

# IC Application Note No. U-9


**Title:** Iodate, chlorite, bromate and nitrite by suppressed ion chromatography applying post column reaction (PCR) and UV/VIS detection

**Summary:** Determination of iodate, chlorite, bromate and nitrite using suppressed anion chromatography with UV/VIS detection after post column reaction.

**Sample:** Standard solution


**Sample Preparation:** –

**Column:** 6.1020.030 Metrosep A Supp 10 – 250  
**Wavelength:** 352 nm  
**Eluent:** 20 mmol/L sodium hydrogencarbonate  
**PCR reagent 1:** 0.2 mmol/L ammonium molybdate in 1 mmol/L sulfuric acid  
**PCR reagent 2:** 40 g/L potassium iodide  
**PCR flow:** 0.5 mL/min  
**Flow:** 1.0 mL/min  
**Injection Volume:** 1000 µL



<b>Results:</b>	Iodate µg/L	Chlorite µg/L	Bromate µg/L	Nitrite µg/L
	1	1	1	1

**Setup of the PCR unit:**



It is important to replace the PTFE mixing coil of the Post Column Reactor (6.2836.000) by a PEEK mixing coil (6.2758.000). Iodine can diffuse through PTFE and this will lead to severe peak tailing.