

IC Application Note No. U-46

Title: Aluminum in an acid extract using UV/VIS detection after post-column reaction with Tiron

Summary: Determination of aluminum in an acidic extract containing metals (e.g. alkali, alkaline earth, iron, chromium, molybdenum, etc.) using cation chromatography with UV detection after post-column reaction with Tiron.

Sample: Acidic extract (4 mol/L nitric acid)

Sample Preparation: 0.1 mL extract diluted in 250 mL eluent, 100 mg ascorbic acid is added to reduce iron(III) to iron(II).

Column: 6.1007.000 Nucleosil 5SA


Wavelength: 312 nm

Eluent: 60 mmol/L potassium sulfate, pH = 2 (sulfuric acid)

Flow: 1.0 mL/min **PCR Flow:** 1.0 mL/min

PCR Reagent: 0.3 mmol/L Tiron in 3 mol/L ammonium acetate

Injection Volume: 250 μ L



Results:	Aluminum mg/L
Sample	1002
Expected value	1000