
IC Application Note U-60

Corrosion inhibitors in cooling water



In industrial cooling water systems, copper and its alloys are widely used because of their superior heat transfer properties. These materials are, however, susceptible to corrosion. Azoles are commonly used to protect copper and its alloys from corrosion. These corrosion inhibitors are quantified by ion chromatography with UV/VIS detection.

Results

	Concentration spiked [mg/L]	Concentration found [mg/L]	RSD (n = 6) [%]
Benzotriazole	1.00	1.07	2.04
Tolyltriazole	1.00	1.04	3.33
Mercaptobenzothiazole	1.00	1.05	3.46



Sample

Artificial cooling water sample spiked with corrosion inhibitors

Sample preparation

Inline Ultrafiltration

Columns

ProntoSil 120-5-C18 AQ - 150/4.0	6.1008.100
Prontosil 120-3-C18 AQ Guard/4.0	

Solutions

Eluent	0.5% phosphoric acid 25% acetonitrile
Simulated cooling water	500 mg/L calcium 250 mg/L magnesium 354 mg/L chloride 240 mg/L sulfate pH = 7.2

Parameters

0.8 mL/min
20 μL
120 MPa
20 min
40 °C
Deuterium lamp
214 nm

Analysis

UV detection

Instrumentation

881 Compact IC pro – Cation	2.881.0010
858 Professional Sample Processor – Pump	2.858.0020
887 Professional UV/VIS Detector	2.887.0010

