

VA Application Note No. V- 3

Title:	Nickel, cobalt and iron in a polyterephthalic acid (PTA) solution
---------------	--

Summary:	Determination of Ni, Co and Fe in a PTA solution containing HCl.
-----------------	--


Sample:	PTA in HCl
Sample Preparation:	none

Nickel and Cobalt:	
Electrolyte:	Dimethylglyoxime in ethanol, NH ₄ Cl / NH ₃ buffer
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-50 mV), HMDE $U_{\text{meas}} = -700 \text{ mV (30s)}$, $U_{\text{start}} = -700 \text{ mV}$, $U_{\text{end}} = -1200 \text{ mV}$ $E_p(\text{Ni}) = -935 \text{ mV}$, $E_p(\text{Co}) = -1050 \text{ mV}$

Iron:	
Electrolyte:	Catechol, Pipes buffer, pH = 7.0
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-50 mV), HMDE $U_{\text{meas}} = -100 \text{ mV (60s)}$, $U_{\text{start}} = -100 \text{ mV}$, $U_{\text{end}} = -700 \text{ mV}$ $E_p(\text{Fe}) = -405 \text{ mV}$

Results:	Ni μg/L	Co μg/L	Fe μg/L
	15.5	61	519.9

Determination of nickel and cobalt



Determination of iron

