

VA Application Note No. V- 16

Title:	Nickel, iron and copper in a silver plating bath
---------------	---

Summary:	Determination of Ni, Fe and Cu in a silver plating bath
-----------------	---


Sample:	Strongly alkaline silver bath containing 2.5 g/L CN^- and 60 g/L Ag
Sample Preparation:	Wet digestion with HNO_3

Nickel:	
Electrolyte:	Dimethylglyoxime, NH_4Cl buffer
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-50 mV), HMDE $U_{\text{meas}} = -700 \text{ mV}$ (30s), $U_{\text{start}} = -700 \text{ mV}$, $U_{\text{end}} = -1200 \text{ mV}$ $E_p(\text{Ni}) = -985 \text{ mV}$

Iron and copper:	
Electrolyte:	Catechol, Pipes buffer pH = 7.0 with NH_3
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-50 mV), HMDE $U_{\text{meas}} = -100 \text{ mV}$ (60s), $U_{\text{start}} = 0 \text{ mV}$, $U_{\text{end}} = -600 \text{ mV}$ $E_p(\text{Cu}) = -175 \text{ mV}$, $E_p(\text{Fe}) = -325 \text{ mV}$

Results:	Ni mg/L	Fe mg/L	Cu mg/L
	25.5	4.1	31.6

Determination of nickel



Determination of iron

Determination of copper and iron

