

# VA Application Note No. V- 30

<b>Title:</b>	<b>Zinc, cadmium, lead, nickel and cobalt in FeCl<sub>3</sub> solution 40%</b>
---------------	--

<b>Summary:</b>	Determination of Zn, Cd, Pb, Ni and Co in FeCl <sub>3</sub> solution 40%
-----------------	--


<b>Sample:</b>	FeCl <sub>3</sub> solution 40%
<b>Sample Preparation:</b>	Dilution 1/100 for nickel and cobalt

<b>Zinc, cadmium and lead:</b>	
<b>Electrolyte:</b>	NH <sub>4</sub> Ac buffer, ascorbic acid
<b>AE:</b>	Pt
<b>RE:</b>	Ag/AgCl/KCl 3M
<b>Parameters:</b>	DPASV (+50 mV), HMDE $U_{\text{meas}} = -1150 \text{ mV (30s)}$ , $U_{\text{start}} = -1150 \text{ mV}$ , $U_{\text{end}} = -200 \text{ mV}$ $E_p(\text{Zn}) = -960 \text{ mV}$ , $E_p(\text{Cd}) = -580 \text{ mV}$ $E_p(\text{Pb}) = -380 \text{ mV}$

<b>Nickel and cobalt:</b>	
<b>Electrolyte:</b>	Dimethylglyoxime in triethanolamine and NH <sub>4</sub> Cl buffer
<b>AE:</b>	Pt
<b>RE:</b>	Ag/AgCl/KCl 3M
<b>Parameters:</b>	DPCSV (-75 mV), HMDE $U_{\text{meas}} = -700 \text{ mV (30s)}$ , $U_{\text{meas}} = -800 \text{ mV (10s)}$ , $U_{\text{start}} = -800 \text{ mV}$ , $U_{\text{end}} = -1000 \text{ mV}$ $E_p(\text{Ni}) = -940 \text{ mV}$ , $E_p(\text{Co}) = -1050 \text{ mV}$

<b>Results:</b>	<b>Zn</b> mg/L	<b>Cd</b> µg/L	<b>Pb</b> mg/L	<b>Ni</b> mg/L	<b>Co</b> mg/L
	<b>212</b>	<b>550</b>	<b>36</b>	<b>141</b>	<b>19</b>

## Determination of zinc (left), cadmium (center) and lead (right)



## Determination of nickel (top) and cobalt (bottom)

