

VA Application Note No. V- 32

Title: Zinc, cadmium, lead, copper, iron, nickel and cobalt in freeze-dried hop

Summary: Determination of Zn, Cd, Pb, Cu, Ni, Co and Fe in freeze-dried hop after a wet digestion.

Sample: Freeze-dried hop samples
Sample Preparation: Wet digestion with H₂SO₄ and H₂O₂


Iron:
Electrolyte: Catechol and Pipes buffer pH = 7.0 ± 0.1 with NH₃
AE: Pt
RE: Ag/AgCl/KCl 3M
Parameters: DPCSV (-50 mV), HMDE
 U_{meas} = -100 mV (20s), U_{start} = -100 mV, U_{end} = -600 mV
 Ep (Fe) = -310 mV

Nickel and cobalt:
Electrolyte: Dimethylglyoxime in ethanol, NH₄Cl buffer
AE: Pt
RE: Ag/AgCl/KCl 3M
Parameters: DPCSV (-50 mV), HMDE
 U_{meas} = -700 mV (30s), U_{start} = -700 mV, U_{end} = -1200 mV
 Ep (Ni) = -970 mV, Ep (Co) = -1110 mV


Zinc, cadmium, lead and copper:
Electrolyte: NH₄Ac buffer, pH = 4.6
AE: Pt
RE: Ag/AgCl/KCl 3M
Parameters: DPASV (+50 mV), HMDE
 U_{meas} = -1150 mV (90s), U_{start} = -1150 mV, U_{end} = +200 mV
 Ep (Zn) = -960 mV, Ep (Cd) = -550 mV, Ep (Pb) = -395 mV, Ep (Cu) = +55 mV

Results:	Zn mg/L	Cd mg/L	Pb mg/L	Cu mg/L	Ni mg/L	Co mg/L	Fe mg/L
	27.2	1.85	3.08	159.3	9.78	0.151	381.5

Determination of iron



Determination of nickel and cobalt



Determination of zinc, cadmium, lead and copper

