

VA Application Note No. V- 41

Title:	Cadmium, lead, copper, nickel and cobalt in soybean oil after digestion
---------------	--

Summary:	Determination of Cd, Pb, Cu, Ni, Co in soybean oil after extraction by boiling with HCl under reflux
-----------------	--


Sample:	Soybean Oil
Sample Preparation:	Extraction with HCl by boiling in a reflux apparatus. Washing of the oil phase with hot water. Separation of the oil phase in a separation funnel. The aqueous phase is analysed by voltammetry.

Cadmium, lead and copper:	
Electrolyte:	Merck reagent cocktail «Suprapure B-931814» KCl / NaAc NH ₄ Ac buffer Suprapure pH = 4.64
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPASV (+50 mV), HMDE U _{meas} = -800 mV (120s), U _{start} = -800 mV, U _{end} = 0 mV Ep (Cd) = -650 mV, Ep (Pb) = -450 mV Ep (Cu) = -135 mV


Nickel and cobalt:	
Electrolyte:	Dimethylglyoxime in ethanol, NH ₄ Cl buffer pH = 8.9 of the solution to analyse
AE:	Pt
RE:	Ag/AgCl/KCl 3M
Parameters:	DPCSV (-75 mV), HMDE U _{meas} = -800 mV (30s), U _{start} = -800 mV, U _{end} = -1250 mV Ep (Ni) = -980 mV, Ep (Co) = -1080 mV

Results:	Cd μg/L	Pb mg/L	Cu mg/L	Ni mg/L	Co μg/L
	/	7.1	1.2	1.1	/


Determination of cadmium




Determination of lead



Determination of copper



Determination of nickel



Determination of cobalt

