


# Uranium in drinking water



Uranium can be determined in drinking water by adsorptive stripping voltammetry (AdSV) at the hanging mercury drop electrode (HMDE). Chloranilic acid is used as complexing agent.

## Results

U in drinking water

0.40  $\mu\text{g/L}$

# Method description

## Sample

Drinking water

## Instruments

797 VA Computrace



## Parameters

Working electrode	HMDE
Drop size	4
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	0.1 V
Deposition time	60 s
Equilibration time	5 s
Start potential	0.05 V
End potential	-0.3 V
Pulse amplitude	0.05 V
Pulse time	0.04 s
Voltage step	0.004 V
Voltage step time	0.1 s
Sweep rate	0.04 V/s
Peak potential U	-0.1 V

## Electrodes

Multi-Mode Electrode pro	6.1246.120
Silanized capillaries	6.1226.050
Ag/AgCl/KCl (3 mol/L) reference electrode. Bridge electrolyte c(KCl) = 3 mol/L	6.0728.020 6.1245.010
Separate Pt rod electrode	6.0343.000

## Reagents

NaNO <sub>3</sub>	Sodium nitrate, 99%
HNO <sub>3</sub>	Nitric acid, for trace analysis*, w(HNO <sub>3</sub> ) = 65%
Chloranilic acid	Chloranilic acid, 99%

\*e.g., Merck suprapur®, Sigma-Aldrich TraceSelect® or equivalent.

## Solutions

Sodium nitrate solution	c(NaNO <sub>3</sub> ) = 1 mol/L
Diluted nitric acid	c(HNO <sub>3</sub> ) = 0.5 mol/L
Chloranilic acid solution	c(Chloranilic acid) = 5 mmol/L

## Analysis

Measuring solution	10 mL sample + 1 mL sodium nitrate solution + 150 µL diluted nitric acid + 110 µL chloranilic acid solution
--------------------	--

[www.metrohm.com](http://www.metrohm.com)

 **Metrohm**