


Zinc, cadmium, lead, and copper in waste water after UV digestion



Zinc, cadmium, lead, and copper can be determined in waste water samples after UV digestion by anodic stripping voltammetry according to DIN 38406 part 16.

Results

Zn	198 µg/L
Cd	1.9 µg/L
Pb	7.5 µg/L
Cu	41.2 µg/L

Method description

Sample

Waste water

*e.g., Merck suprapur®, Sigma-Aldrich TraceSelect® or equivalent.

Instruments

797 VA Computrace & 909 UV Digester



Sample preparation

10 mL waste water sample, 10 µL HCl, and 100 µL H₂O₂ are pipetted into the 12 mL quartz sample vessels. The sample holder with the 12 quartz sample vessels is placed in the 909 UV Digester. The samples are irradiated at 90 °C for 60 min.

Parameters 909 UV Digester

Temperature	90 °C
Irradiation time	60 min

Electrodes

Multi-Mode Electrode pro	6.1246.120
Silanized capillaries	6.1226.050
Ag/AgCl/KCl (3 mol/L) reference electrode. Bridge electrolyte c(KCl) = 3 mol/L	6.0728.020 6.1245.010
Separate Pt rod electrode	6.0343.000

Reagents

HCl	Hydrochloric acid, for trace analysis*, w(HCl) = 30%
H ₂ O ₂	Hydrogen peroxide solution, for trace analysis*, w(H ₂ O ₂) = 30%
CH ₃ COOH	Acetic acid, for trace analysis*
NH ₃	Ammonia solution, for trace analysis*, w(NH ₃) = 25%

Solutions

Supporting electrolyte	Ammonium acetate buffer pH 4.6 c(CH ₃ COOH) = 2 mol/L c(NH ₃) = 1 mol/L
------------------------	---

Analysis

Measuring solution	10.11 mL digested sample solution + 1 mL supporting electrolyte
--------------------	--

Parameters 797 VA Computrace

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	-1.15 V
Deposition time	90 s
Equilibration time	10 s
Start potential	-1.15 V
End potential	0.1 V
Pulse amplitude	0.05 V
Pulse time	0.04 s
Voltage step	0.006 V
Voltage step time	0.1 s
Sweep rate	0.06 V/s
Peak potential Zn	-1.0 V
Peak potential Cd	-0.58 V
Peak potential Pb	-0.42 V
Peak potential Cu	0 V

www.metrohm.com

 **Metrohm**