


Nickel in white wine after UV digestion



For the determination of nickel in white wine, UV digestion is required to mineralize the sample. The determination is done by adsorptive stripping voltammetry (AdSV) at the HMDE in ammonia buffer with dimethylglyoxime (DMG).

Results

Ni in white wine

660 $\mu\text{g/L}$

Method description

Sample

White wine

Instruments

797 VA Computrace & 909 UV Digester



Sample preparation

3 mL ultrapure water, 0.3 mL white wine, 10 µL HCl, and 2 mL H₂O₂ are pipetted into the 12 mL quartz sample vessels. The sample holder with the 12 quartz sample vessels is placed in the 909 UV Digester. The samples are irradiated at 90 °C for 120 min. After 60 min another 1 mL H₂O₂ is added. For the voltammetric determination of nickel 4 mL ultrapure water are added to the digested solution.

Parameters 909 UV Digester

Temperature	90 °C
Irradiation time	120 min

Electrodes

Multi-Mode Electrode pro Non-silanized capillaries	6.1246.120 6.1226.030
Ag/AgCl/KCl (3 mol/L) reference electrode. Bridge electrolyte c(KCl) = 3 mol/L	6.0728.020 6.1245.010
Separate Pt rod electrode	6.0343.000

Reagents

HCl	Hydrochloric acid, for trace analysis*, w(HCl) = 30%
H ₂ O ₂	Hydrogen peroxide solution, for trace analysis*, w(H ₂ O ₂) = 30%
NH ₃	Ammonia solution, for trace

analysis*, w(NH₃) = 25%

DMG Dimethylglyoxime disodium salt
octahydrate, 97%

*e.g., Merck suprapur®, Sigma-Aldrich TraceSelect® or
equivalent.

Solutions

Ammonia buffer c(NH₃) = 3 mol/L
+ c(HCl) = 1 mol/L

DMG solution c(DMG) = 0.1 mol/L in water

Analysis

Measuring
solution 1 mL digested sample solution
+ 10 mL ultrapure water
+ 0.5 mL ammonia buffer
+ 0.1 mL DMG solution
→ pH adjusted to 9.3 with HCl

Parameters 797 VA Computrace

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	-0.7 V
Deposition time	30 s
Equilibration time	5 s
Pulse amplitude	0.05 V
Start potential	-0.7 V
End potential	-1.25 V
Voltage step	0.004 V
Voltage step time	0.3 s
Sweep rate	0.0133 V/s
Peak potential Ni	-0.94 V

www.metrohm.com

 **Metrohm**