

# VA Application Note No. V-99

**Title:** Cadmium, lead and copper in triphosphate

**Summary:** Cadmium, lead and copper are determined by anodic stripping voltammetry (ASV) at the HMDE using aqueous nitric acid as supporting electrolyte.

**Sample:** Pentasodium triphosphate

**Sample preparation:** none

## Determination of cadmium, lead and copper

**Electrolyte** Nitric acid:  
w(HNO<sub>3</sub>) = 65%, suprapur

**Measuring solution** 10 mL water  
+ 250 mg sample  
+ 250 µL nitric acid

**Working electrode (WE)** **MME** (Multi-Mode Electrode) 6.1246.020

**Auxiliary electrode (AE)** **Pt** 6.0343.000

**Reference electrode (RE)** Ag/AgCl/KCl (3 mol/L): 6.0728.020 + 6.1245.010

<b>Parameters</b>	
Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	-800 mV
Deposition time	60 s
Equilibration time	10 s
Pulse amplitude	50 mV
Start potential	-750 mV
End potential	+150 mV
Voltage step	6 mV
Voltage step time	0.1 s
Sweep rate	60 mV/s
Peak potential Cd	-600 mV
Peak potential Pb	-400 mV
Peak potential Cu	+10 mV

<b>Results:</b>	Cd	Pb	Cu
	2.3 ng/g	102 ng/g	1.8 µg/g

## Determination of Cd, Pb and Cu

Segment: asv VR(\*\*)

