

# VA Application Note No. V - 102

**Title:** Manganese in triphosphate

**Summary:** Mn is determined by anodic stripping voltammetry (ASV) at the HMDE in alkaline solution.

**Sample:** pentasodium triphosphate

**Sample preparation:** 2.5 g sample  
 + approx. 30 mL ultrapure water  
 + 1.25 mL c(H<sub>2</sub>SO<sub>4</sub>) = 96%, suprapur  
 The solution is boiled for 30 min. Water is added to keep a constant volume. After cooling to room temperature the solution is made up to 50 mL.

## Analysis of Mn

**Electrolyte** tartrate solution: c(disodium tartrate) = 0.2 mol/L  
 Zn solution:  $\beta(\text{Zn}^{2+}) = 10 \text{ mg/L}$

**Measuring solution** 10 mL ultrapure water  
 + 0.1 mL digestion solution  
 + 1 mL tartrate solution  
 + 100  $\mu\text{L}$  Zn solution  
 adjust the pH to 10 ... 10.5

**Working electrode (WE)** MME (Multi Mode Electrode) 6.1246.020

**Auxiliary electrode (AE)** Pt: 6.0343.000

**Reference electrode (RE)** Ag/AgCl/KCl (3 mol/L): 6.0728.020 + 6.1245.010

### Parameters

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	400 s
Deposition potential	-1700 mV
Deposition time	90 s
Equilibration time	10 s
Pulse amplitude	<b>-75 mV</b>
Start potential	-1650 mV
End potential	-1300 mV

Voltage step	6 mV
Voltage step time	0.3 s
Sweep rate	20 mV/s
Peak potential Mn	-1530 mV

<b>Results:</b>	Mn
	10 µg/g

### Determination of Mn

