

# VA Application Note No. V - 104

|               |                                                |
|---------------|------------------------------------------------|
| <b>Title:</b> | <b>Formaldehyde in metalworking lubricants</b> |
|---------------|------------------------------------------------|

|                 |                                                                               |
|-----------------|-------------------------------------------------------------------------------|
| <b>Summary:</b> | Formaldehyde is determined polarographically at the DME in alkaline solution. |
|-----------------|-------------------------------------------------------------------------------|

|                            |                        |
|----------------------------|------------------------|
| <b>Sample:</b>             | metalworking lubricant |
| <b>Sample preparation:</b> | none                   |

## Analysis of formaldehyde

|                                 |                                                                                                                                                      |
|---------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Electrolyte</b>              | $c(\text{LiOH}) = 0.1 \text{ mol/L} + c(\text{EDTA}) = 0.02 \text{ mol/L}$                                                                           |
| <b>Measuring solution</b>       | 10 mL electrolyte are degassed manually for 240 s. 100 $\mu\text{L}$ sample are added. After additional 60 s degassing the determination is started. |
| <b>Auxiliary electrode (AE)</b> | Pt                                                                                                                                                   |
| <b>Reference electrode (RE)</b> | Ag/AgCl/KCl (3 mol/L)                                                                                                                                |

|                   |                             |          |
|-------------------|-----------------------------|----------|
| <b>Parameters</b> | Working electrode           | DME      |
|                   | Stirrer speed               | 2000 rpm |
|                   | Mode                        | DP       |
|                   | Purge time                  | 60 s     |
|                   | Equilibration time          | 10 s     |
|                   | Pulse amplitude             | 50 mV    |
|                   | Start potential             | -1400 mV |
|                   | End potential               | -1850 mV |
|                   | Voltage step                | 8 mV     |
|                   | Voltage step time           | 0.5 s    |
|                   | Sweep rate                  | 16 mV/s  |
|                   | Peak potential formaldehyde | -1640 mV |

|                 |              |
|-----------------|--------------|
| <b>Results:</b> | formaldehyde |
|                 | 172 mg/L     |

## Determination of formaldehyde

