# VA Application Note V–110

# Total chromium in waste water after UV digestion (polarographic method with ethylene diamine)



Cr(VI) is determined by polarography at the SMDE in acetate solution containing ethylene diamine to mask interfering copper ions.

Only Cr(VI) is electrochemically active. Therefore, all chromium has to be oxidized to Cr(VI) prior to analysis. UV irradiation at pH > 4 is used to oxidize chromium.

# Results

Cr(total) in waste water

8.2 µg/L



# Method description

### Sample

Waste water

#### Instruments

797 VA Computrace & 909 UV Digester



#### Sample preparation

For UV digestion 10 mL waste water, 20  $\mu$ L HCl, and 50  $\mu$ L H<sub>2</sub>O<sub>2</sub> are pipetted into the 12 mL quartz sample vessels. The sample holder with the 12 quartz sample vessels is placed in the 909 UV Digester. The samples are irradiated at 90 °C for 60 min. To oxidize Cr(III) to Cr(VI), the pH of the solutions is adjusted to 4–6 with NaOH and 50  $\mu$ L H<sub>2</sub>O<sub>2</sub> are added. The samples are irradiated at 90 °C for 20 min.

#### Parameters 909 UV Digester

Temperature	90 °C
Irradiation time (UV digestion)	60 min
Irradiation time (Oxidation)	20 min

#### Electrodes

Multi-Mode Electrode pro Silanized capillaries	6.1246.120 6.1226.050
Ag/AgCl/KCl (3 mol/L) reference electrode. Bridge electrolyte c(KCl) = 3 mol/L	6.0728.020 6.1245.010
Separate Pt rod electrode	6.0343.000

## Reagents

HCI	Hydrochloric acid, for trace analysis*, w(HCl) = 30%
$H_2O_2$	Hydrogen peroxide solution, for trace analysis*, $w(H_2O_2) = 30\%$

NaOH	Sodium hydroxide solution, for trace analysis*, w(NaOH) = 30%
$H_2N(CH_2)_2NH_2$	Ethylene diamine, 99.5%
NH <sub>3</sub>	Ammonia solution, for trace analysis <sup>*</sup> , $w(NH_3) = 25\%$
CH₃COOH	Acetic acid, for trace analysis*

\*e.g., Merck suprapur^®, Sigma-Aldrich  $\mathsf{TraceSelect}^{\$}$  or equivalent.

## Analysis

Measuring solution	5 mL digested, oxidized waste water + 15 mL ultrapure water + 10 $\mu$ L ethylene diamine + 200 $\mu$ L NH <sub>3</sub> + 150 $\mu$ L acetic acid $\rightarrow$ pH adjusted to 6.8 ±
	0.1 with NaOH

## Parameters 797 VA Computrace

Working electrode	SMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	600 s
Equilibration time	10 s
Start potential	0.1 V
End potential	-0.15 V
Pulse amplitude	0.05 V
Voltage step	0.004 V
Voltage step time	0.4 s
Sweep rate	0.01 V/s
Peak potential Cr(VI)	-0.02 V

www.metrohm.com

