

VA Application Note No. V - 131

Title:	Determination of nickel and cobalt in sulfuric acid
---------------	--


Summary:	The concentration of Ni and Co is determined by adsorptive stripping voltammetry at the HMDE with dimethylglyoxime (DMG) as complexing agent.
-----------------	---

Sample:	Sulfuric acid
Sample preparation:	Dilution in water 1:50

Analysis of Ni ²⁺ and Co ²⁺		
Ammonia buffer pH 9.8	c(NH ₄ Cl) = 1 mol/L c(NH ₃) = 1 mol/L	
DMG solution	c(DMG) = 0.1 mol/L in water DMG: Dimethylglyoxime disodium salt	
Measuring solution	5 mL water + 5 mL diluted sample + 500 µL ammonia buffer + 500 µL DMG solution → adjust pH to 9.5 with NH ₃	
Working electrode (WE)	MME (Multi Mode Electrode) 6.1246.020	
Auxiliary electrode (AE)	Pt 6.0343.000	
Reference electrode (RE)	Ag/AgCl/KCl (3 mol/L): 6.0728.020 + 6.1245.010	
Parameters	Working electrode	HMDE
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Deposition potential	-0.7
	Deposition time	30 s
	Equilibration time	10 s
	Pulse amplitude	50 mV
	Start potential	-0.8 V
	End potential	-1.25 V
	Voltage step	4 mV
	Voltage step time	0.3 s
	Sweep rate	13.3 mV/s
	Peak potential Ni	-950 mV
	Peak potential Co	-1250 mV

Results:	Ni	Co
	913 µg/L	47.3 µg/L

Determination of Ni and Co



Ni
 $c = 920.515 \text{ µg/L}$
 $\pm 3.321 \text{ µg/L (0.36\%)}$

Co
 $c = 47.666 \text{ µg/L}$
 $\pm 0.461 \text{ µg/L (0.97\%)}$

