

# VA Application Note No. V - 150

<b>Title:</b>	<b>Copper in a nickel plating bath</b>
---------------	--


<b>Summary:</b>	The concentration of Cu in a Ni plating bath is determined by polarography in chloride containing acetate buffer at pH 4.7.
-----------------	---

<b>Sample:</b>	Ni electroplating bath
<b>Sample preparation:</b>	None

Analysis of Cu																									
<b>Acetate buffer pH 4.7</b>	c(Na acetate) = 1.5 mol/L c(CH <sub>3</sub> COOH) = 1.5 mol/L																								
<b>KCl solution</b>	c(KCl) = 3 mol/L																								
<b>Measuring solution</b>	10 mL ultrapure H <sub>2</sub> O + 1 mL acetate buffer pH 4.7 + 1 mL KCl solution + 0.5 mL Ni plating bath																								
<b>Working electrode (WE)</b>	<b>MME</b> (Multi Mode Electrode) 6.1246.020																								
<b>Auxiliary electrode (AE)</b>	<b>Pt</b> 6.0343.000																								
<b>Reference electrode (RE)</b>	Reference system: Ag/AgCl/KCl (3 mol/L) 6.0728.020 Intermediate electrolyte: c(KCl) = 3 mol/L 6.1245.010																								
<b>Parameters</b>	<table border="1"> <tr> <td>Working electrode</td> <td>DME</td> </tr> <tr> <td>Stirrer speed</td> <td>2000 rpm</td> </tr> <tr> <td>Mode</td> <td>DP</td> </tr> <tr> <td>Purge time</td> <td>300 s</td> </tr> <tr> <td>Equilibration time</td> <td>10 s</td> </tr> <tr> <td>Pulse amplitude</td> <td>0.05 V</td> </tr> <tr> <td>Start potential</td> <td>0.05 V</td> </tr> <tr> <td>End potential</td> <td>-0.2 V</td> </tr> <tr> <td>Voltage step</td> <td>0.006 V</td> </tr> <tr> <td>Voltage step time</td> <td>0.6 s</td> </tr> <tr> <td>Sweep rate</td> <td>0.01 V/s</td> </tr> <tr> <td>Peak potential Cu</td> <td>-0.1 V</td> </tr> </table>	Working electrode	DME	Stirrer speed	2000 rpm	Mode	DP	Purge time	300 s	Equilibration time	10 s	Pulse amplitude	0.05 V	Start potential	0.05 V	End potential	-0.2 V	Voltage step	0.006 V	Voltage step time	0.6 s	Sweep rate	0.01 V/s	Peak potential Cu	-0.1 V
Working electrode	DME																								
Stirrer speed	2000 rpm																								
Mode	DP																								
Purge time	300 s																								
Equilibration time	10 s																								
Pulse amplitude	0.05 V																								
Start potential	0.05 V																								
End potential	-0.2 V																								
Voltage step	0.006 V																								
Voltage step time	0.6 s																								
Sweep rate	0.01 V/s																								
Peak potential Cu	-0.1 V																								

<b>Results:</b>	Cu
	22.1 mg/L

**Determination of Cu**



Cu  
c = 22.139 mg/L  
+/- 0.079 mg/L (0.35%)

