

VA Application Note No. V - 151

Title:	Antimony(III) and total antimony in electroless nickel bath
---------------	--

Summary:	The concentration of Sb(III) and Sb(total) in an electroless nickel bath is determined by anodic stripping voltammetry (ASV). In $c(\text{HCl}) = 0.6 \text{ mol/L}$ only Sb(III) shows a signal. In $w(\text{HCl}) = 10\%$ the Sb(total) content is determined.
-----------------	--

Sample:	Electroless Ni bath
----------------	---------------------

Sample preparation:	None
----------------------------	------

Analysis of Sb(III)

Electrolyte for Sb(III)	$c(\text{HCl}) = 0.6 \text{ mol/L}$	
Measuring solution	15 mL electrolyte for Sb(III) + 50 μL Ni plating bath	
Working electrode (WE)	MME (Multi Mode Electrode)	6.1246.020
Auxiliary electrode (AE)	Pt	6.0343.000
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Intermediate electrolyte: $c(\text{KCl}) = 3 \text{ mol/L}$	6.1245.010

Parameters

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	-0.45 V
Deposition time	30 s
Equilibration time	10 s
Pulse amplitude	0.02 V
Start potential	-0.3 V
End potential	0 V
Voltage step	0.004 V
Voltage step time	0.3 s
Sweep rate	0.013 V/s
Peak potential Sb(III)	-0.12 V

Analysis of Sb(total)



Electrolyte for Sb(total)	w(HCl) = 10%	
Measuring solution	15 mL electrolyte for Sb(total) + 20 µL Ni plating bath	
Working electrode (WE)	MME (Multi Mode Electrode)	6.1246.020
Auxiliary electrode (AE)	Pt	6.0343.000
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Outer system: c(KCl) = 3 mol/L	6.1245.010

Parameters

Working electrode	HMDE
Stirrer speed	2000 rpm
Mode	DP
Purge time	300 s
Deposition potential	-0.45 V
Deposition time	30 s
Equilibration time	10 s
Pulse amplitude	0.02 V
Start potential	-0.3 V
End potential	-0.1 V
Voltage step	0.004 V
Voltage step time	0.3 s
Sweep rate	0.013 V/s
Peak potential Sb(total)	-0.19 V

Results:	Sb(III)	Sb(total)
	5.6 mg/L	47.3 mg/L

Determination of Sb(III)



Determination of Sb(total)

