

VA Application Note No. V - 152

Title: Thallium in a cyanidic gold bath


Summary: The concentration of Tl in a cyanidic Au bath is determined by anodic stripping voltammetry (ASV) without further addition of electrolyte.

Sample: Electroless cyanidic Au bath
Sample preparation: None

Analysis of Tl		
Electrolyte	None	
Measuring solution	10 mL ultrapure water + 0.1 mL Au bath	
Working electrode (WE)	MME (Multi Mode Electrode)	6.1246.020
Auxiliary electrode (AE)	Pt	6.0343.000
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Intermediate electrolyte: c(KCl) = 3 mol/L	6.1245.010
Parameters	Working electrode	HMDE
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Deposition potential	-0.8 V
	Deposition time	30 s
	Equilibration time	5 s
	Pulse amplitude	0.05 V
	Start potential	-0.8 V
	End potential	-0.2 V
	Voltage step	0.006 V
	Voltage step time	0.4 s
	Sweep rate	0.015 V/s
	Peak potential Tl	-0.45 V

Results:	Tl
	5.3 mg/L

Determination of Tl



Thallium
c = 5.246 mg/L
+/- 0.057 mg/L (1.08%)

