

VA Application Note No. V - 167

Title:	Cadmium in a phosphatation bath
---------------	--


Summary:	The concentration of Cd in a Zn phosphatation bath is determined by polarography in HCl electrolyte.
-----------------	--

Sample:	Zn phosphatation bath
Sample preparation:	None

Analysis of Cd																									
HCl solution	w(HCl) = 30%																								
Measuring solution	10 mL H ₂ O + 1 mL HCl solution + 1 mL phosphatation bath																								
Working electrode (WE)	MME (Multi Mode Electrode) 6.1246.020																								
Auxiliary electrode (AE)	Pt 6.0343.000																								
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L) 6.0728.020 Intermediate electrolyte: c(KCl) = 3 mol/L 6.1245.010																								
Parameters	<table border="1"> <tr> <td>Working electrode</td> <td>DME</td> </tr> <tr> <td>Stirrer speed</td> <td>2000 rpm</td> </tr> <tr> <td>Mode</td> <td>DP</td> </tr> <tr> <td>Purge time</td> <td>300 s</td> </tr> <tr> <td>Equilibration time</td> <td>5 s</td> </tr> <tr> <td>Pulse amplitude</td> <td>0.05 V</td> </tr> <tr> <td>Start potential</td> <td>-0.48 V</td> </tr> <tr> <td>End potential</td> <td>-0.7 V</td> </tr> <tr> <td>Voltage step</td> <td>0.006 V</td> </tr> <tr> <td>Voltage step time</td> <td>0.6 s</td> </tr> <tr> <td>Sweep rate</td> <td>0.01 V/s</td> </tr> <tr> <td>Peak potential Cd</td> <td>-0.58 V</td> </tr> </table>	Working electrode	DME	Stirrer speed	2000 rpm	Mode	DP	Purge time	300 s	Equilibration time	5 s	Pulse amplitude	0.05 V	Start potential	-0.48 V	End potential	-0.7 V	Voltage step	0.006 V	Voltage step time	0.6 s	Sweep rate	0.01 V/s	Peak potential Cd	-0.58 V
Working electrode	DME																								
Stirrer speed	2000 rpm																								
Mode	DP																								
Purge time	300 s																								
Equilibration time	5 s																								
Pulse amplitude	0.05 V																								
Start potential	-0.48 V																								
End potential	-0.7 V																								
Voltage step	0.006 V																								
Voltage step time	0.6 s																								
Sweep rate	0.01 V/s																								
Peak potential Cd	-0.58 V																								

Results:	Cd
	2.6 mg/L

Determination of Cd



Cd
 c = 2.577 mg/L
 +/- 0.005 mg/L (0.20%)

