

VA Application Note No. V - 169

Title: Lead in tin soldering contacts

Summary: The concentration of Pb in Sn soldering contacts is determined by anodic stripping voltammetry (ASV) in an electrolyte containing citrate, oxalic acid, HCl and cetyl trimethyl ammonium bromide.

Sample: Sn soldering contact

Sample preparation: 2 mL w(HNO₃) = 65% and 0.1 mL w(HCl) = 30% are added onto a sample of soldering contact (approx. 50 mg). The mixture is heated on a heater to 150°C. When the sample is dissolved completely the solution is cooled down to room temperature and made up to 10 mL with ultrapure water.

Analysis of Pb																			
Supporting electrolyte pH 2.5	c(sodium citrate) = 0.1 mol/L c(oxalic acid) = 0.1 mol/L c(HCl) = 0.2 mol/L																		
CTAB solution	c(CTAB) = 0.005 mol/L CTAB: cetyl trimethylammonium bromide																		
Measuring solution	5 mL H ₂ O + 5 mL supporting electrolyte pH 2.5 + 20 µL prepared sample solution + 0.05 mL CTAB solution																		
Working electrode (WE)	MME (Multi Mode Electrode) 6.1246.020																		
Auxiliary electrode (AE)	Pt 6.0343.000																		
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L) 6.0728.020 Intermediate electrolyte: c(KCl) = 3 mol/L 6.1245.010																		
Parameters	<table border="1"> <tbody> <tr> <td>Working electrode</td> <td>HMDE</td> </tr> <tr> <td>Stirrer speed</td> <td>2000 rpm</td> </tr> <tr> <td>Mode</td> <td>DP</td> </tr> <tr> <td>Purge time</td> <td>300 s</td> </tr> <tr> <td>Deposition potential</td> <td>-0.48 V</td> </tr> <tr> <td>Deposition time</td> <td>90 s</td> </tr> <tr> <td>Equilibration time</td> <td>20 s</td> </tr> <tr> <td>Pulse amplitude</td> <td>0.05 V</td> </tr> <tr> <td>Start potential</td> <td>-0.53 V</td> </tr> </tbody> </table>	Working electrode	HMDE	Stirrer speed	2000 rpm	Mode	DP	Purge time	300 s	Deposition potential	-0.48 V	Deposition time	90 s	Equilibration time	20 s	Pulse amplitude	0.05 V	Start potential	-0.53 V
Working electrode	HMDE																		
Stirrer speed	2000 rpm																		
Mode	DP																		
Purge time	300 s																		
Deposition potential	-0.48 V																		
Deposition time	90 s																		
Equilibration time	20 s																		
Pulse amplitude	0.05 V																		
Start potential	-0.53 V																		

End potential	-0.25 V
Voltage step	0.004 V
Voltage step time	0.2 s
Sweep rate	0.02 V/s
Peak potential Pb	-0.4 V

Results:	Pb
	1.4 mg/g

Determination of Pb

