

VA Application Note No. V - 171

Title:	Tellurium(IV) in zinc plant electrolyte
---------------	--

Summary:	The concentration of Te(IV) in Zn plant electrolyte is determined by cathodic stripping voltammetry (CSV) in ammonium sulfate electrolyte containing EDTA and Cu. To get a proper complexation of the interfering Zn a high amount of EDTA at is necessary pH 3.4.
-----------------	--

Sample:	Zinc plant electrolyte (concentrated ZnSO ₄ solution)
Sample preparation:	None

Analysis of Te		
Supporting electrolyte	c((NH ₄) ₂ SO ₄) = 2.5 mol/L c(Na ₂ EDTA) = 0.1 mol/L <i>Na₂EDTA: Ethylenediaminetetraacetic acid disodium salt</i>	
Cu solution	β(Cu) = 1 g/L	
Measuring solution	10 mL supporting electrolyte + 0.2 mL Cu solution + 50 μL zinc plant electrolyte pH 3.4 adjusted with H ₂ SO ₄	
Working electrode (WE)	MME (Multi Mode Electrode)	6.1246.020
Auxiliary electrode (AE)	Pt	6.0343.000
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Intermediate electrolyte: c(KCl) = 3 mol/L	6.1245.010
Parameters	Working electrode	HMDE
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Deposition potential	-0.8 V
	Deposition time	120 s
	Equilibration time	10 s
	Pulse amplitude	0.08 V
	Start potential	-0.9 V
	End potential	-1.15 V
	Voltage step	0.004 V
	Voltage step time	0.1 s
	Sweep rate	0.04 V/s
	Peak potential Te	-1.0 V

Results:	Te
	27.3 µg/L

Determination of Te

