

VA Application Note No. V - 177

Title:	Total iron in a chromium bath (triethanolamine-bromate-method)
---------------	---


Summary:	The concentration of Fe(total) is determined by polarography in alkaline electrolyte containing triethanolamine (TEA) and KBrO_3 . All reagents typically contain Fe impurities. Therefore a subtraction of the reagent blank is recommended.
-----------------	--

Sample:	Cr(III/VI) electroplating bath
Sample preparation:	Reduction of Cr(VI) 250 μL Cr bath is filled up to 25 mL with $c(\text{oxalic acid}) = 0.5 \text{ mol/L}$. The solution is gently heated while stirring until the color changes.

Analysis of Fe(total)																									
Supporting electrolyte	$c(\text{NaOH}) = 0.3 \text{ mol/L}$ $c(\text{KBrO}_3) = 0.1 \text{ mol/L}$ $c(\text{TEA}) = 0.05 \text{ mol/L}$ TEA: triethanolamine																								
Measuring solution	15 mL H_2O + 5 mL supporting electrolyte + 50 μL reduced Cr bath																								
Working electrode (WE)	MME (Multi Mode Electrode) 6.1246.020																								
Auxiliary electrode (AE)	Pt 6.0343.000																								
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L) 6.0728.020 Intermediate electrolyte: $c(\text{KCl}) = 3 \text{ mol/L}$ 6.1245.010																								
Parameters	<table border="1"> <tbody> <tr> <td>Working electrode</td> <td>DME</td> </tr> <tr> <td>Stirrer speed</td> <td>2000 rpm</td> </tr> <tr> <td>Mode</td> <td>DP</td> </tr> <tr> <td>Purge time</td> <td>300 s</td> </tr> <tr> <td>Equilibration time</td> <td>5 s</td> </tr> <tr> <td>Pulse amplitude</td> <td>0.05 V</td> </tr> <tr> <td>Start potential</td> <td>-0.7 V</td> </tr> <tr> <td>End potential</td> <td>-1.3 V</td> </tr> <tr> <td>Voltage step</td> <td>0.005 V</td> </tr> <tr> <td>Voltage step time</td> <td>0.8 s</td> </tr> <tr> <td>Sweep rate</td> <td>0.0062 V/s</td> </tr> <tr> <td>Peak potential Fe</td> <td>-0.95 V</td> </tr> </tbody> </table>	Working electrode	DME	Stirrer speed	2000 rpm	Mode	DP	Purge time	300 s	Equilibration time	5 s	Pulse amplitude	0.05 V	Start potential	-0.7 V	End potential	-1.3 V	Voltage step	0.005 V	Voltage step time	0.8 s	Sweep rate	0.0062 V/s	Peak potential Fe	-0.95 V
Working electrode	DME																								
Stirrer speed	2000 rpm																								
Mode	DP																								
Purge time	300 s																								
Equilibration time	5 s																								
Pulse amplitude	0.05 V																								
Start potential	-0.7 V																								
End potential	-1.3 V																								
Voltage step	0.005 V																								
Voltage step time	0.8 s																								
Sweep rate	0.0062 V/s																								
Peak potential Fe	-0.95 V																								

Results:	Fe (blank subtracted)
	1.5 g/L

Determination of Fe(total)



Fe
c = 2.604 g/L
+/- 0.022 g/L (0.83%)

