

# VA Application Note No. V - 186

**Title:** Chromium(VI) in electronic components as part of electrotechnical products

**Summary:** The EU directive on «Restriction of Hazardous Substances» (RoHS) requires the testing of four regulated heavy metals (Pb, Hg, Cd, Cr(VI)) in electrotechnical products. After sample preparation according to IEC 62321 the determination of chromium(VI) in electronic components can be carried out by polarography in ammonia buffer pH 9.6.

**Sample:** Electronic components

**Sample preparation:** Approx. 5 g of the ground sample is extracted with an alkaline sodium carbonate solution as described in IEC 62321.


## Analysis of Cr(VI)

<b>Electrolyte</b>	Ammonia buffer pH 9.6 c(NH <sub>3</sub> ) = 2 mol/L c(NH <sub>4</sub> Cl) = 1 mol/L	
<b>Measuring solution</b>	5 mL      ultrapure water + 5 mL      ammonia buffer pH 9.6 + 0.1 mL   extraction solution (equals approx. 5 mg of ground sample)	
<b>Working electrode (WE)</b>	<b>MME</b> (Multi Mode Electrode)	6.1246.020
	With silanized capillary	6.1226.050
<b>Auxiliary electrode (AE)</b>	<b>Pt</b>	6.0343.000
<b>Reference electrode (RE)</b>	Reference system: Ag/AgCl/KCl (3 mol/L)	6.0728.020
	Intermediate electrolyte: c(KCl) = 3 mol/L	6.1245.010
<b>Parameters</b>	Working electrode	DME
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Equilibration time	5 s
	Pulse amplitude	0.05 V
	Start potential	-0.1 V
	End potential	-0.45 V
	Voltage step	0.006 V

	Voltage step time	0.6 s
	Sweep rate	0.01 V/s
	Peak potential Cr(VI)	-0.25 V

<b>Results:</b>	Cr(VI)
	934.8 mg/kg

### Determination of Cr(VI)



Cr  
 $c = 934.791 \text{ mg/kg}$   
 $\pm 14.742 \text{ mg/kg (1.58\%)}$

