

VA Application Note No. V - 190

Title:	Mercury in polymer materials as part of electrotechnical products
---------------	--

Summary:	The EU directive on «Restriction of Hazardous Substances» (RoHS) requires the testing of four regulated heavy metals (Pb, Hg, Cd, Cr(VI)) in electrotechnical products. After sample preparation according to IEC 62321 the determination of mercury in polymer materials can be carried out by anodic stripping voltammetry (ASV) at a gold rotating disk electrode (Au-RDE).
-----------------	--


Sample:	Polymer materials
Sample preparation:	Approx. 0.1 g of sample is mineralized by microwave digestion with a mixture of nitric acid, fluoroboric acid and hydrogen peroxide as described in IEC 62321.

Analysis of Hg													
Electrolyte	Perchloric acid electrolyte $c(\text{HClO}_4) = 0.22 \text{ mol/L}$ $c(\text{EDTA}) = 4 \text{ mmol/L}$ $c(\text{NaCl}) = 6 \text{ mmol/L}$												
Measuring solution	10 mL ultrapure water + 10 mL perchloric acid electrolyte + 0.05 mL digested sample solution (equals approx. 0.2 mg of sample)												
Working electrode (WE)	Au-RDE Gold electrode tip 6.1204.140 Driving axle 6.1204.210												
Auxiliary electrode (AE)	GC Glassy carbon rod 6.1247.000 Electrode holder 6.1241.020												
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L) 6.0728.020 Intermediate electrolyte: $c(\text{NaCl}) = 3 \text{ mol/L}$ 6.1245.010												
Parameters	<table border="1"> <tr> <td>Working electrode</td> <td>RDE</td> </tr> <tr> <td>Stirrer speed</td> <td>2000 rpm</td> </tr> <tr> <td>Mode</td> <td>DP</td> </tr> <tr> <td>Purge time</td> <td>300 s</td> </tr> <tr> <td>Cleaning potential</td> <td>+1.3 V</td> </tr> <tr> <td>Cleaning time</td> <td>30 s</td> </tr> </table>	Working electrode	RDE	Stirrer speed	2000 rpm	Mode	DP	Purge time	300 s	Cleaning potential	+1.3 V	Cleaning time	30 s
Working electrode	RDE												
Stirrer speed	2000 rpm												
Mode	DP												
Purge time	300 s												
Cleaning potential	+1.3 V												
Cleaning time	30 s												

Deposition potential	+0.37 V
Deposition time	30 s
Equilibration time	10 s
Pulse amplitude	0.05 V
Start potential	+0.4 V
End potential	+0.8 V
Voltage step	0.002 V
Voltage step time	0.1 s
Sweep rate	0.02 V/s
Peak potential Hg	+0.64 V

Results:	Hg
	1014.0 mg/kg

Determination of Hg



Hg
 $c = 1013.997 \text{ mg/kg}$
 $\pm 50.683 \text{ mg/kg (5.00\%)}$

