

VA Application Note No. V - 192

Title: Chromium(VI) in chromate coating on metallic materials as part of electrotechnical products

Summary: The EU directive on «Restriction of Hazardous Substances» (RoHS) requires the testing of four regulated heavy metals (Pb, Hg, Cd, Cr(VI)) in electrotechnical products. After sample preparation according to IEC 62321 the determination of chromium(VI) in chromate coating on metallic materials can be carried out by adsorptive stripping voltammetry (AdSV) using DTPA (diethylenetriamine pentaacetic acid) as complexing agent.

Sample: Colorless and colored chromate coating on metallic materials

Sample preparation: Boiling-water-extraction on a total surface area of 50 ± 5 cm² as described in IEC 62321.


Analysis of Cr(VI)

Electrolyte	<i>DTPA electrolyte</i> c(sodium acetate) = 0.2 mol/L c(DTPA) = 0.05 mol/L c(NaNO ₃) = 2.5 mol/L pH 6.20 ± 0.1 adjusted	
Measuring solution	10 mL ultrapure water + 2 mL DTPA electrolyte + 0.25 mL extraction solution	
Working electrode (WE)	MME (Multi Mode Electrode) With silanized capillary	6.1246.020 6.1226.050
Auxiliary electrode (AE)	Pt	6.0343.000
Reference electrode (RE)	Reference system: Ag/AgCl/KCl (3 mol/L) Intermediate electrolyte: c(KCl) = 3 mol/L	6.0728.020 6.1245.010
Parameters	Working electrode	HMDE
	Stirrer speed	2000 rpm
	Mode	DP
	Purge time	300 s
	Deposition potential	-0.9 V
	Deposition time	30 s
	Equilibration time	10 s

Pulse amplitude	0.05 V
Start potential	-1.0 V
End potential	-1.5 V
Voltage step	0.01 V
Voltage step time	0.3 s
Sweep rate	0.033 V/s
Peak potential Cr(VI)	-1.3 V

Results:	Cr(VI)
	0.020 mg/kg

Determination of Cr(VI)



Cr(VI)
 c = 0.020 mg/kg
 +/- 0.001 mg/kg (3.86%)

