


# Aluminum in drinking water by adsorptive stripping voltammetry using alizarin red S (DASA) as complexing agent



Aluminum can be determined in drinking water by adsorptive stripping voltammetry at the HMDE using alizarin red S (DASA) as complexing agent. The method is linear up to 35  $\mu\text{g/L}$ . The detection limit for this method is  $\beta(\text{Al}) = 1 \mu\text{g/L}$ , the limit of quantification is  $\beta(\text{Al}) = 3 \mu\text{g/L}$ . The sensitivity of the method cannot be increased by deposition.

## Results

Al in drinking water

4.0  $\mu\text{g/L}$

# Method description

## Sample

- Tap water

## Instrument

797 VA Computrace



## Electrodes

Working electrode (WE)	6.1246.020 MME (Multi-Mode Electrode) with 6.1226.050 silanized glass capillary
Reference electrode (RE)	6.0728.000 reference electrode (Ag/AgCl/ c(KCl) = 3 mol/L) with 6.1245.000 glass electrolyte vessel filled with intermediate electrolyte c(KCl) = 3 mol/L
Auxiliary electrode (AE)	6.0343.000 platinum electrode

## Reagents

HNO <sub>3</sub>	Nitric acid, suprapur, w(HNO <sub>3</sub> ) = 65 %
Ca solution	Commercial Ca standard stock solution, β(Ca) = 1 g/L
NH <sub>4</sub> Cl	Ammonium chloride, suprapur, NH <sub>4</sub> Cl
NaOH	Sodium hydroxide solution, suprapur, w(NaOH) = 30 %
KBrO <sub>3</sub>	Potassium bromate, analytical grade
DASA	Alizarin red S sodium salt (DASA, 3,4-Dihydroxy-9,10-dioxo-2-anthracenesulfonic acid sodium salt, C <sub>14</sub> H <sub>7</sub> NaO <sub>7</sub> S, CAS 130-22-3)
Al standard stock	β(Al <sup>3+</sup> ) = 1 g/L

solution

## Solutions

Ammonia buffer	c(NH <sub>4</sub> Cl) = 1.0 mol/L, adjusted to pH= 9.2 with NaOH
KBrO <sub>3</sub> solution	c(KBrO <sub>3</sub> ) = 0.1 mol/L
DASA solution	c(DASA) = 5.26 mmol/L in water
Al standard addition solution	β(Al <sup>3+</sup> ) = 1 mg/L, acidified with 1 mL/L HNO <sub>3</sub>

## Analysis

Measuring solution	10 mL sample + 0.1 mL Ca solution + 1 mL ammonia buffer + 1 mL KBrO <sub>3</sub> solution + 0.1 mL DASA solution
--------------------	--

## Parameters

Working electrode	HMDE
Drop size	4
Stirrer speed	2000 rpm
Mode	DP
Initial purge time	300 s
Addition purge time	120 s
Deposition potential	0 V
Deposition time	0 s
Equilibration time	5 s
Pulse amplitude	0.05 V
Start potential	-0.9 V
End potential	-1.2 V
Voltage step	0.005 V
Voltage step time	0.2 s
Sweep rate	0.025 V/s
Peak potential Al	-1.07 V