


# Determination of copper in electrolyte solutions for production of copper indium gallium selenide (CIGS) solar cells



Copper can be determined in electrolyte solutions used in production of copper indium gallium selenide solar cells – also called CIGS cells – for the electrodeposition of the CIGS absorber layer.

The determination of copper is carried out by polarography after dilution of the sample in sulfuric acid supporting electrolyte.

## Results

Cu in CIGS electrolyte solution

160 mg/L

# Method description

## Sample

Electrolyte solution containing sulfuric acid

## Sample preparation

No sample preparation necessary.

Pulse amplitude	0.05 V
Potential step	0.005 V
Potential step time	0.8 s
Sweep rate	0.006 V/s
Peak potential Cu	0.07 V

## Electrodes

Multi-Mode Electrode pro	6.1246.120
Non-silanized glass capillaries	6.1226.030
Ag/AgCl/ KCl (3 mol/L) reference electrode. Bridge electrolyte KNO <sub>3</sub> sat.	6.0728.020 6.1245.010
Platinum rod electrode	6.0343.000

## Reagents

H <sub>2</sub> SO <sub>4</sub>	Sulfuric acid, for trace analysis*, w(HCl) = 96%
Cu standard stock solution	β(Cu <sup>2+</sup> ) = 1 g/L, commercially available
Ultrapure water	Resistivity >18 MΩ·cm (25 °C), type I grade (ASTM D1193)

\*e.g., Merck suprapur®, Sigma-Aldrich TraceSelect® or equivalent.

## Solutions

Supporting electrolyte	c(H <sub>2</sub> SO <sub>4</sub> ) = 0.1 mol/L made by dilution of concentrated sulfuric acid.
Cu standard solution	β(Cu <sup>2+</sup> ) = 100 mg/L

## Analysis

Measuring solution	10 mL supporting electrolyte + 50 µL undiluted sample
--------------------	---

## Parameters

Working electrode	DME
Stirring speed	2000 min <sup>-1</sup>
Mode	DP
Purge time	300 s
Equilibration time	5 s
Start potential	0.2 V
End potential	-0.2 V

[www.metrohm.com](http://www.metrohm.com)

