

# Cation suppression in ion chromatography



Straightforward cation analysis at trace levels

## Suppression in cation analysis – benefits and applications

02

Trace analysis of cations, amines, and transition metals is possible with or without suppression. However, some applications require particularly high sensitivity of the analysis. This can only be achieved by sequential suppression, as suppression considerably lowers the detection limits for the analytes. Such analyses are common, for instance, in power plant applications or in pharmaceutical applications. Moreover, there are a number of norms and standards that request cation analysis with suppression.

In short, suppression reduces background conductivity to a minimum and decreases baseline noise. Both effects together improve the signal-to-noise ratio and increase the sensitivity of the measuring system. Thus, whenever the quantification of very low concentrations of cations is required, analysis with sequential suppression is the method of choice.

### The Metrohm Suppressor Module

The patented Metrohm Suppressor Module (MSM) is a very effective and highly robust solution when suppression is needed. The heart of the MSM is a small rotor that contains three cartridges, filled with anion exchange resin. While the first cartridge is being used for suppression, the second one is regenerated. At the same time, the third cartridge is rinsed to remove any residual regen-

erant. As the cartridges are rotated one step before each new injection, there is always a freshly regenerated suppressor cartridge ready for the next injection.


### STREAM – the green way of suppression

STREAM (Suppressor Treatment Re-using Eluent After Measurement) adds even more benefits to the suppression technique described. STREAM refers to a process whereby, instead of simply disposing of the suppressed eluent after detection, it is used to rinse the regenerated suppressor unit. Moreover, the use of the STREAM technique results in considerably lower consumption of regenerant.

#### The advantages of STREAM

- **Less waste**
- **Less use of consumables**
- **No additional rinsing solutions required**
- **Longer system run times due to lower regenerant consumption**

STREAM enables users to perform sustainable, «green chemistry» while saving on operating costs. Moreover, longer system run times frees up more time for the user to perform other tasks in the laboratory.





# Determination of cations, amines, and transition metals after sequential suppression

## Some typical examples are:

- Traces and ultratraces of sodium in the presence of monoethanolamine at high concentrations (typical of sample matrices in nuclear power plants).
- Trace and ultratrace concentrations of alkali and alkaline earth metals such as lithium, sodium, ammonium, potassium, magnesium or calcium in ultrapure water.
- Traces of transition metals, e.g., cobalt, nickel, zinc, manganese and cadmium in various types of water samples.
- Aliphatic and aromatic amines in pharmaceuticals, e.g., piperazine in cetirizine•HCl, tetrabutylammonium in atorvastatin, dimethylamine in meropenem, dimethylamine in imatinib mesylate, meglumine in meglumine salts.
- etc.

## Advantages

- Ultratrace range can be quantified
- Lower detection limits
- Higher measurement sensitivity due to improved signal-to-noise ratios
- Stable baselines with minimal noise, < 0.1 nS/cm
- Rapidly eluting peaks are better evaluated
- Option of using gradients
- 100% pressure stability
- 100% resistance to solvents
- Short conditioning time
- STREAM – the green way of suppression
- Robust analyses



## Lower limits of detection through improved signal-to-noise ratios

Differences between non-suppressed (a) and suppressed (b) cation chromatography. Sensitivity is increased by factor of 10 to 20. Determination of 10 µg/L standard cations in ultrapure water; column: Metrosep C Supp 1 - 150/4.0; eluent: 5 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 20 µL; flow rate: 1.0 mL/min, conductivity detection with sequential suppression




### Trace analysis of cations with the new Metrosep C Supp 1 column

The Metrosep C Supp 1 column is ideally suited for cation chromatography with suppression. The column excels by its outstanding separating efficiency, short retention times and high stability. Applications can be found in all major industries where cations, amines, and transition metals have to be determined at trace levels. Examples to be mentioned here include the pharmaceutical industry


(e.g., piperazine in cetirizine•HCl) and the energy industry (e.g., sodium in the cooling water of nuclear power plants).

Alkali and alkaline earth metals can also be detected with small injection volumes (20 µL) down to the low µg/L concentration range. A linear calibration function can be applied over a wide concentration range.




#### Determination of 1 µg/L standard cations in ultrapure water

Column: Metrosep C Supp 1 - 150/4.0; eluent: 5 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 20 µL; flow rate: 1.0 mL/min, conductivity detection with sequential suppression



#### Linear calibrations in the trace range 1–50 µg/L lithium, sodium, ammonium, potassium, magnesium, calcium

Column: Metrosep C Supp 1 - 150/4.0; eluent: 5 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 20 µL; flow rate: 1.0 mL/min, conductivity detection with sequential suppression



#### Ammonium, linear calibration of 1–50 µg/L, 3 injections per level

Column: Metrosep C Supp 1 - 150/4.0; eluent: 5 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 20 µL; flow rate: 1.0 mL/min, conductivity detection with sequential suppression

# Applications

06

The Metrosep C Supp 1 column in combination with the IC cation suppressor is recommended, for routine analyses as well as for research applications. The free choice of


column length ensures that the chromatography can be flexibly adapted to any application requirement.

## Typical applications

- Cations in the ultratrace range
- Aliphatic and aromatic amines in pharmaceuticals
- Ammonium determination in difficult matrices
- Transition metals in aqueous extracts


### Determination of 1 µg/L alkali, alkaline earth, and transition metals in the cooling water of a nuclear reactor, with fully automatic removal of the matrix


Column: Metrosep C Supp 1 - 250/4.0; eluent: 4 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 100 µL; sample preparation: Inline Matrix Elimination and Inline Preconcentration with Metrosep C PCC 1 HC/4.0; flow rate: 1.0 mL/min, conductivity detection with sequential suppression



### Determination of 1 µg/L lithium (1) and sodium (2) in a power plant sample with 4 mg/L monoethanolamine (3) (red). For comparison (black): Standard cations 0.1 µg/L


Column: Metrosep C Supp 1 - 250/4.0; eluent: 4 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 2000 µL; sample preparation: Inline Matrix Elimination and Inline Preconcentration with Metrosep C PCC 1 HC/4.0; flow rate: 1.0 mL/min, conductivity detection with sequential suppression






### Determination of choline in baby food analogous to standard AOAC 2012.20

Calculated for the original sample weight, the choline content is 82 mg/100 g of milk powder. Column: Metrosep C Supp 1 - 250/4.0; eluent: 4 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 20 µL; flow rate: 1.0 mL/min, conductivity detection with sequential suppression



### Analysis of zinc in a sports beverage

The relative standard deviation (n = 36) makes the precision of the measurement clear: Sodium 0.11%, potassium 1.60%, magnesium 0.31%, zinc 1.16% and calcium 2.01%. Column: Metrosep C Supp 1 - 250/4.0; eluent: 5 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 20 µL; flow rate: 1.0 mL/min, conductivity detection with sequential suppression



### Determination of aliphatic amines (1 mg/L each)

Column: Metrosep C Supp 1 - 250/4.0; eluent: 2.5 mmol/L HNO<sub>3</sub> + 7.5% (v/v) acetonitrile + 50 µg/L Rb<sup>+</sup>; column temperature: 40 °C; sample volume: 20 µL; flow rate: 1.0 mL/min, conductivity detection with sequential suppression

## Technical information

Rotor	MSM-HC C
<b>Design</b>	Micro packed-bed suppressor, robust, chemically resistant
<b>Capacity</b>	> 50 min* (0.25 meq)
<b>Dead volume</b>	< 250 µL
<b>Regeneration</b>	Chemical regeneration with STREAM
<b>Maximum flow rate</b>	Extremely robust: irreversible damage not possible
<b>Backpressure limitations</b>	Extremely robust: irreversible damage not possible, no specific backpressure required
<b>Solvent stability</b>	100% solvent stability
<b>Temperature range</b>	No limit

Separation column	Metrosep C Supp 1
<b>Substrate</b>	Polyvinyl alcohol with carboxyl groups
<b>Particle size</b>	5 µm
<b>Standard eluent</b>	5 mmol/L HNO <sub>3</sub> + 50 µg/L Rb <sup>+</sup>
<b>Standard flow rate</b>	1.0 mL/min
<b>Maximum flow rate</b>	1.5 mL/min
<b>Maximum pressure</b>	15 MPa
<b>Standard temperature</b>	40 °C
<b>Temperature range</b>	20–40 °C
<b>Capacity</b>	30 µmol K <sup>+</sup> (250 mm column)
<b>pH range</b>	1–12
<b>Organic modifier (in eluent)</b>	0–30% (no methanol)
<b>Organic modifier (in sample)</b>	0–100% acetone, acetonitrile and methanol
<b>Storage</b>	At 2–4 °C rinsed with ultrapure water

\*Standard conditions: Metrosep C Supp 1 - 250/4.0; eluent: 5 mmol/L HNO<sub>3</sub> + 50 µg/L Rb<sup>+</sup>, flow rate: 1.0 mL/min

## Ordering information

### Suppressor rotor

6.2842.200 MSM-HC C

### Accessories

6.2835.010 Connecting piece for MSM-HC and SPM

### Dosino Regeneration

2.800.0010 800 Dosino

6.5330.190 IC equipment: Dosino Regeneration

### Separation Columns

6.1052.410 Metrosep C Supp 1 - 100/4.0

6.1052.420 Metrosep C Supp 1 - 150/4.0

6.1052.430 Metrosep C Supp 1 - 250/4.0

6.1052.500 Metrosep C Supp 1 Guard/4.0

6.1052.510 Metrosep C Supp 1 S-Guard/4.0

[www.metrohm.com](http://www.metrohm.com)

