

The whole world of ion chromatography



Metrohm – the comprehensive solution



Metrohm has become a synonym for ion chromatography. For more than 30 years, Metrohm has been offering innovative and creative solutions in the area of IC separation columns and IC devices.



Welcome to the world of ion chromatography



With high-tech, long years of applications know-how, Swiss quality standards and a reasonable pricing policy, for procurement as well as for maintenance, Metrohm guarantees the optimum solution in ion chromatography.

Welcome to Metrohm!

Table of contents

Separation columns from Metrohm	10
iColumn	11
Which column for which application?	12
Preselection	13
A) Anions without chemical suppression	13
B) Anions with chemical suppression	14
C) Oxidizable anions	16
D) Cations without chemical suppression	17
E) Cations with chemical suppression	18
F) Organic acids	18
G) Carbohydrates	19
H) Amino acids	19
Capacity of the separation columns	20
Position of the system peak	21
«MCS» Metrohm CO ₂ Suppressor	22
Flexibility in application thanks to	
free selection of the pH value	23
Standards	24
ABC of practical work	26
Tips for eluent preparation	28
Inline Eluent Preparation	29
	29
Separation columns	30
IC anion-separation columns for analyses	
without chemical suppression	30
Hamilton PRP-X100 - 125/4.0 (6.1005.000)	32
Hamilton PRP-X100 - 250/4.0 (6.1005.010)	34
Super-Sep - 100/4.6 (6.1009.000)	36
IC anion-separation columns for analyses	
with or without chemical suppression	38
Metrosep Anion Dual 2 - 75/4.6 (6.1006.100)	40
Metrosep Anion Dual 3 - 100/4.0 (6.1006.120)	42
Metrosep Dual 4 - 100/4.6 (6.1016.030)	44

Metrosep A Supp 4 - 250/4.0 (6.1006.430)

IC anion-separation columns for analyses	
with chemical suppression	48
Phenomenex Star-Ion™ A 300 - 100/4.6	
(6.1005.100)	50
Metrosep A Supp 1 - 250/4.6 (6.1005.300)	52
Metrosep A Supp 1 HS - 50/4.6 (6.1005.350)	54
Metrosep A Supp 3 - 250/4.6 (6.1005.320)	56
Metrosep A Supp 5 - 50/4.0 (6.1006.550)	58
Metrosep A Supp 5 - 100/4.0 (6.1006.510)	60
Metrosep A Supp 5 - 150/4.0 (6.1006.520)	62
Metrosep A Supp 5 - 250/4.0 (6.1006.530)	64
Metrosep A Supp 7 - 150/4.0 (6.1006.620)	66
Metrosep A Supp 7 - 250/4.0 (6.1006.630)	68
Metrosep A Supp 10 - 50/4.0 (6.1020.050)	70
Metrosep A Supp 10 - 75/4.0 (6.1020.070)	72
Metrosep A Supp 10 - 100/4.0 (6.1020.010)	74
Metrosep A Supp 10 - 250/4.0 (6.1020.030)	76
Metrosep A Supp 16 - 100/4.0 (6.1031.410)	78
Metrosep A Supp 16 - 150/4.0 (6.1031.420)	80
Metrosep A Supp 16 - 250/4.0 (6.1031.430)	82
Metrosep A Supp 17 - 100/4.0 (6.01032.410)	84
Metrosep A Supp 17 - 150/4.0 (6.01032.420)	86
Metrosep A Supp 17 - 250/4.0 (6.01032.430)	88

Microbore IC anion-separation columns for lower

28	eluent consumption and greater sensitivity	90
29	Metrosep A Supp 4 - 250/2.0 (6.01021.230)	92
	Metrosep A Supp 5 - 150/2.0 (6.1006.220)	94
30	Metrosep A Supp 5 - 250/2.0 (6.1006.230)	96
	Metrosep A Supp 7 - 150/2.0 (6.1006.640)	98
30	Metrosep A Supp 7 - 250/2.0 (6.1006.650)	100
32	Metrosep A Supp 10 - 50/2.0 (6.1020.250)	102
34	Metrosep A Supp 10 - 75/2.0 (6.1020.270)	104
36	Metrosep A Supp 10 - 100/2.0 (6.1020.210)	106
	Metrosep A Supp 10 - 150/2.0 (6.1020.220)	108
	Metrosep A Supp 10 - 250/2.0 (6.1020.230)	110
38	Metrosep A Supp 16 - 100/2.0 (6.1031.210)	112
40	Metrosep A Supp 16 - 150/2.0 (6.1031.220)	114
42	Metrosep A Supp 16 - 250/2.0 (6.1031.230)	116
44		

IC separation columns for the determination of organic acids –			
ion-exclusion chromatography	118		
Hamilton PRP-X300 - 250/4.0 (6.1005.030)	120		
Metrosep Organic Acids - 100/7.8 (6.1005.210)	122		
Metrosep Organic Acids - 250/7.8 (6.1005.200)	124		

IC carbohydrate-separation columns – anionexchange chromatography applying pulsed amperometric detection (PAD)

•	
Metrosep Carb 2 - 100/4.0 (6.1090.410)	128
Metrosep Carb 2 - 150/4.0 (6.1090.420)	130
Metrosep Carb 2 - 250/4.0 (6.1090.430)	132
Hamilton RCX-30 - 150/4.6 (6.1018.010)	134
Hamilton RCX-30 - 250/4.6 (6.1018.000)	136

Microbore IC carbohydrate-separation columns for lower eluent consumption and greater sensitivity 138

Sensitivity	150
Metrosep Carb 2 - 100/2.0 (6.01090.210)	140
Metrosep Carb 2 - 150/2.0 (6.01090.220)	142
Metrosep Carb 2 - 250/2.0 (6.01090.230)	144

IC amino acid-separation column with

optical detection after post-column reaction	146
Metrosep Amino Acids 1 - 100/4.0 (6.4001.410)	148

weuosep	Amino	Acius	1 -	100/4.0	(0.4001.4	10)

IC cation-separation columns for analyses without chemical suppression

Nucleosil 5SA - 125/4.0 (6.1007.000)
Metrosep C 3 - 100/4.0 (6.1010.410)
Metrosep C 3 - 150/4.0 (6.1010.420)
Metrosep C 3 - 250/4.0 (6.1010.430)
Metrosep C 4 - 50/4.0 (6.1050.450)
Metrosep C 4 - 100/4.0 (6.1050.410)
Metrosep C 4 - 150/4.0 (6.1050.420)
Metrosep C 4 - 250/4.0 (6.1050.430)
Metrosep C 5 - 150/4.6 (6.4000.320)
Metrosep C 6 - 100/4.0 (6.1051.410)
Metrosep C 6 - 150/4.0 (6.1051.420)
Metrosep C 6 - 250/4.0 (6.1051.430)

Microbore IC cation-separation columns for lower eluent consumption and greater sensitivity 176 Metrosep C 4 - 100/2.0 (6.1050.210) 178 Metrosep C 4 - 150/2.0 (6.1050.220) 180 Metrosep C 4 - 250/2.0 (6.1050.230) 182 Metrosep C 6 - 100/2.0 (6.01051.210) 184

Metrosep C 4 - 250/2.0 (6.1050.230)	182
Metrosep C 6 - 100/2.0 (6.01051.210)	184
Metrosep C 6 - 150/2.0 (6.01051.220)	186
Metrosep C 6 - 250/2.0 (6.01051.230)	188

126 IC cation-separation columns for analyses with chemical suppression 190 30 Metrosep C Supp 1 - 100/4.0 (6.1052.410) 192 Metrosep C Supp 1 - 150/4.0 (6.1052.420) 194 32 Metrosep C Supp 1 - 250/4.0 (6.1052.430) 196 Metrosep C Supp 2 - 100/4.0 (6.01053.410) 6 198 Metrosep C Supp 2 - 150/4.0 (6.01053.420) 200

Separation column	n for the	determination
-------------------	-----------	---------------

Metrosep C Supp 2 - 250/4.0 (6.01053.430)

of organic substances	204
MetroSil RP 3 - 150/4.0 (6.01070.420)	206

144		
	IC guard columns (precolumns)	208
	IC guard column cartridge for Hamilton PRP-X100	
146	(6.1005.020)	210
148	Super-Sep Guard/4.6 (6.1009.010)	211
	Metrosep Dual 4 Guard Column Kit (6.1016.500)	212
	Metrosep A Supp 1 Guard/4.6 (6.1005.340)	213
150	Metrosep A Supp 4 Guard/4.0 (6.01021.500)	214
152	Metrosep A Supp 4 S-Guard/4.0 (6.01021.510)	214
154	Metrosep A Supp 4 Guard/2.0 (6.01021.600)	215
156	Metrosep A Supp 4 S-Guard/2.0 (6.01021.610)	215
158	Metrosep A Supp 5 Guard/4.0 (6.1006.500)	216
160	Metrosep A Supp 5 S-Guard/4.0 (6.1006.540)	216
162	Metrosep A Supp 5 Guard/2.0 (6.1006.600)	217
164	Metrosep A Supp 5 S-Guard/2.0 (6.1006.610)	217
166	Metrosep A Supp 10 Guard/4.0 (6.1020.500)	218
168	Metrosep A Supp 10 S-Guard/4.0 (6.1020.510)	218
170	Metrosep A Supp 10 Guard HC/4.0 (6.1020.520)	218
172	Metrosep A Supp 10 Guard/2.0 (6.1020.600)	219
174	Metrosep A Supp 10 S-Guard/2.0 (6.1020.610)	219

7

Metrosep A Supp 16 Guard/4.0 (6.1031.500)	220
Metrosep A Supp 16 S-Guard/4.0 (6.1031.510)	220
Metrosep A Supp 16 Guard/2.0 (6.1031.600)	221
Metrosep A Supp 16 S-Guard/2.0 (6.1031.610)	221
Metrosep A Supp 17 Guard/4.0 (6.01032.500)	222
Metrosep A Supp 17 S-Guard/4.0 (6.01032.510)	222
Metrosep A Supp 17 S-Guard - 50/4.0 (6.01032.530)	222
Metrosep Organic Acids Guard/4.6 (6.1005.250)	223
Metrosep Carb 2 Guard/4.0 (6.1090.500)	224
Metrosep Carb 2 S-Guard/4.0 (6.1090.510)	224
Metrosep Carb 2 Guard/2.0 (6.01090.600)	225
Metrosep Carb 2 S-Guard/2.0 (6.01090.610)	225
Nucleosil 5SA 2 Guard Cartridge/4.0 (6.1007.110)	226
Metrosep C 3 Guard/4.0 (6.1010.450)	227
Metrosep C 3 S-Guard/4.0 (6.1010.460)	227
Metrosep C 4 Guard/4.0 (6.1050.500)	228
Metrosep C 4 S-Guard/4.0 (6.1050.510)	228
Metrosep C 4 S-Guard - 50/4.0 (6.1050.530)	228
Metrosep C 4 Guard/2.0 (6.1050.600)	229
Metrosep C 4 S-Guard/2.0 (6.1050.610)	229
Metrosep C 6 Guard/4.0 (6.1051.500)	230

Metrosep C 6 S-Guard/4.0 (6.1051.510)

Metrosep C 6 Guard/2.0 (6.01051.600)

Metrosep C 6 S-Guard/2.0 (6.01051.610)

Metrosep C PCC 1 VHC/4.0 (6.1010.320)

Metrosep Chel PCC 1 VHC/4.0 (6.01010.350)

Metrosep C Supp 1 Guard/4.0 (6.1052.500)

Metrosep C Supp 1 S-Guard/4.0 (6.1052.510)

Metrosep C Supp 2 Guard/4.0 (6.01053.500)	233
Metrosep RP 2 Guard/3.5 (6.1011.030)	234
Metrosep RP 3 Guard HC/4.0 (6.1011.040)	235
MetroSil RP 3 Guard/4.0 (6.01070.500)	236
Metrosep BP 1 Guard/2.0 (6.1015.100)	237
Preconcentration columns	238
Metrosep A PCC 2/4.0 (6.1006.330)	240
Metrosep A PCC 2 HC/4.0 (6.1006.340)	240
Metrosep A PCC 2 VHC/4.0 (6.1006.350)	240
Metrosep C PCC 1/4.0 (6.1010.300)	241
Metrosep C PCC 1 HC/4.0 (6.1010.310)	241

IC trap columns Metrosep A Trap 1 - 100/4.0 (6.1014.000) Metrosep C Trap 1 - 100/4.0 (6.1015.000) Metrosep C Trap 1 - 30/4.0 (6.01015.030) Metrosep RP Trap 1 - 50/4.0 (6.1014.100) Metrosep RP Trap 2 - 100/4.0 (6.1014.150) Metrosep | Trap 1 - 100/4.0 (6.1014.200)

Metrosep BO₃³⁻ Trap 1 - 100/4.0 (6.1015.200)

Metrosep CO ₃ ²⁻ Trap 1 - 100/4.0 (6.1015.300)	253
IC sample-preparation cartridges	254
IC sample-preparation cartridge IC-RP (6.1012.X00)	256
IC sample-preparation cartridge IC-H (6.1012.X10)	256
IC sample-preparation cartridge IC-Ag (6.1012.X20)	256
IC sample-preparation cartridge IC-OH (6.1012.X30)	257
IC sample-preparation cartridge IC-Na (6.1012.X40)	257
IC sample-preparation cartridge IC-C18 (6.1012.X50)	257

IC accessories parts	258
PEEK inline filter (6.2821.120)	259
Coupling safety olive with PEEK inline filter	
(6.2744.180)	259

Separation columns from Metrohm

Separation columns from Metrohm are the backbone of high-performance analytics in ion chromatography. The combination of Metrohm IC separation columns and Metrohm IC systems guarantees:

- High separating efficiency
- Short analysis times
- Excellent reproducibility
- Long lifetime
- Low costs

Metrohm offers the right separation columns for all disciplines in ion chromatography:

- Anions with or without suppression
- Organic acids
- Cations with or without suppression
- Transition metals
- Carbohydrates
- Amines
- Amino acids

A small number of separation columns are sufficient for solving the majority of application problems. Metrohm offers the correct column material, both for standard applications and for complex separation tasks:

- Polystyrene divinylbenzene
- Polyvinyl alcohol
- Polymethacrylate
- Silica gel
- Monolith

This flexibility guarantees top performance. The 2 mm columns enable reduced eluent consumption, low detection limits, and universal use.

Ion chromatography with Metrohm systems can be carried out with or without chemical suppression. Therefore, the application determines the pH value and not the other way around. A large selection of anion eluents is thus available:

- Phthalate
- Benzoate
- Borate
- Hydrogen carbonate
- Carbonate
- Hydroxide

Metrohm is the comprehensive solution for applications in ion chromatography. For more than 30 years, Metrohm has been offering innovative and creative solutions in the area of IC separation columns and IC instruments. A comprehensive network of application laboratories is available worldwide to provide the correct answers to analytical questions: on-site, quickly, and focused on the customer. The Internet (www.metrohm.com) provides an extensive database with information and problem solutions for all areas of ion analysis.

Metrohm – the whole world of ion chromatography

iColumn

The world's first intelligent column generation in ion chromatography – just click and start

- All column data available immediately
- Active monitoring of all important column functions
- Uninterrupted retraceability of all column parameters

Intelligent columns, built into all Metrohm IC instruments - the MagIC Net software registers immediately which separation columns are available to the IC system. One click, and the software detects typical standard conditions for columns such as standard eluents and flow rate, the permissible maximum values for pressure and flow, and the correct guard column. This data and other information are applied in the method if desired. As a result, it is no longer necessary to enter all of this data into the method. The iColumn knows how many hours it has already been in operation and how many samples it has already analyzed. MagIC Net monitors the separating efficiency of the column and switches the system off in the event that the permissible pressure is exceeded. If the performance of the column falls below a previously defined value, then MagIC Net can even have the column reordered automatically by E-mail.

If the Metrosep separation column is used in different systems, then it will take along all information stored on its memory chip to the next IC system. This allows for uninterrupted traceability and GLP monitoring, even on different IC devices, for all columns used. System validation is simplified significantly.

The iColumn concept is flexible and encompasses the following data types:

- Freely definable data, e.g.
 - Column name
 - Comment field in which, for example, the name of the application can be entered
- Data permanently linked with the column, e.g.
 - Column type (e.g. Metrosep A Supp 10 100/4.0)
 - Order number
 - Serial number
 - Standard flow
 - Standard eluent
 - Standard injection volume
- Standard temperature
- Length x inner diameter
- Particle size
- pH range
- Maximum permitted pressure
- Maximum permitted flow rate
- Data entered by the IC system and the MagIC Net software, e.g.
 - Operating hours
 - Number of injections
 - Maximum pressure
 - Maximum flow used

With the exception of the Metrosep Anion Dual 2 - 75/4.6 (6.1006.100), all Metrosep separation columns are available exclusively as iColumns. Excluded from this are the other column types, i.e. guard, preconcentration, and trap columns as well as separation columns which do not have the Metrosep designation.







Which column for which application?

Find the column finder online under https://www.metrohm.com/Column-Finder

Preselection

Anions	Anions	High concentration Simple setup HBO_3^{2-} $H_2SiO_4^{2-}$ no HPO_4^{2-}	without suppression	→ A
		Entire concentration spectrum	with suppression	→ B
	Oxidizable anions	Amperometric detection		\rightarrow C
Cations			without suppression with suppression	\rightarrow D \rightarrow E
Additional analyt- es	Organic acids		with or without suppression	→ F
	Carbohydrates			\rightarrow G
	Amino acids	Post-column reaction with	ninhydrin	→H

This symbol indicates the respective standard column

A) Anions without chemical suppression

Requirements or application	Column	Page
No F	IC anion column	
Simple separation problems	Metrosep A Supp 4 - 250/x.0	
Simple matrices	6.1006.430 (250/4.0)	46
Rapid separation	6.01021.230 (250/2.0)	92
F	IC anion column	
Difficult separation problems	Metrosep Anion Dual 2 - 75/4.6	
Difficult matrices	6.1006.100	40
No F⁻	IC anion column	,
Difficult separation problems	Metrosep Anion Dual 3 - 100/4.0	\checkmark
Difficult matrices	6.1006.120	42
Biological samples		
Cl ⁻ , NO ₃ ⁻ , SO ₄ ²⁻	IC anion columns	
Difficult matrices (e.g. dyes)	Hamilton PRP-X100	
HBO ₃ ²⁻ , H ₂ SiO ₄ ²⁻	6.1005.000 (125/4.0)	32
	6.1005.010 (250/4.0)	34
No F⁻, acetate	IC anion column	
Difficult matrices	Super-Sep - 100/4.6	
Special applications (e.g. BF_4^{-})	6.1009.000	36

B) Anions with chemical suppression

Requirements or application	Column	Page
F ⁻ , formate, acetate, Cl ⁻ SO ₃ ²⁻ , SO ₄ ²⁻ Simple matrices Standard anions ClO_2^{-} , ClO_3^{-} , ClO_4^{-} , BrO_3^{-} Special separation problems	IC anion column Metrosep Anion Dual 2 - 75/4.6 6.1006.100	40
F ⁻ , formate, acetate, Cl ⁻ SO ₃ ²⁻ , SO ₄ ²⁻ Simple matrices Standard anions ClO_2^- , ClO_3^- , ClO_4^- , BrO_3^- Special separation problems Biological samples	IC anion column Metrosep Anion Dual 3 - 100/4.0 6.1006.120	42
Perchlorate in difficult matrices, EPA 314 Very high ionic strength	IC anion column (Monolith) Metrosep Dual 4 - 100/4.6 6.1016.030	44
Great differences in concentration High ionic strength ClO_{2}^{-} , ClO_{3}^{-} , ClO_{4}^{-} , BrO_{3}^{-}	IC anion column Metrosep A Supp 1 - 250/4.6 6.1005.300	52
SCN ⁻ , SO ₃ ²⁻ , SO ₄ ²⁻ , S ₂ O ₃ ²⁻ Polyphosphates	IC anion column Metrosep A Supp 3 - 250/4.6 6.1005.320	56
Standard anions Difficult matrices Critical samples I ⁻	IC anion column Metrosep A Supp 4 - 250/x.0 6.1006.430 (250/4.0) 6.01021.230 (250/2.0)	√ 46 92
Standard anions F^- , Cl^- , Br^- , l^- ClO_2^- , ClO_3^- , ClO_4^- , BrO_3^- BrO_3^- at high ionic strength $Cr(VI)$ (CrO_4^{-2-}) l^- (not with 250 mm) Method development	IC anion columns Metrosep A Supp 5 6.1006.550 (50/4.0) 6.1006.510 (100/4.0) 6.1006.520 (150/4.0)	58 60 √ 62
Universal applications Difficult matrices Difficult separation problems Rapid separation (with 50 and 100 mm) PO ₄ ³⁻ in soft drinks (with 100 mm) IC-MS coupling Applications with gradient	6.1006.530 (250/4.0) 6.1006.220 (150/2.0) 6.1006.230 (250/2.0)	64 94 96
Standard anions Oxohalides, EPA 300 (with 250 mm) Isocratic separation of glycolate and acetate Difficult separations	IC anion columns Metrosep A Supp 7 6.1006.620 (150/4.0) 6.1006.630 (250/4.0)	66 √ 68
Bayer liquors Applications with gradient	6.1006.640 (150/2.0) 6.1006.650 (250/2.0)	68 98 100

Requirements or application	Column	Page
Anions in salt solutions	Metrosep Carb 2 - 100/x.0	
	6.1090.410 (100/4.0)	128
	6.01090.210 (100/2.0)	140
PO_4^{3-} in soft drinks with cyclamate	IC anion columns	
Standard anions (no F⁻)	Metrosep A Supp 10	
SCN ⁻ , SO ₃ ²⁻ , SO ₄ ²⁻ , S ₂ O ₃ ²⁻	6.1020.050 (50/4.0)	70
Separation SO_3^{2-} , SO_4^{2-}	6.1020.070 (75/4.0)	72
Aerosols with PILS/MARGA (75 mm)	6.1020.010 (100/4.0)	74
Air analytics	6.1020.030 (250/4.0)	76
IC-MS coupling	6.1020.250 (50/2.0)	102
Aggressive matrices	6.1020.270 (75/2.0)	104
	6.1020.210 (100/2.0)	106
	6.1020.220 (150/2.0)	108
	6.1020.230 (250/2.0)	110
Standard anions	IC anion columns	
Universal applications	Metrosep A Supp 16 - 100/x.0	
Non-critical matrices	6.1031.410 (100/4.0)	78
Bro ₃ " (EPA 326, DIN EN ISO 11206)	6.1031.210 (100/2.0)	112
IC-MS coupling		
Standard anions	IC anion columns	1
Universal applications	Metrosep A Supp 16 - 150/x.0	\checkmark
Complex matrices	6.1031.420 (150/4.0)	80
IC-MS coupling	6.1031.220 (150/2.0)	114
Standard anions	IC anion columns	
Universal applications	Metrosep A Supp 16 - 250/x.0	
Oligosaccharides and polysaccharides	6.1031.430 (250/4.0)	82
Cl^{-} , SO_4^{2-} in electroplating baths	6.1031.230 (250/2.0)	116
Silicate in addition to standard anions (4 mm column)		
Quality monitoring of high-purity chemicals		
(e.g. conc. acids)		
Complex separation problems		
Difficult matrices		
IC-MS coupling		
Standard anions	IC anion columns	
Water analysis	Metrosep A Supp 17 - xx0/4.0	
	6.01032.410 (100/4.0)	84
	6.01032.420 (150/4.0)	\checkmark
		86
	6.01032.430 (250/4.0)	88
Short retention times	IC anion column	
No Cl ⁻	Phenomenex Star-Ion™ A 300 - 100/4.6	
	6.1005.100	50

C) Oxidizable anions

Requirements or application	Column	Page
CN ⁻ S ²⁻	IC anion column Metrosep A Supp 1 - 250/4.6 6.1005.300	52
CN [−] S ^{2−}	IC anion column Metrosep A Supp 10 - 100/x.0 6.1020.010 (100/4.0) 6.1020.210 (100/2.0)	√ 74 106
ClO ₂ ⁻ , NO ₂ ⁻ , S ₂ O ₃ ²⁻ , SCN ⁻ , I ⁻	IC anion column Super-Sep - 100/4.6 6.1009.000	36
Br⁻, I⁻	IC anion column Metrosep Anion Dual 2 - 75/4.6 6.1006.100	40
NO ₂ ⁻ , ClO ₂ ⁻ S ₂ O ₃ ²⁻ , SCN ⁻ , I ⁻	IC anion column Metrosep A Supp 5 - 100/4.0 6.1006.510	60

D) Cations without chemical suppression

Requirements or application	Column	Page
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ Very rapid separations Simple matrices	IC cation column Metrosep C 4 - 50/4.0 6.1050.450	160
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ Lipophilic amines with short retention times Rapid separations	IC cation columns Metrosep C 4 - 100/x.0 6.1050.410 (100/4.0) 6.1050.210 (100/2.0)	162 178
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ Transition metals Amines	IC cation columns Metrosep C 4 - 150/x.0 6.1050.420 (150/4.0) 6.1050.220 (150/2.0)	164 180
Li [*] , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ , Mn ²⁺ , Co ²⁺ , Ni ²⁺ , Zn ²⁺ , Cd ²⁺ , Pb ²⁺ , amines NH ₄ ⁺ , ethanolamines Na ⁺ /NH ₄ ⁺ separation NH ₄ ⁺ , methylamines, and ethylamines Transition metals Difficult separation problems Great differences in concentration	IC cation columns Metrosep C 4 - 250/x.0 6.1050.430 (250/4.0) 6.1050.230 (250/2.0)	166 182
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ Lipophilic amines with short retention times Rapid separations	IC cation column Metrosep C 6 - 100/x.0 6.1051.410 (100/4.0) 6.01051.210 (100/2.0)	170 184
Amines Transition metals Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺	IC cation column Metrosep C 6 - 150/x.0 6.1051.420 (150/4.0)	√ 172
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ , Co ²⁺ , Ni ²⁺ , Zn ²⁺ , Cd ²⁺ , Pb ²⁺ , Amine Very good Na ⁺ /NH ₄ ⁺ separation NH ₄ ⁺ , (CH ₃)NH ₃ ⁺ , (CH ₃) ₂ NH ₂ ⁺ , (CH ₃) ₃ NH ⁺ , (CH ₃) ₄ N ⁺ , and the respective ethanolamines Difficult separation problems Great differences in concentration Transition metals	6.01051.220 (150/2.0) IC cation columns Metrosep C 6 - 250/x.0 6.1051.430 (250/4.0) 6.01051.230 (250/2.0)	186 174 188
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ , Co ²⁺ , Ni ²⁺ , Zn ²⁺ , Cd ²⁺ , Pb ²⁺ , amines NH ₄ ⁺ , monoethanolamine Transition metals Na ⁺ /NH ₄ ⁺ separation Matrices with high pH	IC cation columns Metrosep C 3 6.1010.410 (100/4.0) 6.1010.420 (150/4.0) 6.1010.430 (250/4.0)	154 156 158
Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , Fe ²⁺ , Co ²⁺ , Ni ²⁺ , Cd ²⁺ , Zn ²⁺ , Mn ²⁺ Mg ²⁺ , Ca ²⁺ in addition to large amounts of Na ⁺	IC cation column Nucleosil 5SA - 125/4.0 6.1007.000	152
Transition metals	IC cation column Metrosep C 5 - 150/4.6 6.4000.320	168
Transition metals, U, and Pu	see footnote on next page	

E) Cations with chemical suppression

Requirements or application	Column	Page
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ Lipophilic amines with short retention times Rapid separations	IC cation columns Metrosep C Supp 1 - 100/4.0 6.1052.410	192
Trace analysis	Metrosep C Supp 2 - 100/4.0 6.01053.410	198
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ Transition metals Amines Trace analysis	IC cation columns Metrosep C Supp 1 - 150/4.0 6.1052.420 Metrosep C Supp 2 - 150/4.0 6.01053.420	194 √ 200
Li ⁺ , Na ⁺ , K ⁺ , Rb ⁺ , Cs ⁺ , Mg ²⁺ , Ca ²⁺ , Sr ²⁺ , Ba ²⁺ , NH ₄ ⁺ , Mn ²⁺ , Co ²⁺ , Ni ²⁺ , Zn ²⁺ , Cd ²⁺ , Pb ²⁺ , amines NH_4^{+} , ethanolamines	IC cation columns Metrosep C Supp 1 - 250/4.0 6.1052.430	196
Na ⁺ /NH ₄ ⁺ separation NH ₄ ⁺ , methylamines, and ethylamines Transition metals Difficult separation problems Great differences in concentration Trace analysis	Metrosep C Supp 2 - 250/4.0 6.01053.430	202
Transition metals, U, and Pu	*	

F) Organic acids

Requirements or application	Column	Page
Organic acids: Citrate, tartrate, malate, ascorbate,	IC exclusion columns	
succinate, etc.	Metrosep Organic Acids	
Short-chain fatty acids: Formate, acetate, propionate,	6.1005.210 (100/7.8)	122
butyrate, etc.	6.1005.200 (250/7.8)	124
Simple matrices (100 mm)		1
Difficult matrices (250 mm)		v
Simple separation problems (100 mm)		
Difficult separation problems (250 mm)		
Glycolic acid, monochloroacetic acid	IC exclusion column	
Simple matrices	Hamilton PRP-X300 - 250/4.0	
Simple separation problems	6.1005.030	120
Formate determination		

* Reliable determinations of transition metals as well as uranium and plutonium can be made into the ultra trace range using voltammetry and polarography.

G) Carbohydrates

Requirements or application	Column	Page
Monosaccharides	IC carbohydrate column	
Disaccharides	Metrosep Carb 2 - 100/x.0	
Sugar alcohols	6.1090.410 (100/4.0)	128
Oligosaccharides	6.01090.210 (100/2.0)	140
Simple separation problems		
Very rapid separations		
Monosaccharides	IC carbohydrate column	1
Disaccharides	Metrosep Carb 2 - 150/x.0	\checkmark
Sugar alcohols	6.1090.420 (150/4.0)	130
Anhydrosugars	6.01090.220 (150/2.0)	142
Oligosaccharides		
Rapid separations		
Monosaccharides	IC carbohydrate column	
Disaccharides	Metrosep Carb 2 - 250/x.0	
Sugar alcohols	6.1090.430 (250/4.0)	132
Anhydrosugars	6.01090.230 (250/2.0)	144
Complex separations		
Monosaccharides	IC carbohydrate column	
Disaccharides	Hamilton RCX-30 - 150/4.6	
Oligosaccharides	6.1018.010	134
Sugar alcohols		
Simple separation problems Rapid separation		
Monosaccharides	IC carbohydrate column	
Disaccharides	Hamilton RCX-30 - 250/4.6	
Sugar alcohols	6.1018.000	136
Difficult separation problems		
Difficult matrices		

H) Amino acids

Requirements or application	Column	Page
Amino acids	IC amino acid column	
	Metrosep Amino Acids 1 - 100/4.0	
	6.4001.410	148

Capacity of the separation columns

The capacity of a separation column is determined by the type of the stationary phase used. The capacity has no direct influence on selectivity, whereas the column material does.

In addition, the capacity of a separation column changes in proportion to the quantity of packaging material used. This means that the capacity of a separation column also increases as the column length and diameter increase.


The capacity of separation columns can be determined with a wide variety of methods, all of which can be justified scientifically. The capacities specified here have to do with chloride or potassium exchange capacities, respectively, which are calculated by means of static charging. The specifications of other manufacturers are based to some extent on proton exchange and neutralization methods. The latter results in disproportionately higher numerical values. The capacity specifications of a column manufacturer can be used to compare different columns of this manufacturer. Capacity values from different manufacturers that apply different determination methods are not comparable.

Which capacity is right? The following rules apply:

- Simple separation tasks, weakly ionic matrix
 → Small capacity and therefore rapid separation of
 the analytes
- Complex separation tasks, strongly ionic matrix
 → High capacity and therefore long retention times for the analytes

The practical solution is often to be found somewhere in between. Most separations can thus be solved with just a limited number of column types:


- Determination of anions: Metrosep A Supp 4 -250/4.0, Metrosep A Supp 5 - 150/4.0, Metrosep A Supp 16 - 150/4.0
- Determination of cations: Metrosep C 4 150/4.0



Position of the system peak

When work is performed with carbonate eluent, a system peak is always present. Its size and position are determined by various factors. The position of the system peak, however, corresponds in the first approximation to the retention time of carbonate. For this reason this peak is also called the system peak. To ensure that the measured values are not skewed, it is important to know the position of the system peak in the chromatogram. With columns based on polystyrene in particular the system peak often lies directly beneath the chloride peak. The position of the system peak with the respective standard carbonate eluents is shown below.

If the CO₂ suppressor «MCS» is used, then the influence of carbonate can be virtually excluded. Both system peak and injection peak are minimized and do not interfere with the integration.




Position of the system peaks on Metrohm anion columns using the respective standard eluent.

«MCS» Metrohm CO₂ Suppressor

The «MCS» eliminates both carbonate from the sample and CO_2 that develops during the suppression reaction. As a result, the injection peak is practically nonexistent and the peak areas of the analyte ions are significantly larger. In addition, the system peak is effectively eliminated. The MCS is based on the gas permeability of a fluoropolymer membrane. The special system setup with an integrated vacuum cell, a fluoropolymer membrane, and a CO₂ adsorber cartridge is controlled by MagIC Net.

Chromatography the way it should always be




Chromatography with (red) and without sequential suppression (blue)

A unique feature of Metrohm IC is the combination of chemical suppression with the «MSM» (Metrohm Suppressor Module) and CO_2 -suppression using the «MCS» (Metrohm CO_2 Suppressor) for sequential suppression. This technique achieves lowest background conductivities for anion and cation suppression.

- Extremely low background conductivity
- Very small injection peak (1)
- Larger peak areas
- Lower detection limits (3)
- Minimized carbonate influence (2)
- No system peak (2)

Flexibility in application thanks to free selection of the pH value

Metrohm enables successful working with a large number of eluents. The correct choice of eluent has a decisive influence on the analytics. With Metrohm, there are no limitations and you can use virtually the entire pH range to obtain an optimum separation.



pH ranges of various eluents for anion chromatography

Standards

24

The world of standards is also changing. Today it is the result that matters when it comes to compliance with standards. Which instrument was used to calculate the result is of secondary importance. This was not always the case. The fact that today you have a free choice of analytical instruments is in part due to the commitment of Metrohm application chemists, who have demonstrated to standards committees that new methods can also produce the correct result and therefore meet requirements.

The following standards deal with the determination of anions and cations in water using ion chromatography. The analytical requirements of these standards can be met with Metrohm IC systems. The separation columns with which the standard can be fulfilled are indicated in brackets.

EPA 300.1 Part A and Part B

Determination of inorganic anions in drinking water by ion chromatography. (Metrosep A Supp 7 - 250/4.0; 6.1006.630; see page 68)

EPA 314.0

Determination of perchlorate in drinking water by ion chromatography. (Metrosep Dual 4 - 100/4.6; 6.1016.030; see page 44)

EPA 218.7

Determination of dissolved hexavalent chromium by means of ion chromatography (post-column reaction and VIS detection). (Metrosep A Supp 5 - 150/4.0; 6.1006.520; see page 62)

EPA 317.0

Determination of inorganic oxyhalide disinfection byproducts in drinking water using ion chromatography. Trace analysis of bromate by means of post-column reaction – «o-dianisidine method». (Metrosep A Supp 5 -250/4.0; 6.1006.530; see page 64)

EPA 326, DIN EN ISO 11206

Determination of inorganic oxidation by-products in drinking water using ion chromatography. Trace analysis of bromate by means of post-column reaction and UV detection – «triiodide» method. Improved method without acidification using a suppressor. (Metrosep A Supp 16 - 100/4.0; 6.1031.410; see page 78)

EPA 332.0

Determination of perchlorate in drinking water by ion chromatography with chemical suppression as well as conductivity detection and ESI/MS detection. (Metrosep A Supp 5 - 100/4.0; 6.1006.510; see page 60)

DIN 38413-8

Determination of the dissolved complexing agents nitrilotriacetic acid (NTA), ethylenediaminetetraacetic acid (EDTA), and diethylenetriaminepentaacetic acid (DTPA) by liquid chromatography (LC). (2 x MetroSil RP 3 - 150/4.0; 6.01070.420; see page 206)

DIN EN ISO 14911

Water quality – Determination of dissolved lithium, sodium, ammonium, potassium, manganese(II), calcium, magnesium, strontium, and barium using ion chromatography – Method for water and wastewater. (Metrosep C 4 - 150/4.0; 6.1050.420; see page 164)

DIN EN 13368-1

Determination of chelating agents in fertilizers by ion chromatography, Part 1: HEDTA, EDTA, DTPA. (Metrosep A Supp 3 - 250/4.6; 6.1005.320; see page 56)

DIN 38405-7

Determination of cyanides in slightly polluted water by ion chromatography or potentiometric titration (Metrosep A Supp 10 - 100/2.0; 6.1020.210; see page 106)

DIN EN ISO 10304-1

Water quality – Determination of the dissolved anions fluoride, chloride, nitrite, orthophosphate, bromide, nitrate, and sulfate by liquid chromatography – Part 1: Method for slightly polluted wastewater. (Metrosep A Supp 17 - 250/4.0; 6.01032.430; Metrosep A Supp 16 - 150/4.0; 6.1031.420 or Metrosep A Supp 5 - 250/4.0; 6.1006.530)

DIN EN ISO 10304-3

Water quality – determination of dissolved anions by liquid chromatography (LC) – Part 3: Determination of chromate, iodide, sulfite, thiocyanate, and thiososulfate in wastewater. (Metrosep A Supp 17 - 150/4.0; 6.01032.420 or Metrosep A Supp 5 - 150/4.0; 6.1006.520)

DIN EN ISO 10304-4

Water quality – Determination of dissolved anions by liquid chromatography (LC) – Part 4: Determination of chlorate, chloride, and chlorite in slightly polluted waste-water. (Metrosep Anion Dual 2 - 75/4.6; 6.1006.100 or Metrosep A Supp 5 - 250; 6.1006.530)

DIN EN ISO 15061

Determination of dissolved bromate in water. (Metrosep A Supp 5 - 250/4.0; 6.1006.530 or Metrosep A Supp 7 - 250/4.0; 6.1006.630)

ABC of practical work

Bacterial growth

Bacterial growth has a significant negative effect on chromatography and destroys the analytical columns. A large number of chromatographic problems can be traced back to the growth of algae, bacteria, and molds. In order to prevent bacterial growth, eluents, rinsing, and regeneration solutions should always be prepared fresh and not reused after prolonged periods. We recommend that all vessels be thoroughly rinsed with ultrapure and UV-treated water and then rinsed with methanol/water or acetone/water and finally again with water before being refilled. If bacteria or algae should form despite this treatment, then 5% methanol or acetone can be added to the eluent. This is not possible when using membrane suppressors, because these could be destroyed by organic solvents. The «MSM», «MSM-HC», and «MSM-LC» Metrohm Suppressor Modules are 100% solvent-resistant. Methanol should not be used with some cation columns.

Cation analyses

For all analyses we recommend that the samples be acidified with nitric acid (approximately 100 μ L 2 mol/L HNO₃ per 100 mL of sample) (pH 2.5–3.5), otherwise divalent cation results may be overestimated depending on the age of the injection loop.

Chemical stress

Although many separation phases cover a wide pH range in terms of specification, this does not mean that they are chemically inert. Separation columns achieve their longest service life under constant chemical conditions. A column must never be allowed to dry out and must always be kept well-sealed.

\mathbf{CO}_2

Carbon dioxide from air affect alkaline eluents. To avoid this, the eluent bottle should always be furnished with a CO_2 adsorber material («soda lime»). Eluents with a weak buffer capacity must also be protected against CO_2 .

Degassing the eluent

In order to prevent bubble formation, we recommended to use the Eluent Degasser in the IC instrument. Alternatively this is done by applying a vacuum created by a water-jet pump or vacuum pump for approximately 10 minutes or by means of an ultrasonic bath.

Eluent bottles

Eluents are positioned in special eluent bottles, usually directly on the IC system. To prevent moisture and carbon dioxide from being absorbed by the eluent, the bottles are equipped with a drying tube which normally has a molecular sieve and is filled with soda lime (as a weak CO_2 adsorber material) for sodium hydroxide and carbonate eluents.

Environmental protection

A great advantage of ion chromatography is that most work is carried out with aqueous media. The chemicals used in ion chromatography are therefore as non-toxic as possible and do not pollute the environment. Nevertheless, when work is carried out with acids, bases, organic solvents, or heavy metal standards, they must be disposed of properly after use.

Filter

If problems occur with IC systems, they are usually due to particles introduced by bacterial growth, unfiltered eluents, by the sample or by rinsing and regeneration solutions. This risk can be reduced to an absolute minimum by using an aspiration filter (6.2821.090), inline filter (6.2821.120), and guard columns (starting on page 208). The filters are part of the basic equipment of the Metrohm ion chromatographs and are included in the scope of delivery. We strongly recommend their use. Care should be taken to ensure that the filters are replaced regularly.

Filtration of the eluent

All eluents should be microfiltered (0.45 $\mu\text{m})$ immediately before being used.

Fun

Ion chromatography should be fun and not get on your nerves. Metrohm does everything it can to ensure that your IC systems work reliably with a minimum of upkeep, maintenance, and cost. Metrosep separation columns stand for quality, long lifetime, and outstanding results.

Guard columns (precolumns)

Guard columns (starting on page 208) are used to protect the valuable separation columns. We strongly recommend their use. As a rule they contain the same stationary phase as the separation column, although in a considerably smaller quantity to avoid influencing the chromatography. Guard columns eliminate critical contaminations which might react with the column material and they effectively eliminate particles and bacterial contamination. Guard columns need to be replaced if

- the backpressure in the system rises
- the chromatograohy gets worse

It is recommended to use 3–4 guard columns during the lifetime of an analytical column. Guard columns are available for all Metrosep separation columns.

Long-term storage of the ion chromatograph

If the ion chromatograph will not be used for a prolonged period (>1 week), then the separation column should be removed and sealed with the stoppers provided. The ion chromatograph should be rinsed with methanol/water (1:4). Care should be taken to ensure that all three chambers of the suppressor are rinsed during this process. The separation column should be stored in the medium listed on the column data sheet, optimally between 4 and 8 °C. When the instrument is restarted, rinse the system with fresh eluent before installing the separation column and bring it up to room temperature.

Particles

All solutions, samples, regeneration solutions, the water and the eluents should be free of particles because they may clog the separation columns over time (increase in column pressure). This must be taken into account particularly when eluents are being produced, because eluents flow continuously through the column (500–1000 mL per working day in contrast to approximately 0.5 mL of sample solution). The sample can be filtered or dialyzed fully automatically with the «MISP» Metrohm Inline Sample Preparation systems.

Pulsation absorber

We recommend the use of a pulsation absorber (6.2620.150). In particular, polymethacrylate and polyvinyl alcohol columns should be protected against brief pressure surges which inevitably occur when the valves are switched. This protection is ensured when a pulsation absorber is used.

Quality of chemicals

All chemicals should be at least of p.a. or puriss. quality. The standards must be specially suited to ion chromatography.

Regeneration of separation columns

As a rule, if separation columns are operated with clean eluents and charged with particle-free samples, then a very long lifetime is guaranteed. A regeneration of the column is then not necessary and is also no longer possible after a large number of injections. Nevertheless, if the pressure in the column should rise unexpectedly or the separating efficiency decrease, then the regeneration steps which are indicated for each separation column can be carried out. In general, it must be noted that the regeneration takes place outside the analytical line. This means that the separation column is connected directly to the pump and the regeneration solution feeds through the column directly into the waste vessel. Before the separation column is reinstalled, it should be rinsed sufficiently - for 30 minutes at standard flow - with fresh eluent.

Sample-preparation cartridges

Sample-preparation cartridges are used for the preparation of critical samples which cannot be injected directly on the separation columns. Thus, for example, samplepreparation cartridges remove organic contamination or neutralize strongly alkaline or acidic samples. Samplepreparation cartridges are consumable materials which, as a rule, cannot be regenerated. They do not replace the guard column (precolumn), which should always be used with each separation column. «MISP» (Metrohm Inline Sample Preparation) offers an alternative to sample cartridges, e.g. for the fully automated neutralization of alkaline samples.

Water quality

Ion chromatography primarily involves work in aqueous media. Water quality is therefore of decisive importance for obtaining good chromatographic results. If the water quality is unsatisfactory, then the results will certainly be unsatisfactory as well. In addition, there is the risk of damaging instruments and separation columns due to insufficient water quality. The ultrapure water used should have a specific resistance greater than 18 M Ω ·cm and be particle-free. It is therefore recommended that the water be filtered through a 0.45 µm filter and treated with UV. Modern ultrapure water plants for laboratory use guarantee this water quality (Type I).

Tips for eluent preparation

Please note that the eluents must be degassed once in order to avoid bubble formation during the measurements. Degassing can be carried out fully automatically by the eluent degasser in the IC instrument. Alternatively, the ultrapure water used can already be degassed before the reagents are added.

Excellent water quality (high resistance, absence of particles, and bacteria) is crucial for good ion chromatography (see also the chapter «ABCs of practical work», starting on page 26). The exact concentration specifications of the recommended standard eluents are listed in the chapter «Separation columns», starting on page 30.




Inline Eluent Preparation

Inline Eluent Preparation means that eluent is refilled fully automatically while the ion chromatograph is in operation. An eluent concentrate is diluted in portions with ultra pure water for the purpose of producing the required eluent.

For automatic Inline Eluent Preparation, the 940 Professional IC Vario or 930 Compact IC Flex need only be expanded to include an 941 Eluent Production Module.

Tests with repeated injections of $250 \mu g/L$ standard solutions over a time period of approximately 20 days have demonstrated outstanding stability with respect to reten-

tion times. After more than 800 sample injections, the relative standard deviations for a series of anions and cations were less than 0.55 and 0.41 percent, respectively. During a test sequence over a 24-hour period, the precision of the retention times for anions and cations were better than 0.09 and 0.08 percent, respectively. In short, this increases the reproducibility of retention times, thus permitting the exact analysis of anions and cations over extended periods, and does so without manual eluent production.



Superimposition of 200 sequential cation chromatograms (250 µg/L of the standard cation)



Separation columns



IC anion-separation columns for analyses without chemical suppression

Hamilton PRP-X100 - 125/4.0 (6.1005.000)

The Hamilton PRP-X100 - 125/4.0 IC anion column is a robust separation column based on a polystyrene/divinylbenzene copolymer. It is especially suited for the separation of chloride, nitrate, and sulfate without chemical suppression. Fluoride can also be determined if the cations are first removed with an H⁺ cartridge. The Hamilton PRP-X100 - 125/4.0 is also the separation column of choice for the determination of silicate. The column is characterized by a very good price-performance ratio.

Applications

- Cl⁻, NO₃⁻, SO₄²⁻
- Difficult matrices, e.g. dyes
- HBO₃²⁻, H₂SiO₄²⁻

Technical information

	rechnical information		
Substrate		Polystyrene/divinylbenzene	
		copolymer with quaternary	
		ammonium groups	
	Column dimensions	125 x 4.0 mm	
	Column body	Stainless steel	
	Standard flow	2.0 mL/min	
	Maximum flow	8.0 mL/min	
	Maximum pressure	34 MPa	
	Particle size	10 µm	
	Organic modifier	0-100%	
	pH range	1–13 (T > 30 °C: 1–8)	
	Capacity	62 µmol (Cl⁻)	

Eluents

Phthalic acid eluent (standard eluent)	Phthalic acid Acetone	665 mg/2 L 152 mL/2 L or 200 mL/2 L	2.0 mmol/L 7.6% or 10%
	NaOH	200 11122 1	pH = 5
Silicate eluent	Sodium hydroxide (c = 10 mol/L) Sodium carbonate	0.64 mL/2 L 106 mg/2 L	3.2 mmol/L 0.5 mmol/L


Care

Regeneration

Rinse the column with 0.5 mol/L tartaric acid or with For short periods (days) in the eluent, for longer periods 60 mmol/L HNO₃ in methanol at a flow rate of (weeks) in methanol/water (1:4) 0.5 mL/min for 2 h.

Storage

In case of contamination with iron: Rinse the column overnight with 0.1 mmol/L Na_2H_2EDTA at a flow rate of 0.5 mL/min.



Ordering information

Hamilton PRP-X100 - 125/4.0 Guard column cartridge for Hamilton PRP-X100 Guard cartridge holder, 20 mm 6.1005.000 6.1005.020 6.02821.000

Hamilton PRP-X100 - 250/4.0 (6.1005.010)

The Hamilton PRP-X100 - 250/4.0 IC anion column is a 34 robust separation column based on a polystyrene/divinylbenzene copolymer. It is mainly used with difficult matrices, e.g. dyes.

Applications

Conductivity detection

- Cl⁻, NO₃⁻, SO₄²⁻
- Difficult matrices, e.g. dyes

Technical information	
Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	250 x 4.0 mm
Column body	Stainless steel
Standard flow	2.0 mL/min
Maximum flow	8.0 mL/min
Maximum pressure	34 MPa
Particle size	10 µm
Organic modifier	0-100%
pH range	1–13 (T > 30 °C: 1–8)
Capacity	120 µmol (Cl⁻)

Eluent

Phthalic acid eluent (standard eluent)	Phthalic acid Acetone	665 mg/2 L 152 mL/2 L or 200 mL/2 L	2.0 mmol/L 7.6% or 10%
	NaOH	200 m22 2	pH = 5

Care

Regeneration

Rinse the column with 0.5 mol/L tartaric acid or with For short periods (days) in the eluent, for longer periods 60 mmol/L HNO₃ in methanol at a flow rate of (weeks) in methanol/water (1:4) 0.5 mL/min for 2 h.

In case of contamination with iron: Rinse the column overnight with 0.1 mmol/L Na₂H₂EDTA at a flow rate of 0.5 mL/min.

Storage



Ordering information

Hamilton PRP-X100 - 250/4.0 Guard column cartridge for Hamilton PRP-X100 Guard cartridge holder, 20 mm 6.1005.010 6.1005.020 6.02821.000

Super-Sep - 100/4.6 (6.1009.000)

In addition to the analysis of standard anions without chemical suppression, this column can be used for a variety of special applications. The Super-Sep - 100/4.6 IC anion column can be used for successful phosphate detection with alkaline eluent. Formate, acetate, and fluoride can be separated with suitable eluents. Overall, it is a column with very good separation performance.

Applications

- Conductivity detection
- F[−], acetate
- Difficult matrices
- Special applications, e.g. BF_4^-

Amperometric detection

ClO₂⁻, NO₂⁻, S₂O₃²⁻, SCN⁻, I⁻

Technical information

Substrate	Polymethacrylate
Column dimensions	100 x 4.6 mm
Column body	Stainless steel
Standard flow	1.5 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	2.5 MPa
Particle size	12 µm
Organic modifier	0–20%
pH range	1–13
Temperature range	20–50 °C
Capacity	32 µmol (Cl⁻)

Liuent				
Phthalic acid eluent	Phthalic acid	831 mg/2 L	2.5 mmol/L	
(standard eluent)	Acetonitrile	100 mL/2 L	5.0%	
	TRIS		pH = 4.0	

Care

Fluont


Regeneration

Rinse the column with 20% acetonitrile in 0.1 mol/L nitric acid; flow rate 0.3 mL/min for approx. 24 h. If insufficient:

- Metal contaminants: 0.1 mol/L sodium tartrate
- Protein contaminants: 0.1 mol/L sodium hydroxide or 20% acetic acid
- Organic contaminants: 20% acetonitrile in ultrapure water

Storage In the eluent





| 5

6

7

5.00

5.00

10.00

Nitrate

Sulfate

System peak

Conc. (mg/L)

10.00


_

Phthalic acid eluent, standard

1 Fluoride 5.00 Chloride 2

3 Nitrite

4 Bromide



Conc. (mg/L) Phthalic acid eluent, drinking water sample Fluoride 0.03 | 4 5.35 1 Sulfate 10.00 2 Chloride 6.43 5 System peak _ 7.83 3 Nitrate

Ordering information

Super-Sep - 100/4.6 Super-Sep Guard/4.6 (no holder required)

6.1009.000 6.1009.010



Separation columns



IC anion-separation columns for analyses with or without chemical suppression

Metrosep Anion Dual 2 - 75/4.6 (6.1006.100)

The Metrosep Anion Dual 2 - 75/4.6 column is based on an acrylate polymer. Without chemical suppression, it can be used with conventional phthalic acid eluents and a range of other eluent systems. The benefits of working without chemical suppression (simple analytical method and linear calibration) can thus be taken advantage of.

When used with chemical suppression, the column excels in its outstanding separating efficiency, particularly with respect to early eluting ions (fluoride, acetate, formate). The separation between fluoride and the «water dip» as well as between chloride and nitrite is remarkable.

In order to prolong the lifetime of the column even further, we recommend the use of the Metrosep RP 2 Guard/3.5 (6.1011.030) or the Metrosep RP 3 Guard HC/4.0 (6.1011.040).

Applications

Conductivity detection without chemical suppression

- Difficult separation problems
- Difficult matrices
- F
- Conductivity detection with chemical suppression
- Standard anions
- F⁻, formate, acetate, Cl⁻, SO_3^{2-} , SO_4^{2-}
- Simple matrices
- ClO₂⁻, ClO₃⁻, ClO₄⁻, BrO₃⁻

• Special separation problems Amperometric detection

• Br⁻, I⁻

Technical information

Polymethacrylate with qua-	
ternary ammonium groups	
75 x 4.6 mm	
Stainless steel	
0.8 mL/min	
1.2 mL/min	
7 MPa	
6 µm	
0–20%	
1-12	
17 µmol (Cl⁻)	

Eluents

Without chemical suppression Phthalic acid eluent Phthalic acid 1660 mg/2 L 5.0 mmol/L (standard eluent) Acetonitrile 40 mL/2 L 2% NaOH pH = 4.5With chemical suppression Carbonate eluent Sodium hydrogen carbonate 336 mg/2 L 2.0 mmol/L (standard eluent) Sodium carbonate 276 mg/2 L 1.3 mmol/L

Care

Regeneration

Rinse with 0.1 mol/L HNO_3 at 0.3 mL/min for approx. 2 h, then switch to standard eluent.


Rinse with an eluent containing 100 mmol/L of the respective buffer salt.

Organic contaminants:

Rinse the column with eluent containing organic solvent (max. 20%).

Transition metals:

When using carbonate eluent: Add 0.1 mmol/L dipicolinic acid to the standard carbonate eluent; rinse in the opposite flow direction with this eluent for 3 hours at 0.5 mL/min; then rinse the column again under standard conditions in the flow direction for at least 2 hours.



Metrosep Anion Dual 2 - 75/4.6 Metrosep RP 2 Guard/3.5 Replacement filters for RP 2 Guard/3.5 (10 pcs.) Metrosep RP 3 Guard HC/4.0 6.1006.100 6.1011.030 6.1011.130 6.1011.040

Metrosep Anion Dual 3 - 100/4.0 (6.1006.120)

The Metrosep Anion Dual 3 - 100/4.0 made of PEEK is particularly suitable for biological samples. It can be used to solve separation problems both with and without chemical suppression. Its performance profile is identical to that of the Metrosep Anion Dual 2 - 75/4.6. This means that early eluting ions are separated very well and that oxohalides can be determined very efficiently. The Metrosep Anion Dual 3 - 100/4.0 can also be used for samples with high chloride but low nitrite contents.

When working with amperometric detection, the Metrosep Anion Dual 3 - 100/4.0 can be used for the determination of bromide and iodide.

In order to prolong the lifetime of the column, we recommend the use of the Metrosep RP 2 Guard/3.5 (6.1011.030) or the Metrosep RP 3 Guard HC/4.0 (6.1011.040).

Applications

Conductivity detection without chemical suppression

- Difficult separation problems
- Difficult matrices
- F⁻
- Conductivity detection with chemical suppression
- Standard anions
- F⁻, formate, acetate ,Cl⁻, SO_3^{2-} , SO_4^{2-}
- Simple matrices
- CIO₂⁻, CIO₃⁻, CIO₄⁻, BrO₃⁻

• Special separation problems Amperometric detection

• Brī, lī

Technical information

Substrate	Polymethacrylate with qua-
	ternary ammonium groups
Column dimensions	100 x 4.0 mm
Column body	PEEK
Standard flow	0.8 mL/min
Maximum flow	1.2 mL/min
Maximum pressure	7 MPa
Particle size	6 µm
Organic modifier	0–20%
pH range	1-12
Capacity	14 µmol (Cl⁻)

Eluents

Without chemical suppression

without chemical suppressio	n		
Phthalic acid eluent	Phthalic acid	1660 mg/2 L	5.0 mmol/L
(standard eluent)	Acetonitrile	40 mL/2 L	2%
	NaOH		pH = 4.5
With chemical suppression			
Carbonate eluent	Sodium hydrogen carbonate	336 mg/2 L	2.0 mmol/L
(standard eluent)	Sodium carbonate	276 mg/2 L	1.3 mmol/L

Care

Regeneration

Rinse with 0.1 mol/L HNO_3 at 0.3 mL/min for approx. 2 h, then switch to standard eluent.

Rinse with an eluent containing 100 mmol/L of the respective buffer salt.


Organic contaminants:

Rinse the column with eluent containing organic solvent (max. 20%).

Transition metals:

When using carbonate eluent: Add 0.1 mmol/L dipicolinic acid to the standard carbonate eluent; rinse in the opposite flow direction with this eluent for 3 hours at 0.5 mL/min; then rinse the column again under standard conditions in the flow direction for at least 2 hours.





Ordering information	
Metrosep Anion Dual 3 - 100/4.0	6.1006.120
Metrosep RP 2 Guard/3.5	6.1011.030
Replacement filters for RP 2 Guard/3.5 (10 pcs.)	6.1011.130
Metrosep RP 3 Guard HC/4.0	6.1011.040

43

Metrosep Dual 4 - 100/4.6 (6.1016.030)

The Metrosep Dual 4 separation columns are based on a functionalized Monolith based on silica gel. The Monolith permits an eluent flow of up to 5 mL/min. Despite the high flow, the column is characterized by low backpressure. In contrast to traditional materials, the Monolith has a much larger surface due to its structure of macropores and mesopores. This contributes to the high column capacity with simultaneously very low dead volume.

The Metrosep Dual 4 - 100/4.6 is suitable for a great number of applications. All standard anions can thus be separated in less than nine minutes. The high column capacity makes it largely insensitive to matrix influences. Even in a matrix of 3 g/L chloride, carbonate and sulfate, 0.5 μ g/L perchlorate can be detected. The column is therefore used in perchlorate analysis in accordance with EPA standard 314.

The column can be used with or without chemical suppression. When p-cyanophenol is used as the eluent, it is recommended that the Metrosep RP Trap 1 - 50/4.0 (6.1014.100) be installed between pulsation absorber and injection valve.

Applications

- Rapid separations
- Complex sample matrices
- EPA 314
- Detection of perchlorate

Technical information


Substrate	Monolithic silica gel
Column dimensions	100 x 4.6 mm
Column body	PEEK
Standard flow	2.0 mL/min
Maximum flow	5.0 mL/min
Maximum pressure	20 MPa
Particle size	Monolith with 2 μm macro-
	pores and 13 nm mesopores
Organic modifier	0–5% (methanol or
	acetonitrile only)
pH range	2–8
Capacity	74 µmol (Cl⁻)

Eluent			
p-cyanophenol eluent	p-cyanophenol	2859 mg/2 L	12.0 mmol/L
(standard eluent)	КОН		$pH = 7.4 \pm 0.1$

Care

Preparation Rinse the column with eluent for 0.5–1 h. Storage In the eluent

Regeneration Rinse the column with a maximum of 5% acetonitrile at a flow rate of 0.5 mL/min for 30 min.



Ordering information	
Metrosep Dual 4 - 100/4.6	6.1016.030
Guard column kit for the Metrosep Dual 4, comprised of three guard column cartridges and	
one guard column cartridge holder	6.1016.500
Guard column cartridges for the Metrosep Dual 4 (3 pcs.)	6.1016.510

Metrosep A Supp 4 - 250/4.0 (6.1006.430)

The Metrosep A Supp 4 - 250/4.0 is an extremely robust column with very good separation properties. The separation phase is comprised of polyvinyl alcohol particles with quaternary ammonium groups and a diameter of 9 μ m. This structure guarantees great stability and a greater tolerance to very small particles which could pass through the integrated filter plate. The Metrosep A Supp 4 - 250/4.0 has a medium ion-exchange capacity; sulfate elutes after 12.5 minutes. The number of plates which can be achieved with this separation column is higher than those on the Metrosep Anion Dual 2 - 75/4.6. The A Supp 4 - 250/4.0 is particularly suitable for all routine tasks in water analysis.

To protect the IC separation column – even though it is not particularly sensitive to contaminants – we recommend the use of the Metrosep A Supp 4 Guard/4.0 or the A Supp 4 S-Guard/4.0.

Applications

- Standard anions
- Water analysis
- Difficult matrices
- Critical samples
- Iodide

Technical information

Substrate	Polyvinyl alcohol with
	quaternary
	ammonium groups
Column dimensions	250 x 4.0 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	12 MPa
Particle size	9 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Temperature range	20–60 °C
Capacity	37 µmol (Cl⁻)

Eluent

Without chemical suppression				
Phthalic acid eluent	Phthalic acid	1660 mg/2 L	5.0 mmol/L	
(standard eluent)	Acetone	40 mL/2 L	2.0%	
	TRIS		pH = 4.4	
With chemical suppression				
Carbonate eluent	Sodium hydrogen carbonate	286 mg/2 L	1.7 mmol/L	
(standard eluent)	Sodium carbonate	382 mg/2 L	1.8 mmol/L	
Carbonate eluent, mod.	Sodium hydrogen carbonate	672 mg/2 L	4.0 mmol/L	
	Sodium carbonate	212 mg/2 L	1.0 mmol/L	

Care

Regeneration

Contamination with hydrophilic ions:

a) Rinse with ultrapure water (15 min at 0.5 mL/min)

b) Rinse with 10x concentrated eluent


(60 min at 0.5 mL/min)

c) Rinse with ultrapure water (15 min at 0.5 mL/min)


d) Rinse with eluent (60 min at 0.5 mL/min)

Contamination with lipophilic ions:

a) Rinse with ultrapure water (15 min at 0.5 mL/min)
b) Rinse with 5% acetonitrile (10 min at 0.5 mL/min)
c) Rinse with 100% acetonitrile (60 min at 0.5 mL/min)
d) Rinse with 50% acetonitrile (10 min at 0.5 mL/min)
e) Rinse with ultrapure water (30 min at 0.5 mL/min)
f) Rinse with eluent (60 min at 0.5 mL/min)



Ordering information	
Metrosep A Supp 4 - 250/4.0	6.1006.430
Metrosep A Supp 4 Guard/4.0	6.01021.500
Metrosep A Supp 4 S-Guard/4.0	6.01021.510



Separation columns



IC anion-separation columns for analyses with chemical suppression

Phenomenex Star-Ion[™] A 300 - 100/4.6 (6.1005.100)

The Star-Ion[™] A 300 - 100/4.6 IC anion column's separation properties are based on the use of styrene-divinylbenzene resin. The separation between the «water dip» and chloride is not as sharp as with acrylate-based and vinyl alcohol-based columns. The Star-Ion[™] A 300 -100/4.6 features lower separation efficiency in the fluoride range. This column nevertheless excels in its extremely short analysis times.

Standard anions can be determined in less than 8 minutes. Moreover, this column can be used to determine fluoride in the ppm range in the presence of very large amounts of lactate (Application Note S-37).

The system peak and high carbonate concentrations interfere with the integration of the chloride peak. It is therefore recommended that the «MCS» (Metrohm $\rm CO_2$ Suppressor) be used.

Applications

• Rapid separation of standard anions

- Fluoride in addition to an excess of lactate
- Chromate

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	100 x 4.6 mm
Column body	PEEK
Standard flow	1.5 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	7 MPa
Particle size	7 µm
Organic modifier	0%
pH range	1–12
Capacity	21 µmol (Cl⁻)

Eluents

Carbonate eluent	Sodium hydrogen carbonate	286 mg/2 L	1.7 mmol/L
(standard eluent)	Sodium carbonate	382 mg/2 L	1.8 mmol/L
Carbonate eluent	Sodium hydrogen carbonate	841 mg/2 L	5.0 mmol/L
(modified)	Sodium carbonate	743 mg/2 L	3.5 mmol/L


Care

RegenerationStoraRinse for 30 min at 1.0 mL/min with a solution con-
taining 18 mmol/L Na_2CO_3 (1908 mg) and 17 mmol/LIn the
NaHCO_3 (1428 mg).

Storage

In the eluent





| 5

6

7

0.50

0.50

1.00

Carbonate eluent, standard


- 1 Fluoride 0.20
- 2 Chloride 3 Nitrite

4 Bromide

Conc. (mg/L) Nitrate Phosphate Sulfate

1.00 1.00

1.00



Carbonate eluent, modified, dye, Dil. 1:100 Chromate (VI) 3.24 1

Conc. (mg/L)

Ordering information

Phenomenex Star-Ion™ A 300 - 100/4.6 Metrosep RP 2 Guard/3.5 Replacement filters for RP 2 Guard/3.5 (10 pcs.) Metrosep RP 3 Guard HC/4.0

6.1005.100 6.1011.030 6.1011.130 6.1011.040

Metrosep A Supp 1 - 250/4.6 (6.1005.300)

The Metrosep A Supp 1 - 250/4.6 is a universal anion column which is characterized by medium capacity and special selectivity. With this column it is possible to process samples with great differences in concentration. For example, 4 µg/L of sulfate can be determined in a solution containing 150 g/L sodium chloride. An additional advantage is that bromide elutes after nitrate. Particularly in the area of oxohalide analysis, the A Supp 1 - 250/4.6 excels in its outstanding separation properties. Pressure fluctuations, constantly changing eluents, and large sample through-put do not influence the separating efficiency of this column, even after very long periods. It is the «workhorse» for development and routine laboratories.

Applications

- Conductivity detection
- Great differences in concentration
- High ionic strength
- ClO₂⁻, ClO₃⁻, ClO₄⁻
- Amperometric detection
- CN⁻
 S²⁻
- 5

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	250 x 4.6 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.5 mL/min
Maximum pressure	15 MPa
Particle size	7 µm
Organic modifier	0–100%
pH range	1–13
Capacity	12 µmol (Cl⁻)

Eluents				
Carbonate eluent	Sodium carbonate	636 mg/2 L	3.0 mmol/L	
(standard eluent)				
Sodium hydroxide eluent	Sodium hydroxide (c = 10 mol/L)	20 mL/2 L	100 mmol/L	

Care

Regeneration


Rinse with 50 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.5 mL/min. Then rinse with 0.1 mol/L NaOH at 0.5 mL/min for 1 h.

Organic contaminants: Rinse with 70% methanol at 1.0 mL/min for 12 h. The addition of 1% acetic acid may be useful.

Storage

In the eluent. For a longer period (weeks), store the column in a refrigerator at minimum +4 $^{\circ}$ C.





Carbonate eluent, standard

Fluoride
 Chlorite


3

4

5

Bromate Chloride Nitrite

Conc. (mg/L) 10.00 2.00 6 Chlorate 5.00 10.00 7 Nitrate 5.00 Bromide 10.00 8 5.00 9 Phosphate 10.00 5.00 10.00 10 Sulfate



Sodium hydroxide eluent, standard, amperometric detection 1 Cyanide 4.0

Conc. (µg/L)

Ordering information

Metrosep A Supp 1 - 250/4.6 Metrosep A Supp 1 Guard/4.6 6.1005.300 6.1005.340

Metrosep A Supp 1 HS - 50/4.6 (6.1005.350)

The Metrosep A Supp 1 HS - 50/4.6 permits the separation of standard anions in a very short time. The Metrosep A Supp 1 HS - 50/4.6 is also the column of choice for the determination of only a few anions in an uncomplicated sample matrix. For example, the analysis of phosphate as well as chloride and sulfate in cola beverages can be carried out in less than three minutes.

Applications

- Cl^{-} , PO_{a}^{3-} , SO_{a}^{2-} in cola beverages
- Very rapid separation
- Standard anions in uncomplicated sample matrices

Technical information	
Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	50 x 4.6 mm
Column body	PEEK
Standard flow	1.3 mL/min
Maximum flow	2.5 mL/min
Maximum pressure	4.0 MPa
Particle size	7 µm
Organic modifier	0-100%
pH range	1–13
Capacity	3.1 µmol (Cl⁻)

Eluents			
Carbonate eluent	Sodium carbonate	636 mg/2 L	3.0 mmol/L
(standard eluent)			

Care

Regeneration

Rinse with 15 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.25 mL/min. Then rinse with 0.1 mol/L NaOH at 0.25 mL/min for 1 h.

Organic contaminants: Rinse with 70% methanol at 0.4 mL/min for 12 h. The addition of 1% acetic acid may be useful.


Storage

In the eluent. For a longer period (weeks), store the column in a refrigerator at minimum +4 $^{\circ}\mathrm{C}.$




2

Nitrate



20.0 4

Sulfate



Larbo	onate eluent, cola pe	verage		Con	c. (mg/L)
1	Chloride	5.0	3	Phosphate	496.3
2	Nitrate	8.2	4	Sulfate	10.4

Ordering information	
Metrosep A Supp 1 HS - 50/4.6	6.1005.350
Metrosep RP 2 Guard/3.5	6.1011.030
Replacement filters for RP 2 Guard/3.5 (10 pcs.)	6.1011.130
Metrosep RP 3 Guard HC/4.0	6.1011.040

55

Metrosep A Supp 3 - 250/4.6 (6.1005.320)

The Metrosep A Supp 3 - 250/4.6 solves separation problems in aqueous and organic media. It can be used reliably with a wide range of eluents – even those with high proportions of organic solvents. With the Metrosep A Supp 3 - 250/4.6, highly demanding samples can be analyzed in routine operation, for example the measurement of biological samples or the determination of inorganic anions in organic matrices. In isocratic operation, the column is also suitable for the separation of sulfite, sulfate, and thiosulfate in less than 20 minutes.

Applications

- SCN⁻, SO₃²⁻, SO₄²⁻, S₂O₃²⁻
- Polyphosphates
- Organic matrices

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	250 x 4.6 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	1.5 mL/min
Maximum pressure	15 MPa
Particle size	9 µm
Organic modifier	0-100%
pH range	1–13
Capacity	35 µmol (Cl⁻)

EI	uont
E	uent

Carbonate eluent	Sodium hydrogen carbonate	286 mg/2 L	1.7 mmol/L
(standard eluent)	Sodium carbonate	382 mg/2 L	1.8 mmol/L

Care

Regeneration

Rinse with 50 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.5 mL/min. Then rinse with 0.1 mol/L NaOH at 0.5 mL/min for 1 h.


Organic contaminants:

Rinse with 70% methanol at 1.0 mL/min for 12 h. The addition of 1% acetic acid may be useful.

Storage

In the eluent. For a longer period (weeks), store the column in a refrigerator at minimum +4 $^{\circ}$ C.





_

2.00

5.00

5.00

µS/cm 26 24 22 20 18 2 1 16 ż 3 4 ÷ 8 ó 10 11 12 13 14 15 16 5 Carbo er)

4

5

6

Carbonate eluent, standard

1 Injection peak 2 Fluoride Chloride

Nitrite

3

4

5 Bromide 6 Nitrate 7 Phosphate 8 Sulfate

10.00 10.00 10.00 10.00

Conc. (mg/L)


1

2

3

onate eluent, drinking	wate
Injection peak	-
Fluoride	0.04
Chloride	5.24

Conc.	(mg/L)
System peak	-
Nitrate	10.25
Sulfate	6.92



Ordering information Metrosep A Supp 3 - 250/4.6 Metrosep RP 2 Guard/3.5 Replacement filters for RP 2 Guard/3.5 (10 pcs.) Metrosep RP 3 Guard HC/4.0

6.1005.320 6.1011.030 6.1011.130 6.1011.040

Metrosep A Supp 5 - 50/4.0 (6.1006.550)

The Metrosep A Supp 5 - 50/4.0 separates the seven standard anions in less than six minutes. Even fluoride is still separated from the injection peak and can be integrated perfectly. Like all columns in the A-Supp-5 product range, the column, which is based on a polyvinyl alcohol polymer, is characterized by high plate numbers and therefore by outstanding separating efficiency. The Metrosep A Supp 5 - 50/4.0 is the column of choice when simple separation tasks must be solved in a short time – and that without having to sacrifice very low detection limits.

Applications

- Rapid separation of standard anions
- Simple sample matrices
- Method development

Technical information

Substrate	Polyvinyl alcohol with qua-
	ternary ammonium groups
Column dimensions	50 x 4.0 mm
Column body	PEEK
Standard flow	0.7 mL/min
Maximum flow	0.8 mL/min
Maximum pressure	15 MPa
Particle size	5 µm
Organic modifier	0–100%, (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Temperature range	20-60 °C
Capacity	18 µmol (Cl⁻)

Deliver district a landle a localitation

Eluent

Carbonate eluent	Sodium hydrogen carbonate	168 mg/2 L	1.0 mmol/L
(standard eluent)	Sodium carbonate	678 mg/2 L	3.2 mmol/L
Perchlorate eluent	Sodium carbonate	3.178 g/2 L	15.0 mmol/L
	Acetone	200 mL/2 L	10%

Care

Regeneration

Contamination with hydrophilic ions:

- a) Rinse with ultrapure water (25 min at 0.3 mL/min)
- b) Rinse with 10x concentrated eluent (100 min at 0.3 mL/min)
- c) Rinse with ultrapure water (25 min at 0.3 mL/min) d) Rinse with eluent (100 min at 0.3 mL/min)

Contamination with lipophilic ions:

- a) Rinse with ultrapure water (25 min at 0.3 mL/min)
 b) Rinse with 5% acetonitrile (20 min at 0.3 mL/min)
 c) Rinse with 100% acetonitrile (60 min at 0.3 mL/min)
 d) Rinse with 50% acetonitrile (10 min at 0.3 mL/min)
- e) Rinse with ultrapure water (50 min at 0.3 mL/min)
- f) Rinse with eluent (100 min at 0.3 mL/min)


With shifted system peak (regeneration method with column oven):


- a) Rinse with concentrated eluent 1 mol/L $\rm Na_2CO_3$ (25 min at 0.4 mL/min)
- b) Maintain for 10-12 h at 45-50 °C (without rinsing)
- c) Rinse with standard eluent (at least 40 min at 0.4 mL/min)

Storage

In the eluent







| 5

6

7

5.00

5.00

10.00


Nitrate

Sulfate

Phosphate

Carbonate eluent, standard 2.00

- Fluoride 1
- 2 Chloride 3 Nitrite
- 4 Bromide



Carbonate eluent, drinking water					Conc. (mg/L)
1	Fluoride	0.04	4	Nitrate	10.36
2	Chloride	5.25	5	Sulfate	6.90
3	System peak	-			

Ordering information	
Metrosep A Supp 5 - 50/4.0	6.1006.550
Metrosep A Supp 5 Guard/4.0	6.1006.500
Metrosep A Supp 5 S-Guard/4.0	6.1006.540

10.00

10.00

10.00

Metrosep A Supp 5 - 100/4.0 (6.1006.510)

The Metrosep A Supp 5 - 100/4.0 allows highly efficient, rapid separations. This property makes the Metrosep A Supp 5 - 100/4.0 the standard column for short analysis times and the determination of late eluting anions (e.g. perchlorate).

Applications

- Conductivity detection
- Standard anions
- F⁻, Cl⁻, Br⁻, l⁻, ClO₂⁻, ClO₃⁻, ClO₄⁻, BrO₃⁻
- ClO₄⁻
- Cr(VI) (CrO₄²⁻), I⁻
- Method development
- Universal applications
- Determination of phosphate in cola beverages
- Rapid separation
- Amperometric detection
- NO₂⁻, ClO₂⁻
- S₂O₃²⁻, SCN⁻, I⁻

Technical information		Maximum pressure	15 MPa
Substrate	Polyvinyl alcohol with qua-	Particle size	5 µm
	ternary ammonium groups	Organic modifier	0–100%, (particularly
Column dimensions	100 x 4.0 mm		acetone, acetonitrile,
Column body	PEEK		methanol)
Standard flow	0.7 mL/min	pH range	3–12
Maximum flow	0.8 mL/min	Temperature range	20-60 °C
		Capacity	32 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium hydrogen carbonate	168 mg/2 L	1.0 mmol/L
(standard eluent)	Sodium carbonate	678 mg/2 L	3.2 mmol/L
Cola eluent	Sodium hydrogen carbonate	504 mg/2 L	3.0 mmol/L
	Sodium carbonate	1484 mg/2 L	7.0 mmol/L

Care

Regeneration

Contamination with hydrophilic ions:

a) Rinse with ultrapure water (25 min at 0.3 mL/min)

- b) Rinse with 10x concentrated eluent (100 min at 0.3 mL/min)
- c) Rinse with ultrapure water (25 min at 0.3 mL/min) d) Rinse with eluent (100 min at 0.3 mL/min)

Contamination with lipophilic ions:

a) Rinse with ultrapure water (25 min at 0.3 mL/min)
b) Rinse with 5% acetonitrile (20 min at 0.3 mL/min)
c) Rinse with 100% acetonitrile (60 min at 0.3 mL/min)
d) Rinse with 50% acetonitrile (10 min at 0.3 mL/min)

e) Rinse with ultrapure water (50 min at 0.3 mL/min) f) Rinse with eluent (100 min at 0.3 mL/min)


With shifted system peak (regeneration method with column oven):

- a) Rinse with concentrated eluent 1 mol/L $\rm Na_2CO_3$ (25 min at 0.4 mL/min)
- b) Maintain for 10-12 h at 45-50 °C (without rinsing)
- c) Rinse with standard eluent (at least 40 min at 0.4 mL/min)




4


Bromide



Cola eluent, cola beverage Conc. (mg/L) 1 Phosphate 532.53 | 2 Sulfate 36.63




10.00



RSD retention time < 0.1% RSD concentration < 0.2%

Number of analyses n = 400



Ordering information	
Metrosep A Supp 5 - 100/4.0	6.1006.510
Metrosep A Supp 5 Guard/4.0	6.1006.500
Metrosep A Supp 5 S-Guard/4.0	6.1006.540
Metrosep A Supp 10 Guard HC/4.0	6.1020.520

Metrosep A Supp 5 - 150/4.0 (6.1006.520)

The 150 mm version of the Metrosep A Supp 5 is characterized by its very good separation properties. High plate numbers and excellent peak symmetries simplify working in the lower μ g/L range. The particle size of 5 μ m makes a decisive contribution to the separating efficiency of this column. The Metrosep A Supp 5 -150/4.0 offers the optimum combination of selectivity and capacity, with which even complex separation tasks can be solved within a short time. This characteristic makes the Metrosep A Supp 5 - 150/4.0 one of the best universally applicable standard IC columns.

Applications

- Standard anions
- F⁻, Cl⁻, Br⁻, l⁻
- ClO₂⁻, ClO₃⁻, ClO₄⁻, BrO₃⁻
- Cr(VI) (CrO4²⁻)
- Method development
- Difficult matrices
- Difficult separation problems

Technical information

Substrate	Polyvinyl alcohol with qua-
	ternary ammonium groups
Column dimensions	150 x 4.0 mm
Column body	PEEK
Standard flow	0.7 mL/min
Maximum flow	0.8 mL/min
Maximum pressure	15 MPa
Particle size	5 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Temperature range	20-60 °C
Capacity	52 µmol (Cl⁻)

Eluents

Carbonate eluent	Sodium hydrogen carbonate	168 mg/2 L	1.0 mmol/L
(standard eluent)	Sodium carbonate	678 mg/2 L	3.2 mmol/L
Chromate eluent	Sodium hydrogen carbonate	672 mg/2 L	4.0 mmol/L
	Sodium carbonate	2714 mg/2 L	12.8 mmol/L

Care

Regeneration

Contamination with hydrophilic ions:

a) Rinse with ultrapure water (25 min at 0.3 mL/min)

- b) Rinse with 10x concentrated eluent (100 min at 0.3 mL/min)
- c) Rinse with ultrapure water (25 min at 0.3 mL/min) d) Rinse with eluent (100 min at 0.3 mL/min)


Contamination with lipophilic ions:

a) Rinse with ultrapure water (25 min at 0.3 mL/min)
b) Rinse with 5% acetonitrile (20 min at 0.3 mL/min)
c) Rinse with 100% acetonitrile (60 min at 0.3 mL/min)
d) Rinse with 50% acetonitrile (10 min at 0.3 mL/min)


e) Rinse with ultrapure water (50 min at 0.3 mL/min) f) Rinse with eluent (100 min at 0.3 mL/min)


- With shifted system peak (regeneration method with column oven):
- a) Rinse with concentrated eluent 1 mol/L Na_2CO_3 (25 min at 0.4 mL/min)
- b) Maintain for 10-12 h at 45-50 °C (without rinsing)
- c) Rinse with standard eluent (at least 40 min at 0.4 mL/min)





Carbo	nate eluent, standa	Conc.	(mg/L)		
1	Fluoride	2.00	5	Nitrate	10.00
2	Chloride	5.00	6	Phosphate	10.00
3	Nitrite	5.00	7	Sulfate	10.00
4	Bromide	10.00			





Conc. (µg/L)

Ordering information	
Metrosep A Supp 5 - 150/4.0	6.1006.520
Metrosep A Supp 5 Guard/4.0	6.1006.500
Metrosep A Supp 5 S-Guard/4.0	6.1006.540

Metrosep A Supp 5 - 250/4.0 (6.1006.530)

The high-performance separation column from Metrohm with an extremely high number of plates for the most demanding separation tasks. Even complex separation problems can be solved easily and reproducibly with the Metrosep A Supp 5 - 250/4.0. The high capacity of the column allows, for example, the detection of 1 μ g/L bromate along with 150 mg/L chloride without sample preparation. The range of applications possible with this column far exceeds the detection of standard anions. The Metrosep A Supp 5 - 250/4.0 is the column of choice when it comes to reliable monitoring of the high purity standards in the semiconductor industry or of the boiler feed water of power plants.

Applications

- Standard anions
- F[−], Cl[−], Br[−], l[−]
- CIO₂⁻, CIO₃⁻, CIO₄⁻, BrO₃⁻
- ClO_4^- at high ionic strength
- BrO₃⁻ at high ionic strength
- Method development
- Universal applications
- Difficult matrices
- Difficult separation problems
- Applications with gradient

feed water of power plants		Technical information	
		Substrate	Polyvinyl alcohol with qua-
			ternary ammonium groups
		Column dimensions	250 x 4.0 mm
		Column body	PEEK
		Standard flow	0.7 mL/min
		Maximum flow	0.8 mL/min
		Maximum pressure	15 MPa
		Particle size	5 µm
		Organic modifier	0–100% (particularly
			acetone, acetonitrile,
			methanol)
		pH range	3–12
		Temperature range	20–60 °C
		Capacity	85 µmol (Cl⁻)
Eluent			
Carbonate eluent	Sodium hydrogen carbonate	168 mg/2 L	1.0 mmol/L
(standard eluent)	Sodium carbonate	678 mg/2 L	3.2 mmol/L

Care

Regeneration

Contamination with hydrophilic ions:

a) Rinse with ultrapure water (25 min at 0.3 mL/min)

b) Rinse with 10x concentrated eluent (100 min at 0.3 mL/min)

c) Rinse with ultrapure water (25 min at 0.3 mL/min) d) Rinse with eluent (100 min at 0.3 mL/min)

Contamination with lipophilic ions:


a) Rinse with ultrapure water (25 min at 0.3 mL/min)
b) Rinse with 5% acetonitrile (20 min at 0.3 mL/min)
c) Rinse with 100% acetonitrile (60 min at 0.3 mL/min)
d) Rinse with 50% acetonitrile (10 min at 0.3 mL/min)

e) Rinse with ultrapure water (50 min at 0.3 mL/min) f) Rinse with eluent (100 min at 0.3 mL/min)

With shifted system peak (regeneration method with column oven):

- a) Rinse with concentrated eluent 1 mol/L $\rm Na_2CO_3$ (25 min at 0.4 mL/min)
- b) Maintain for 10-12 h at 45-50 °C (without rinsing)


c) Rinse with standard eluent (at least 40 min at 0.4 mL/min)




Training and

Curr	vonace cracine, i	not and last q	C 510	an loan o,	conc. (mg/L)
215	0 injections				
1	Fluoride	1.00	5	Nitrate	10.00
2	Chlorido	2 00	G	Dhocphata	10.00

2	Chloride	2.00	6	Phosphate	10.00
3	Nitrite	5.00	7	Sulfate	10.00
4	Bromide	10.00			



Carbonate eluent, fluoride in HCl (32%) dil. 1:500 Conc. (µg/L) 1 Fluoride 20.0



Carbonate eluent, VIS detection (λ = 450 nm), EPA 317.0 Conc. (µg/L) 1 Bromate 10.0

5 10 11 12

1.95

Ordering information Metrosep A Supp 5 - 250/4.0

Metrosep A Supp 5 Guard/4.0 Metrosep A Supp 5 S-Guard/4.0 6.1006.530 6.1006.500 6.1006.540

min

Metrosep A Supp 7 - 150/4.0 (6.1006.620)

The Metrosep A Supp 7 - 150/4.0 is the shorter A Supp 7 column. It allows similarly complex separation tasks to be solved the same way as with the 250 mm version, with no significant loss in separating efficiency. Chlorite and bromate can thus be easily separated from standard anions with this separation column. With the Metrosep A Supp 7 - 150/4.0, these ions are determined with certainty and precision down to the lower μ g/L range. The high detection sensitivity is achieved through the use of the 5 μ m polyvinyl alcohol polymer, with which extremely high plate numbers and thus outstanding separation and detection properties are achieved. In addition, the separation can be adapted to the specific requirements of the application by modifying the temperature.

Applications

- Standard anions
- Determination of standard anions and ClO $_{\rm 2}^-$, ClO $_{\rm 3}^-$, BrO $_{\rm 3}^-$
- Complex separation tasks
- Applications with gradient

Technical information

:	Substrate	Polyvinyl alcohol with qua- ternary ammonium groups
I	Column dimensions	150 x 4.0 mm
•	Column body	PEEK
;	Standard flow	0.7 mL/min
	Maximum flow	1.0 mL/min
	Maximum pressure	15 MPa
	Particle size	5 µm
	Organic modifier	0–100% (particularly
		acetone, acetonitrile,
		methanol)
	pH range	3–12
	Temperature range	20–60 °C
	Capacity	76 µmol (Cl⁻)
	763 ma/2 L	3.6 mmol/L

Eluent			
Carbonate eluent	Sodium carbonate	763 mg/2 L	3.6 mmol/L
(standard eluent)		Column temperature 45 °C	

Care

Regeneration

Contamination with low-valency hydrophilic ions:

a) Rinse with ultrapure water (25 min at 0.3 mL/min)

b) Rinse with 10x concentrated eluent (100 min at 0.3 mL/min)


c) Rinse with ultrapure water (25 min at 0.3 mL/min) d) Rinse with eluent (100 min at 0.3 mL/min) Contaminations with high-valency hydrophobic ions and organic contaminations:

- a) Rinse with ultrapure water (25 min at 0.3 mL/min) b) Rinse with 100% acetonitrile (20 min at 0.3 mL/min)
- c) Rinse with ultrapure water (25 min at 0.3 mL/min)
- d) Rinse with 10x concentrated eluent (100 min at 0.3 mL/min)
- e) Rinse with ultrapure water (25 min at 0.3 mL/min) f) Rinse with eluent (100 min at 0.3 mL/min)

Storage In the eluent at max. 8 °C



Chromatograms



Carbonate eluent, standard, 45 °C,

Fluoride
 Chloride
 Nitrite

Bromide

4


2.00 5 5.00 6 5.00 7

10.00

Sulfate

Nitrate

Phosphate



Gradient: Carbonate eluent 1–6 mmol/L,

stand	ard, 45 °C			Conc. (mg/L)
1	Fluoride	5.00	11	Nitrate	5.00
2	Glycolate	5.00	12	Dibromoacetate	5.00
3	Propionate	5.00	13	Phosphate	5.00
4	Butyrate	5.00	14	Sulfate	5.00
5	Methacrylate	5.00	15	Tartrate	5.00
6	Monochloroacetate	5.00	16	Selenate	5.00
7	Chloride	5.00	17	Arsenate	5.00
8	Nitrite	5.00	18	lodide	5.00
9	Bromide	5.00	19	Thiosulfate	5.00
10	Dichloroacetate	5.00			

Ordering inf	ormation
--------------	----------

Metrosep A Supp 7 - 150/4.0	6.1006.620
Metrosep A Supp 5 Guard/4.0	6.1006.500
Metrosep A Supp 5 S-Guard/4.0	6.1006.540
Metrosep A Supp 16 Guard/4.0	6.1031.500
Metrosep A Supp 16 S-Guard/4.0	6.1031.510
Metrosep RP 2 Guard/3.5	6.1011.030
Replacement filters for RP 2 Guard/3.5 (10 pcs.)	6.1011.130

Conc. (mg/L)

10.00

10.00

10.00

Metrosep A Supp 7 - 250/4.0 (6.1006.630)

Disinfection by products from water purification are suspected carcinogens. Oxohalides have therefore become the subject of many investigations and standards (e.g. EPA 300.1 Part B, EPA 317.0, EPA 326, DIN EN ISO 11206). Of primary concern is bromate, which forms from bromide during the ozonization of drinking water. The Metrosep A Supp 7 - 250/4.0 is a high-performance separation column for the parallel determination of standard anions, oxohalides, and dichloroacetic acid. With this column, these ions are determined with certainty and precision down to the lower μ g/L range. The high detection sensitivity is achieved through the use of the 5 µm polyvinyl alcohol polymer, with which extremely high plate numbers and thus outstanding separation and detection properties are achieved. In addition, the separation can be adapted to the specific requirements of the application by modifying the temperature.

Applications

- Standard anions
- EPA Method 300.1 Part B, simultaneous determination of standard anions and CIO₂⁻, CIO₃⁻, BrO₃⁻ and DCA (dichloroacetic acid)

- Isocratic separation of glycolate, acetate, and formate
- Complex separation tasks
- Applications with gradient

Technical information

Substrate	Polyvinyl alcohol with qua-
	ternary ammonium groups
Column dimensions	250 x 4.0 mm
Column body	PEEK
Standard flow	0.7 mL/min
Maximum flow	1.0 mL/min
Maximum pressure	15 MPa
Particle size	5 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Temperature range	20–60 °C
Capacity	110 µmol (Cl⁻)

Eluent				
Carbonate eluent	Sodium carbonate	763 mg/2 L	3.6 mmol/L	
(standard eluent)		Column temperature 4	15 °C	
Carbonate eluent	Sodium carbonate	763 mg/2 L	3.6 mmol/L	
(modified)	Acetone	40 mL/2 L	2%	
		Column temperature 4	15 °C	

Care

Regeneration

Contamination with low-valency hydrophilic ions:

a) Rinse with ultrapure water (25 min at 0.3 mL/min) b) Rinse with 10x concentrated eluent (100 min at

0.3 mL/min)

c) Rinse with ultrapure water (25 min at 0.3 mL/min) d) Rinse with eluent (100 min at 0.3 mL/min) Contaminations with high-valency hydrophobic ions and organic contaminations:

a) Rinse with ultrapure water (25 min at 0.3 mL/min)


b) Rinse with 100% acetonitrile (20 min at 0.3 mL/min)

- c) Rinse with ultrapure water (25 min at 0.3 mL/min) d) Rinse with 10x concentrated eluent (100 min at
- 0.3 mL/min) e) Rinse with ultrapure water (25 min at 0.3 mL/min)

f) Rinse with eluent (100 min at 0.3 mL/min)


Storage

In the eluent at max. 8 °C




Carbonate eluent, standard, 45 °C, EPA 300.1,

Par	ts 1 and 2 (without	MCS)		Cor	nc. (mg/L)		
1	Fluoride	2.00	7	Bromide	10.00		
2	Chlorite	10.00	8	Chlorate	20.00		
3	Bromate	20.00	9	DCA	5.00		
4	Chloride	3.00	10	Nitrate	10.00		
5	System peak	-	11	Phosphate	20.00		
6	Nitrite	10.00	12	Sulfate	15.00		




Carbonate eluent, nuclear power plant,

C
te
t




Carbonate eluent, mod. 1, «Bayer liquor»

after	inline neutralization		35 °C	Conc	. (mg/L)
1	Fluoride	0.35	6	Sulfate	22.59
3	Acetate	9.13	7	Malonate	3.87
4	Formate	3.71	8	Succinate	1.93
5	Chloride	5.66	9	Oxalate	2.07
5	Chloride	5.66	9	Oxalate	2.07



Carb	oonate eluent, c	lrinking water,	45 °(C, EPA 300.1	,
Part	s 1 and 2 (with	MCS)			Conc. (mg/L)
1	Fluoride	0.099	6	Nitrite	0.002
2	Acetate	n.q.	7	Bromide	0.008
3	Formate	n.q.	8	Nitrate	4.378
4	Bromate	0.002	9	Sulfate	35.62
5	Chloride	6.94			



Gradient: Carbonate eluent 1–5 mmol/L,

stanc	lard, 45 °C			Conc	:. (mg/L)
1	Galacturonate	10.00	12	Bromide	10.00
2	Fluoride	2.00	13	DCA	10.00
3	Glycolate	10.00	14	Nitrate	5.00
4	Acetate	10.00	15	Phosphite	10.00
5	Propionate	10.00	16	Phosphate	10.00
6	Formate	5.00	17	Sulfate	5.00
7	Pyruvate	10.00	18	Tartrate	10.00
8	Methacrylate	10.00	19	Selenate	5.00
9	Monochloroacetate	e 5.00	20	Oxalate	5.00
10	Chloride	5.00	21	Arsenate	10.00
11	Nitrite	5.00	22	n. ident.	-

Ordering information	
Metrosep A Supp 7 - 250/4.0	6.1006.630
Metrosep A Supp 5 Guard/4.0	6.1006.500
Metrosep A Supp 5 S-Guard/4.0	6.1006.540
Metrosep A Supp 16 Guard/4.0	6.1031.500
Metrosep A Supp 16 S-Guard/4.0	6.1031.510
Metrosep RP 2 Guard/3.5	6.1011.030
Replacement filters for RP 2 Guard/3.5 (10 pcs.)	6.1011.130

Conc. (µg/L)

2.26 2.06

2.12

1.91

2.18

Metrosep A Supp 10 - 50/4.0 (6.1020.050)

The Metrosep A Supp 10 - 50/4.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 μ m. This proven column concept optimized by Metrohm is characterized by its robust construction, high selectivity, and outstanding separating efficiency. High plate numbers and the favorable position of the system peak between fluoride and chloride complete its properties. Temperature, flow, and eluent composition can be used to modify the properties of the column to accommodate current applications directly.

The short length in conjunction with the relatively low overall capacity of this 50 mm column enable rapid separations of standard anions. They can be determined in less than nine minutes at a flow rate of 1.0 mL/min. The Metrosep A Supp 10 - 50/4.0 is well-suited to simple separation problems and uncomplicated matrices.

Applications

- Standard anions
- Separation of azide and nitrate
- Simple separation problems
- Uncomplicated matrices
- Short analysis times

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary
	ammonium groups
Column dimensions	50 x 4.0 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0–100%
pH range	0–14
Temperature range	10–70 °C
Capacity	17 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium hydrogen carbonate	840 mg/2 L	5.0 mmol/L
(standard eluent)	Sodium carbonate	1060 mg/2 L	5.0 mmol/L
		Column temperature 45 °C	
Hydroxide eluent	Sodium hydroxide (30%)	20 mL/2 L	100 mmol/L


Care


RegenerationSRinse with 50 mL of a 0.05 mol/L solution of Na $_4$ EDTA atIa flow rate of 0.5 mL/min. Then rinse with 0.1 mol/LNaOH at 0.5 mL/min for 1 h.

Organic contaminants:

Rinse with 70% methanol at 1.0 mL/min for 12 h. The addition of 1% acetic acid may be useful.







5.00 | 4


10.00 6

5.00

5

Conc. (mg/L)

Bromide 10.00 Nitrate 10.00 Sulfate 10.00



		5,			, J,
1	System peak	-	3	Nitrate	9.64
2	Chloride	10.05	4	Sulfate	5.19

Ordering information

Chloride

Phosphate

Nitrite

1

2

3

Metrosep A Supp 10 - 50/4.0 Metrosep A Supp 10 Guard/4.0 Metrosep A Supp 10 S-Guard/4.0 6.1020.050 6.1020.500 6.1020.510

Metrosep A Supp 10 - 75/4.0 (6.1020.070)

The Metrosep A Supp 10 - 75/4.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 µm. Under standard conditions, phosphate elutes between nitrite and bromide. Applications can be optimized by modifying temperature, composition of the eluent and flow.

The capacity of the A Supp 10 - 75/4.0 has been optimized with respect to two aspects: matrix and speed. Baseline separation is achieved in samples with high ionic strength, e.g. for phosphate in cola beverages. Even in the presence of large quantities of nitrate and sulfate, the analysis time remains less than 7.5 minutes. High sample through-put is also of great importance in air analytics.

Applications

- Standard anions
- Air monitoring
- Aerosols with PILS/MARGA
- Separation of sulfite and sulfate
- Phosphate in addition to cyclamate in cola beverages

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary
	ammonium groups
Column dimensions	75 x 4.0 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0-100%
pH range	0-14
Temperature range	10–70 °C
Capacity	22 µmol (Cl⁻)

Eluent

Carbonate eluent (standard eluent)	Sodium hydrogen carbonate Sodium carbonate	840 mg/2 L 1060 mg/2 L Column temperature 45 °C	5.0 mmol/L 5.0 mmol/L
Cola eluent	Sodium hydrogen carbonate Sodium carbonate	67 mg/2 L 1695 mg/2 L	0.4 mmol/L 8.0 mmol/L
		5	
		Column temperature 30 °C	
Carbonate eluent	Sodium hydrogen carbonate		4.0 mmol/L
Carbonate eluent (modified)	Sodium hydrogen carbonate Sodium carbonate		4.0 mmol/L 6.0 mmol/L
	, ,	672 mg/2 L	

Care


Regeneration

Rinse with 50 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.5 mL/min. Then rinse with 0.1 mol/L NaOH at 0.5 mL/min for 1 h.

Organic contaminants:

Rinse with 70% methanol at 1.0 mL/min for 12 h. The addition of 1% acetic acid may be useful.






S/cm	+																
2.2	+									8	1						
2.0	+										2						
1.8	+										1						
1.6	+																
1.4	+																
1.2	+																
1.0	+																
0.8	+																
0.6	+										1						
0.4	+			~		•					1						
0.2	+	-	_	-	-			-	_	+	1	-			-	/	-
0.0	누	-	-		-		-		-			-			-	-	_
	0.0	0.5	1.0	1.5	2.0	2.5	3.0	3.5	4.0	4.5	5.0	5.5	6.0	6.5	7.0	7.5	min

Cola eluent, cola beverage, sample volume 250 nL, 30 °C 1 Phosphate 587.3

Conc. (mg/L)



4 Unknown –

Ordering information	
Metrosep A Supp 10 - 75/4.0	6.1020.070
Metrosep A Supp 10 Guard/4.0	6.1020.500
Metrosep A Supp 10 S-Guard/4.0	6.1020.510
Metrosep A Supp 10 Guard HC/4.0	6.1020.520

Metrosep A Supp 10 - 100/4.0 (6.1020.010)

The Metrosep A Supp 10 - 100/4.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 μ m. This column is characterized by high plate numbers and high selectivity. Sulfite and sulfate thus can be reliably separated in the eluent without the addition of organic modifiers. These characteristics are completed by great flexibility with respect to column temperature, flow, and the composition of the eluent.

The Metrosep A Supp 10 - 100/4.0 is the column of choice for routine applications. Its robust construction, excellent price-performance ratio, and very good separating efficiency, in conjunction with moderate chromatography times, make the Metrosep A Supp 10 - 100/4.0 a highly universal anion-separation column.

Applications

- Standard anions
- Separation of sulfite and sulfate
- Simple separation problems
- Uncomplicated matrices

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with guaternary
	ammonium groups
	5 1
Column dimensions	100 x 4.0 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0–100%
pH range	0–14
Temperature range	10–70 °C
Capacity	37 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium hydrogen carbonate	840 mg/2 L	5.0 mmol/L
(standard eluent)	Sodium carbonate	1060 mg/2 L	5.0 mmol/L
		Column temperature 45 °C	
Carbonate eluent	Sodium hydrogen carbonate	672 mg/2 L	4.0 mmol/L
(modified)	Sodium carbonate	1272 mg/2 L	6.0 mmol/L
	Sodium perchlorate	1.2 mg/2 L	5.0 µmol/L
		Room temperature	


Care

Regeneration


Rinse with 50 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.5 mL/min. Then rinse with 0.1 mol/L NaOH at 0.5 mL/min for 1 h.

Organic contaminants: Rinse with 70% methanol at 1.0 mL/min for 12 h. The addition of 1% acetic acid may be useful.












Carbonate eluent, top (black), temperature 45 °C, without 5 μ mol/L ClO₄⁻; bottom (red), temperature 25 °C, with 5 µmol/L ClO₄⁻

tem	perature 25 °C, wit	h 5 µmol/L C	10_4^{-}	Col	nc. (mg/L)
1	System peak	-	5	Bromide	10.00
2	Chloride	2.00	6	Nitrate	10.00
3	Nitrite	5.00	7	Sulfate	10.00

4 Phosphate 10.00








Carbonate eluent, mod., standard, room

tempe	erature				Conc. (mg/L)
1	Chloride	50.00	4	Bromide	10.00
2	Phosphate	10.00	5	Nitrate	10.00
3	Sulfite	10.00	6	Sulfate	50.00



Carbonate eluent, temperature 30...70 ° C in 5 °C increments (from top to bottom), with CO₂ suppressor Conc. (ma/L)

(110111	top to bottom),	with CO ₂ su	phie	:5501	CONC. (ING/L)
1	Fluoride	2.00	5	Bromide	10.00
2	Chloride	2.00	6	Nitrate	10.00
3	Nitrite	5.00	7	Sulfate	10.00
4	Phosphate	10.00			

6.1020.010 6.1020.500 6.1020.510 6.1020.520

Ordering information

Metrosep A Supp 10 - 100/4.0 Metrosep A Supp 10 Guard/4.0 Metrosep A Supp 10 S-Guard/4.0 Metrosep A Supp 10 Guard HC/4.0

Metrosep A Supp 10 - 250/4.0 (6.1020.030)

The Metrosep A Supp 10 - 250/4.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 μ m. The longest column of the A-Supp-10 product range offers the greatest selectivity and flexibility. Utilization of the MSM-HC is particularly recommended with longer chromatogram duration. Changes in temperature, flow, and composition of the eluent also enable a wide variety of separations of anions on this separation column.

Metrosep A Supp 10 - 250/4.0 has a very high capacity. It is suitable for samples with high ionic strength, for complex separation tasks and for analyzing samples in which great differences in concentration between the individual components are present.

Applications

- Standard anions
- Complex separation problems
- Difficult matrices
- Anions in concentrated acids
- Aggressive matrices

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary
	ammonium groups
Column dimensions	250 x 4.0 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0–100%
pH range	0–14
Temperature range	10-70 °C
Capacity	100 µmol (Cl⁻)

Eluent

Carbonate eluent (standard eluent)	Sodium hydrogen carbonate Sodium carbonate	1060 mg/2 L	5.0 mmol/L 5.0 mmol/L
		Column temperature 45 °C	
Carbonate eluent	Sodium hydrogen carbonate	840 mg/2 L	5.0 mmol/L
(modified 1)	Sodium carbonate	1060 mg/2 L	5.0 mmol/L
	Acetone	50 mL	2.5%
		Column temperature 50 °C	
Carbonate eluent	Sodium hydrogen carbonate	672 mg/2 L	4.0 mmol/L
(modified 2)	Sodium carbonate	1272 mg/2 L	6.0 mmol/L
		Column temperature 45 °C	

Care

Regeneration

Column purification:


Rinse with 50 mL of a 0.05 mol/L solution of Na_4EDTA at a flow rate of 0.5 mL/min. Then rinse with 0.1 mol/L NaOH at 0.5 mL/min for 1 h.

Organic contaminants:


Rinse with 70% methanol at 1.0 mL/min for 12 h. The addition of 1% acetic acid may be useful.








			· ·		····
1	Chloride	5.00	4	Bromide	10.00
2	Nitrite	5.00	5	Nitrate	10.00
3	Phosphate	10.00	6	Sulfate	10.00



1 HZ MARKED

 $\begin{array}{c} \mbox{Carbonate eluent, modified 1, artificial} \\ \mbox{seawater 50°C, UV detection } (\lambda = 218 \mbox{ nm}) & \mbox{Conc. (mg/L)} \\ 1 & \mbox{Nitrite} & 2.1 & \mbox{3 Nitrate} & 51.2 \\ 2 & \mbox{Bromide} & 4.4 & \mbox{} \end{array}$



Carbonate eluent, modified 2, ointment based on glycol, 45 °C Conc. (mg/kg) 1 Sulfite 1028° 2 Sulfat n.q. (from metabisulfite) °calculated as metabisulfite

Ordering information

Metrosep A Supp 10 - 250/4.0 Metrosep A Supp 10 Guard/4.0 Metrosep A Supp 10 S-Guard/4.0 Metrosep A Supp 10 Guard HC/4.0

6.1020.030
6.1020.500
6.1020.510
6.1020.520

Metrosep A Supp 16 - 100/4.0 (6.1031.410)

The Metrosep A Supp 16 is a high-capacity column based on a surface-functionalized polystyrene/divinylbenzene copolymer. The functional groups are bonded covalently. The morphology of the anion exchanger results in unique selectivity. In addition, this column type is noteworthy for its high mechanical and chemical resistance.

The column is well-suited to applications with a high ionic load but which require only relatively low resolution. Determination of bromate in water by means of the triiodide method (EPA 326, DIN EN ISO 11206) is another of the numerous applications of the Metrosep A Supp 16 - 100/4.0.

Applications

- Standard anions
- Universal applications
- Bromate (EPA 326, DIN EN ISO 11206)

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	100 x 4.0 mm
Column body	PEEK
Standard flow	0.8 mL/min
Maximum flow	1.2 mL/min
Maximum pressure	20 MPa
Particle size	4.6 µm
Organic modifier	0-10%
pH range	0-14
Temperature range	10-70 °C
Capacity	80 µmol (Cl⁻)

Eluent

Carbonate/hydroxide eluent	Sodium carbonate	1590 mg/2 L	7.5 mmol/L
(standard eluent)	Sodium hydroxide	6.0 mL/2 L	0.75 mmol/L
	(c = 0.25 mol/L)	Column temperature 45 °C	
Sulfuric acid eluent	Sulfuric acid (c = 1 mol/L)	200 mL/2 L	100 mmol/L
	Ammonium heptamolybdate	19.3 mL/2 L	19.3 µmol/L
	(c = 2 mmol/L)	Column temperature 45 °C	
PCR reagent	Potassium iodide	90 g/2 L	0.27 mol/L

Care


Regeneration

Rinse the column overnight (12 h) with standard eluent at a low flow rate (0.4 mL/min).

Rinse the column with one third of the standard flow in the opposite direction for 2 h with 15 mmol/L Na_2CO_3 and then for 2 h with ultrapure water.

Eluent change

When installing or changing to eluents which have an organic modifier to avoid high backpressure, adjust the flow in small increments from 0.4 mL/min to match standard conditions within one hour while maintaining the direction of flow.



Carbonate/hydroxide eluent, drinking water, 45 °CConc. (mg/L)1Fluoriden.q.4Nitrate9.72Chloride9.25Sulfate10.2

_

3 System peak

Sulfuric acid eluent, triiodide method with UV/VIS detection drinking water, 45 °C 1 Bromate 0.6

Ordering information
Metrosep A Supp 16 - 100/4.0

Metrosep A Supp 16 - 100/4.0 Metrosep A Supp 16 Guard/4.0 Metrosep A Supp 16 S-Guard/4.0 6.1031.410 6.1031.500 6.1031.510

Metrosep A Supp 16 - 150/4.0 (6.1031.420)

The Metrosep A Supp 16 is ideal for high-capacity separation problems and excels in its outstanding resolution. The Metrosep A Supp 16 - 150/4.0 is based on a surfacefunctionalized polystyrene/divinylbenzene copolymer. The functional groups are bonded covalently.

The Metrosep A Supp 16 - 150/4.0 is characterized by very good resolution and solves difficult separation problems. The column is well-suited to applications with a high ionic load but which do not require the highest resolution. It is one of the standard columns in anion chromatography.

Applications

- Standard anions
- Universal applications
- Azide/nitrate separation
- Matrices with high ionic strength
- Applications with gradient

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	150 x 4.0 mm
Column body	PEEK
Standard flow	0.8 mL/min
Maximum flow	1.2 mL/min
Maximum pressure	20 MPa
Particle size	4.6 µm
Organic modifier	0-10%
pH range	0–14
Temperature range	10–70 °C
Capacity	125 µmol (Cl⁻)

Eluent

Carbonate/hydroxide eluent	Sodium carbonate	1590 mg/2 L	7.5 mmol/L
(standard eluent)	Sodium hydroxide	6.0 mL/2 L	0.75 mmol/L
	(c = 0.25 mol/L)	Column temperature 45 °C	

Care


Regeneration

Rinse the column overnight (12 h) with standard eluent at a low flow rate (0.4 mL/min).

Rinse the column with one third of the standard flow in the opposite direction for 2 h with 15 mmol/L Na_2CO_3 and then for 2 h with ultrapure water.

Eluent change

When installing or changing to eluents which have an organic modifier to avoid high backpressure, adjust the flow in small increments from 0.4 mL/min to match standard conditions within one hour while maintaining the direction of flow.



Carb	onate/hydroxide e	luent, stand	dard,	45 °C (Conc. (mg/L)	Carbo
1	Fluoride	2.00	5	Bromide	10.00	1
2	Chloride	2.00	6	Nitrate	10.00	2
3	Nitrite	5.00	7	Sulfate	10.00	3
4	System peak	-	8	Phosphate	10.00	

Carbonate/hydroxide eluent, drinking water, 45 °C					Conc. (mg/L)
1	Fluoride	n.q.	4	Nitrate	9.7
2	Chloride	9.2	5	Sulfate	10.2
3	System peak	-			

Ordering information

Metrosep A Supp 16 - 150/4.0 Metrosep A Supp 16 Guard/4.0 Metrosep A Supp 16 S-Guard/4.0 6.1031.420 6.1031.500 6.1031.510

Metrosep A Supp 16 - 250/4.0 (6.1031.430)

The Metrosep A Supp 16 is ideal for high-capacity separation problems and distinguishes itself with its outstanding resolution, even in complex separation problems. The Metrosep A Supp 16 separation column is based on a surface-functionalized polystyrene/divinylbenzene copolymer. The functional groups are bonded covalently. This and the surface structure of the anion exchanger results in unique selectivity. The high-capacity Metrosep A Supp 16 is used for solving complex problems.

The Metrosep A Supp 16 - 250/4.0 is characterized by outstanding resolution and solves the most difficult separation problems. The column is very well-suited to monitoring electroplating baths. Traces of anions can be determined in concentrated acids. Utilization in food analysis for the determination of maltose derivatives is only one more of the numerous applications of the high-capacity Metrosep A Supp 16 - 250/4.0.

Applications

- Standard anions
- Universal applications
- Oligosaccharides and polysaccharides
- Separation of organic acids
- Cl^{-} , SO_4^{2-} in electroplating baths
- Quality monitoring of high-purity chemicals (e.g. conc. acids)
- Complex separation problems
- Difficult matrices

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary
	ammonium groups
Column dimensions	250 x 4.0 mm
Column body	PEEK
Standard flow	0.8 mL/min
Maximum flow	1.2 mL/min
Maximum pressure	20 MPa
Particle size	4.6 µm
Organic modifier	0–10%
pH range	0–14
Temperature range	10-70 °C
Capacity	195 µmol (Cl⁻)

Eluent

Carbonate/hydroxide eluent	Sodium carbonate	1590 mg/2 L	7.5 mmol/L
(standard eluent)	Sodium hydroxide	6.0 mL/2 L	0.75 mmol/L
	(c = 0.25 mol/L)	Column temperature 45 °C	
Hydroxide eluent	Sodium hydroxide	4.0 mL/2 L	20 mmol/L
	(c = 10 mol/L)	Column temperature 32 °C	
Carbonate eluent	Sodium hydrogen carbonate	420 mg/2 L	2.5 mmol/L
	Sodium carbonate	1166 mg/2 L	5.5 mmol/L
		Column temperature 45 °C	

Care

Regeneration

Rinse the column overnight (12 h) with standard eluent at a low flow rate (0.4 mL/min).

Rinse the column with one third of the standard flow in the opposite direction for 2 h with 15 mmol/L Na_2CO_3 and then for 2 h with ultrapure water.

Eluent change


When installing or changing to eluents which have an organic modifier to avoid high backpressure, adjust the flow in small increments from 0.4 mL/min to match standard conditions within one hour while maintaining the direction of flow.




83



Car	bonate/hydroxid	e eluent, stanc	dard,	45 °C	Conc. (mg/L)
1	Fluoride	2.00	5	Nitrate	10.00
2	Chloride	2.00	6	Sulfate	10.00
3	Nitrite	5.00	7	Phospha	te 10.00
4	Bromide	10.00			



Hydi	roxide eluent, malte	ose derivati	ve, 32	2 °C; PAD Conc. (mg/L)
1	Glucose	5.00	5	Maltopentaose 20.00
2	Maltose	5.00	6	Maltohexaose 30.00
3	Maltotriose	10.00	7	Maltoheptaose 40.00
4	Maltotetraose	10.00		



Carbonate eluent, diluted sulfuric acid

after r	neutralization			Conc.	(mg/L)
1	Fluoride	0.50	6	Bromide	0.50
2	Formate	n.q.	7	Nitrate	0.50
3	Acetate	n.q.	8	Phosphate	0.50
4	Chloride	0.50	9	Sulfate	n.q.
5	Nitrite	0.50			

Ordering information	
Metrosep A Supp 16 - 250/4.0	6.1031.430
Metrosep A Supp 16 Guard/4.0	6.1031.500
Metrosep A Supp 16 S-Guard/4.0	6.1031.510

Metrosep A Supp 17 - 100/4.0 (6.01032.410)

The Metrosep A-Supp-17 columns are anion separation columns for use at room temperature with a very good price-performance ratio. The Metrosep A Supp 17 -100/4.0 enables rapid separation of the standard anions. Thanks to its great flexibility with respect to flow rates (up to 1.8 mL/min), very short analysis times can be achieved, depending on the separation problem.

Applications

- Anion determinations at room temperature
- Simple water analysis

Technical information	
Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	100 x 4.0 mm
Column body	PEEK
Standard flow	0.6 mL/min
Maximum flow	1.8 mL/min
Maximum pressure	18 MPa
Particle size	5.0 µm
Organic modifier	0–100% methanol
	0–40% acetone or
	acetonitrile
pH range	0–14
Temperature range	10-70 °C
Standard temperature	25 °C
Capacity	43 µmol (Cl⁻)

Carbonate eluent	Sodium hydrogen carbonate	33.6 mg/2 L	0.2 mmol/L
(standard eluent)	Sodium carbonate	1060 mg/2 L	5.0 mmol/L

Care

Eluent

Preparation Rinse the column with eluent for 2-3 h.

Regeneration


Inorganic contamination


- 1. Rinse with ultrapure water (20 min at 0.3 mL/min)
- 2. Rinse with 10x concentrated standard eluent (120 min at 0.3 mL/min)
- 3. Rinse with ultrapure water (20 min at 0.3 mL/min)
- 4. Rinse with standard eluent (120 min at 0.3 mL/min)

Organic contamination

- 1. Rinse with 70% methanol (16 h at 0.3 mL/min)
- 2. Rinse with standard eluent (120 min at 0.3 mL/min)







Ordering information	
Metrosep A Supp 17 - 100/4.0	
Metrosep A Supp 17 Guard/4.0	
Metrosep A Supp 17 S-Guard/4.0	

6.01032.410 6.01032.500 6.01032.510 85

Metrosep A Supp 17 - 150/4.0 (6.01032.420)

The separation column Metrosep A Supp 17 - 150/4.0 is the column of choice for anion determinations that require good separating efficiency and short separation times at room temperature. The maximum flow rate of 1.4 mL/min offers in this connection the possibility of optimizing the determination. The Metrosep A-Supp-17 columns convince with their good price-performance ratio.

Applications

- Anion determinations at room temperature
- Water analysis

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	150 x 4.0 mm
Column body	PEEK
Standard flow	0.6 mL/min
Maximum flow	1.4 mL/min
Maximum pressure	18 MPa
Particle size	5.0 µm
Organic modifier	0–100% methanol
	0–40% acetone or
	acetonitrile
pH range	0–14
Temperature range	10-70 °C
Capacity	65 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium hydrogen carbonate	33.6 mg/2 L	0.2 mmol/L
(standard eluent)	Sodium carbonate	1060 mg/2 L	5.0 mmol/L

Care

Preparation Rinse the column with eluent for 2–3 h.

Regeneration

Inorganic contamination

- 1. Rinse with ultrapure water (20 min at 0.3 mL/min)
- 2. Rinse with 10x concentrated standard eluent (120 min at 0.3 mL/min)
- 3. Rinse with ultrapure water (20 min at 0.3 mL/min)
- 4. Rinse with standard eluent (120 min at 0.3 mL/min)

Organic contamination

1. Rinse with 70% methanol (16 h at 0.3 mL/min)

2. Rinse with standard eluent (120 min at 0.3 mL/min)




3

4

Nitrite

Bromide




1.25


1.25

7

Phosphate

1.25





Ordering information		
Metrosep A Supp 17 - 150/4.0	6.01032.420	
Metrosep A Supp 17 Guard/4.0	6.01032.500	
Metrosep A Supp 17 S-Guard/4.0	6.01032.510	

Metrosep A Supp 17 - 250/4.0 (6.01032.430)

The Metrosep A Supp 17 - 250/4.0 combines high separating efficiency with a good price-performance ratio without requiring the use of a column oven. The polystyrene/divinylbenzene basic material used guarantees a long service life for the column. Complex separation tasks can be solved on this column.

Applications

- Anion determinations at room temperature
- Water analysis
- Wastewater analysis
- Difficult matrices

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	250 x 4.0 mm
Column body	PEEK
Standard flow	0.6 mL/min
Maximum flow	0.9 mL/min
Maximum pressure	18 MPa
Particle size	5.0 µm
Organic modifier	0–100% methanol
	0–40% acetone or
	acetonitrile
pH range	0–14
Temperature range	10-70 °C
Capacity	109 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium hydrogen carbonate	33.6 mg/2 L	0.2 mmol/L
(standard eluent)	Sodium carbonate	1060 mg/2 L	5.0 mmol/L
Hydroxide eluent	Sodium hydroxide (c = 10 mol/L)	10 mL/2 L	50 mmol/L

Care

Preparation Rinse the column with eluent for 2-3 h.

Regeneration


Inorganic contamination

- 1. Rinse with ultrapure water (20 min at 0.3 mL/min)
- 2. Rinse with 10x concentrated standard eluent (120 min at 0.3 mL/min)
- 3. Rinse with ultrapure water (20 min at 0.3 mL/min)
- 4. Rinse with standard eluent (120 min at 0.3 mL/min)


Organic contamination


- 1. Rinse with 70% methanol (16 h at 0.3 mL/min)
- 2. Rinse with standard eluent (120 min at 0.3 mL/min)





10-5	0 mmol/L, standard				
1	Fluoride	2.0	5	Nitrate	10.0
2	Chloride	2.0	6	Sulfate	10.0
3	Nitrite	10.0	7	Phosphate	10.0
4	Bromide	10.0			





Ordering information			
Metrosep A Supp 17 - 250/4.0	6.01032.430		
Metrosep A Supp 17 Guard/4.0	6.01032.500		
Metrosep A Supp 17 S-Guard/4.0	6.01032.510		
Metrosep A Supp 17 S-Guard - 50/4.0	6.01032.530		



Separation columns



Microbore IC anion-separation columns for lower eluent consumption and greater sensitivity

Metrosep A Supp 4 - 250/2.0 (6.01021.230)

The microbore column Metrosep A Supp 4 - 250/2.0 is an extremely robust column with very good separation properties. The separation phase is comprised of polyvinyl alcohol particles with quaternary ammonium groups and a diameter of 9 µm. This structure guarantees great stability and a greater tolerance to very small particles that could pass through the integrated filter pad. The Metrosep A Supp 4 - 250/2.0 has a medium ion exchange capacity and is particularly suitable for all routine tasks in water analysis.

To protect the IC separation column – even though it is not particularly sensitive to contaminants – we recommend the use of the Metrosep A Supp 4 Guard/2.0 or the A Supp 4 S-Guard/2.0.

Applications

- Standard anions
- Water analysis
- Difficult matrices
- Critical samples
- lodide
- IC-MS coupling

Technical information

Substrate	Polyvinyl alcohol with
	quaternary
	ammonium groups
Column dimensions	250 x 2.0 mm
Column body	PEEK
Standard flow	0.25 mL/min
Maximum flow	0.7 mL/min
Maximum pressure	15 MPa
Particle size	9 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Temperature range	20-60 °C
Capacity	11 µmol (Cl⁻)

Eluent

With chemical suppression			
Carbonate eluent	Sodium hydrogen carbonate	286 mg/2 L	1.7 mmol/L
(standard eluent)	Sodium carbonate	382 mg/2 L	1.8 mmol/L

Care

Regeneration

Contamination with hydrophilic ions:

- a) Rinse with ultrapure water (15 min at 0.1 mL/min)
- b) Rinse with 10x concentrated eluent


(60 min at 0.1 mL/min)

c) Rinse with ultrapure water (15 min at 0.1 mL/min)

d) Rinse with eluent (60 min at 0.1 mL/min)

Contamination with lipophilic ions:

a) Rinse with ultrapure water (15 min at 0.1 mL/min)
b) Rinse with 5% acetonitrile (10 min at 0.1 mL/min)
c) Rinse with 100% acetonitrile (60 min at 0.1 mL/min)
d) Rinse with 50% acetonitrile (10 min at 0.1 mL/min)
e) Rinse with ultrapure water (30 min at 0.1 mL/min)
f) Rinse with eluent (60 min at 0.1 mL/min)



4.0			3	4 1 1 1 1		10.0 - 5.0 - 0.0 -			12	l,	4 ^	,	5	6 ^		+
	0.0 2.0 4.0	6.0 8.0 10	0	12.0 14.0 16.0	18.0 min		0.0 2.0	4.0	6.0		10.0	12.0	14.0	16.0	18.0	min
Carb	onate eluent, s	tandard		Co	nc. (mg/L)	Carbo	onate elu	ent, 1	treated	wastev	wate	r		Cor	nc. (m	ng/L)
1	Fluoride	0.008	4	Phosphate	0.080	1	Fluoride	5		0.16	4	Ν	itrate		2	6.75
2	Chloride	4.000	5	Sulfate	1.200	2	Unknov	vn		-	5	Pl	nosph	ate		0.99
3	Nitrate	0.800				3	Chlorid	e	11	3.72	6	Si	ulfate		3	0.66

Ordering information	
Metrosep A Supp 4 - 250/2.0	6.01021.230
Metrosep A Supp 4 Guard/2.0	6.01021.600
Metrosep A Supp 4 S-Guard/2.0	6.01021.610

Metrosep A Supp 5 - 150/2.0 (6.1006.220)

The Metrosep A Supp 5 - 150/2.0 in the microbore version is distinguished for its excellent separation properties. The particle size of 5 μ m makes a decisive contribution to the separating efficiency of this column. The Metrosep A Supp 5 - 150/2.0 offers the optimum combination of selectivity and capacity, with which even complex separation tasks can be solved within a short time. The 2 mm Metrosep A Supp 5 separation columns are packed with the same material as the corresponding 4 mm separation columns. The 150 mm version of this column type is used for universal applications at low eluent consumption.

With its low eluent flow, this column is particularly suitable for IC/MS coupling.

Applications

- Standard anions
- F^{*}, Cl^{*}, Br^{*}, l^{*}
- CIO₂, CIO₃, CIO₄, BrO₃
- Cr (VI) (CrO₄²⁻)
- Method development
- IC-MS coupling

Technical information

Substrate	Polyvinyl alcohol with
	quarternary ammonium
	groups
Column dimensions	150 x 2.0 mm
Column body	PEEK
Standard flow	0.18 mL/min
Maximum flow	0.21 mL/min
Maximum pressure	20 MPa
Particle size	5 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Temperature range	20–60 °C
Capacity	12 µmol (Cl⁻)

Eluent

94

Carbonate eluent	Sodium hydrogen carbonate	168 mg/2 L	1.0 mmol/L
(standard eluent)	Sodium carbonate	678 mg/2 L	3.2 mmol/L

Care

Regeneration

Contamination with low-valence hydrophilic ions


- 1. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 2. Rinse with 10x concentrated eluent (100 min at 0.1 mL/min)
- 3. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 4. Rinse with eluent (100 min at 0.1 mL/min)


Contamination with high-valence hydrophobic ions or organic contaminations

- 1. Rinse with eluent (100 min at 0.1 mL/min)
- 2. Rinse with 5% acetonitrile (20 min at 0.1 mL/min)
- 3. Rinse with 100% acetonitrile (60 min at 0.1 mL/min)
- 4. Rinse with 50% acetonitrile (10 min at 0.1 mL/min)
- 5. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 6. Rinse with eluent (100 min at 0.1 mL/min)

Shifted system peak

- 1. Regeneration method with column oven
- 2. Rinse with concentrated eluent of 1 mol/L Na_2CO_3 (25 min at 0.1 mL/min)
- 3. Maintain for 10–12 hours at 45–50 °C (without rinsing)
- 4. Rinse with the normal eluent (at least 40 min at 0.1 mL/min)





Carbo	nate eluent, polyeti	nyiene pe	ellets r	'y	
Comb	oustion IC (CIC)			Co	onc. (mg/kg)
1	Chloride	94.2	3	Sulfate	74.7
2	Bromide	84.0			

Ordering information

Metrosep A Supp 5 - 150/2.0 Metrosep A Supp 5 Guard/2.0 Metrosep A Supp 5 S-Guard/2.0 6.1006.220 6.1006.600 6.1006.610

Metrosep A Supp 5 - 250/2.0 (6.1006.230)

The Metrosep A Supp 5 - 250/2.0 is the microbore highperformance separation column with which even complex separation problems can be solved easily and reproducibly. The range of applications possible with this column far exceeds the detection of standard anions. The Metrosep A Supp 5 - 250/2.0 is used wherever maximum separating efficiency must be combined with both the lowest of detection limits and low eluent consumption.

With its low eluent flow, this column is particularly suitable for IC/MS coupling.

Applications

- Standard anions
- F^{*}, Cl^{*}, Br^{*}, l^{*}
- CIO₂, CIO₃, CIO₄, BrO₃
- ClO_4^- at high ionic strength
- BrO₃⁻ at high ionic strength
- Method development
- Universal applications
- Difficult matrices
- Complex separation problems
- Applications with gradient
- IC-MS coupling

Technical information

Substrate	Polyvinyl alcohol with
	quarternary ammonium
	groups
Column dimensions	250 x 2.0 mm
Column body	PEEK
Standard flow	0.18 mL/min
Maximum flow	0.21 mL/min
Maximum pressure	20 MPa
Particle size	5 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Temperature range	20–60 °C
Capacity	21 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium hydrogen carbonate	168 mg/2 L	1.0 mmol/L
(standard eluent)	Sodium carbonate	678 mg/2 L	3.2 mmol/L

Care

Regeneration

Contamination with low-valence hydrophilic ions

- 1. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 2. Rinse with 10x concentrated eluent (100 min at 0.1 mL/min)
- 3. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 4. Rinse with eluent (100 min at 0.1 mL/min)

Contamination with high-valence hydrophobic ions or organic contaminations


- 1. Rinse with eluent (100 min at 0.1 mL/min)
- 2. Rinse with 5% acetonitrile (20 min at 0.1 mL/min)
- 3. Rinse with 100% acetonitrile (60 min at 0.1 mL/min)
- 4. Rinse with 50% acetonitrile (10 min at 0.1 mL/min)


- 5. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 6. Rinse with eluent (100 min at 0.1 mL/min)


Shifted system peak

- 1. Regeneration method with column oven
- 2. Rinse with concentrated eluent of 1 mol/L $\rm Na_2CO_3$ (25 min at 0.1 mL/min)
- 3. Maintain for 10–12 hours at 45–50 $^{\circ}\mathrm{C}$ (without rinsing)
- 4. Rinse with the normal eluent (at least 40 min at 0.1 mL/min)









Ordering information

Metrosep A Supp 5 - 250/2.0 Metrosep A Supp 5 Guard/2.0 Metrosep A Supp 5 S-Guard/2.0 6.1006.230 6.1006.600 6.1006.610 97

Metrosep A Supp 7 - 150/2.0 (6.1006.640)

The Metrosep A Supp 7 - 150/2.0 is the shorter of the two A Supp 7 columns in the microbore version. It allows similarly complex separation tasks to be solved the same way as with the corresponding 250 mm version, with no significant loss in separating efficiency.

With the Metrosep A Supp 7 - 150/2.0, these ions are determined with certainty and precision down to the lower μ g/L range. High detection sensitivity is achieved by using the 5 μ m polyvinyl alcohol polymer, which allows extremely high plate numbers and therefore outstanding separation and detection properties can be achieved. In addition, the separation can be adapted to the specific requirements of the application by modifying the temperature.

This microbore column is particularly suitable for use with an MS detector.

Applications

- Standard anions
- Fast analysis (high flow rate)
- Applications with gradient
- IC-MS

Technical information

Substrate	Polyvinyl alcohol with
	quarternary ammonium
	groups
Column dimensions	150 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	0.6 mL/min
Maximum pressure	20 MPa
Particle size	5 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile and
	methanol)
pH range	3–12
Temperature range	20–60 °C
Capacity	18 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium carbonate	763 mg/2 L	3.6 mmol/L
(standard eluent)	(column temperature: 45 °C)		
Carbonate eluent	Sodium carbonate	878 mg/2 L	4.0 mmol/L
(modified)	(column temperature: 55 °C)		

Care

Regeneration procedure for contamination with lowvalency hydrophilic ions:


- 1. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 2. Rinse with 10x concentrated eluent (100 min at 0.1 mL/min)
- 3. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 4. Rinse with eluent (100 min at 0.1 mL/min)


Regeneration procedure for contamination with highvalency hydrophobic ions and organic contaminations:

- 1. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 2. Rinse with 100% acetonitrile (20 min at 0.1 mL/min)
- 3. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 4. Rinse with 10x concentrated eluent (100 min at 0.1 mL/min)
- 5. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 6. Rinse with eluent (100 min at 0.1 mL/min)

Storage In the eluent at maximum +8 °C.







Ordering information

Metrosep A Supp 7 - 150/2.0 Metrosep A Supp 5 Guard/2.0 Metrosep A Supp 5 S-Guard/2.0 Metrosep A Supp 16 Guard/2.0 Metrosep A Supp 16 S-Guard/2.0 6.1006.640 6.1006.600 6.1006.610 6.1031.600 6.1031.610

Metrosep A Supp 7 - 250/2.0 (6.1006.650)

Disinfection byproducts from water treatment are suspected not only of being health hazards but also of being carcinogenic. Oxyhalides have therefore become the subject of many investigations and standards (e.g., EPA 300.1 Part B, EPA 317.0, EPA 326.0). Of primary concern is bromate, which forms from bromide during the ozonization of drinking water.

The microbore version of the Metrosep A Supp 7 - 250/4.0 is a high-performance separation column for the parallel determination of standard anions, oxyhalides and dichloro-acetic acid. With this column, these ions are determined with certainty and precision down to the lower μ g/L range. High detection sensitivity is achieved by using the 5 μ m polyvinyl alcohol polymer, which allows extremely high plate numbers and therefore outstanding separation and detection properties can be achieved. In addition, the separation can be adapted to the specific requirements of the application by modifying the temperature.

This microbore column is particularly suitable for use with an MS detector.

Applications

- Standard anions
- EPA 300.1 Part B, simultaneous determination of standard anions and ClO₂⁻, ClO₃⁻, BrO₃⁻ and DCA (dichloroacetic acid)
- Isocratic separation of glycolate, acetate and formate
- Complex separation tasks
- Applications with gradient
- IC-MS

Technical information	
Substrate	Polyvinyl alcohol with
	quarternary ammonium
	groups
Column dimensions	250 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	0.4 mL/min
Maximum pressure	20 MPa
Particle size	5 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile and
	methanol)
pH range	3–12
Temperature range	20–60 °C
Capacity	27 µmol (Cl⁻)

Eluent

Carbonate eluent	Sodium carbonate	763 mg/2 L	3.6 mmol/L
(standard eluent)	(column temperature: 45 °C)		

Care

Regeneration procedure for contamination with low-valency hydrophilic ions:

- 1. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 2. Rinse with 10x concentrated eluent (100 min at 0.1 mL/min)
- 3. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 4. Rinse with eluent (100 min at 0.1 mL/min)

Regeneration procedure for contamination with highvalency hydrophobic ions and organic contaminations:

- 1. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 2. Rinse with 100% acetonitrile (20 min at 0.1 mL/min)
- 3. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 4. Rinse with 10x concentrated eluent (100 min at 0.1 mL/min)
- 5. Rinse with ultrapure water (25 min at 0.1 mL/min)
- 6. Rinse with eluent (100 min at 0.1 mL/min)

Storage

In the eluent at maximum +8 °C.



µS/cm 1.24

1.22 -

1.20 -

1.0

1.0

30.0


15.0

40.0

1

1 .

Chromatogram



1.1	8 -	1 7				
1.1	6 -	3	5			
1.1	4 -	2	Ĩ.	7		
1.1	2 -		6		9	12
1.1	0 -	IW I		٨	8 10 11	^
1.0	har	1144	A.A.	Λ	~N_AA	1
1.0		8.0 12.0	16.0 20.0	24.0	28.0 32.0 36.0	40.0
	0.0 4.0	8.0 12.0	16.0 20.0	24.0	28.0 32.0 36.0	40.0 min
Car	bonate elue	nt, stan	dard, 45 °0	2	Co	onc. (mg/L)
1	Fluoride		0.1	7	Nitrate	0.1
1 2	Fluoride Acetate		0.1 0.1	7 8	Nitrate Phosphate	0.1 0.1
•				1.		
2	Acetate		0.1	8	Phosphate	0.1
2 3	Acetate Fromate		0.1 0.1	8 9	Phosphate Sulfate	0.1 0.1

2.0 1 -Fluoride 1

Chlorite 2 3 Bromate 4 Chloride

Nitrite 5

6 Bromide

Ordering information Metrosep A Supp 7 - 250/2.0 Metrosep A Supp 5 Guard/2.0 Metrosep A Supp 5 S-Guard/2.0 Metrosep A Supp 16 Guard/2.0 Metrosep A Supp 16 S-Guard/2.0

,		,
2.0	7	Chlorate
1.0	8	Dichloroacetate
1.0	9	Nitrate
10.0	10	Phosphate
5.0	11	Sulfate
1.0		

6.1006.650
6.1006.600
6.1006.610
6.1031.600
6.1031.610

101

Metrosep A Supp 10 - 50/2.0 (6.1020.250)

102

The Metrosep A Supp 10 - 50/2.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 μ m. Temperature, flow and eluent composition can be used to modify the properties of the column to accommodate the current application directly. The 2 mm Metrosep A Supp 10 separation columns are packed with the same material as the corresponding 4 mm separation columns. The short length and associated relatively low overall capacity of this 50 mm column enable very rapid separations of standard anions.

The Metrosep A Supp 10 - 50/2.0 is well-suited to simple separation problems and uncomplicated matrices. Thanks to its low flow, this microbore separation column is ideal for IC-MS applications.

Applications

- Standard anions
- Simple separation problemsUncomplicated matrices
- Short analysis times
- IC-MS coupling

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	50 x 2.0 mm
Column body	PEEK
Standard flow	0.25 mL/min
Maximum flow	1.3 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0-100%
pH range	0-14
Temperature range	10–70 °C
Capacity	4.6 µmol (Cl⁻)

Eluent

Carbonate eluentSodium hydrogen carbonate840 mg/2 L(standard eluent)Sodium carbonate1060 mg/2

nate 840 mg/2 L 1060 mg/2 L Column temperature 45 °C

5.0 mmol/L 5.0 mmol/L

Care

Regeneration


Rinse with 50 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.12 mL/min. Then rinse with 0.1 mol/L NaOH at 0.12 mL/min for 1 h.


Organic contaminants:

Rinse with 70% methanol at 0.12 mL/min for 12 h. The addition of 1% acetic acid may be useful.










Nitrate 8.98

2



2 Nitrate 2.15

Ordering information

Metrosep A Supp 10 - 50/2.0 Metrosep A Supp 10 Guard/2.0 Metrosep A Supp 10 S-Guard/2.0 6.1020.250 6.1020.600 6.1020.610 103

Metrosep A Supp 10 - 75/2.0 (6.1020.270)

104

The Metrosep A Supp 10 - 75/2.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 μ m. This proven column concept optimized by Metrohm is characterized by its robust construction, great selectivity and outstanding separating efficiency. The 2 mm Metrosep A Supp 10 separation columns are packed with the same material as the corresponding 4 mm separation columns. Temperature, flow and eluent composition can be used to modify the properties of the column to accommodate the current application directly.

The capacity of the Metrosep A Supp 10 - 75/2.0 has been optimized with respect to two aspects: matrix and speed. A rapid baseline separation of the standard anions can also be achieved in samples of high ionic strength. Particularly suitable for IC-MS applications.

Applications

- Standard anions
- IC-MS coupling
- Separation of sulfite and sulfate
- Fermentation samples

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary
	ammonium groups
Column dimensions	75 x 2.0 mm
Column body	PEEK
Standard flow	0.25 mL/min
Maximum flow	1.1 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0–100%
pH range	0–14
Temperature range	10-70 °C
Capacity	6.9 µmol (Cl⁻)

Eluent	

 Carbonate eluent
 Sodium hydrogen carbonate
 840 mg/2 L
 5.0 mmol/L

 (standard eluent)
 Sodium carbonate
 1060 mg/2 L
 5.0 mmol/L

 Column temperature 45 °C


Care

Regeneration

Rinse with 50 mL of a 0.05 mol/L solution of Na_4EDTA at a flow rate of 0.12 mL/min. Then rinse with 0.1 mol/L NaOH at 0.12 mL/min for 1 h.


Organic contaminants: Rinse with 70% methanol at 0.12 mL/min for 12 h. The addition of 1% acetic acid may be useful.





Conc. (g/L) Carbonate eluent, fermentation broth Dilution 1 : 100, 45 °C 1.249

Chloride 0.025 | 3 Phosphate Nitrite 6.461



1

2

6.1020.270
6.1020.600
6.1020.610

Metrosep A Supp 10 - 100/2.0 (6.1020.210)

The Metrosep A Supp 10 - 100/2.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 µm. Temperature, flow and eluent composition can be used to modify the properties of the column to accommodate the current application directly. The 2 mm Metrosep A Supp 10 separation columns are packed with the same material as

106

The Metrosep A Supp 10 - 100/2.0 is the microbore column of choice for routine applications. Thanks to the high flow and pressure stability of this separation column, very rapid chromatograms with good separation of the ions can be achieved. The standard anions can thus be separated within less than three minutes. Particularly suitable for IC-MS applications.

the corresponding 4 mm separation columns.

Applications

- Standard anions
- Simple separation problems
- Traces of cyanide and sulfide with PAD
- Uncomplicated matrices
- IC-MS coupling

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary
	ammonium groups
Column dimensions	100 x 2.0 mm
Column body	PEEK
Standard flow	0.25 mL/min
Maximum flow	0.9 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0-100%
pH range	0–14
Temperature range	10–70 °C
Capacity	8.7 µmol (Cl⁻)

Eluent

Carbonate eluent (standard eluent)	Sodium hydrogen carbonate Sodium carbonate	840 mg/2 L 1060 mg/2 L Column temperature 45 °C	5.0 mmol/L 5.0 mmol/L
Hydroxide/EDTA eluent	Sodium hydroxide (c = 20 mol/L)	10 mL/2 L	100 mmol/L
	EDTA	2.0 mg/2 L Column temperature 35 °C	0.007 mmol/L


Care


Regeneration

Rinse with 50 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.12 mL/min. Then rinse with 0.1 mol/L NaOH at 0.12 mL/min for 1 h.

Organic contaminants: Rinse with 70% methanol at 0.12 mL/min for 12 h. The addition of 1% acetic acid may be useful.







Ordering information

Metrosep A Supp 10 - 100/2.0 Metrosep A Supp 10 Guard/2.0 Metrosep A Supp 10 S-Guard/2.0 6.1020.210 6.1020.600 6.1020.610

Metrosep A Supp 10 - 150/2.0 (6.1020.220)

The Metrosep A Supp 10 - 150/2.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 µm. This proven column concept optimized by Metrohm is characterized by its robust construction, great selectivity and outstanding separating efficiency. The 2 mm Metrosep A Supp 10 separation columns are packed with the same material as the corresponding 4 mm separation columns. Temperature, flow and eluent composition can be used to modify the properties of the column to accommodate the current application directly.

The Metrosep A Supp 10 - 150/2.0 separation column is suitable for complex separation tasks with wide differences in concentrations. The microbore version exhibits low eluent consumption and is therefore particularly suitable for IC-MS applications.

Applications

- Standard anions
- Universal applications
- Different matrices
- Transition metal complexes
- Chromium(VI) in toys (EU directive 2009/48/EC)
- IC-MS coupling

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	150 x 2.0 mm
Column body	PEEK
Standard flow	0.25 mL/min
Maximum flow	0.7 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0–100%
pH range	0-14
Temperature range	10-70 °C
Capacity	15 µmol (Cl⁻)

Eluent

Carbonate eluent (standard eluent)	Sodium hydrogen carbonate Sodium carbonate	840 mg/2 L 1060 mg/2 L	5.0 mmol/L 5.0 mmol/L
		Column temperature 45 °C	
Dipicolinic acid eluent	Dipicolinic acid	1.0 g/2 L	3.0 mmol/L
	Sodium sulfate	1.42 g/2 L	10 mmol/L
	Sodium hydroxide (c = 20 mol/L)	6.6 mL/2 L	66 mmol/L
	Formic acid		pH = 4.33

Care


Regeneration

Column purification:

Rinse with 50 mL of a 0.05 mol/L solution of Na_4EDTA at a flow rate of 0.12 mL/min. Then rinse with 0.1 mol/L NaOH at 0.12 mL/min for 1 h.

Organic contaminants:

Rinse with 70% methanol at 0.12 mL/min for 12 h. The addition of 1% acetic acid may be useful.



Ordering information Metrosep A Supp 10 - 150/2.0 Metrosep A Supp 10 Guard/2.0 Metrosep A Supp 10 S-Guard/2.0

6.1020.220 6.1020.600 6.1020.610

Metrosep A Supp 10 - 250/2.0 (6.1020.230)

The Metrosep A Supp 10 - 250/2.0 separation column is based on a high-capacity polystyrene/divinylbenzene copolymer with a particle size of only 4.6 µm. This proven column concept optimized by Metrohm is characterized by its robust construction, great selectivity and outstanding separating efficiency. The 2 mm Metrosep A Supp 10 separation columns are packed with the same material as the corresponding 4 mm separation columns. Temperature, flow and eluent composition can be used to modify the properties of the column to accommodate the current application directly.

The Metrosep A Supp 10 - 250/2.0 has a very high capacity for a microbore column. It is suitable for samples with high ionic strength, for complex separation tasks and for analyses in which great differences in concentration between the individual components are present. Thanks to its low flow, this microbore separation column is ideal for IC-MS applications.

Applications

- Standard anions
- Complex separation problems
- Difficult matrices
- Anions in concentrated acids
- Aggressive matrices
- IC-MS coupling

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	250 x 2.0 mm
Column dimensions	250 x 2.0 mm
Column body	PEEK
Standard flow	0.25 mL/min
Maximum flow	0.7 mL/min
Maximum pressure	25 MPa
Particle size	4.6 µm
Organic modifier	0–100%
pH range	0–14
Temperature range	10-70 °C
Capacity	24 µmol (Cl⁻)

Eluent

Carbonate eluent (standard eluent)	Sodium hydrogen carbonate Sodium carbonate	840 mg/2 L 1060 mg/2 L	5.0 mmol/L 5.0 mmol/L
(, 		Column temperature 45 °C	
Ammonium sulfate eluent	Ammonium sulfate	66.0 g/2 L	250 mmol/L
	Ammonium hydroxide (c = 5.0 mol/L)	40 mL/2 L	100 mmol/L


Care


Regeneration Column purification: Rinse with 50 mL of a 0.05 mol/L solution of Na₄EDTA at a flow rate of 0.12 mL/min. Then rinse with 0.1 mol/L NaOH at 0.12 mL/min for 1 h. Organic contaminants:


Rinse with 70% methanol at 0.12 mL/min for 12 h. The addition of 1% acetic acid may be useful.

Storage In the eluent









Ordering information Metrosep A Supp 10 - 250/2.0 Metrosep A Supp 10 Guard/2.0 Metrosep A Supp 10 S-Guard/2.0

6.1020.230 6.1020.600 6.1020.610

Metrosep A Supp 16 - 100/2.0 (6.1031.210)

112

In the case of the microbore version of the Metrosep A Supp 16 - 100, lower flows are applied through the smaller inner diameter. Eluent consumption is reduced drastically as a result. The dwell time of the ions in the detector becomes longer and the sensitivity or the peak area is increased accordingly (with the same sample amount). Microbore separation columns are used together with the MSM-LC (6.2844.000). The 2 mm Metrosep A Supp 16 separation columns are packed with the same material as the corresponding 4 mm separation columns. The short version of this column type enables extremely rapid separations.

The column is well-suited to applications with a high ionic load but which require only relatively low resolution. With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Applications

- Standard anions
- Universal applications
- Rapid analysis (standard anions in 5 min)
- IC-MS coupling

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	100 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	0.6 mL/min
Maximum pressure	16 MPa
Particle size	4.6 µm
Organic modifier	0-10%
pH range	0–14
Temperature range	10–70 °C
Capacity	22 µmol (Cl⁻)

Eluent

Carbonate/hydroxide eluent	Sodium carbonate	1590 mg/2 L	7.5 mmol/L
(standard eluent)	Sodium hydroxide	6.0 mL/2 L	0.75 mmol/L
	(c = 0.25 mol/L)	Column temperature 45 °C	

Care

Regeneration

Rinse the column overnight (12 h) with standard eluent at low flow rate (0.1 mL/min).


Rinse the column with half of the standard flow in the opposite direction for 2 h with 15 mmol/L $\rm Na_2CO_3$ and then for 2 h with ultrapure water.

Eluent change

When installing or changing to eluents which have an organic modifier to avoid high backpressure, adjust the flow in small increments from 0.1 mL/min within one hour to match standard conditions while maintaining the direction of flow.


Storage In the eluent






Carbonate/hydroxide eluent, standard, 40 °C

Carbo						
flow	rate 0.6 mL/min.			Conc	:. (mg/L)	
1	Fluoride	2.00	5	Nitrate	10.0	
2	Chloride	5.00	6	Sulfate	10.0	
3	Nitrite	5.00	7	Phosphate	10.0	
4	Bromide	10.0				





		4	9 1	0 11 12 13	14 15 min	
onate/hydrox	ide eluent	t, stand	lard, 4	45 °C	Conc. (mg	/L) Carl
Fluoride	:	2.00	5	Bromide	10.0	00 1
Chloride	:	2.00	6	Nitrate	10.0	00 2
Nitrite	!	5.00	7	Sulfate	10.0	00

3 4 System peak

Ordering information Metrosep A Supp 16 - 100/2.0

Metrosep A Supp 16 Guard/2.0

Metrosep A Supp 16 S-Guard/2.0

2

- 5.00 8
- Sulfate Phosphate

10.00

6.1031.210
6.1031.600
6.1031.610

Metrosep A Supp 16 - 150/2.0 (6.1031.220)

114

The microbore version of the Metrosep A Supp 16 - 150 is well-suited to medium-capacity separation problems. Eluent consumption is drastically reduced as a result of the smaller inner diameter of this column type and the correspondingly lower flows. As a result of the lower flows, the dwell time of the anions in the detector, and thus also the peak areas with identical sample amounts, are increased. Microbore separation columns are used together with the MSM-LC (6.2844.000). The 2 mm Metrosep A Supp 16 separation columns are packed with the same material as the corresponding 4 mm separation columns. The medium version of this column type is used for universal applications.

The column is well-suited to applications with a high ionic load but which do not require the highest resolution. With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Applications

- Standard anions
- Universal applications
- Difficult matrices with high ionic strength
- Applications with gradient
- IC-MS coupling

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	150 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	0.3 mL/min
Maximum pressure	16 MPa
Particle size	4.6 µm
Organic modifier	0-10%
pH range	0–14
Temperature range	10-70 °C
Capacity	31 µmol (Cl⁻)

Eluent

Carbonate/hydroxide eluent	Sodium carbonate	1590 mg/2 L	7.5 mmol/L
(standard eluent)	Sodium hydroxide	6.0 mL/2 L	0.75 mmol/L
	(c = 0.25 mol/L)	Column temperature 45 °C	
Carbonate eluent	Sodium carbonate	763 mg/2 L	3.6 mmol/L

Care

Regeneration

Rinse the column overnight (12 h) with standard eluent at low flow rate (0.1 mL/min).

Rinse the column with half of the standard flow in the opposite direction for 2 h with 15 mmol/L $\rm Na_2CO_3$ and then for 2 h with ultrapure water.


Eluent change

When installing or changing to eluents which have an organic modifier to avoid high backpressure, adjust the flow in small increments from 0.1 mL/min to match standard conditions within one hour while maintaining the direction of flow.


Storage In the eluent












Conc. (mg/L) Carbonate/hydroxide eluent, standard, 45 °C Fluoride 1 2.00 | 5 Bromide 10.00 Chloride 2.00 2 6 Nitrate 10.00 3 Nitrite 5.00 7 10.00 Sulfate 4 System peak - 8 10.00 Phosphate



Carb	onate/hydroxide el	uent, drinkii	ng wa	ater, 45 °C	Conc. (mg/L)
1	Fluoride	n.q.	4	Nitrate	9.7
2	Chloride	9.2	5	Sulfate	10.2
3	System peak	-			

Ordering information

Metrosep A Supp 16 - 150/2.0 Metrosep A Supp 16 Guard/2.0 Metrosep A Supp 16 S-Guard/2.0

6.1031.220 6.1031.600 6.1031.610

Metrosep A Supp 16 - 250/2.0 (6.1031.230)

116

The microbore version of the Metrosep A Supp 16 - 250 is well suited to high-capacity separation problems. Lower flows are applied due to the smaller inner diameter of this column type. Eluent consumption is reduced drastically as a result. The dwell time of the ions in the detector becomes longer and the sensitivity or the peak area is increased accordingly (with the same sample amount). Microbore separation columns are used with the MSM-LC (6.2844.000). The 2 mm Metrosep A Supp 16 separation columns are packed with the same material as the corresponding 4 mm separation columns. The separation column is based on a surface-functionalized polystyrene/ divinylbenzene copolymer. The functional groups are bonded covalently. The morphology of the anion exchanger results in unique selectivity. The high-capacity Metrosep A Supp 16 - 250/2.0 is used for solving complex problems.

The Metrosep A Supp 16 - 250/2.0 is characterized by outstanding resolution and solves the most difficult separation problems. With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Applications

- Standard anions
- Universal applications
- Azide/nitrate separation
- Divalent organic acids besides standard anions
- Matrices with high ionic strength
- Applications with gradient
- IC-MS coupling

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	250 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	0.3 mL/min
Maximum pressure	16 MPa
Particle size	4.6 µm
Organic modifier	0-10%
pH range	0–14
Temperature range	10–70 °C
Capacity	51 µmol (Cl⁻)

Eluent

Carbonate/hydroxide eluent	Sodium carbonate	1590 mg/2 L	7.5 mmol/L
(standard eluent)	Sodium hydroxide	6.0 mL/2 L	0.75 mmol/L
	(c = 0.25 mol/L)	Column temperature 45 °C	

Care

Regeneration

Rinse the column overnight (12 h) with standard eluent at low flow rate (0.1 mL/min).


Rinse the column with half of the standard flow in the opposite direction for 2 h with 15 mmol/L $\rm Na_2CO_3$ and then for 2 h with ultrapure water.

Eluent change


When installing or changing to eluents which have an organic modifier to avoid high backpressure, adjust the flow in small increments from 0.1 mL/min to match standard conditions within one hour while maintaining the direction of flow.

Storage In the eluent





Carbonate/hydroxide eluent, standard, 65 °C Conc. (mg/L)			
0.0			
0.0			
0.0			
0.0			
0.0			
0.0			
0.0			
0.0			




4

System peak

Carbonate/hydroxide eluent, standard, 45 $^{\circ}\mathrm{C}$ Conc. (mg/L) 1 Fluoride 2.00 5 Bromide 10.00 Chloride 2.00 10.00 2 6 Nitrate Nitrite 5.00 3 7 10.00 Sulfate

- 8

Phosphate



Carbonate/hydroxide eluent,				Conc. (mg/L)
drinking water, 45 °C				
1 Fluoride	n.q.	4	Nitrate	9.7
2 Chloride	9.2	5	Sulfate	10.2
3 System peak	-			

Ordering information	
Metrosep A Supp 16 - 250/2.0	6.1031.230
Metrosep A Supp 16 Guard/2.0	6.1031.600
Metrosep A Supp 16 S-Guard/2.0	6.1031.610

10.00



Separation columns



IC separation columns for the determination of organic acids – ion-exclusion chromatography

«Inverse suppression» – dissociation desired!

The use of the Metrohm «MSM» suppressor module is recommended to improve sensitivity in the detection of organic acids which are only weakly dissociated. A non-conventional approach is used: The suppressor is charged with lithium ions instead of hydrogen ions. As a result, it is possible to transfer the protonated and thus undissociated acids into their nearly completely dissociated salts. This increases sensitivity in the conductivity detector considerably. The construction is the same as with chemical suppression, except that the suppressor is regenerated with lithium chloride instead of with sulfuric acid. The «MSM» is used as a post-column reactor between the ion-exclusion column and the conductivity detector.

Hamilton PRP-X300 - 250/4.0 (6.1005.030)

120

The Hamilton-PRP-X300 ion-exclusion column is a cationexchanger column with low capacity. The combination of a polystyrene/divinylbenzene copolymer with sulfonic acid groups as ion exchanger is ideal for the solution of simple separation problems. The column is characterized by the possibility of determining the salts of organic acids, in particular the very sensitive determination of formate.

Applications

- Glycolic acid, monochloroacetic acid
- Simple matrices
- Simple separation problems
- Formate determination

Technical information

sulfonic acid groupsColumn dimensions250 x 4.0 mmColumn bodyStainless steelStandard flow1.0 mL/minMaximum flow8.0 mL/min
Column bodyStainless steelStandard flow1.0 mL/min
Standard flow 1.0 mL/min
Maximum flow 8.0 mL/min
Maximum pressure 34 MPa
Particle size 7 µm
Organic modifier 0–100%
pH range 1–13
Temperature range 5–60 °C
Capacity 19 µmol (K ⁺)

Eluents			
Sulfuric acid eluent	Sulfuric acid (c = 0.1 mol/L)	10 mL/2 L	0.5 mmol/L

Care

Regeneration


Divalent cations remain on the column and form complexes with citrate that falsify the citrate peak. Injection of 100 μ L 0.1 mol/L Na₂H₂EDTA.

Rinse the column with 0.01 mol/L H_2SO_4 with 20% methanol at a flow rate of 0.5 mL/min for 6 h.

Storage

For short periods (days) in the eluent, for longer periods (weeks) in methanol/water (1:4)





Ordering information	
Hamilton PRP-X300 - 250/4.0	6.1005.030
Metrosep RP 2 Guard/3.5	6.1011.030
Replacement filters for RP 2 Guard/3.5 (10 pcs.)	6.1011.130
Metrosep RP 3 Guard HC/4.0	6.1011.040

Metrosep Organic Acids - 100/7.8 (6.1005.210)

122 The separation column for the determination of organic acids and weak mineral acids. The low capacity in comparison with the Metrosep Organic Acids - 250/7.8 (6.1005.200) allows the rapid separation of organic acids. This column is suitable primarily for small and medium concentrations in uncomplicated sample matrices.

Applications

- Organic acids: Citrate, tartrate, malate, ascorbate, succinate
- Short-chain fatty acids: Formate, acetate, propionate, butyrate, etc.
- F⁻, CO₃²⁻
- Simple matrices
- Simple separation problems

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with
	sulfonic acid groups
Column dimensions	100 x 7.8 mm
Column body	Stainless steel
Standard flow	0.5 mL/min
Maximum flow	0.6 mL/min
Maximum pressure	7 MPa
Particle size	9 µm
Organic modifier	0-20%
pH range	1–13
Temperature range	5–90 °C

Eluent

Sulfuric acid eluent	Sulfuric acid (c = 2 mol/L)	0.5 mL/2 L	0.5 mmol/L
(standard eluent)	Acetone	300 mL/2 L	15%
Oxalic acid eluent	Oxalic acid	45 mg/2 L	0.25 mmol/L

Care

Regeneration

Column purification: Rinse the column in the opposite direction with 20% acetonitrile in 0.01 mol/L H_2SO_4 at a flow rate of 0.1 mL/min for 4 hours at 65 °C.

Contaminations with metals:

If retention times are shortened: Rinse the column in the opposite direction with approx. 30 mL 0.1 mol/L H_2SO_4 at a flow rate of 0.1 mL/min.


Organic contaminants:

Rinse the column in the opposite direction with approx. 30 mL 0.01 mol/L H_2SO_4 /acetonitrile (80/20) at a flow rate of 0.1 mL/min.

Storage

For short periods (days) in the eluent, for longer periods (weeks) in ultrapure water. The column can be stored in a refrigerator at no colder than +4 °C.





50.0

Sulfuric acid eluent, standard

Tartrate 1 2 Malate

3 Succinate


4 Lactate

Conc. (mg/L) 25.0 5 Formate 50.0 6 Acetate 100.0 7 System peak

20.0

100.0

_



Oxalio	c acid eluent, sta	ndard		Cor	nc. (mg/L)
1	Acetate	10.0	5	Caproate	10.0
2	Propionate	10.0	6	Enantate	10.0
3	Butyrate	10.0	7	Octanate	10.0
4	Valerate	10.0			

Ordering information

Metrosep Organic Acids - 100/7.8 Metrosep Organic Acids Guard/4.6 6.1005.210 6.1005.250

Metrosep Organic Acids - 250/7.8 (6.1005.200)

124 The Metrosep Organic Acids - 250/7.8 is is a polymerbased cation-exchanger column. It is the high-performance column for the determination of organic acids and for the solution of difficult and complex separation problems. In addition, carbonate (with inverse suppression), fluoride (hydrofluoric acid), and phosphate (phosphoric acid) can be determined along with organic acids. In comparison with the Hamilton PRP-X300 - 250/4.0, the Metrosep Organic Acids column - 250/7.8 has greater capacity and enhanced selectivity.

Applications

- Organic acids: Citrate, tartrate, malate, ascorbate, succinate
- Short-chain fatty acids: Formate, acetate, propionate, butyrate, etc.
- F⁻, PO₄³⁻, CO₃²⁻
- Difficult matrices
- Difficult separation problems

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with sulfonic acid groups
Column dimensions	250 x 7.8 mm
Column body	Stainless steel
Standard flow	0.5 mL/min
Maximum flow	0.6 mL/min
Maximum pressure	7 MPa
Particle size	9 µm
Organic modifier	0–20%
pH range	1–13
Temperature range	5–90 °C

Eluent			
Sulfuric acid eluent	Sulfuric acid (c = 2 mol/L)	0.5 mL/2 L	0.5 mmol/L
(standard eluent)	Acetone	300 mL/2 L	15%

Care

Regeneration

Column purification: Rinse the column in the opposite direction with 20% acetonitrile in 0.01 mol/L H_2SO_4 at a flow rate of 0.1 mL/min for 4 hours at 65 °C.

Contaminations with metals:


If retention times are shortened: Rinse the column in the opposite direction with approx. 30 mL 0.1 mol/L $\rm H_2SO_4$ at a flow rate of 0.1 mL/min.

Organic contaminants:

Rinse the column in the opposite direction with approx. 30 mL 0.01 mol/L H_2SO_4 /acetonitrile (80/20) at a flow rate of 0.1 mL/min.

Storage

For short periods (days) in the eluent, for longer periods (weeks) in ultrapure water. The column can be stored in a refrigerator at no colder than +4 $^{\circ}$ C.



Ordering information Metrosep Organic Acids - 250/7.8 Metrosep Organic Acids Guard/4.6

6.1005.200 6.1005.250



Separation columns

IC carbohydrate-separation columns – anionexchange chromatography applying pulsed amperometric detection (PAD)

Metrosep Carb 2 - 100/4.0 (6.1090.410)

128

The Metrosep Carb 2 - 100/4.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The high-capacity anion-exchange column is based on a styrene/divinylbenzene copolymer. It is stable in the range of pH = 0-14 and provides separation of glucose, fructose, sucrose and lactose. It is also suitable for the analysis of some sugar alcohols and oligosaccharides. Short analysis times can be achieved on the 100 mm version of the Metrosep Carb 2 separation column.

Applications

- Monosaccharides
- Disaccharides
- Sugar alcohols
- Oligosaccharides
- Simple separation problems
- Very rapid separations

Technical information

copolymer with quaternary ammonium groups Column dimensions 100 x 4.0 mm Column body PEEK Standard flow 0.8 mL/min Maximum flow 1.6 mL/min Maximum pressure 20 MPa Particle size 5.0 µm	Substrate	Styrene/divinylbenzene
Column dimensions100 x 4.0 mmColumn bodyPEEKStandard flow0.8 mL/minMaximum flow1.6 mL/minMaximum pressure20 MPa		copolymer with quaternary
Column bodyPEEKStandard flow0.8 mL/minMaximum flow1.6 mL/minMaximum pressure20 MPa		ammonium groups
Standard flow0.8 mL/minMaximum flow1.6 mL/minMaximum pressure20 MPa	Column dimensions	100 x 4.0 mm
Maximum flow1.6 mL/minMaximum pressure20 MPa	Column body	PEEK
Maximum pressure 20 MPa	Standard flow	0.8 mL/min
	Maximum flow	1.6 mL/min
Particle size 5.0 µm	Maximum pressure	20 MPa
	Particle size	5.0 µm
Organic modifier In the eluent: 0–50 %	Organic modifier	In the eluent: 0–50 %
acetonitrile or methanol		acetonitrile or methanol
In the sample: 0–100 %		In the sample: 0–100 %
acetone, acetonitrile or		acetone, acetonitrile or
methanol		methanol
pH range 0–14	pH range	0–14
Temperature range 20–60 °C	Temperature range	20-60 °C

Eluent

Hydroxide/acetate eluent	Sodium hydroxide (c = 20 mol/L)	10 mL/2 L	100 mmol/L
(standard eluent)	Sodium acetate	1640.7 mg/2 L	10 mmol/L

Note

- 1. Use a flow ramp to establish the standard flow in the column within 5 min.
- 2. Rinse the column with the desired eluent for 2 h at 30 °C.

Care

Organic contamination:


Rinse the column in the flow direction with standard eluent in 50% acetonitrile at a flow rate of 0.5 mL/min for 3 h.

Inorganic contamination:

Rinse the column in the flow direction with a mixture of 100 mmol/L sodium hydroxide and 500 mmol/L sodium acetate at a flow rate of 0.5 mL/min for at least 3 h.

After regeneration, rinse the column with standard eluent for at least 3 h.

Storage In the standard eluent



6.1090.410
6.1090.500
6.1090.510

Metrosep Carb 2 - 150/4.0 (6.1090.420)

130

The Metrosep Carb 2 - 150/4.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The ainon-exchange column is based on a styrene/ divinylbenzene copolymer. It is stable in the range of pH = 0-14 and provides separation of monosaccharides and disaccharides. It is also suitable for the analysis of sugar alcohols, anhydrous sugars, oligosaccharides, etc. The column capacity has been optimized to enable the combination of rapid separations and excellent separation properties.

Applications

- Monosaccharides
- Disaccharides
- Sugar alcohols
- Anhydrosugars
- Oligosaccharides
- Rapid separations

Technical information

Substrate	Styrene/divinylbenzene
	copolymer with
	quaternary ammonium
	groups
Column dimensions	150 x 4.0 mm
Column body	PEEK
Standard flow	0.5 mL/min
Maximum flow	1.2 mL/min
Maximum pressure	20 MPa
Particle size	5.0 µm
Organic modifier	In the eluent: 0–50 %
	acetonitrile or methanol
	In the sample: 0–100 %
	acetone, acetonitrile or
	methanol
pH range	0–14
Temperature range	20-60 °C

Eluent

Hydroxide/acetate eluent	Sodium hydroxide (c = 20 mol/L)	10 mL/2 L	100 mmol/L
(standard eluent)	Sodium acetate	1640.7 mg/2 L	10 mmol/L
Hydroxide eluent	Sodium hydroxide (c = 20 mol/L)	1.0 mL/2 L	20 mmol/L
Hydroxide/acetate eluent	Sodium hydroxide (c = 20 mol/L)	0.5 mL/2 L	5 mmol/L
(modified)	Sodium acetate	328.1 mg/2 L	2 mmol/L
Hydroxide eluent (modified)	Sodium hydroxide (c = 20 mol/L)	5.0 mL/2 L	100 mmol/L

Note

- 1. Use a flow ramp to establish the standard flow in the column within 5 min.
- 2. Rinse the column for 2 h at 30 $^{\circ}\mathrm{C}$ with the desired eluent.

Care

Organic contamination:

Rinse the column in the flow direction with standard eluent in 50% acetonitrile at a flow rate of 0.5 mL/min for 3 h.


Inorganic contamination:

Rinse the column in the flow direction with a mixture of 100 mmol/L sodium hydroxide and 500 mmol/L sodium acetate at a flow rate of 0.5 mL/min for at least 5 h.

After regeneration, rinse the column with standard eluent for at least 5 h.


Storage In the standard eluent






Hydroxide eluent, (modified), ASTM D 7591, Conc. (mg/kg) free and total glycerin in biodiesel

Total glycerin 98.15 1 Free glycerin 6.52 2




Hydroxide/acetate eluent, standard, 30 °C 1.0 1 Inositol 1.0 5 Xylose 2 Arabitol 1.0 1.0 6 Fructose 3 Sorbitol 1.0 7 Lactose 1.0

4	Glucose	1.0	8



Conc. (mg/L) Hydroxide eluent, standard, anhydrosugars in aerosols, 45 °C

1	Inositol	0.6	7	Mannosan	3.2
2	Erythritol	0.6	8	Galactosan	3.2
3	Arabitol	1.3	9	Rhamnose	3.2
4	Sorbitol	1.3	10	Glucose	3.2
5	Mannitol	1.3	11	Xylose	3.2
6	Levoglucosan	3.2	12	Sucrose	3.2



Sucrose

1.0

Hydroxide/acetate eluent, (mod.), standard, 35 °C

5.0

5.0

Galactose 1

N-acetyl-

nA

200

1600

2

3

Mannose

Ordering information Metrosep Carb 2 - 150/4.0 Metrosep Carb 2 Guard/4.0 Metrosep Carb 2 S-Guard/4.0

- glucosamine 20.0

Conc. (mg/L) 4 N-acetylgalactosamine 20.0

10 12 14 16 18

Hydroxide/acetate eluent, (mod.), lactose-free milk, diluted 1 : 100, Inline Dialysis spiked with 100 mg/L Lactose, 28 °C Conc. (mg/L)

					conc. (mg/L)
1	Galactose	n.q.	3	Lactose	100.0
2	Glucose	n.q.			

6	1090.420
6	1090.500
6	1090.510

Metrosep Carb 2 - 250/4.0 (6.1090.430)

The Metrosep Carb 2 - 250/4.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The high-capacity anion-exchange column is based on a styrene/divinylbenzene copolymer. It is stable in the range of pH = 0-14 and provides separation of mono-saccharides and disaccharides. It is also suitable for the analysis of sugar alcohols, anhydrosugars, amino sugars, etc. The 250 mm version of the Metrosep Carb 2 separation column is optimized for complex separations.

Applications

- Monosaccharides
- Disaccharides
- Sugar alcohols
- Anhydrosugars
- Difficult matrices
- Complex separations

Technical information

Substrate	Styrene/divinylbenzene copolymer with quaternary ammonium groups
Column dimensions	250 x 4.0 mm
Column body	PEEK
Standard flow	0.5 mL/min
Maximum flow	0.8 mL/min
Maximum pressure	20 MPa
Particle size	5.0 µm
Organic modifier	In the eluent: 0–50 %
	acetonitrile or methanol
	In the sample: 0–100 %
	acetone, acetonitrile or
	methanol
pH range	0–14
Temperature range	20–60 °C

Eluent

Hydroxide/acetate eluent	Sodium hydroxide (c = 20 mol/L)	10 mL/2 L	100 mmol/L
(standard eluent)	Sodium acetate	1640.7 mg/2 L	10 mmol/L
Hydroxide/acetate eluent	Sodium hydroxide (c = 20 mol/L)	0.5 mL/2 L	5 mmol/L
(modified)	Sodium acetate	328.1 mg/2 L	2 mmol/L

Note

- 1. Use a flow ramp to establish the standard flow in the column within 5 min.
- 2. Rinse the column with the desired eluent for 2 h at 30 °C.

Care

Organic contamination:


Rinse the column in the flow direction with standard eluent in 50% acetonitrile at a flow rate of 0.5 mL/min for 7 h.

Inorganic contamination:

Rinse the column in the flow direction with a mixture of 100 mmol/L sodium hydroxide and 500 mmol/L sodium acetate at a flow rate of 0.5 mL/min for at least 7 h.


After regeneration, rinse the column with standard eluent for at least 7 h.

Storage In the standard eluent



4

Glucose




5.0


8

Sucrose

15.0



Hydroxide/acetate eluent, (mod.), standard, 40 °C				Cor	ic. (mg/L)
1	Inositol	0.5	9	Fucose	2.0
2	Xylitol	0.5	10	Sucrose	2.0
3	Sorbitol	0.5	11	Galactose	2.0
4	Mannitol	0.5	12	Glucose	2.0
5	Lactitol	0.5	13	Mannose	2.0
6	Levoglucosan	2.0	14	Sorbose	5.0
7	Mannosan	2.0	15	Fructose	5.0
8	Galactosan	2.0	16	Lactose	5.0



6.1090.430 6.1090.500 6.1090.510 133

Ordering information

Metrosep Carb 2 - 250/4.0 Metrosep Carb 2 Guard/4.0 Metrosep Carb 2 S-Guard/4.0

Hamilton RCX-30 - 150/4.6 (6.1018.010)

The Hamilton RCX-30 - 150/4.6 is a column for the separation of monosaccharides, disaccharides, oligosaccharides, and sugar alcohols. It is an anion-exchange column based on polystyrene/divinylbenzene resin. The RCX-30 - 150/4.6 can be used for universal applications.

> The Hamilton RCX-30 - 150/4.6 separation column features an outstanding separation of fructose and lactose. The column also offers the advantage that flows of up to 2 mL/min can be used in order to accelerate the chromatography. The column is used both for the rapid separation of small carbohydrates and for the separation of oligosaccharides.

Applications

- Monosaccharides
- Disaccharides
- Oligosaccharides
- Sugar alcohols
- Simple separation problems
- Rapid separations

Technical information

Substrate	Styrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	150 x 4.6 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	34 MPa
Particle size	7 µm
pH range	1–13 (T>35 °C max. pH 8)
Temperature range	20–60 °C


Eluent

Hydroxide eluent	Sodium hydroxide (c = 20 mol/L)	15 mL/2 L	150 mmol/L
		Column temperature 32 °C	
Hydroxide eluent	Sodium hydroxide (c = 20 mol/L)	20 mL/2 L	200 mmol/L
(modified)		Column temperature 32 °C	

Care

Regeneration Rinse the column with 150 mL 0.1 mol/L NaOH at a flow rate of 1 mL/min. Storage In ultrapure water with 1 mmol/L sodium azide





5.00 6

7

8

10.00


10.00

Hydroxide eluent, standard

1 Inositol 2 Arabitol 3 Sorbitol 4 Fucose

4 Fucose 5 Arabinose

10.00 9 ose 10.00



Hydroxide eluent, modified, glycerin

Conc. (mg/L)

10.00

10.00

10.00

10.00

Glucose

Fructose

Sucrose

Lactose

in f	fermentation solu	tion		Cor	nc. (mg/L)
1	Glycerin	20.31	4	Trehalose	n.q.
2	Arabitol	n.q.	5	Glucose	n.q.
3	Sorbitol	n.q.			

Ordering	information
----------	-------------

Hamilton RCX-30 - 150/4.6 Metrosep RP 2 Guard/3.5 Replacement filters for RP 2 Guard/3.5 (10 pcs.) Metrosep RP 3 Guard HC/4.0 6.1018.010 6.1011.030 6.1011.130 6.1011.040

Hamilton RCX-30 - 250/4.6 (6.1018.000)

The Hamilton RCX-30 - 250/4.6 is a column for the separation of monosaccharides, disaccharides, and sugar alcohols. It is an anion-exchange column based on polystyrene/divinylbenzene resin. The RCX-30 - 250/4.6 can be used for universal applications.

The Hamilton RCX-30 - 250/4.6 separation column excels in an outstanding separation of fructose and lactose. The column also offers the advantage that flows of up to 2 mL/min can be used in order to accelerate the chromatography. The long version of the column (250 mm) is preferred for the determination of small carbohydrates (monosaccharides, disaccharides, and sugar alcohols).

Applications

- Monosaccharides
- Disaccharides
- Sugar alcohols
- Difficult separation problems
- Difficult matrices

Technical information

Substrate	Styrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	250 x 4.6 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	34 MPa
Particle size	7 µm
pH range	1–13 (T>35 °C max. pH 8)
Temperature range	20-60 °C

150 mmol/L

Eluent

Hydroxide eluent (standard eluent) Sodium hydroxide (c = 20 mol/L) 15 mL/2 L Column temperature 32 °C

Care

Regeneration Rinse the column with 150 mL 0.1 mol/L NaOH at a flow rate of 1 mL/min.

Storage In ultrapure water with 1 mmol/L sodium azide



Conc. (mg/L)

10.00

10.00

10.00

10.00


Glucose

Fructose

Lactose

Sucrose

Chromatograms



Hydroxide eluent, standard

1	Inositol	
2	Arabitol	
3	Sorbitol	
	_	


4 Fucose 5 Arabinose

10.00 8 10.00 9 e 10.00

5.00 | 6

7

10.00



1 Glucose 15 3 Sucrose 71 2 Fructose 11

Ordering information

Hamilton RCX-30 - 250/4.6 Metrosep RP 2 Guard/3.5 Replacement filters for RP 2 Guard/3.5 (10 pcs.) Metrosep RP 3 Guard HC/4.0 6.1018.000 6.1011.030 6.1011.130 6.1011.040



Separation columns

Microbore IC carbohydrate-separation columns for lower eluent consumption and greater sensitivity

Metrosep Carb 2 - 100/2.0 (6.01090.210)

The Metrosep Carb 2 - 100/2.0 IC column is the short microbore version of the Metrosep Carb 2 columns and is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The high-capacity anion exchanger column is based on a styrene-divinylbenzene copolymer. It is stable in the range of pH = 0-14 and provides separation of glucose, fructose and sucrose. It is also suitable for the analysis of some sugar alcohols and oligosaccharides. Short analysis times can be achieved on the 100 mm version of the Metrosep Carb 2 separation column.

With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Applications

- Monosaccharides
- Disaccharides
- Sugar alcohols
- Oligosaccharides
- Simple separation problems
- Very rapid separations
- Anions in sea water
- IC-MS

	Technical information	
t-	Substrate	Polystyrene-divinylbenzene
		copolymer with quaternary
		ammonium groups
	Column dimensions	100 x 2.0 mm
	Column body	PEEK
	Standard flow	0.2 mL/min
	Maximum flow	0.7 mL/min
	Maximum pressure	20 MPa
	Particle size	5.0 µm
	Organic modifier	In the eluent: 0–50 %
		acetonitrile or methanol
		In the sample: 0–100 %
		acetone, acetonitrile or
		methanol
	pH range	0–14
	Temperature range	20-60 °C
	10 mL/2 L	100 mmol/L
	IV IIIL/Z L	
	1640.7 mg/2 L	10 mmol/L
	20 a/2 l	10 a/l

Eluents

Eluents			
Hydroxide/acetate eluent (standard eluent)	Sodium hydroxide (c = 20 mol/L)	10 mL/2 L	100 mmol/L
	Sodium acetate	1640.7 mg/2 L	10 mmol/L
Sodium chloride eluent	Sodium chloride	20 g/2 L	10 g/L
Ammonium nitrate eluent	Ammonium nitrate	16.0 g/2 L	100 mmol/L
	Ammonium hydroxide		pH = 9.0

Care

Regeneration Note:

- 1. Use a flow ramp to establish the standard flow in the column within 5 min.
- 2. Rinse the column with the desired eluent for 2 h at 30 °C.

Organic contamination:


Rinse the column in the flow direction with 25 mL of solution (standard eluent in 50% acetonitrile) at a flow rate of 0.13 mL/min.

Inorganic contamination:

Rinse the column in the flow direction with a mixture of 100 mmol/L sodium hydroxide and 500 mmol/L sodium acetate at a flow rate of 0.13 mL/min for at least 7 h.

After regeneration, rinse the column with standard eluent for at least 7 h.

Storage In the standard eluent



Ordering information	
Metrosep Carb 2 - 100/2.0	6.01090.210
Metrosep Carb 2 Guard/2.0	6.01090.600
Metrosep Carb 2 S-Guard/2.0	6.01090.610

Metrosep Carb 2 - 150/2.0 (6.01090.220)

142

The microbore version of the Metrosep Carb 2 - 150/2.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The anion exchanger column is based on a styrene-divinylbenzene copolymer. It is stable in the range of pH = 0-14 and provides separation of monosaccharides and disaccharides. It is also suitable for the analysis of sugar alcohols, anhydrous sugars, oligosaccharides, etc. The column capacity has been optimized to enable the combination of rapid separations and excellent separation properties.

Thanks to its low eluent consumption, it is particularly suitable for IC-MS coupling.

Applications

- Monosaccharides
- Disaccharides
- Sugar alcohols
- Oligosaccharides
- Anhydrous sugars
- Rapid separations
- IC-MS

Technical information

Substrate	Polystyrene-divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	150 x 2.0 mm
Column body	PEEK
Standard flow	0.13 mL/min
Maximum flow	0.45 mL/min
Maximum pressure	20 MPa
Particle size	5.0 µm
Organic modifier	In the eluent: 0–50 %
	acetonitrile or methanol
	In the sample: 0–100 %
	acetone, acetonitrile or
	methanol
pH range	0–14
Temperature range	20-60 °C

Eluents

Eluents				
Hydroxide/acetate eluent	Sodium hydroxide	10 mL/2 L	100 mmol/L	
(standard eluent)	(c = 20 mol/L)			
	Sodium acetate	1640.7 mg/2 L	10 mmol/L	
Hydroxide/acetate eluent	Sodium hydroxide	0.5 mL/2 L	5 mmol/L	
(modified)	(c = 20 mol/L)			
	Sodium acetate	328.1 mg/2 L	2 mmol/L	

Care

Regeneration Note:

- 1. Use a flow ramp to establish the standard flow in the column within 5 min.
- 2. Rinse the column with the desired eluent for 2 h at 30 °C.

Organic contamination:

Rinse the column in the flow direction with 25 mL of solution (standard eluent in 50% acetonitrile) at a flow rate of 0.13 mL/min.


Inorganic contamination:

Rinse the column in the flow direction with a mixture of 100 mmol/L sodium hydroxide and 500 mmol/L sodium acetate at a flow rate of 0.13 mL/min for at least 7 h.

After regeneration, rinse the column with standard eluent for at least 7 h.


Storage In the standard eluent






Hydroxide/acetate eluent, (mod.), yoghurt,

1:1	000 diluted, 40 °C				Conc. (g/L)
1	Sucrose	64.5	4	Fructose	15.2
2	Galactose	4.2	5	Lactose	30.6
3	Glucose	10.0			





1	Sucrose	5.0	3	unknown	-
2	Glucose	26.8	4	Fructose	59.4

Ordering information

Metrosep Carb 2 - 150/2.0 Metrosep Carb 2 Guard/2.0 Metrosep Carb 2 S-Guard/2.0

6.01090.220 6.01090.600 6.01090.610

Metrosep Carb 2 - 250/2.0 (6.01090.230)

The Metrosep Carb 2 - 250/2.0 IC column is particularly suitable for the determination of carbohydrates using alkaline eluents and pulsed amperometric detection. The high-capacity anion exchanger column is based on a styrol-divinylbenzene copolymer. It is stable in the range of pH = 0–14 and provides separation of monosaccharides and disaccharides. It is also suitable for the analysis of sugar alcohols, anhydrous sugars, amino sugars, etc. The 250 mm microbore version of the Metrosep Carb 2 separation column is optimized for complex separations.

Thanks to its low eluent consumption, it is particularly suitable for IC-MS coupling.

Applications

- Monosaccharides
- Disaccharides
- Sugar alcohols
- Anhydrous sugars
- Oligosaccharides
- Difficult matrices
- Complex separations
- IC-MS

Technical information		
Substrate	Polystyrene-divinylbenzene	
	copolymer with quaternary	
	ammonium groups	
Column dimensions	250 x 2.0 mm	
Column body	PEEK	
Standard flow	0.13 mL/min	
Maximum flow	0.30 mL/min	
Maximum pressure	20 MPa	
Particle size	5.0 µm	
Organic modifier	In the eluent: 0–50 %	
	acetonitrile or methanol	
	In the sample: 0–100 %	
	acetone, acetonitrile or	
	methanol	
pH range	0–14	
Temperature range	20–60 °C	
10	100	
10 mL/2 L	100 mmol/L	

Eluents

Eluents			
Hydroxide/acetate eluent (standard eluent)	Sodium hydroxide (c = 20 mol/L)	10 mL/2 L	100 mmol/L
	Sodium acetate	1640.7 mg/2 L	10 mmol/L
Hydroxid eluent	Sodium hydroxide	1.0 mL/2 L	10 mmol/L

Care

Regeneration

Note:

- 1. Use a flow ramp to establish the standard flow in the column within 5 min.
- 2. Rinse the column with the desired eluent for 2 h at 30 °C.

Organic contamination:

Rinse the column in the flow direction with 25 mL of solution (standard eluent in 50% acetonitrile) at a flow rate of 0.13 mL/min.

Inorganic contamination:


Rinse the column in the flow direction with a mixture of 100 mmol/L sodium hydroxide and 500 mmol/L sodium acetate at a flow rate of 0.13 mL/min for at least 7 h.


After regeneration, rinse the column with standard eluent for at least 7 h.

Storage In the standard eluent




Chromatograms





Ordering information

Metrosep Carb 2 - 250/2.0 Metrosep Carb 2 Guard/2.0 Metrosep Carb 2 S-Guard/2.0 6.01090.230 6.01090.600 6.01090.610



Separation columns

IC amino acid-separation column with optical detection (VIS) after post-column reaction

Metrosep Amino Acids 1 - 100/4.0 (6.4001.410)

148

The Metrosep Amino Acids 1 - 100/4.0 is the standard separation column for amino acids. The column is based on a sulfonated polystyrene-divinylbenzene material. The determination of amino acids is accomplished by means of photometric detection following a post-column reaction with ninhydrin.

The Metrosep Amino Acids 1 - 100/4.0 permits the separation of up to 44 amino acids in research and routine applications including all naturally occurring amino acids.

Lithium citrate

Lithium citrate

Lithium chloride

Phenol HCl

Phenol HCl

Ninhydrin Hydrindantin Dimethyl sulfoxide Lithium acetate buffer (2 mol/L, pH = 5.2 with acetic acid)

B:

Applications

• Amino acids

Ammonium

ans		
eac-	Technical information	
pa-	Substrate	Sulfonated polystyrene- divinylbenzene copolymer, lithium form
tine	Column dimensions	100 x 4.0 mm
ids.	Column body	Stainless steel
	Standard flow	0.4 mL/min
	Maximum flow	0.5 mL/min
	Maximum pressure	10 MPa
	Particle size	5 µm
	Organic modifier	10% acetonitrile,
		0–5% other organic
		solvents
	pH range	1–14
	Temperature range	30–90 °C
	Capacity	2.9 mmol (K ⁺)
	17.8 g/2 L	42.6 mmol/L
	2.0 g/2L	10.6 mmol/L
	2.0 g/ 22	pH = 2.8
	17.8 g/2L	42.6 mmol/L
	86.0 g/2L	1.0 mol/L
	2.0 g/2L	10.6 mmol/L
	5	pH = 4.2
	Column temperature 50 °C	
		0.44
	4.0 g/200 mL	0.11 mol/L
	0.16 g/200 mL	2.5 mmol/L 100 mL
		100 mL
	Reactor temperature 120 °C	

Care

Eluents

Gradient

Citrate/phenol eluent A:

Regeneration

PCR reagents Ninhydrin

In the event of temporary loss of column performance:

• Apply fresh eluent, rinse the instrument and column for 1 h at 0.20 mL/min at 65 °C

For minor contaminations:

• 120 min 0.3 mol/L lithium hydroxide with 0.25 g/L EDTA (0.20 mL/min, 90 °C)

With contaminations caused by organic components: Rinse the column with the following solutions in sequence (0.2 mL/min, 65 $^{\circ}$ C):


- 30 min ultrapure water
- 60 min 20% acetonitrile/water
- 60 min ultrapure water to completely remove the acetonitrile

Storage

Short-term: Storage in eluent with 2.5% acetonitrile Long-term: Storage in 0.3 mol/L lithium hydroxide with 5% acetonitrile.



Chromatogram



2	L-serine	2.5	11	L-isoleucine	2.5
3	L-threonine	2.5	12	L-leucine	2.5
4	L-glutamic acid	2.5	13	L-tyrosine	2.5
5	L-proline	2.5	14	L-phenylalanine	2.5
6	Glycine	2.5	15	Ammonium	2.5
7	L-alanine	2.5	16	L-lysine	2.5
8	L-valine	2.5	17	L-histidine	2.5
9	L-cystine	1.25	18	L-arginine	2.5

Ordering information	
Metrosep Amino Acids 1 - 100/4.0	6.4001.410
Metrosep RP 2 Guard/3.5	6.1011.030
Replacement filters for RP 2 Guard/3.5 (10 pcs.)	6.1011.130
Metrosep RP 3 Guard HC/4.0	6.1011.040



Separation columns



IC cation-separation columns for analyses without chemical suppression

Nucleosil 5SA - 125/4.0 (6.1007.000)

152

The Nucleosil 5SA IC cation column uses sulfonic acid groups for the separation of cations. With eluents containing organic acids and ethylenediamine, this column separates divalent cations such as magnesium and calcium as well as some of the transition metal elements (e.g. nickel, zinc, cobalt, manganese). The Nucleosil 5SA -125/4.0 is therefore the inexpensive and robust separation column for the determination of transition metals by direct conductivity measurement without post-column reaction. In addition to high concentrations of alkaline metals, calcium, and magnesium can be determined reliably. The column is only suitable for divalent cations. Monovalent cations elute at almost the same time as the injection peak.

Applications

- Mg^{2+} , Ca^{2+} , Sr^{2+} , Ba^{2+} , Fe^{2+} , Co^{2+} , Ni^{2+} , Cd^{2+} , Zn^{2+} , Mn^{2+}
- Mg^{2+} , Ca^{2+} in addition to a high amount of sodium

Technical information	
Substrate	Spherical silica gel with
	sulfonic acid groups
Column dimensions	125 x 4.0 mm
Column body	Stainless steel
Standard flow	1.5 mL/min
Maximum flow	5.0 mL/min
Maximum pressure	40 MPa
Particle size	5 µm
pH range	2–8
Temperature range	0-40°C
Capacity	95 µmol (K⁺)

Eluents

Tartaric acid/	Tartaric acid	1200 mg/2 L	4.0 mmol/L
citric acid eluent	Citric acid	192 mg/2 L	0.5 mmol/L
(standard eluent)	Ethylenediamine	360 mg/2 L	3.0 mmol/L
	Acetone	100 mL/2 L	5%


Care

Regeneration

Injection of 100 μL Na_2H_2EDTA (0.1 mol/L) – do not use alkaline EDTA solutions – or rinse with 30 mL HNO₃ (0.1 mol/L) at a flow rate of 0.5 mL/min.

Storage

For short periods (days) in the eluent, for longer periods (weeks) in methanol/water (1:4).



Ordering information Nucleosil 5SA - 125/4.0

Nucleosil 55A 2 Guard cartridge/4.0 Holder to Nucleosil 55A 2 Guard Cartridge/4.0 (holder for guard column cartridges 6.1007.110) 6.1007.000 6.1007.110 6.2821.140

Metrosep C 3 - 100/4.0 (6.1010.410)

154 The substrate on a polyvinyl alcohol base increases selectivity for monovalent and divalent cations significantly on this cation column. The peak forms are highly symmetrical.

The shortest separation column of the Metrosep-C-3 product range is particularly suitable for rapid separations of standard cations and for the separation of larger organic amines.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺, Mn²⁺, Co²⁺, Zn²⁺, Ni²⁺
- Larger organic amines
- Low detection limits
- Matrices with high pH

Technical information

Substrate	Polyvinyl alcohol with
	carboxyl groups
Column dimensions	100 x 4.0 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	1.5 mL/min
Maximum pressure	15 MPa
Particle size	5 µm
Organic modifiers	0–50% acetonitrile,
	0–30% acetone,
	no methanol
pH range	2–12
Temperature range	20-40 °C
Capacity	12 µmol (K ⁺)

Eluents

Nitric acid eluent (standard eluent) Nitric acid (c = 1 mol/L)


.) 10 mL/2 L Column temperature 40 °C 5 mmol/L

Care

Regeneration Add 30% acetonitrile to the standard eluent. Storage For 1–3 days in the eluent; in ultrapure water for longer storage. Recommended temperature: 4–8 °C



Chromatograms



Nitric acid eluent, standard

1 Lithium 1.00 5.00

Sodium 2 3 Ammonium


Conc. (mg/L) | 4 Potassium 5 Magnesium 5.00 Calcium 6

10.00

10.00

10.00

2



Sodium	5.83	3	Magnesium	18.91
Potassium	1.45	4	Calcium	87.51

Ordering information Metrosep C 3 - 100/4.0

Metrosep C 3 Guard/4.0 Metrosep C 3 S-Guard/4.0 6.1010.410 6.1010.450 6.1010.460

Metrosep C 3 - 150/4.0 (6.1010.420)

The substrate on a polyvinyl alcohol base increases selectivity for monovalent and divalent cations significantly. The peak forms on this cation column are highly symmetrical.

The middle separation column of the Metrosep-C-3 product range is particularly suitable for rapid separations of standard cations and certain transition metal cations as well as for the separation of mid-sized organic amines.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺, Mn²⁺, Co²⁺, Zn²⁺, Ni²⁺
- Organic amines
- Low detection limits
- Matrices with high pH

Technical information

Substrate	Polyvinyl alcohol with
	carboxyl groups
Column dimensions	150 x 4.0 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	1.5 mL/min
Maximum pressure	15 MPa
Particle size	5 µm
Organic modifiers	0–50% acetonitrile,
	0-30% acetone,
	no methanol
pH range	2-12
Temperature range	20–40 °C
Capacity	18 µmol (K+)


Eluents

Nitric acid eluent	Nitric acid (c = 1 mol/L)	10 mL/2 L	5 mmol/L
(standard eluent)		Column temperature 40 °C	
Nitric acid eluent	Nitric acid (c = 1 mol/L)	5 mL/2 L	2.5 mmol/L
(modified)		Column temperature 40 °C	

Care

Regeneration Add 30% acetonitrile to the standard eluent. Storage For 1–3 days in the eluent; in ultrapure water for longer storage. Recommended temperature: 4–8 °C





5

6


7

Nitric acid eluent, modified,

traces of cations, 40 °C				
1	Lithium	0.050		
2	Sodium	0.050		
3	Ammonium	0.050		
4	Monoethylamine	0.100		

Cond	c. (μg/L)
Potassium	0.050
Diethylamine	0.100
Triethylamine	0.100

Monoethy	amine	0.10



5.00 6

		Cor
1.00	4	Potassium
5.00	5	Magnesium

Calcium


10.00

10.00

10.00

1

2



	5			- () /
Sodium	5.86	3	Magnesium	18.90
Potassium	1.41	4	Calcium	87.48

Ordering information

Lithium

Sodium

Ammonium

1

2

3

Metrosep C 3 - 150/4.0 Metrosep C 3 Guard/4.0 Metrosep C 3 S-Guard/4.0 6.1010.420 6.1010.450 6.1010.460

Metrosep C 3 - 250/4.0 (6.1010.430)

The substrate on a polyvinyl alcohol base increases selectivity for monovalent and divalent cations significantly. A characteristic of this is the number of «theoretical plates per meter». On the Metrosep C 3 - 250/4.0, for example, 42,000 plates are achieved for sodium, 51,000 for ammonium, and 31,000 for barium with its delayed elution. The peak forms on this cation column are highly symmetrical.

The selectivity of the Metrosep C 3 - 250/4.0 also permits the separation of transition metals. Because Metrohm ion chromatographs generally determine the cations without chemical suppression, the transition metals can be analyzed on the Metrosep C 3 - 250/4.0 together with the alkaline and earth alkaline metals.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺, Mn²⁺, Co²⁺, Zn²⁺, Ni²⁺
- Good Na⁺/NH₄⁺ separation
- Low detection limits
- Matrices with high pH

Technical information

Polyvinyl alcohol with carboxyl groups
250 x 4.0 mm
PEEK
1.0 mL/min
1.5 mL/min
15 MPa
5 µm
0–50% acetonitrile,
0–30% acetone,
no methanol
2–12
20-40 °C
30 µmol (K⁺)

Eluents


Nitric acid eluent	Nitric acid (c = 1 mol/L)	10 mL/2 L	5 mmol/L
(standard eluent)		Column temperature 40 °C	
Nitric acid/crown ether	Nitric acid (c = 1 mol/L)	7 mL/2 L	3.5 mmol/L
eluent	Crown ether 18-crown-6	264 mg/2 L	0.5 mmol/L
		Column temperature 40 °C	

Care

Regeneration Add 30% acetonitrile to the standard eluent.

Storage

For 1–3 days in the eluent; in ultrapure water for longer storage. Recommended temperature: 4–8 °C



Ordering information	
Metrosep C 3 - 250/4.0	
Metrosep C 3 Guard/4.0	
Metrosep C 3 S-Guard/4.0	

6.1010.430 6.1010.450 6.1010.460

Metrosep C 4 - 50/4.0 (6.1050.450)

The Metrosep C 4 - 50/4.0 is the shortest separation column in the Metrosep-C-4 product range. With a capacity of 5 µmol (K⁺), it is particularly suitable for very rapid separations. The low capacity makes it possible to quickly analyze the earth alkaline metals with their delayed elution. Thanks to the short retention times, applications that, in terms of analysis duration, were previously possible only with an FIA system (Flow Injection Analysis system) can now be transferred over to ion chromatography.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Alkylamines
- Very rapid separations
- Simple sample matrices

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	50 x 4.0 mm
Column body	PEEK
Standard flow	0.9 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2–7
Temperature range	20–60 °C
Capacity	5 µmol (K⁺)

Eluents

160

Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L
dipicolinic acid eluent	Dipicolinic acid	234 mg/2 L	0.7 mmol/L
(standard eluent)			
Nitric acid eluent	Nitric acid (c = 1 mol/L)	4 mL/2 L	2.0 mmol/L
(modified)			

Care


Regeneration

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.9 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.9 mL/min.

Storage

In the eluent or in ultrapure water



Ordering information

Metrosep C 4 - 50/4.0 Metrosep C 4 Guard/4.0 Metrosep C 4 S-Guard/4.0 6.1050.450 6.1050.500 6.1050.510

Metrosep C 4 - 100/4.0 (6.1050.410)

162

The 100 mm version of the Metrosep C 4 column is intended for rapid determinations of the standard cations. Very short retention times are achieved, for which the elution times of sodium and ammonium nevertheless differ by 25 s. When a special eluent is used, the six cations lithium, ammonium, sodium, calcium, magnesium, and potassium can be determined in less than 5 minutes with the Metrosep C 4 - 100/4.0.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Lipophilic amines with short retention times
- Rapid separations

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	100 x 4.0 mm
Column body	PEEK
Standard flow	0.9 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2–7
Temperature range	20–60 °C
Capacity	10 µmol (K ⁺)

Eluents

Nitric acid/ dipicolinic acid eluent (standard eluent)	Nitric acid (c = 1 mol/L) Dipicolinic acid	3.4 mL/2 L 234 mg/2 L	1.7 mmol/L 0.7 mmol/L
Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L
dipicolinic acid/acetone	Dipicolinic acid	234 mg/2 L	0.7 mmol/L
eluent	Acetone	100 mL/2 L	5%
Nitric acid/	Nitric acid (c = 1 mol/L)	4.0 mL/2 L	2.0 mmol/L
dipicolinic acid eluent (modified)	Dipicolinic acid	401 mg/2 L	1.2 mmol/L

Care

Regeneration

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.9 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.


Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.9 mL/min.

Storage

In the eluent or in ultrapure water

Chromatograms

2




Nitric acid/dipicolinic acid eluent, standard


1 Lithium 2 Sodium


Ammonium

3

- 1.00 4 5 5.00 5.00 6
- 10.00 Potassium Calcium 10.00 Magnesium 10.00







Nitric acid/dipicolinic acid/acetone eluent,

			••,	
tandard			Cor	nc. (mg/L)
Sodium	5.00	4	Guanidine	15.00
Ammonium	5.00	5	Aminoguanidi	ne 15.00
Methylamine	5.00			

6.1050.410
6.1050.500
6.1050.510
6.1050.530

Ordering information

Metrosep C 4 - 100/4.0 Metrosep C 4 Guard/4.0 Metrosep C 4 S-Guard/4.0 Metrosep C 4 S-Guard - 50/4.0

Metrosep C 4 - 150/4.0 (6.1050.420)

164

The Metrosep C 4 - 150/4.0 is the universal standard column in cation analysis. High separating efficiency in a brief time. The Metrosep C 4 - 150/4.0 is the ideal separation column for the analysis of alkaline and earth alkaline metals in aqueous media.

Applications

- Standard column
- Amines
- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Universal applications
- Different matrices

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	150 x 4.0 mm
Column body	PEEK
Standard flow	0.9 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2–7
Temperature range	20-60 °C
Capacity	15 µmol (K+)

Eluents

Nitric acid/ dipicolinic acid eluent (standard eluent)	Nitric acid (c = 1 mol/L) Dipicolinic acid	3.4 mL/2 L 234 mg/2 L	1.7 mmol/L 0.7 mmol/L
Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L
dipicolinic acid/	Dipicolinic acid	234 mg/2 L	0.7 mmol/L
crown ether eluent	18-crown-6	26.4 mg/2 L	0.05 mmol/L

Care


Regeneration

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.9 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L HNO $_3$ + 4 mmol/L dipicolinic acid for 1 h at a flow rate of 0.9 mL/min.

Storage

In the eluent or in ultrapure water



n.q. 3 4

5

Monoethanolamine 30.0

Diethanolamine

9

30.0

1.13

min

10.0

10.0

10.0

Ordering information Metrosep C 4 - 150/4.0 Metrosep C 4 Guard/4.0 Metrosep C 4 S-Guard/4.0 Metrosep C 4 S-Guard - 50/4.0

6.1050.420 6.1050.500 6.1050.510 6.1050.530

Magnesium

Metrosep C 4 - 250/4.0 (6.1050.430)

The Metrosep C 4 - 250/4.0 is the cation column with the greatest capacity in the C 4 series. It is predestined for applications which require the highest separating efficiency. Samples with extreme differences in concentrations can be analyzed reliably with this column. The performance capability of the column is demonstrated, for example, in connection with the analysis of boiler feed water for which the requirement is the perfect quantification of 7 μ g/L sodium in addition to 7 mg/L monoethanolamine (MEA). With the C 4 - 250/4.0, not only amines and transition metals but also alkaline and alkaline earth metals can be determined in a single run.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺, Co²⁺, Ni²⁺, Zn²⁺, Cd²⁺, Pb²⁺, amines
- Very good Na⁺/NH₄⁺ separation
- NH_4^+ , $(CH_3)NH_3^+$, $(CH_3)_2NH_2^+$, $(CH_3)_3NH^+$, $(CH_3)_4N^+$, and the respective ethanolamines
- Difficult separation problems
- Great differences in concentration
- Transition metals

Technical information	
Substrate	Silica gel with
	carboxyl groups
Column dimensions	250 x 4.0 mm
Column body	PEEK
Standard flow	0.9 mL/min
Maximum flow	2.0 mL/min
Maximum pressure	25 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2-7
Temperature range	20–60 °C
Capacity	25 µmol (K⁺)

Eluents

166

Eluents				
Nitric acid/	Nitric acid ($c = 1 \text{ mol/L}$)	3.4 mL/2 L	1.7 mmol/L	
dipicolinic acid eluent	Dipicolinic acid	234 mg/2 L	0.7 mmol/L	
(standard eluent)				
Amine eluent	Nitric acid ($c = 1 \text{ mol/L}$)	3.4 mL/2 L	1.7 mmol/L	
	Dipicolinic acid	234 mg/2 L	0.7 mmol/L	
	18-crown-6	26.4 mg/2 L	0.05 mmol/L	
	Acetone	25 mL/2 L	2.5%	

Care


Regeneration

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.9 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4 \text{ mmol/L}$ dipicolinic acid for 1 h at a flow rate of 0.9 mL/min.

Storage

In the eluent or in ultrapure water



Ordering information	
Metrosep C 4 - 250/4.0	6.1050.430
Metrosep C 4 Guard/4.0	6.1050.500
Metrosep C 4 S-Guard/4.0	6.1050.510
Metrosep C 4 S-Guard - 50/4.0	6.1050.530

Metrosep C 5 - 150/4.6 (6.4000.320)

The Metrosep C 5 - 150/4.6 is based on a sulfonated polystyrene/divinylbenzene polymer. The strongly acidic cation-exchanger groups make it a preferred separation column for the determination of divalent cations, particularly of transition metals. The column is preferably also used with UV/VIS detection after post-column reaction. This column may also be used with sequential suppression.

Applications

• Special column for transition metals

Technical information

Substrate	Sulfonated polystyrene/
	divinylbenzene polymer
Column dimensions	150 x 4.6 mm
Column body	PEEK
Standard flow	1.0 mL/min
Maximum flow	1.5 mL/min
Maximum pressure	4.9 MPa
Particle size	12 µm
Organic modifier	10% acetonitrile,
	0–5 % other organic
	solvents
pH range	1–14 (optimal 2–6)
Temperature range	20–70 °C
Capacity	15 µmol (K ⁺)

Eluents

Oxalic acid/citric acid	Oxalic acid	1.080 g/2 L	6.0 mmol/L
(standard eluent)	Citric acid	1.153 g/2 L	3.0 mmol/L
	КОН		pH = 4.2
PCR reagents			
PAR	PAR (4(2-pyridylazo) resorcinol	64.6 mg/2 L	0.15 mmol/L
	Ammonium hydroxide (c = 1 mol/L)	800 mL/2 L	0.4 mol/L
	Nitric acid (c = 1 mol/L)	160 mL/2 L	80 mmol/L
			pH = 10–11

Care

Regeneration

Slight contamination (e.g. divalent cations): Operate the column with eluent (0.5 mL/min, room temperature) and inject 100 μ L 1 mol/L nitric acid 4–6 times in succession.

More extreme contamination: In the event of heavier contamination, rinse the column with the following solutions in sequence (0.5 mL/min, room temperature):

- 60 min (30 mL) 100 mmol/L tartaric acid
- 60 min (30 mL) 100 mmol/L nitric acid
- Eluent for rinsing the column

Contamination from proteins and nitrogen components: Rinse the column with the following solutions in sequence (0.5 mL/min, room temperature):


- 30 min ultrapure water
- 60 min (30 mL) 100 mmol/L sodium hydroxide
- 60 min (30 mL) 100 mmol/L nitric acid
- Eluent for rinsing the column

Storage

Immediately after use, rinse to 3.0 mmol/L nitric acid.



Chromatogram



Oxalic acid/citric acid eluent, PCR with PAR,

VIS detection ($\lambda = 530$ nm), standard			Cone	c. (mg/L)	
1	Copper	5.00	5	Lead	30.0
2	Nickel	3.00	6	Manganese	4.00
3	Zinc	4.00	7	Cadmium	8.00
4	Cobalt	5.00			

Ordering information	
Metrosep C 5 - 150/4.6	6.4000.320
Metrosep RP 2 Guard/3.5	6.1011.030
Replacement filters for RP 2 Guard/3.5 (10 pcs.)	6.1011.130
Metrosep RP 3 Guard HC/4.0	6.1011.040

169

Metrosep BP 1 Guard/2.0

6.1015.100

Metrosep C 6 - 100/4.0 (6.1051.410)

170

The 100 mm version of the Metrosep C 6 column is intended for the determination of standard cations, e.g., in drinking water. Excellent separation of sodium and ammonium is still achieved, even with the very short retention times applied. The high capacity of the C ${\rm 6}$ mater

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Lipophilic amines with short retention times
- Rapid separations

material permits larger sa	mple volumes.	Technical information	hnical information		
		Substrate	Silica gel with		
			carboxyl groups		
		Column dimensions	100 x 4.0 mm		
		Column body	PEEK		
		Standard flow	0.9 mL/min		
		Maximum flow	3.5 mL/min		
		Maximum pressure	20 MPa		
		Particle size	5 µm		
		Organic modifier	Eluent: 0–100% acetone		
			and acetonitrile (no		
			alcohols)		
			Sample: 0–100%		
			acetone, acetonitrile,		
			and alcohols		
		pH range	2–7		
		Temperature range	20–60 °C		
		Standard temperature	20–30 °C		
		Capacity	20 µmol (K⁺)		
Eluents					
Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L		
dipicolinic acid eluent	Dipicolinic acid	568 mg/2 L	1.7 mmol/L		
(standard eluent)	1	0 -			
Oxalic acid/	Oxalic acid	360 mg/2 L	2.0 mmol/L		
dipicolinic acid/	Dipicolinic acid	668 mg/2 L	2.0 mmol/L		
	•	2			

40 mL/2 L

Care

Regeneration

acetonitril eluent

The column must be rinsed with ultrapure water before and after the regeneration.

Acetonitril

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.9 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60).

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.9 mL/min.

2%


Storage Standard eluent at 10-22 °C



Chromatograms

3

Ammonium




5.00


6

Magnesium

10.00



Nitric acid/dipicolinic acid eluent, drinking waterConc. (mg/L)1Sodium112.123Calcium33.442Potassium0.754Magnesium6.88



Oxalic acid/dipicolinic acid/acetonitril eluent, standard

				Conc.	(mg/L)
1	Sodium	20	6	Calcium	20
2	Ammonium	20	7	Magnesium	20
3	Monoethanolamine	20	8	Strontium	20
4	Potassium	20	9	Barium	40
5	Diethanolamine	20			

Ordering information	
Metrosep C 6 - 100/4.0	6.1051.410
Metrosep C 6 Guard/4.0	6.1051.500
Metrosep C 6 S-Guard/4.0	6.1051.510

Metrosep C 6 - 150/4.0 (6.1051.420)

172

The high-capacity C 6 material makes the Metrosep C 6 - 150/4.0 separation column the optimum solution for the separation of standard cations with high differences in concentration in conjunction with reasonable retention times. Drinking water with low ammonium contents can be determined with this column.

Applications

- Standard column
- Amines
- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Universal applications
- Different matrices
- Transition metals

Technical information

Substrate	Silica gel with carboxyl groups
Column dimensions	150 x 4.0 mm
Column body	PEEK
Standard flow	0.9 mL/min
Maximum flow	2.5 mL/min
Maximum pressure	20 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2–7
Temperature range	20–60 °C
Standard temperature	20–30 °C
Capacity	30 µmol (K⁺)

Eluents

Nitric acid/ dipicolinic acid eluent (standard eluent) Nitric acid (c = 1 mol/L) Dipicolinic acid 3.4 mL/2 L 568 mg/2 L 1.7 mmol/L 1.7 mmol/L

Care


Regeneration


The column must be rinsed with ultrapure water before and after the regeneration.


Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.9 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60).

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L HNO₃ + 4 mmol/L dipicolinic acid for 1 h at a flow rate of 0.9 mL/min.

Storage Standard eluent at 10–22 °C







10.00

10.00

10.00

Ordering information	
Metrosep C 6 - 150/4.0	6.1051.420
Metrosep C 6 Guard/4.0	6.1051.500
Metrosep C 6 S-Guard/4.0	6.1051.510

Metrosep C 6 - 250/4.0 (6.1051.430)

174

The Metrosep C 6 - 250/4.0 is the cation column with the greatest capacity in the C 6 series. It is predestined for applications which require the highest separating efficiency. Samples with extreme differences in concentrations can be analyzed reliably with this column. The separation of sodium and ammonium is particularly outstanding here.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺, Co²⁺, Ni²⁺, Zn²⁺, Cd²⁺, Pb²⁺, amines
- Excellent Na⁺/NH₄⁺ separation
- NH_4^+ , $(CH_3)NH_3^+$, $(CH_3)_2NH_2^+$, $(CH_3)_3NH^+$, $(CH_3)_4N^+$, and the respective ethanolamines
- Difficult separation problems
- Great differences in concentration
- Transition metals

Technical information

	rectificat information	
	Substrate	Silica gel with
		carboxyl groups
	Column dimensions	250 x 4.0 mm
	Column body	PEEK
	Standard flow	0.9 mL/min
	Maximum flow	1.5 mL/min
	Maximum pressure	20 MPa
	Particle size	5 µm
	Organic modifier	Eluent: 0–100% acetone
		and acetonitrile (no
		alcohols)
		Sample: 0–100%
		acetone, acetonitrile,
		and alcohols
	pH range	2–7
	Temperature range	20–60 °C
	Standard temperature	20–30 °C
	Capacity	50 μmol (K⁺)
	cupacity	
cid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L
nic acid	568 mg/2 L	1.7 mmol/L
	5	

Eluents

LIUCIIIS				
Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L	
dipicolinic acid eluent	Dipicolinic acid	568 mg/2 L	1.7 mmol/L	
(standard eluent)				
Nitric acid/	Nitric acid (c = 1 mol/L)	16 mL/2 L	8.0 mmol/L	
dipicolinic acid eluent	Dipicolinic acid	434 mg/2 L	1.3 mmol/L	
(modified)				

Care


Regeneration

The column must be rinsed with ultrapure water before and after the regeneration.

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.9 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60).

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.9 mL/min.

Storage Standard eluent at 10–22 °C



Ordering information Metrosep C 6 - 250/4.0

Metrosep C 6 Guard/4.0 Metrosep C 6 S-Guard/4.0 6.1051.430 6.1051.500 6.1051.510



Separation columns

Microbore IC cation-separation columns for lower eluent consumption and greater sensitivity

Metrosep C 4 - 100/2.0 (6.1050.210)

178

The short version of the Metrosep C 4 column with 2 mm inner diameter is intended for rapid determinations of the standard cations. Very short retention times are achieved, for which the elution times of sodium and ammonium nevertheless differ by 25 s. When a special eluent is used, the six cations lithium, ammonium, sodium, calcium, magnesium, and potassium can be determined in less than 5 minutes with the Metrosep C 4 - 100/2.0. With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Lipophilic amines with short retention times
- High flow rate fast separations
- Fast analysis

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	100 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	1.6 mL/min
Maximum pressure	25 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2–7
Temperature range	20-60 °C
Capacity	3 µmol (K⁺)

Eluents

Nitric acid/ dipicolinic acid eluent (standard eluent)	Nitric acid (c = 1 mol/L) Dipicolinic acid	3.4 mL/2 L 234 mg/2 L	1.7 mmol/L 0.7 mmol/L
Nitric acid/ dipicolinic acid eluent (modified)	Nitric acid (c = 1 mol/L) Dipicolinic acid	4.0 mL/2 L 401 mg/2 L	2.0 mmol/L 1.2 mmol/L

Care


Regeneration

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.2 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.2 mL/min.

Storage

In the eluent or in ultrapure water



Ordering information

Metrosep C 4 - 100/2.0 Metrosep C 4 Guard/2.0 Metrosep C 4 S-Guard/2.0 6.1050.210 6.1050.600 6.1050.610

Metrosep C 4 - 150/2.0 (6.1050.220)

180

The Metrosep C 4 - 150/2.0 is the universal standard column in cation analysis of microbore separation columns. High separating efficiency in a brief time. The Metrosep C 4 - 150/2.0 is the ideal separation column for the analysis of alkaline and earth alkaline metals in aqueous media. With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Amines
- Transition metals

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	150 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	1.1 mL/min
Maximum pressure	25 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2-7
Temperature range	20–60 °C
Capacity	4 µmol (K⁺)

Eluents

Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L
dipicolinic acid eluent (standard eluent)	Dipicolinic acid	234 mg/2 L	0.7 mmol/L
Nitric acid/	Nitric acid (c = 1 mol/L)	4.0 mL/2 L	2.0 mmol/L
dipicolinic acid eluent (modified)	Dipicolinic acid	43.6 mg/2 L	0.13 mmol/L
Nitric acid eluent	Nitric acid (c = 1 mol/L)	4.0 mL/2 L	2.0 mmol/L


Care

Regeneration


Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.2 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water. Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.2 mL/min.

Storage


In the eluent or in ultrapure water



Chromatograms



Nitric	acid/dipicolinic acid	eluent, st	andard	Conc.	(mg/L)
1	Lithium	1.00	4	Potassium	10.00
2	Sodium	5.00	5	Calcium	10.00
3	Ammonium	5.00	6	Magnesium	10.00




Nitric acid eluent, traces of cations (MiPCT), 40 °C Lithium 1 С 2 Sodium С

Ammonium

3

		Con	c. (µg/L)
0.50	4	Potassium	0.50
0.50	5	Magnesium	0.50
0.50	6	Calcium	0.50



Nitric	acid/dipicolinic acid e	eluent (m	od.), s [.]	tandard C	Conc. (mg/L)
1	Zinc	2.50	7	Lead	2.50
2	Lithium	0.25	8	Cesium	2.50
3	Cobalt	2.50	9	Manganese	2.50
4	Sodium	1.25	10	Magnesium	n 2.50
5	Ammonium	1.25	11	Calcium	2.50
6	Potassium	2.50			

Ordering information	
Metrosep C 4 - 150/2.0	
Metrosep C 4 Guard/2.0	
Metrosep C 4 S-Guard/2.0	

6.1050.220 6.1050.600 6.1050.610

Metrosep C 4 - 250/2.0 (6.1050.230)

The Metrosep C 4 - 250/2.0 is the cation column with the greatest capacity in the C 4 series with 2 mm inner diameter. It is predestined for applications which require the highest separating efficiency. Samples with high differences in concentrations can be analyzed reliably with this column. The performance capability of the column is demonstrated, for example, in connection with the analysis of sodium traces in addition to monoethanolamine (MEA). With the C 4 - 250/2.0, not only amines and transition metals but also alkaline and alkaline earth metals can be determined in a single run. With its low eluent flow, this column is particularly suitable for IC-MS coupling.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺, Co²⁺, Ni²⁺, Zn²⁺, Cd²⁺, Pb²⁺, amines
- Good Na⁺/NH₄⁺ separation
- NH₄⁺, (CH₃)NH₃⁺, (CH₃)₂NH₂⁺, (CH₃)₃NH⁺, (CH₃)₄N⁺, and the respective ethanolamines
- Difficult separation problems
- High differences in concentration
- Transition metals

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	250 x 2.0 mm
Column body	PEEK
Standard flow	0.2 mL/min
Maximum flow	0.8 mL/min
Maximum pressure	25 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2–7
Temperature range	20–60 °C
Capacity	6 µmol (K⁺)

Eluents

LIUEIIIS				
Nitric acid/	Nitric acid ($c = 1 \text{ mol/L}$)	3.4 mL/2 L	1.7 mmol/L	
dipicolinic acid eluent	Dipicolinic acid	234 mg/2 L	0.7 mmol/L	
(standard eluent)				
Nitric acid/	Nitric acid ($c = 1 \text{ mol/L}$)	5.0 mL/2 L	2.5 mmol/L	
oxalic acid eluent	Oxalic acid	90 mg/2 L	0.5 mmol/L	

Care


Regeneration

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.2 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.2 mL/min.

Storage

In the eluent or in ultrapure water



I	Vitric a	acid/dipicolinic aci	d eluent,	drinking	water	Conc.	(mg/L)
	1	Sodium	3.90	4	Calcium		82.81
	2	Ammonium	n.q.	5	Magnesiu	um	18.76
	3	Potassium	1.12				

Ordering information Metrosep C 4 - 250/2.0 Metrosep C 4 Guard/2.0 Metrosep C 4 S-Guard/2.0

1.00 | 4

5

6

5.00

5.00

Potassium

Calcium Magnesium 10.00

10.00

1

2

3

Sodium

Ammonium

6.1050.230 6.1050.600 6.1050.610

Metrosep C 6 - 100/2.0 (6.01051.210)

184 The 100 mm version of the microbore Metrosep C 6 column is intended for the determination of standard cations, for example in drinking water. Short retention times are attained with a relatively good sodium/ammonium separation. The high capacity of the C 6 material permits larger sample volumes.

The column is suitable for use in IC-MS coupling.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Lipophilic amines with short retention times
- Rapid separations
- IC-MS coupling

Technical information

Substrate	Silica gel with carboxyl groups
Column dimensions	100 x 2.0 mm
Column body	PEEK
Standard flow	0.25 mL/min
Maximum flow	1.0 mL/min
Maximum pressure	20 MPa
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone
	and acetonitrile (no
	alcohols)
	Sample: 0–100%
	acetone, acetonitrile,
	and alcohols
pH range	2–7
Temperature range	20–60 °C
Standard temperature	20–30 °C
Capacity	5 µmol (K⁺)

Eluent

Nitric acid/ dipicolinic acid eluent (standard eluent) Nitric acid (c = 1 mol/L) Dipicolinic acid 3.4 mL/2 L 568 mg/2 L 1.7 mmol/L 1.7 mmol/L

Care

Regeneration

The column must be rinsed with ultrapure water before and after the regeneration.

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.25 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.


Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L HNO₃ + 4 mmol/L dipicolinic acid for 1 h at a flow rate of 0.25 mL/min.


Storage

Standard eluent at 10–22 °C



Chromatograms





Ordering information	
Metrosep C 6 - 100/2.0	6.01051.210
Metrosep C 6 Guard/2.0	6.01051.600
Metrosep C 6 S-Guard/2.0	6.01051.610

Metrosep C 6 - 150/2.0 (6.01051.220)

186

The high-capacity C 6 material makes the microbore version of the Metrosep C 6 - 150/2.0 column the optimum solution for the separation of standard cations with high differences in concentration with reasonable retention times. Drinking water with low ammonium contents can be determined with this column.

The column is suitable for use in IC-MS coupling.

Applications

- Standard column
- Amines
- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Universal applications
- Different matrices
- Transition metals
- IC-MS coupling

Technical information

Substrate Column dimensions Column body Standard flow Maximum flow Maximum pressure Particle size Organic modifier pH range Temperature range Standard temperature Capacity	Silica gel with carboxyl groups 150 x 2.0 mm PEEK 0.25 mL/min 0.7 mL/min 20 MPa 5 µm Eluent: 0–100% acetone and acetonitrile (no alcohols) Sample: 0–100% acetone, acetonitrile, and alcohols 2–7 20–60 °C 20–30 °C 8 µmol (K*)
 3.4 mL/2 L 568 mg/2 L	1.7 mmol/L 1.7 mmol/L

Eluents

Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L
dipicolinic acid eluent	Dipicolinic acid	568 mg/2 L	1.7 mmol/L
(standard eluent)			

Care

Regeneration

The column must be rinsed with ultrapure water before and after the regeneration.


Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.25 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.

Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.25 mL/min.

Storage Standard eluent at 10–22 °C



Chromatograms



Ordering information
Metrosep C 6 - 150/2.0
Metrosep C 6 Guard/2.0
Metrosep C 6 S-Guard/2.0

6.01051.220 6.01051.600 6.01051.610

Metrosep C 6 - 250/2.0 (6.01051.230)

188

The Metrosep C 6 - 250/2.0 is the microbore cation column with the greatest capacity in the C 6 series. It is predestined for applications which require the highest separating efficiency. Samples with extreme differences in concentrations can be analyzed reliably with this column. The separation of sodium and ammonium is particularly outstanding here. The column is suitable for use in IC-MS coupling.

Applications

- Li⁺, Na⁺, K⁺, Rb⁺, Cs⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺, Co²⁺, Ni²⁺, Zn²⁺, Cd²⁺, Pb²⁺, amines
- Excellent Na⁺/NH₄⁺ separation
- NH₄⁺, (CH₃)NH₃⁺, (CH₃)₂NH₂⁺, (CH₃)₃NH⁺, (CH₃)₄N⁺, and the respective ethanolamines
- Difficult separation problems
- Great differences in concentration
- IC-MS coupling

		Technical information	
		Substrate	Silica gel with
			carboxyl groups
		Column dimensions	250 x 2.0 mm
		Column body	PEEK
		Standard flow	0.25 mL/min
		Maximum flow	0.4 mL/min
		Maximum pressure	20 MPa
		Particle size	5 µm
		Organic modifier pH range Temperature range	Eluent: 0–100% acetone and acetonitrile (no alcohols) Sample: 0–100% acetone, acetonitrile, and alcohols 2–7 20–60 °C
		Standard temperature	20-30 °C
		Capacity	13 µmol (K⁺)
/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L
acid eluent	Dipicolinic acid	568 mg/2 L	1.7 mmol/L
eluent)		5	-
eluent	Nitric acid (c = 1 mol/L)	13.5 mL/2 L	6.75 mmol/L

Eluents

Nitric acid/	Nitric acid (c = 1 mol/L)	3.4 mL/2 L	1.7 mmol/L	
dipicolinic acid eluent	Dipicolinic acid	568 mg/2 L	1.7 mmol/L	
(standard eluent)				
Nitric acid eluent	Nitric acid (c = 1 mol/L)	13.5 mL/2 L	6.75 mmol/L	

Care


Regeneration

The column must be rinsed with ultrapure water before and after the regeneration.

Organic contamination: Rinse the column in the opposite flow direction at a flow rate of 0.25 mL/min for 1 h with ultrapure water, then for 1 h with acetonitrile/water (40/60), and finally for 1 h with ultrapure water.


Inorganic contamination: Rinse the column in the opposite flow direction with 10 mmol/L $HNO_3 + 4$ mmol/L dipicolinic acid for 1 h at a flow rate of 0.25 mL/min.

Storage Standard eluent at 10-22 °C



6.01051.230 6.01051.600 6.01051.610

Metrosep C 6 - 250/2.0 Metrosep C 6 Guard/2.0 Metrosep C 6 S-Guard/2.0



Separation columns



IC cation-separation columns for analyses with chemical suppression

Metrosep C Supp 1 - 100/4.0 (6.1052.410)

192

The short version of the Metrosep C Supp 1 is used for the rapid determination of cations in the µg/L range with conductivity detection following sequential suppression.

The baseline noise in cation analysis is improved by the suppression. This results in lower detection limits for the cations to be determined.

Applications

- Li⁺, Na⁺, K⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Samples with low concentrations
- Larger amines
- Low limits of detection
- Fast analysis
- Excellent peak shape
- Matrix with high pH

		Technical information	
		Substrate	Polyvinyl alcohol with
			carboxyl groups
		Column dimensions	100 x 4.0 mm
		Column body	PEEK
		Standard flow	1.0 mL/min
		Maximum flow	1.5 mL/min
		Maximum pressure	15 MPa
		Particle size	5 µm
		Organic modifier	0–50% acetonitrile,
			0–30% acetone,
			no methanol
		pH range	1–12
		Temperature range	20–40 °C
		Standard temperature	40 °C
		Capacity	12 µmol (K⁺)
Eluents			
Nitric acid eluent	Nitric acid (c = 1 mol/L)	10 mL/2 L	5.0 mmol/L
(standard eluent)	Rubidium	172.5 µg/2 L (RbNO₃)	50 µg/L Rb⁺

Care

Note:

Ensure that the maximum pressure is never exceeded during regeneration. If the pressure becomes too high, reduce the flow rate.

Regeneration:

- Disconnect the column outlet from the downstream function units such as suppressor or detector and collect the flow of liquid in a beaker instead.
- 2. Rinse the column with ultrapure water before and after regeneration.

Depending on the type of contamination, proceed in accordance with one of the following instructions:

Organic contaminations:

Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min with the following solutions in succession:

- 1. 1 h with ultrapure water
- 2. 1 h with acetonitrile-water mixture (30:70)
- 3.1 h with ultrapure water

Inorganic contaminations:


- 1. Add 30% acetonitrile to the standard eluent.
- 2. Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min for 1 h.

Storage

Store the column in ultrapure water at 4–8 °C. Do not store the column below 0 °C.



Chromatogram



Ordering information		
Metrosep C Supp 1 - 100/4.0		
Metrosep C Supp 1 Guard/4.0		
Metrosep C Supp 1 S-Guard/4.0		

6.1052.410 6.1052.500 6.1052.510

Metrosep C Supp 1 - 150/4.0 (6.1052.420)

194 The Metrosep C Supp 1 - 150/4.0 separation column is the column of choice for the determination of low concentrations of standard cations.

Detection limits below one μ g/L are achieved through low baseline noise after sequential suppression.

Applications

- Li⁺, Na⁺, K⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Samples with low concentrations
- Organic amines
- Low limits of detection
- Transition metals
- Fast analysis
- Excellent peak shape
- Matrix with high pH

		Technical information	
		Substrate	Polyvinyl alcohol with carboxyl groups
		Column dimensions	150 x 4.0 mm
		Column body	PEEK
		Standard flow	1.0 mL/min
		Maximum flow	1.5 mL/min
		Maximum pressure	15 MPa
		Particle size	5 µm
		Organic modifier	0–50% acetonitrile,
			0–30% acetone,
			no methanol
		pH range	1-12
		Temperature range	20–40 °C
		Standard temperature	40 °C
		Capacity	18 µmol (K⁺)
Eluents			
Nitric acid eluent	Nitric acid ($c = 1 \text{ mol/L}$)	10 mL/2 L	5.0 mmol/L
(standard eluent)	Rubidium	172.5 µg/2 L (RbNO₃)	50 µg/L Rb⁺

Care

Note:

Ensure that the maximum pressure is never exceeded during regeneration. If the pressure becomes too high, reduce the flow rate.

Regeneration:

- Disconnect the column outlet from the downstream function units such as suppressor or detector and collect the flow of liquid in a beaker instead.
- 2. Rinse the column with ultrapure water before and after regeneration.

Depending on the type of contamination, proceed in accordance with one of the following instructions:

Organic contaminations:

Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min with the following solutions in succession:

- 1.1 h with ultrapure water
- 2. 1 h with acetonitrile-water mixture (30:70)
- 3.1 h with ultrapure water

Inorganic contaminations:


- 1. Add 30% acetonitrile to the standard eluent.
- 2. Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min for 1 h.


Storage

Store the column in ultrapure water at 4–8 °C. Do not store the column below 0 °C




Chromatograms





The second second



Ordering information	
Metrosep C Supp 1 - 150/4.0	6.1052.420
Metrosep C Supp 1 Guard/4.0	6.1052.500
Metrosep C Supp 1 S-Guard/4.0	6.1052.510

Metrosep C Supp 1 - 250/4.0 (6.1052.430)

196

The Metrosep C Supp 1 - 250/4.0 separation column is used for difficult separations of standard cations, some transition cations and amines in the low concentration range. Conductivity detection after sequential suppression enables low detection limits thanks to low baseline noise.

Applications

- Li⁺, Na⁺, K⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Samples with low concentrations
- Good Na⁺/NH₄⁺ separation
- Low limits of detection
- Matrix with high pH

Technical information

Polyvinyl alcohol with
carboxyl groups
250 x 4.0 mm
PEEK
1.0 mL/min
1.5 mL/min
15 MPa
5 µm
0–50% acetonitrile,
0–30% acetone,
no methanol
1-12
20-40 °C
40 °C
30 µmol (K ⁺)

Eluents

Nitric acid eluent	Nitric acid (c = 1 mol/L)	10 mL/2 L	5.0 mmol/L
(standard eluent)	Rubidium	172.5 μg/2 L (RbNO₃)	50 µg/L Rb⁺

Care

Note:

Ensure that the maximum pressure is never exceeded during regeneration. If the pressure becomes too high, reduce the flow rate.

Regeneration:

- Disconnect the column outlet from the downstream function units such as suppressor or detector and collect the flow of liquid in a beaker instead.
- 2. Rinse the column with ultrapure water before and after regeneration.

Depending on the type of contamination, proceed in accordance with one of the following instructions:

Organic contaminations:

Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min with the following solutions in succession:

- 1.1 h with ultrapure water
- 2. 1 h with acetonitrile-water mixture (30:70)
- 3. 1 h with ultrapure water

Inorganic contaminations:

- 1. Add 30% acetonitrile to the standard eluent.
- 2. Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min for 1 h.

Storage

Store the column in ultrapure water at 4–8 °C. Do not store the column below 0 °C

Chromatograms µS/cm µS/cn 1.10 0.34 3 1.00 0.32 0.90 0.30 0.80 0.28 0.70 0.26 0.60 0.24 0.50 0.22 0.40 0.20 0.30 0.18 0.20 0.16 4.0 2.0 3.0 5.0 8.0 4.0 6.0 8.0 1.0 6.0 7.0 9.0 10.0 11.0 min 0.0 2.0 10.0 12.0 14.0 16.0 min Nitric acid eluent, geological leachate, 40 °C Conc. (mg/L) Nitric acid eluent, power plant sample, Conc. (µg/L) 1 Lithium 164.6 4 Magnesium 0.35 MiPCT-ME, 2000 µL, 40 °C 5 12.9 2 Sodium Calcium 1.02 1 Lithium Monoethanol- 4000 1.0 3 3 0.81 Potassium 2 Sodium 1.0 amine (MEA) 3.2 4.5 2.8 4.0 3.5 2/ 3.0 2.0 2.5 1.6 2.0 1.5 1.2 1.0 0. 0.0 0.4 10.0 12.0 14.0 16.0 18.0 10.0 12.0 0.0 2.0 4.0 6.0 8.0 min Nitric acid eluent, standard, 40 °C Conc. (mg/L) Ni

2

1

2

3

Lithium

Sodium

Ammonium

1.00 | 4

5.00

5.00

5

6

Potassium

Calcium

Magnesium

	0.0 2.0	4.0	6.0	8.0	10.0	12.0	14.0	10.0	19/0	min
Nitri	c acid elue	ent, m	agnesi	ium sp	ort dri	ink, 40) °C	Co	nc. (n	ng/L)
1	Sodium			227	4	Zi	inc			6.0
2	Potassiu	um		202	5	С	alcium			0.6
3	Magnes	sium		165						
	5				I					

Ordering information	
Metrosep C Supp 1 - 250/4.0	6.1052.430
Metrosep C Supp 1 Guard/4.0	6.1052.500
Metrosep C Supp 1 S-Guard/4.0	6.1052.510

10.00

10.00

10.00

Metrosep C Supp 2 - 100/4.0 (6.01053.410)

198

The Metrosep C Supp 2 separation material is based on a polystyrene/divinylbenzene copolymer with carboxyl groups. It is suitable for the separation and determination of monovalent and divalent cations. The Metrosep C Supp 2 - 100/4.0 column is the shortest separation column in the Metrosep C Supp 2 product range. It is especially suitable for trace analysis of standard cations. Limits of guantification below the µg/L range are achieved thanks to the extremely low baseline noise following sequential suppression.

Applications

- Li⁺, Na⁺, K⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Samples with low concentrations
- Larger amines
- Low limits of detection
- Fast analysis
- Excellent peak shape
- Matrix with high pH

Technical information

	Substrate	Polystyrene-divinyl- benzene copolymer with carboxyl groups
	Column dimensions	100 x 4.0 mm
	Column body	PEEK
	Standard flow	1.0 mL/min
	Maximum flow	3.8 mL/min
	Maximum pressure	25 MPa
	Particle size	5 μm
	Organic modifier	Eluent: 0–100% acetone
	organie mounter	and acetonitrile (no alcohol)
		Sample: 0–100% acetone,
		acetonitrile and alcohols
	pH range	Eluent: 0–12
	P	Sample: 0–14
	Temperature range	10–60 °C
	Standard temperature	40 °C
	Capacity	23 µmol (K⁺)
		•
Nitric acid (c = 1 mol/L)	10 mL/2 L	5.0 mmol/L
Rubidium	172.5 μg/2 L (RbNO ₃)	50 μg/L Rb⁺
Nitric acid (c = 1 mol/L)	14 mL/2 L	7.0 mmol/L
Rubidium	172.5 μg/2 L (RbNO₃)	50 μg/L Rb⁺
	·· -· · · · · · · · · · · · · · · · · ·	50 Mg/ L 110

Eluents

Nitric acid eluent (standard eluent) Nitric acid eluent (modified)

Care

Note:

Ensure that the maximum pressure is never exceeded during regeneration. If the pressure becomes too high, reduce the flow rate.

Preparation:

Rinse the column with eluent for 3 h.

Regeneration:

- 1. Disconnect the column outlet from the downstream function units such as suppressor or detector and collect the flow of liquid in a beaker instead.
- 2. Depending on the type of contamination, proceed in accordance with one of the following instructions:

a. Organic contaminations:

Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min with the following solutions in succession:

- 1. 1 h with ultrapure water
- 2. 1 h with acetonitrile-water mixture (40:60)
- 3.1 h with ultrapure water

b. Inorganic contaminations:


Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min for 1 h with 50 mmol/L nitric acid.

Storage

Store the column in standard eluent at ambient temperature.



Chromatogram



Ordering information

Metrosep C Supp 2 - 100/4.0 Metrosep C Supp 2 Guard/4.0 6.01053.410 6.01053.500

Metrosep C Supp 2 - 150/4.0 (6.01053.420)

200

The Metrosep C Supp 2 - 150/4.0 column is the standard separation column of the Metrosep C Supp 2 product range. It is suitable for the separation and determination of monovalent and divalent cations with an excellent sodium/ammonium separation. The Metrosep C Supp 2 separation material is based on a polystyrene/divinylbenzene copolymer with carboxyl groups. The column is used with sequential suppression. It is accordingly particularly suitable for the determination of concentrations in the middle µg/L range and below.

Applications

- Li⁺, Na⁺, K⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Samples with low concentrations
- Organic amines
- Low limits of detection
- Transition metals
- Fast analysis
- Excellent peak shape
- Matrix with high pH

		Technical information	
		Substrate	Polystyrene-divinylbenzene copolymer with carboxyl groups
		Column dimensions	150 x 4.0 mm
		Column body	PEEK
		Standard flow	1.0 mL/min
		Maximum flow	3.1 mL/min
		Maximum pressure	25 MPa
		Particle size	5 µm
		Organic modifier	Eluent: 0–100% acetone and acetonitrile (no alcohol) Sample: 0–100% acetone, acetonitrile and alcohols
		pH range	Eluent: 0–12 Sample: 0–14
		Temperature range	10–60 °C
		Standard temperature	40 °C
		Capacity	35 µmol (K⁺)
eluent	Nitric acid (c = 1 mol/L)	10 mL/2 L	5.0 mmol/L
eluent)	Rubidium	172.5 μg/2 L (RbNO ₃)	50 µg/L Rb⁺

Nitric acid ((standard e

Eluents

Care

Note:

Ensure that the maximum pressure is never exceeded during regeneration. If the pressure becomes too high, reduce the flow rate.

Preparation:

Rinse the column with eluent for 3 h.

Regeneration:

- Disconnect the column outlet from the downstream function units such as suppressor or detector and collect the flow of liquid in a beaker instead.
- 2. Depending on the type of contamination, proceed in accordance with one of the following instructions:

a. Organic contaminations:

Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min with the following solutions in succession:

- 1.1 h with ultrapure water
- 2. 1 h with acetonitrile-water mixture (40:60)
- 3. 1 h with ultrapure water

b. Inorganic contaminations:


Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min for 1 h with 50 mmol/L nitric acid.

Storage


Store the column in standard eluent at ambient temperature.




Chromatograms



Nitric acid eluent, standard, 40 °C				Con	c. (mg/L)
1	Lithium	1.00	4	Potassium	10.00
2	Sodium	5.00	5	Magnesium	10.00
3	Ammonium	5.00	6	Calcium	10.00



Marine 100 2



Ordering information

Metrosep C Supp 2 - 150/4.0 Metrosep C Supp 2 Guard/4.0 6.01053.420 6.01053.500

Metrosep C Supp 2 - 250/4.0 (6.01053.430)

The longest separation column in the Metrosep C Supp 2 product range is the Metrosep C Supp 2 - 250/4.0. The Metrosep C Supp 2 separation material is based on a polystyrene/divinylbenzene copolymer with carboxyl groups. Thanks to the optimized sodium/ammonium separation of this separation material, this column is perfectly suitable for determination of the smallest concentrations of ammonium in addition to a large amount of sodium. The column is used with sequential suppression. It is accordingly particularly suitable for the determination of concentrations in the middle below.

Applications

- Li⁺, Na⁺, K⁺, Mg²⁺, Ca²⁺, Sr²⁺, Ba²⁺, NH₄⁺
- Samples with low concentrations
- Good Na⁺/NH⁺₄ separation
- Low limits of detection
- Transition metals
- Excellent peak shape
- Matrix with high pH

Technical information

e μg/L range and	Substrate	Polystyrene-divinyl- benzene copolymer
		with carboxyl groups
	Column dimensions	250 x 4.0 mm
	Column body	PEEK
	Standard flow	1.0 mL/min
	Maximum flow	2.0 mL/min
	Maximum pressure	25 MPa
	Particle size	5 µm
	Organic modifier	Eluent: 0–100% acetone and acetonitrile (no alcohol)
		Sample: 0–100% acetone,
		acetonitrile and alcohols
	pH range	Eluent: 0–12
		Sample: 0–14
	Temperature range	10-60 °C
	Standard temperature	40 °C
	Capacity	58 µmol (K⁺)
	14 mal (2)	7.0 mmol/L
d (c = 1 mol/L)	14 mL/2 L	
1	172.5 μg/2 L (RbNO ₃)	50 µg/L Rb⁺
rile	40 mL/2 L	2%
d (c = 1 mol/L)	2.5 mL/2 L	1.25 mmol/L
1	43.1 μg/2 L (RbNO ₃)	12.5 µg/L Rb⁺
d (c = 1 mol/L)	25 mL/2 L	12.5 mmol/L
1	431 μg/2 L (RbNO ₃)	125 µg/L Rb⁺

Eluents			
Nitric acid eluent	Nitric acid (c = 1 mol/L)	14 mL/2 L	7.0 mmol/L
(modified)	Rubidium	172.5 µg/2 L (RbNO₃)	50 µg/L Rb⁺
	Acetonitrile	40 mL/2 L	2%
Gradient: Eluent A	Nitric acid (c = 1 mol/L)	2.5 mL/2 L	1.25 mmol/L
	Rubidium	43.1 μg/2 L (RbNO ₃)	12.5 µg/L Rb⁺
Eluent B	Nitric acid (c = 1 mol/L)	25 mL/2 L	12.5 mmol/L
	Rubidium	431 µg/2 L (RbNO₃)	125 µg/L Rb⁺

Care

Note:

Ensure that the maximum pressure is never exceeded during regeneration. If the pressure becomes too high, reduce the flow rate.

Preparation:

Rinse the column with eluent for 3 h.

Regeneration:

- 1. Disconnect the column outlet from the downstream function units such as suppressor or detector and collect the flow of liquid in a beaker instead.
- 2. Depending on the type of contamination, proceed in accordance with one of the following instructions:

a. Organic contaminations:

Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min with the following solutions in succession:


- 1. 1 h with ultrapure water
- 2.1 h with acetonitrile-water mixture (40:60)
- 3. 1 h with ultrapure water

b. Inorganic contaminations:

Regenerate the column in the direction opposite to the flow at a flow rate of 1.0 mL/min for 1 h with 50 mmol/L nitric acid.

Storage

Store the column in standard eluent at ambient temperature.



Ordering information Metrosep C Supp 2 - 250/4.0 Metrosep C Supp 2 Guard/4.0

To a

6.01053.430 6.01053.500



Separation columns



Separation column for the determination of organic substances

MetroSil RP 3 - 150/4.0 (6.01070.420)

The MetroSil RP 3 - 150/4.0 reversed phase column with medium capacity which can be used universally with aqueous eluents as well as with aqueous samples. This characteristic is especially important for applications in ion chromatography. The MetroSil RP 3 - 150/4.0 can be used to solve application problems which lie in the boundary range between HPLC and ion chromatography. The MetroSil RP 3 material is an "endcapped" C18 silica gel with a pore width of 120 angstroms.

Applications

- Determination of organic substances with low polarity and low charge
- Caffeine
- Determination of pharmaceutical products
- NTA, EDTA, DTPA (with UV/VIS detection)

Technical information

Substrate	Silica gel C ₁₈
Column dimensions	150 x 4.0 mm
Column body	Stainless steel
Standard flow	0.7 mL/min
Maximum flow	5.0 mL/min
Maximum pressure	40 MPa
Particle size	5 µm
Organic modifier	0-100%
pH range	2–9

Eluents

Acetonitrile/water	Acetonitrile	300 mL/2 L	15%
(standard eluent)	Water	1700 mL/2 L	85%
Nitrate/sulfuric acid/	Potassium nitrate	44.044 g/2 L	20 mmol/L
Methanol (phenol eluent)	Sulfuric acid (c = 1 mol/L)	1.0 mL/2 L	0.5 mmol/L
	Methanol	1000 mL/2 L	50%
Acetonitrile/water/sulfuric acid	Acetonitrile	300 mL/2 L	15%
(paracetamol eluent)	Water	1660 mL/2 L	85%
	Sulfuric acid (c = 1 mol/L)	40 mL/2 L	20 mmol/L

Care

Regeneration


Rinse for 15 min at 1.0 mL/min with each of the following: 100% water, followed by acetonitrile, isopropanol, hexane, isopropanol, and back to acetonitrile.


Storage

For short periods (< 48 h): Rinse the column for 30 min at 0.5 mL/min with acetonitrile/water 50:50 (v:v). For prolonged periods (> 48 h): Rinse the column with water for 30 min at 0.5 mL/min. Rinse the column with acetonitrile for 30 min at 0.5 mL/min.



Chromatograms






20

| 3

 $(\lambda = 256 \text{ nm})$

1 Paracetamol 2 Caffeine

Conc. (mg/L) 2,6 Hydroxybenzoic acid 20



100

Ordering information MetroSil RP 3 - 150/4.0 MetroSil RP 3 - Guard/4.0 Cartridge holder for MetroSil RP 3 Guard/4.0

2

Cresol

6.01070.420 6.01070.500 6.02821.010



Guard columns



IC guard columns (precolumns)

Optimum protection for the separation columns, minimal dead volume, the same phase, and therefore nearly no influence on the chromatography are the characteristics of the Metrohm IC guard columns. They are extremely efficient, easy to handle and yet economical.

IC guard column cartridge for Hamilton PRP-X100 (6.1005.020)

210

For the preservation of PRP-X100 analytical separation columns. The cartridge effectively removes contaminations in the form of particles, as might arise, for example, from inadvertent bacteria and algae growth.

Applications

• Anions

Technical information

20 x 4.0 mm
Stainless steel
10 µm
Cartridge



Ordering information

Guard column cartridge for Hamilton PRP-X100 Guard cartridge holder, 20 mm

For use with Hamilton PRP-X100 - 100/4.0 Hamilton PRP-X100 - 250/4.0 6.1005.020 6.02821.000

> 6.1005.000 6.1005.010

Super-Sep Guard/4.6 (6.1009.010)

For the protection of the Super-Sep - 100/4.6 analytical separation column

Applications

• Anions

Technical information

Column dimensions	12 x 4.6 mm
Column body	Stainless steel
Particle size	12 µm
Туре	Column



Ordering information

Super-Sep Guard/4.6

For use with Super-Sep - 100/4.6 6.1009.010

6.1009.000

Metrosep Dual 4 Guard Column kit (6.1016.500)

Even if the Dual-4 columns based on monolithic silica gel are very durable, the use of the Dual 4 guard column is advised in order to increase the safety of the analytical separation column even more. The Dual 4 guard column is a PEEK cartridge which is also filled with monolithic silica gel. This cartridge is easy to replace and is screwed directly onto the analytical column in an aluminum holder. The proven «On Column Guard System» is simple to use and additionally offers the advantage of minimal dead volume.

Applications

• Anions

Technical information

Substrate	Monolithic silica gel
Column dimensions	5 x 4.6 mm
Column body	PEEK cartridge in an
	aluminum cartridge holder
	(replaceable)
Particle size	Monolith with 2 μ m
	Macropores and 13 nm
	Mesopores
Organic modifier	0–5% methanol or
	acetonitrile only
pH range	0–8
Туре	Cartridge



Ordering information	
Guard column kit for the Metrosep Dual 4, comprised of three guard column cartridges and	
one guard column cartridge holder	6.1016.500
Guard column cartridges for the Metrosep Dual 4 (3 pcs.)	6.1016.510
For use with	
Metrosep Dual 4 - 100/4.6	6.1016.030

Metrosep A Supp 1 Guard/4.6 (6.1005.340)

The Metrosep A Supp 1 Guard/4.6 protects the Metrosep A Supp 1 - 250/4.6 separation column against contamination from particles and bacteria.

Applications

- Anions
- Oxhalogenides

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	50 x 4.6 mm
Column body	PEEK
Particle size	7 µm
Organic modifier	0–100%
pH range	0–13
Туре	Column



Ordering information

Metrosep A Supp 1 Guard/4.6

For use with Metrosep A Supp 1 - 250/4.6 6.1005.340

6.1005.300

Metrosep A Supp 4 Guard/4.0 (6.01021.500) Metrosep A Supp 4 S-Guard/4.0 (6.01021.510)

214

The Metrosep A Supp 4 Guard/4.0 reliably protects the Metrosep A Supp 4 IC anion column against contamination from the sample or eluent. It contains the same separation material as the Metrosep A Supp 4, is also made of PEEK, and is screwed directly onto the separation column with nearly no dead volume («On Column Guard System»). The guard column prolongs the lifetime of the analytical column, with practically no influence on its chromatographic separation performance. The economical price and simple handling make using the Metrosep A Supp 4 Guard/4.0 highly recommended.

Applications

• Anions

Technical information

Substrate	Polyvinyl alcohol with
	quaternary
	ammonium groups
Column dimensions	5 x 4.0 mm
Column body	PEEK
Particle size	9 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 4 S-Guard/4.0 (6.01021.510) must be used instead of the Metrosep A Supp 4 Guard/4.0 (6.01021.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information

Metrosep A Supp 4 Guard/4.0 Metrosep A Supp 4 S-Guard/4.0

For use with Metrosep A Supp 4 - 250/4.0 6.01021.500 6.01021.510

6.1006.430

Metrosep A Supp 4 Guard/2.0 (6.01021.600) Metrosep A Supp 4 S-Guard/2.0 (6.01021.610)

The Metrosep A Supp 4 Guard/2.0 reliably protects the Metrosep A Supp 4 microbore anion column against contamination from the sample or eluent. It contains the same separation material as the Metrosep A Supp 4, is also made of PEEK, and is screwed directly onto the separation column with nearly no dead volume («On Column Guard System»). The guard column prolongs the lifetime of the analytical column, with practically no influence on its chromatographic separation performance. The economical price and simple handling make using the Metrosep A Supp 4 Guard/2.0 highly recommended.

Applications

• Anions

Technical information

Substrate	Polyvinyl alcohol with
	quaternary
	ammonium groups
Column dimensions	5 x 2.0 mm
Column body	PEEK
Particle size	9 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 4 S-Guard/2.0 (6.01021.610) must be used instead of the Metrosep A Supp 4 Guard/2.0 (6.01021.600). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information

Metrosep A Supp 4 Guard/2.0 Metrosep A Supp 4 S-Guard/2.0

For use with Metrosep A Supp 4 - 250/2.0 6.01021.600 6.01021.610

6.01021.230

Metrosep A Supp 5 Guard/4.0 (6.1006.500) Metrosep A Supp 5 S-Guard/4.0 (6.1006.540)

216

The Metrosep A Supp 5 Guard/4.0 reliably protects the Metrosep A Supp 5 and 7 IC anion columns against contamination from the sample or eluent. It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the separation column with nearly no dead volume («On Column Guard System»). The guard column prolongs the lifetime of the analytical column, with practically no influence on its chromatographic separation performance. The economical price and simple handling make using the Metrosep A Supp 5 Guard/4.0 highly recommended.

Applications

• Anions

Technical information

Substrate	Polyvinyl alcohol with
	quaternary
	ammonium groups
Column dimensions	5 x 4.0 mm
Column body	PEEK
Particle size	5 µm
Organic modifier	0–100% (particularly
	acetone, acetonitrile,
	methanol)
pH range	3–12
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 5 S-Guard/4.0 (6.1006.540) must be used instead of the Metrosep A Supp 5 Guard/4.0 (6.1006.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep A Supp 5 Guard/4.0	6.1006.500
Metrosep A Supp 5 S-Guard/4.0	6.1006.540
For use with	
Metrosep A Supp 5 - 50/4.0	6.1006.550
Metrosep A Supp 5 - 100/4.0	6.1006.510
Metrosep A Supp 5 - 150/4.0	6.1006.520
Metrosep A Supp 5 - 250/4.0	6.1006.530
Metrosep A Supp 7 - 150/4.0	6.1006.620
Metrosep A Supp 7 - 250/4.0	6.1006.630

Metrosep A Supp 5 Guard/2.0 (6.1006.600) Metrosep A Supp 5 S-Guard/2.0 (6.1006.610)

The Metrosep A Supp 5 Guard/2.0 reliably protects the Metrosep A Supp 5 and 7 microbore anion columns against contamination from the sample or eluent. It contains the same separation material as the Metrosep A Supp 5, is also made of PEEK, and is screwed directly onto the separation column with nearly no dead volume («On Column Guard System»). The guard column prolongs the lifetime of the analytical column, with practically no influence on its chromatographic separation performance. The economical price and simple handling make using the Metrosep A Supp 5 Guard/2.0 highly recommended.

Applications

• Anions

Technical information

Substrate	Polyvinyl alcohol with	
	quaternary	
	ammonium groups	
Column dimensions	5 x 2.0 mm	
Column body	PEEK	
Particle size	5 µm	
Organic modifier	0–100% (particularly	
	acetone, acetonitrile,	
	methanol)	
pH range	3–12	
Туре	Column	



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 5 S-Guard/2.0 (6.1006.610) must be used instead of the Metrosep A Supp 5 Guard/2.0 (6.1006.600). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Metrosep A Supp 5 Guard/2.0	6.1006.600
Metrosep A Supp 5 S-Guard/2.0	6.1006.610
For use with	
Metrosep A Supp 5 - 150/2.0	6.1006.220
Metrosep A Supp 5 - 250/2.0	6.1006.230
Metrosep A Supp 7 - 150/2.0	6.1006.640
Metrosep A Supp 7 - 250/2.0	6.1006.650

Metrosep A Supp 10 Guard/4.0 (6.1020.500) Metrosep A Supp 10 S-Guard/4.0 (6.1020.510) Metrosep A Supp 10 Guard HC/4.0 (6.1020.520)

The Metrosep A Supp 10 Guard/4.0 reliably protects the A Supp 10 analytical separation columns against contamination. Thanks to the «On Column Guard System», the guard column is very easy to handle. The guard column screws easily and directly onto the analytical column without tools.



Applications

Anions

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with quaternary
	ammonium groups
Column dimensions	6.1020.500: 5 x 4.0 mm
	6.1020.510: 5 x 4.0 mm
	6.1020.520: 12.5 x 4.0 mm
Column body	PEEK
Particle size	4.6 µm
Organic modifier	0-100%
pH range	0-14
Туре	Column

Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 10 S-Guard/4.0 (6.1020.510) must be used instead of the Metrosep A Supp 10 Guard/4.0 (6.1020.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

The Metrosep A Supp 10 Guard HC/4.0 is the high-capacity variant of the Metrosep A Supp 10 Guard/4.0.

The separation of cyclamate and phosphate is significantly improved when the Metrosep A Supp 5 - 100/4.0 (6.1006.510) is combined with the Metrosep A Supp 10 Guard HC/4.0.

Ordering information	
Metrosep A Supp 10 Guard/4.0	6.1020.500
Metrosep A Supp 10 S-Guard/4.0	6.1020.510
Metrosep A Supp 10 Guard HC/4.0	6.1020.520
For use with	
Metrosep A Supp 5 - 100/4.0 (with Metrosep A Supp 10 Guard HC/4.0; 6.1020.520)	6.1006.510
Metrosep A Supp 10 - 50/4.0	6.1020.050
Metrosep A Supp 10 - 75/4.0	6.1020.070
Metrosep A Supp 10 - 100/4.0	6.1020.010
Metrosep A Supp 10 - 250/4.0	6.1020.030

Metrosep A Supp 10 Guard/2.0 (6.1020.600) Metrosep A Supp 10 S-Guard/2.0 (6.1020.610)

The Metrosep A Supp 10 Guard/2.0 column reliably protects the Metrosep A Supp 10 microbore separation columns against contamination. Thanks to the «On Column Guard System», the guard column is very easy to handle. The guard column screws easily and directly onto the analytical column without tools.



Applications

• Anions

Technical information

Substrate	Polystyrene/divinylbenzene	
	copolymer with quaternary	
	ammonium groups	
Column dimensions	5 x 2.0 mm	
Column body	PEEK	
Particle size	4.6 µm	
Organic modifier	0–100%	
pH range	0–14	
Туре	Column	

Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 10 S-Guard/2.0 (6.1020.610) must be used instead of the Metrosep A Supp 10 Guard/2.0 (6.1020.600). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep A Supp 10 Guard/2.0	6.1020.600
Metrosep A Supp 10 S-Guard/2.0	6.1020.610
For use with	
Metrosep A Supp 10 - 50/2.0	6.1020.250
Metrosep A Supp 10 - 75/2.0	6.1020.270
Metrosep A Supp 10 - 100/2.0	6.1020.210
Metrosep A Supp 10 - 150/2.0	6.1020.220
Metrosep A Supp 10 - 250/2.0	6.1020.230

Metrosep A Supp 16 Guard/4.0 (6.1031.500) Metrosep A Supp 16 S-Guard/4.0 (6.1031.510)

The Metrosep A Supp 16 Guard/4.0 reliably protects the A Supp 16 analytical separation columns against contamination. Thanks to the «On Column Guard System», the guard column is very easy to handle. The guard column screws easily onto the analytical column without tools.

Applications

Anions

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	5 x 4.0 mm
Column body	PEEK
Particle size	4.6 µm
Organic modifier	0-10%
pH range	0-14
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 16 S-Guard/4.0 (6.1031.510) must be used instead of the Metrosep A Supp 16 Guard/4.0 (6.1031.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information		
Metrosep A Supp 16 Guard/4.0	6.1031.500	
Metrosep A Supp 16 S-Guard/4.0	6.1031.510	
For use with		
Metrosep A Supp 16 - 100/4.0	6.1031.410	
Metrosep A Supp 16 - 150/4.0	6.1031.420	
Metrosep A Supp 16 - 250/4.0	6.1031.430	
Metrosep A Supp 7 - 150/4.0	6.1006.620	
Metrosep A Supp 7 - 250/4.0	6.1006.630	

Metrosep A Supp 16 Guard/2.0 (6.1031.600) Metrosep A Supp 16 S-Guard/2.0 (6.1031.610)

The Metrosep A Supp 16 Guard/2.0 reliably protects the A Supp 16 analytical separation columns with 2 mm inner diameter against contamination. Thanks to the «On Column Guard System», the guard column is very easy to handle. The guard column screws easily onto the analytical column without tools.

Applications

• For anions

Technical information

Substrate	Polystyrene/divinylbenzene	
	copolymer with quaternary	
	ammonium groups	
Column dimensions	5 x 2.0 mm	
Column body	PEEK	
Particle size	4.6 µm	
Organic modifier	0–10%	
pH range	0–14	
Туре	Column	



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 16 S-Guard/2.0 (6.1031.610) must be used instead of the Metrosep A Supp 16 Guard/2.0 (6.1031.600). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information		
Metrosep A Supp 16 Guard/2.0	6.1031.600	
Metrosep A Supp 16 S-Guard/2.0	6.1031.610	
For use with		
Metrosep A Supp 16 - 100/2.0	6.1031.210	
Metrosep A Supp 16 - 150/2.0	6.1031.220	
Metrosep A Supp 16 - 250/2.0	6.1031.230	
Metrosep A Supp 7 - 150/2.0	6.1006.640	
Metrosep A Supp 7 - 250/2.0	6.1006.650	

Metrosep A Supp 17 Guard/4.0 (6.01032.500) Metrosep A Supp 17 S-Guard/4.0 (6.01032.510) Metrosep A Supp 17 S-Guard - 50/4.0 (6.01032.530)

The Metrosep A Supp 17 Guard/4.0 reliably protects the 222 A Supp 17 analytical separation columns against contamination. Thanks to the «On Column Guard System», the guard column is very easy to handle. The guard column screws easily onto the analytical column without tools.

Applications

• Anions

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	5 x 4.0 mm, and
	50 x 4.0 mm respectively
Column body	PEEK
Particle size	5.0 µm
Organic modifier	0–100% methanol
	0–40% acetone or
	acetonitrile
pH range	0-14
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep A Supp 17 S-Guard/4.0 (6.01032.510) must be used instead of the Metrosep A Supp 17 Guard/4.0 (6.01032.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep A Supp 17 Guard/4.0	6.01032.500
Metrosep A Supp 17 S-Guard/4.0	6.01032.510
Metrosep A Supp 17 S-Guard - 50/4.0	6.01032.530
For use with	
Metrosep A Supp 17 - 100/4.0	6.01032.410
Metrosep A Supp 17 - 150/4.0	6.01032.420
Metrosep A Supp 17 - 250/4.0	6.01032.430

Metrosep Organic Acids Guard/4.6 (6.1005.250)

The Metrosep Organic Acids Guard/4.6 effectively removes contamination, thus protecting the analytical separation column.

Applications

• Organic acids

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with
	sulfonic acid groups
Column dimensions	50 x 4.6 mm
Column body	Stainless steel
Particle size	9 µm
Organic modifier	0–20%
pH range	1–13
Туре	Column



Orc	lering	information	

Metrosep Organic Acids Guard/4.6

For use with Metrosep Organic Acids - 100/7.8 Metrosep Organic Acids - 250/7.8 6.1005.250

6.1005.210 6.1005.200

Metrosep Carb 2 Guard/4.0 (6.1090.500) Metrosep Carb 2 S-Guard/4.0 (6.1090.510)

The Metrosep Carb 2 Guard/4.0 and the the Metrosep Carb 2 S-Guard/4.0 effectively removes contaminations, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on chromatographic separation can be ignored.

Applications

• Carbohydrates

Technical information

Substrate	Polystyrene-divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	5 x 4.0 mm
Column body	PEEK
Particle size	5 µm
Organic modifier	0–50% acetonitrile or
	methanol (eluent)
	0–100% acetone,
	acetonitrile or methanol
	(sample)
pH range	0–14
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep Carb 2 S-Guard/4.0 (6.1090.510) must be used instead of the Metrosep Carb 2 Guard/4.0 (6.1090.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep Carb 2 Guard/4.0	6.1090.500
Metrosep Carb 2 S-Guard/4.0	6.1090.510
For use with	
Metrosep Carb 2 - 100/4.0	6.1090.410
Metrosep Carb 2 - 150/4.0	6.1090.420
Metrosep Carb 2 - 250/4.0	6.1090.430

Metrosep Carb 2 Guard/2.0 (6.01090.600) Metrosep Carb 2 S-Guard/2.0 (6.01090.610)

The microbore guard column, Metrosep Carb 2 Guard/2.0 effectively removes contamination, thus protecting the analytical separation column. The design of the guard column has been selected in such a way that its influence on the chromatographic separation can be ignored.

Applications

• Carbohydrates

Technical information

Styrene/divinylbenzene
copolymer with quaternary
ammonium groups
5 x 2.0 mm
PEEK
5 µm
0–50% acetonitrile or
methanol (eluent)
0–100% acetone,
acetonitrile or methanol
(sample)
0–14
Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep Carb 2 S-Guard/2.0 (6.01090.610) must be used instead of the Metrosep Carb 2 Guard/2.0 (6.01090.600). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep Carb 2 Guard/2.0	6.01090.600
Metrosep Carb 2 S-Guard/2.0	6.01090.610
For use with	
Metrosep Carb 2 - 100/2.0	6.01090.210
Metrosep Carb 2 - 150/2.0	6.01090.220
Metrosep Carb 2 - 250/2.0	6.01090.230

Nucleosil 5SA 2 Guard Cartridge/4.0 (6.1007.110)

226 For the protection of the Nucleosil 5SA - 125/4.0 analytical separation column.

Applications

• Cations

Technical information

Substrate	S
	S
Column dimensions	2
Column body	S
Particle size	5
Туре	C

Spherical silica gel with sulfonic acid groups 20 x 4.0 mm Stainless steel 5 µm Cartridge



Ordering information

Nucleosil 5SA 2 Guard Cartridge/4.0 Holder to Nucleosil 5SA 2 Guard Cartridge/4.0

For use with IC Cation Column Nucleosil 5SA - 125/4.0 6.1007.110 6.2821.140

6.1007.000

Metrosep C 3 Guard/4.0 (6.1010.450) Metrosep C 3 S-Guard/4.0 (6.1010.460)

The Metrosep C 3 Guard/4.0 contains the C 3 column material and is used to protect Metrosep C 3 cation columns. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be extended considerably. The Metrosep C 3 Guard/4.0 also functions according to the «On Column Guard System» and is screwed directly onto the separation column with nearly no dead volume.



Technical information

Substrate	Polyvinyl alcohol with
	carboxyl groups
Column dimensions	5 x 4.0 mm
Column body	PEEK
Particle size	5 µm
Organic modifier	50% acetonitrile or 30%
	acetone (no methanol)
pH range	2–12
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep C 3 S-Guard/4.0 (6.1010.460) must be used instead of the Metrosep C 3 Guard/4.0 (6.1010.450). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep C 3 Guard/4.0	6.1010.450
Metrosep C 3 S-Guard/4.0	6.1010.460
For use with	
Metrosep C 3 - 100/4.0	6.1010.410
Metrosep C 3 - 150/4.0	6.1010.420
Metrosep C 3 - 250/4.0	6.1010.430

Metrosep C 4 Guard/4.0 (6.1050.500) Metrosep C 4 S-Guard/4.0 (6.1050.510) Metrosep C 4 S-Guard - 50/4.0 (6.1050.530)

228

The Metrosep C 4 Guard/4.0 contains the C 4 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be extended considerably. The economical price is an additional plus. The Metrosep C 4 Guard/4.0 also functions according to the «On Column Guard System» and is screwed directly onto the separation column with nearly no dead volume. In comparison with the Standard C 4 guard columns, the Metrosep C 4 S-Guard - 50/4.0 exhibits greater capacity and therefore an even longer service life.

Applications

• Cations

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	5 x 4.0 mm, and
	50 x 4.0 mm respectively
Column body	PEEK
Particle size	5 µm
Organic modifier	0–100% (no methanol)
pH range	2–7
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep C 4 S-Guard/4.0 (6.1050.510) must be used instead of the Metrosep C 4 Guard/4.0 (6.1050.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep C 4 Guard/4.0	6.1050.500
Metrosep C 4 S-Guard/4.0	6.1050.510
Metrosep C 4 S-Guard - 50/4.0	6.1050.530
For use with	
Metrosep C 4 - 50/4.0	6.1050.450
Metrosep C 4 - 100/4.0	6.1050.410
Metrosep C 4 - 150/4.0	6.1050.420
Metrosep C 4 - 250/4.0	6.1050.430

Metrosep C 4 Guard/2.0 (6.1050.600) Metrosep C 4 S-Guard/2.0 (6.1050.610)

The Metrosep C 4 Guard/2.0 contains the C 4 column material and is used to protect all Metrosep cation columns with 2 mm inner diameter which have a substrate based on silica gel. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be extended considerably. The economical price is an additional plus. The Metrosep C 4 Guard/2.0 also functions according to the «On Column Guard System» and is screwed directly onto the separation column with nearly no dead volume.

Applications

• Cations

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	5 x 2.0 mm
Column body	PEEK
Particle size	5 µm
Organic modifier	0–100% (no methanol)
pH range	2–7
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep C 4 S-Guard/2.0 (6.1050.610) must be used instead of the Metrosep C 4 Guard/2.0 (6.1050.600). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep C 4 Guard/2.0	6.1050.600
Metrosep C 4 S-Guard/2.0	6.1050.610
For use with	
Metrosep C 4 - 100/2.0	6.1050.210
Metrosep C 4 - 150/2.0	6.1050.220
Metrosep C 4 - 250/2.0	6.1050.230

Metrosep C 6 Guard/4.0 (6.1051.500) Metrosep C 6 S-Guard/4.0 (6.1051.510)

230

The Metrosep C 6 Guard/4.0 contains the C 6 column material and is used to protect all Metrosep cation columns that have a substrate based on silica gel. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be extended considerably. The economical price is an additional plus. The Metrosep C 6 Guard/4.0 also functions according to the «On Column Guard System» and is screwed directly onto the separation column with nearly no dead volume.

Applications

• Cations

Technical information

Silica gel with
carboxyl groups
5 x 4.0 mm
PEEK
5 µm
0–100% (no alcohol)
2–7
Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep C 6 S-Guard/4.0 (6.1051.510) must be used instead of the Metrosep C 6 Guard/4.0 (6.1051.500). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep C 6 Guard/4.0	6.1051.500
Metrosep C 6 S-Guard/4.0	6.1051.510
For use with	
Metrosep C 6 - 100/4.0	6.1051.410
Metrosep C 6 - 150/4.0	6.1051.420
Metrosep C 6 - 250/4.0	6.1051.430

Metrosep C 6 Guard/2.0 (6.01051.600) Metrosep C 6 S-Guard/2.0 (6.01051.610)

The Metrosep C 6 Guard/2.0 contains the C 6 column material and is used to protect all Metrosep cation columns with 2 mm inner diameter which have a substrate based on silica gel. Particles and contaminations are reliably retained, allowing the lifetime of the analytical separation column to be extended considerably. The economical price is an additional plus. The Metrosep C 6 Guard/2.0 also functions according to the «On Column Guard System» and is screwed directly onto the separation column with nearly no dead volume.

Applications

• Cations

Technical information

Substrate	Silica gel with
	carboxyl groups
Column dimensions	5 x 2.0 mm
Column body	PEEK
Particle size	5 µm
Organic modifier	0–100% (no alcohol)
pH range	2–7
Туре	Column



Important note

If separation columns with a length of 250 mm are used in a column oven of the «Advanced» instrument generation, then the Metrosep C 6 S-Guard/2.0 (6.01051.610) must be used instead of the Metrosep C 6 Guard/2.0 (6.01051.600). It is connected to the 250 mm column by means of capillary connection and can thus also be placed in the column oven.

Ordering information	
Metrosep C 6 Guard/2.0	6.01051.600
Metrosep C 6 S-Guard/2.0	6.01051.610
For use with	
Metrosep C 6 - 100/2.0	6.01051.210
Metrosep C 6 - 150/2.0	6.01051.220
Metrosep C 6 - 250/2.0	6.01051.230

Metrosep C Supp 1 Guard/4.0 (6.1052.500) Metrosep C Supp 1 S-Guard/4.0 (6.1052.510)

232

The Metrosep C Supp 1 Guard/4.0 contains the C Supp 1 column material and is used to protect Metrosep C Supp 1 cation columns. Particles and contaminations are reliably retained, allowing the service life of the analytical separation column to be prolonged considerably. The Metrosep C Supp 1 Guard/4.0 also functions according to the "On Column Guard System" and is screwed directly onto the separation column with nearly no dead volume.

Applications

• Cations

Technical information

Substrate	Polyvinyl alcohol with
	carboxyl groups
Column dimensions	5 x 4.0 mm
Column body	PEEK
Particle size	5 µm
Organic modifier	50 % Acetonitril or
	30 % Aceton
pH range	1-12
Туре	Column



Ordering information	
Metrosep C Supp 1 Guard/4.0	6.1052.500
Metrosep C Supp 1 S-Guard/4.0	6.1052.510
For use with	
Metrosep C Supp 1 - 100/4.0	6.1052.410
Metrosep C Supp 1 - 150/4.0	6.1052.420
Metrosep C Supp 1 - 250/4.0	6.1052.430

Metrosep C Supp 2 Guard/4.0 (6.01053.500)

The Metrosep C Supp 2 Guard/4.0 contains the C Supp 2 column material and is used to protect Metrosep C Supp 2 cation columns. Particles and contaminations are reliably retained, considerably prolonging the service life of the analytical separation column. The Metrosep C Supp 2 Guard/4.0 also functions according to the "On Column Guard System" and is screwed directly onto the respective separation column with virtually no dead volume.



Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with
	carboxyl groups
Column dimensions	5 x 4.0 mm
Column body	PEEK
Particle size	5 µm
Organic modifier	Eluent: 0–100% acetone,
	acetonitrile (no alcohol)
	Sample: 0–100% aceton,
	acetonitrile, alcohols
pH range	Eluent: 0–12
	Sample: 0–14
Туре	Column



Ord	lering	inf	ormation

Metrosep C Supp 2 Guard/4.0

For use with Metrosep C Supp 2 - 100/4.0 Metrosep C Supp 2 - 150/4.0 Metrosep C Supp 2 - 250/4.0

6.01053.500

6.01053.410 6.01053.420 6.01053.430

Metrosep RP 2 Guard/3.5 (6.1011.030)

234 The Metrosep RP 2 Guard/3.5 is a guard column for universal use. It reliably protects the analytical separation column against contamination, removing the smallest particles, traces of iron oxide, and bacteria. The Metrosep RP 2 Guard/3.5 helps to reduce costs; its filter disk can be replaced very easily.

Applications

• Universal guard column

Technical information

Substrate	Polymer
Column dimensions	1.0 x 3.5 mm
Column body	PEEK
Pore size	0.2 µm
Organic modifier	0–100%
pH range	1–13
Туре	Column



|--|

Metrosep RP 2 Guard/3.5 Replacement filters for RP 2 Guard/3.5 (10 pcs.)

For use with Phenomenex Star-IonTM A 300 - 100/4.6 Metrosep Anion Dual 2 - 75/4.6 Metrosep A Supp 1 HS - 50/4.6 Metrosep A Supp 3 - 250/4.6 Hamilton PRP-X300 - 250/4.6 Hamilton RCX-30 - 150/4.6 Metrosep Amino Acids 1 - 100/4.0 Metrosep C 5 - 150/4.6 6.1011.030 6.1011.130 6.1005.100 6.1006.100 6.1005.350 6.1005.320 6.1005.030 6.1018.010 6.1018.000 6.4001.410

6.4000.320

Metrosep RP 3 Guard HC/4.0 (6.1011.040)

The Metrosep RP 3 Guard HC/4.0 is a guard column for universal use. It reliably protects the analytical separation column against contamination, securely removing lipophilic organic contamination, the smallest particles, traces of iron oxide and bacteria. The guard column is based on a polymer material and thanks to its larger pack volume, has a significantly higher capacity that the Metrosep RP 2 Guard/3.5. It can be used in the entire pH range.

Applications

• Universal guard column

Technical information

Substrate	Styrene/divinylbenzene
	copolymer
Column dimensions	5 x 4.0 mm
Column body	PEEK
Organic modifier	0–100%
pH range	1–14
Туре	Column



Ordering information

Metrosep RP 3 Guard HC/4.0

For use with

Phenomenex Star-IonTM A 300 - 100/4.6 Metrosep Anion Dual 2 - 75/4.6 Metrosep Anion Dual 3 - 100/4.0 Metrosep A Supp 1 HS - 50/4.6 Metrosep A Supp 3 - 250/4.6 Hamilton PRP-X300 - 250/4.6 Hamilton RCX-30 - 150/4.6 Metrosep Amino Acids 1 - 100/4.0 Metrosep C 5 - 150/4.6 6.1005.100 6.1006.120 6.1005.350 6.1005.320 6.1005.030 6.1018.010 6.1018.000 6.4001.410 6.4000.320

6.1011.040

MetroSil RP 3 Guard/4.0 (6.01070.500)

236

The MetroSil RP 3 Guard/4.0 is used to protect the MetroSil RP 3 - 150/4.0 against contamination from particles and bacteria.

Applications

• Organic substances

Technical information

Substrate	Silica gel C ₁₈
Column dimensions	14 x 4.0 mm
Column body	Stainless steel
Particle size	5 µm
Organic modifier	0–100%
pH range	2–9
Туре	Cartridge



Ordering information

MetroSil RP 3 Guard/4.0 Cartridge holder for MetroSil RP 3 Guard/4.0

For use with MetroSil RP 3 - 150/4.0 6.01070.500 6.02821.010

6.01070.420

Metrosep BP 1 Guard/2.0 (6.1015.100)

The Metrosep BP 1 Guard/2.0 is used to generate a sufficiently high working pressure in the flow path of postcolumn reagents. They are used in combination with a high-pressure pump for conveying the post-column reagent in the Professional Reactor.

Applications

• Backpressure column

Technical information

Substrate	Diamond
Column dimensions	50 x 2.0 mm
Column body	PEEK
Particle size	6–10 µm
Organic modifier	0-100%
Туре	Column



237

Ordering information Metrosep BP 1 Guard/2.0

6.1015.100





Preconcentration columns

Metrosep A PCC 2/4.0 (6.1006.330), Metrosep A PCC 2 HC/4.0 (6.1006.340), and Metrosep A PCC 2 VHC/4.0 (6.1006.350)

The Metrosep A PCC 2/4.0 is used for the preconcentration of anions from small sample volumes. The small dead volume of the column guarantees an excellent peak shape.

> The Metrosep A PCC 2 HC/4.0 and the Metrosep A PCC 2 VHC/4.0, on the other hand, are high-capacity preconcentration columns for anions. They are used primarily where large sample volumes with very low anion concentrations must be preconcentrated. The high capacity prevents premature elution of the anions by the matrix itself (in most cases water). Reliable determinations can now be made using these high-capacity columns. All preconcentration columns are made of PEEK.

Applications

Substrato

• Preconcentration of anions

Technical information

Substrate
Column dimensions
Column body Maximum flow Maximum pressure Particle size Organic modifier

pH range Туре

Polymethacrylate with quaternary ammonium groups 6.1006.330: 1.0 x 4.0 mm 6.1006.340: 13.0 x 4.0 mm 6.1006.350: 30.0 x 4.0 mm PFFK 5.0 mL/min 20 MPa 65 µm Eluent: 0–10% (acetone, acetonitrile, methanol, isopropanol) Sample: 0-100% (acetone, acetonitrile, methanol, isopropanol) 2-12 Column 6.1006.330: 0.5 µmol (Cl⁻) 6.1006.340: 5 µmol (Cl⁻) 6.1006.350: 10 µmol (Cl⁻)



Ordering information

Metrosep A PCC 2/4.0 Metrosep A PCC 2 HC/4.0 Metrosep A PCC 2 VHC/4.0 6.1006.330 6.1006.340 6.1006.350

Metrosep C PCC 1/4.0 (6.1010.300), Metrosep C PCC 1 HC/4.0 (6.1010.310), and Metrosep C PCC 1 VHC/4.0 (6.1010.320)

The Metrosep C PCC 1/4.0 in the various versions are suitable for the preconcentration of monovalent and divalent cations. They are used primarily where large sample volumes with very low cation concentrations must be preconcentrated. In addition, they fulfill the function of a trap column when working with matrix elimination, i.e. the cations to be determined are held back and allow the removal of the disruptive matrix before the sample is fed to the IC system.

The greater the capacity of the column (in the first approximation, this is proportional to the length of the packing bed) the larger the sample volume which can be preconcentrated. This allows detection limits into the lower ppt range. On the other hand, the packing bed increases the dead volume of the preconcentration column; with increasing size the injection peak in the chromatogram increases in peak area. Three different capacities provide the needed flexibility for all preconcentration tasks.

The preconcentration columns are distinguished by very low noise and very low backpressure. They are suitable for preconcentration using a peristaltic pump or Metrohm Dosino technology.

Applications

• Preconcentration of cations

Technical information

Spherical polymethacrylate
with carboxyl groups
6.1010.300: 8.5 x 4.0 mm
6.1010.310: 16.5 x 4.0 mm
6.1010.320: 30.0 x 4.0 mm
PEEK
15 MPa
35 µm
0–20% methanol, ethanol,
isopropanol or acetonitrile
1-14
Column
6.1010.300: 20 mL
6.1010.310: 60 mL
6.1010.320: 90 mL

* A solution was preconcentrated with Li⁺ = 2 μ g/L, Na⁺, NH₄⁺ = 10 μ g/L and K⁺ = 20 μ g/L. The maximum preconcentration volume is determined by the fact that the peak area of the lithium does not continue to increase. This means that at greater volumes the lithium is already eluting again from the column.



Care Storage In the eluent

Ordering information

Metrosep C PCC 1/4.0 Metrosep C PCC 1 HC/4.0 Metrosep C PCC 1 VHC/4.0 6.1010.300 6.1010.310 6.1010.320

Metrosep Chel PCC 1 VHC/4.0 (6.01010.350)

The Metrosep Chel PCC 1 VHC/4.0 is a preconcentration column for the preconcentration of all kinds of cations. Due to its chelating exchanger group, it is particularly well suited for the preconcentration of polyvalent cations. Earthalkaline metals and transition metals are retained much stronger on the preconcentration column than the alkaline metals. Using matrix elimination, it is possible to remove the disturbing monovalent cations from the preconcentration column prior to injection onto the separation column.

> The chelating exchange group allows the determination of transition metals in the ppb range in matrices with high salt contents.

Applications

 Preconcentration of transition metal ions and polyvalent cations.

Technical information

Substrate	Polystyrene-divinylbenzene
	copolymer with chelating
	exchanger groups
Column dimensions	30 x 4.0 mm
Column body	PEEK
Maximum pressure	20 MPa
Particle size	30 µm
Maximum flow	2.0 mL/min
Organic modifier	0–50% acetonitrile,
	acetone, methanol
pH range	1–13
Туре	Column



Care Storage: In 5 mmol/L nitric acid.

Ordering information Metrosep Chel PCC 1 VHC/4.0

6.01010.350





IC trap columns

Metrosep A Trap 1 - 100/4.0 (6.1014.000)

246 The Metrosep A Trap 1 - 100/4.0 is a high capacity anion column, which is used to purify the eluent flow. Even reagents of the highest quality, e.g. «ultrapure» or «puriss.» can still contain minimal anionic contaminants. These are reliably held back by the Metrosep A Trap 1 - 100/4.0. This column is primarily used with gradient applications.

South Street Party

Applications

• Purification of the anion eluent flow

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with quaternary
	ammonium groups
Column dimensions	100 x 4.0 mm
Column body	PEEK
Maximum pressure	25 MPa
Particle size	570 µm
Organic modifier	0–20%
pH range	1-14
Туре	Column

Care

Regeneration

a) First rinse with 30 mL 0.5 mol/L Na_2CO_3 at a flow rate of 1.0 mL/min.

b) Then rinse with 30 mL ultrapure water at a flow of 1.0 mL/min.

Storage In the eluent

Ordering information

Metrosep A Trap 1 - 100/4.0

Metrosep C Trap 1 - 100/4.0 (6.1015.000)

This is a high capacity cation column, which is used to purify the eluent flow. Even reagents of the highest quality, e.g. «ultrapure» or «puriss.» can still contain minimal cationic contaminants. These are reliably held back by the Metrosep C Trap 1 - 100/4.0.

Participant Providence Providence

Applications

• Purification of the cation eluent flow

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with
	sulfonic acid groups
Column dimensions	100 x 4.0 mm
Column body	PEEK
Maximum pressure	25 MPa
Particle size	37–74 µm
Organic modifier	0-20%
pH range	1–14
Туре	Column

Care

Regeneration

a) First rinse with 20 mL of 5% $\rm H_2SO_4$ at a flow of 1.0 mL/min.

b) Then rinse with 15 mL ultrapure water at a flow of 1.0 mL/min.

Storage In the eluent

Ordering information Metrosep C Trap 1 - 100/4.0

6.1015.000

Metrosep C Trap 1 - 30/4.0 (6.01015.030)

This is a cation column, which is used to purify the eluent flow. Even reagents of the highest quality, e.g. «ultrapure» or «puriss.» can still contain minimal cationic contaminants. These are reliably held back by the Metrosep C Trap 1 - 30/4.0. The small dimension of the trap column reduces the dead volume of the instrument.

Applications

• Purification of the cation eluent flow

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer with
	sulfonic acid groups
Column dimensions	30 x 4.0 mm
Column body	PEEK
Maximum pressure	25 MPa
Particle size	37–74 µm
Organic modifier	0-20%
pH range	1-14
Туре	Column

Care

Regeneration a) First rinse with 20 mL of 5% H_2SO_4 at a flow of 1.0 mL/min. b) Then rinse with 15 mL ultrapure water at a flow of

1.0 mL/min.

Storage In the eluent

Ordering information Metrosep C Trap 1 - 30/4.0

6.01015.030

Metrosep RP Trap 1 - 50/4.0 (6.1014.100)

The Metrosep RP Trap 1 - 50/4.0 column is used to eliminate organic contaminants from the eluent. The Metrosep RP Trap 1 - 50/4.0 column helps avoid eluent-related interference at the baseline, especially with gradient systems. Its use is also recommended for the purification of the p-cyanophenol eluent of the Metrosep Dual 4 separation columns.

Applications

• Purification of the eluent flow

Technical information

Substrate	Silica gel
Column dimensions	50 x 4.0 mm
Column body	PEEK
Maximum pressure	25 MPa
Particle size	10 µm
pH range	1–9
Туре	Column



Care

Regeneration

- a) Rinse with 10 mL 80% acetonitrile/water at a flow rate of 2.0 mL/min.
- b) Rinse with 20 mL 100% acetonitrile at a flow rate of 2.0 mL/min.
- c) Rinse with 10 mL 80% acetonitrile/water at a flow rate of 2.0 mL/min.

Note

If the Metrosep RP Trap 1 - 50/4.0 is used with the Metrosep Dual 4 (6.1016.0X0), then it must be rinsed with 40 mL water at a flow rate of 2.0 mL/min after the regeneration.

Storage In the eluent

Ordering information

Metrosep RP Trap 1 - 50/4.0

6.1014.100

249

For use with Metrosep Dual 4 - 100/4.6

6.1016.030

Metrosep RP Trap 2 - 100/4.0 (6.1014.150)

The Metrosep RP Trap 2 - 100/4.0 column is used to eliminate organic contaminants from the eluent. The Metrosep RP Trap 2 - 100/4.0 column helps avoid eluentrelated interference at the baseline, especially with gradient systems. It is based on a polymer material. Its presence enables the use of the Metrosep RP Trap 2 -100/4.0 in the acidic as well as in the alkaline pH range.

Applications

• Elimination of organic contamination from the eluent.

Technical information

Polystyrene/divinylbenzene
copolymer
100 x 4.0 mm
PEEK
25 MPa
1–14
Column



Metrosep I Trap 1 - 100/4.0 (6.1014.200)

The Metrosep I Trap 1 - 100/4.0 column is used to eliminate ionic, i.e. cationic and anionic, contaminants from aqueous solutions. Its use is especially recommended for the purification of the transfer water in combination with «MISP» (Metrohm Inline Sample Preparation). Using the Metrosep I Trap 1 - 100/4.0 column can significantly reduce the influence of the transfer water on the system blank.

Hone Hones Inde

Applications

• For the elimination of traces of anionic and cationic contaminants from ultrapure water

Technical information

Substrate	Polystyrene/divinylbenzene copolymer with anionic and cationic ion exchangers
Column dimensions	100 x 4.0 mm
Column body	PEEK
Maximum pressure	25 MPa
Particle size	300–840 µm
Organic modifier	0–100%
pH range	0–14
Туре	Column

Care

Regeneration not possible

Storage in ultrapure water

Ordering information Metrosep I Trap 1 - 100/4.0

Metrosep BO₃³⁻ Trap 1 - 100/4.0 (6.1015.200)

252 Trap column for the removal of borate contaminants from the eluent. The Metrosep BO₃³⁺ Trap 1 - 100/4.0 is mainly used in carbohydrate analysis with hydroxide eluents. The removal of borate from the eluent improves the peak shape of sorbitol.

Applications

• Elimination of borate traces from hydroxide eluents.

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer
Column dimensions	100 x 4.0 mm
Column body	PEEK
Maximum pressure	25 MPa
pH range	0–14
Туре	Column



Care

Conditioning

Rinse the column with hydroxide eluent for 90 min at a maximum flow of 0.5 mL/min.

Regeneration

Rinse the column with the following solutions in succession in the direction against the flow:

- during 30 min with 0.1 mol/L hydrochloric acid at a flow rate of 0.3 mL/min
- during 30 min with 1 mol/L sodium chloride solution at a flow rate of 0.3 mL/min
- during 30 min with ultrapure water at a flow rate of 0.5 mL/min
- during 90 min with hydroxide eluent at a flow rate of maximum 0.5 mL/min

Storage in ultrapure water

Ordering information

Metrosep BO₃³⁻ Trap 1 - 100/4.0

Metrosep CO₃²⁻ Trap 1 - 100/4.0 (6.1015.300)

Trap column for the removal of carbonate traces in hydroxide eluents. The Metrosep CO_3^{2-} Trap 1 - 100/4.0 is used in carbohydrate analysis with hydroxide eluents.

Applications

• Elimination of carbonate contamination from hydroxide eluents

Technical information

Substrate	Polystyrene/divinylbenzene
	copolymer
Column dimensions	100 x 4.0 mm
Column body	PEEK
Maximum pressure	25 MPa
pH range	0-14
Туре	Column

The second secon Recommended eluent concentration 5-40 mmol/L hydroxide eluent

Regeneration

Rinse the column during 840 min with 0.3 mol/L sodium hydroxide at a flow rate of 1 mL/min.

Storage Store the column in 0.3 mol/L sodium hydroxide.

Ordering information

Metrosep CO₃²⁻ Trap 1 - 100/4.0



IC sample-preparation cartridges

IC-RP sample-preparation cartridge (6.1012.X00)

Material	RP	
Application	For the non-polar solid-phase	e extraction. The cartridge removes organic substances.
Quantity	50	10
Bed volume	0.5 mL	0.5 mL
Connection	Luer	Luer
Order number	6.1012.000	6.1012.100

IC-H sample-preparation	cartridge (6.1012.X10)		
Material	Cation exchanger in acid for	m	
Application	For the elimination of interfering cations. The cartridge can also be used for the neutralization of alkaline samples.		
Quantity	50	10	25
Bed volume	0.5 mL	0.5 mL	1.5 mL
Capacity	0.8 mmol	0.8 mmol	2.0 mmol
Connection	Luer	Luer	Luer
Order number	6.1012.010	6.1012.110	6.1012.210

IC-Ag sample-preparation cartridge (6.1012.X20)			
Material	Cation exchanger in silver form		
Application	For the removal of halides.		
Quantity	50	10	25
Bed volume	0.5 mL	0.5 mL	1.5 mL
Capacity	0.8 mmol	0.8 mmol	2.0 mmol
Connection	Luer	Luer	Luer
Order number	6.1012.020	6.1012.120	6.1012.220

IC-OH sample-preparation cartridge (6.1012.X30)		
Material	Anion exchanger in hydroxid	de form
Application	For the neutralization of ext	remely acidic samples.
Quantity	50	10
Bed volume	0.5 mL	0.5 mL
Capacity	0.6 mmol	0.6 mmol
Connection	Luer	Luer
Order number	6.1012.030	6.1012.130

IC-Na sample-preparation cartridge (6.1012.X40)		
Material	Cation exchanger in sodium form	
Application	For the elimination of cations.	
Quantity	50	
Bed volume	0.5 mL	
Capacity	0.8 mmol	
Connection	Luer	
Order number	6.1012.040	

IC-C18 sample-preparation cartridge (6.1012.X50)		
Material	C18	
Application	For the removal of polar substances; not suitable for F^- determination.	
Quantity	50	
Bed volume	0.5 mL	
Connection	Luer	
Order number	6.1012.050	



IC accessory parts

PEEK inline filter (6.2821.120)

The inline filter in the PEEK housing not only removes all particles of mineral origin, but also algae and bacteria. The exclusion diameter of 2 μ m ensures that no contamination can damage the column or the suppressor.



Ordering information PEEK inline filter Replacement filters (10 pcs.)

6.2821.1206.2821.130

Coupling safety olive with PEEK inline filter (6.2744.180)

The coupling safety olive with PEEK inline filter connects the tube of the peristaltic pump with the following system, e.g. with the suppressor or with the post-column reactor (PCR). On the one hand, this prevents the tube of the peristaltic pump from detaching unintentionally, while on the other hand all particles with a diameter of greater than 2 μ m are effectively removed from the flow of liquid.



Ordering information Coupling safety olive with PEEK inline filter Replacement filters (10 pcs.)

6.2744.1806.2821.130

Subject to change without notice. Design Edmaner-Schoch ASW, Printing Metrohm AG, CH-9100 Herisau, Switzerland 8.000.5347EN – 2021-05 www.metrohm.com

