

Säulenhandbuch



Metrosep A Supp 20 (6.01035.4x0)

Manual

8.0107.8018EN / v1 / 2026-03-27



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

Säulenhandbuch

Metrosep A Supp 20 (6.01035.4x0)

Manual

Technical Communication
Metrohm AG
CH-9100 Herisau

This documentation is protected by copyright. All rights reserved.

This documentation is an original document.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g., conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	General information	1
1.1	Ordering information	1
1.2	Technical specifications	1
2	Key aspects of working with separation columns	3
3	Eluent production	7
3.1	Production of standard eluent	7
4	Start-up	8
4.1	Connecting and rinsing the guard column	8
4.2	Connecting and rinsing the separation column	10
4.3	Conditioning	13
5	Applications	15
5.1	Standard chromatogram	15
5.2	Temperature influence	16
5.3	Variation of the eluent flow rate	19
5.4	Variation of the eluent	21
5.4.1	Constant Na ₂ CO ₃ /NaHCO ₃ ratio under EPA 300.1, Part A and Part B	21
5.4.2	Constant Na ₂ CO ₃ /NaHCO ₃ ratio under ISO 10304, Part 1 and Part 4 – Conditions	22
5.4.3	NaHCO ₃ variation at constant Na ₂ CO ₃	24
5.4.4	Na ₂ CO ₃ variation at constant NaHCO ₃	25
5.4.5	NaOH variation at constant Na ₂ CO ₃	27
5.5	Variation with organic modifier	28
5.5.1	Variation of the acetone concentration	28
5.5.2	Variation of the methanol concentration	30
5.5.3	Variation of the acetonitrile concentration	32
5.6	Determination of standard anions and oxyhalides (chlorite, bromate, chlorate) in mineral water samples according to ISO 10304, Part 1 and Part 4.	33
5.6.1	Fast version	33
5.6.2	Improved matrix compatibility	34
5.7	Determination of standard anions, oxyhalides (chlorite, bromate, chlorate), and dichloroacetates in mineral water samples according to US EPA 300.1, part A and part B	35
5.7.1	Fast version	35



5.7.2	Improved matrix compatibility	36
6	Troubleshooting	38
6.1	Regeneration	38
6.2	Decreasing resolution and asymmetrical peaks	39
6.3	Unstable retention times	40
6.4	Unknown peaks	40
6.5	Increasing backpressure	41
7	Literature	42
	Index	43

1 General information

The Metrosep A Supp 20 is a high-performance separation column and is particularly suitable for the determination of inorganic anions and disinfection byproducts such as chlorite, bromate, and chlorate with chemical and sequential suppression. Due to their high capacity, the Metrosep A Supp 20 can process samples with high ionic strength and large fluctuations in concentration without any problems. The outstanding peak symmetry and the high number of theoretical bases enable universal use in ion chromatography.

1.1 Ordering information

Table 1 Separation columns

Order number	Designation
6.01035.420	Metrosep A Supp 20 - 150/4.0
6.01035.430	Metrosep A Supp 20 - 250/4.0

Table 2 Guard column

Order number	Designation
6.01035.500	Metrosep A Supp 20 Guard/4.0

1.2 Technical specifications

Column material Hydrophilic polystyrene/divinylbenzene copolymer with quaternary ammonium groups

Particle size 4.6 µm

<i>Measurements</i>	Order number	Measurements
	6.01035.420	150 x 4.0 mm
	6.01035.430	250 x 4.0 mm

pH range eluent 0–14

ph range sample 0–14

Temperature range 10–70 °C

Recommended standard temperature 30 °C



Maximum pressure

Order number	Maximum pressure
6.01035.420	25 MPa (250 bar)
6.01035.430	25 MPa (250 bar)

Flow rate

Order number	Recommended flow rate	Maximum flow rate
6.01035.420	0.85 mL/min	1.50 mL/min
6.01035.430	0.75 mL/min	1.20 mL/min

Standard eluent

5.6 mmol/L sodium carbonate (Na_2CO_3) and 3.1 mmol/L sodiumhydrogencarbonate (NaHCO_3)

Permitted organic additives

0–100% acetonitrile, acetone and methanol

Capacity

Order number	Capacity
6.01035.420	187 μmol (Cl^-)
6.01035.430	312 μmol (Cl^-)

Preparation

Use a flow gradient to set the column to the standard flow within 5 minutes. Then wait until the baseline is given.

Storage

Store the column in the standard eluent and at 4–30 °C.

Typical pressure

For columns with a guard column under standard conditions with chemical suppression:

Order number	Typical pressure
6.01035.420	15 MPa
6.01035.430	18 MPa

Column housing

Smart column with a chip, called an iColumn, made of PEEK

Application

Determination of inorganic anions and disinfection byproducts such as chlorite, bromate, and chlorate with chemical and sequential suppression



2 Key aspects of working with separation columns

- Storage* Once the backpressure in the ion chromatograph has dissipated, remove the column at ambient temperature. Seal the column at both ends using the original stoppers (6.2744.060). Store the column in the standard eluent and at 4 to 30 °C.
- Bacterial growth* Bacterial growth significantly affects chromatography and ruins separation columns. A vast array of problems in chromatography are caused by the growth of algae, bacteria and fungi.
- In order to prevent bacterial growth, always use fresh eluents, rinsing solutions and regeneration solutions. Do not use any solutions that have not been used for a prolonged period. Metrohm recommends cleaning all vessels as follows before filling them:
1. Thoroughly rinse with ultrapure, UV-treated water (> 18.2 MΩ).
 2. Swirl an acetonitrile-water mixture around in the vessel.
 3. Rinse again with ultrapure water.
- If you notice the growth of bacteria or algae despite these precautionary measures, add 5% methanol, acetonitrile or acetone to the eluent. This is only possible if you are *not using membrane suppressors*. Organic solvents can destroy membrane suppressors. The Metrohm Suppressor Modules (MSM, MSM-HC and MSM-LC) are 100% resistant to solvents.
- Chemical quality* All chemicals must have at least a quality of p.a. or puriss. Standard solutions must be intended specifically for ion chromatography.
- Chemical stress* Even though specifications may indicate that separation phases do cover a large pH range, this does not mean they are chemically inert. Separation columns last longest when subjected to constant chemical conditions. Never allow a column to dry out and ensure the column is sealed well at all times.
- CO₂* Carbon dioxide from the air affects the carbonate / hydrogen carbonate balance in the eluent. The eluent becomes weaker over time. In order to prevent this, always outfit the eluent bottle with a CO₂ adsorber with the adsorber material soda lime.
- Degassing the eluent* In order to prevent bubbles from forming, degas the produced eluent before using it in the IC system. To do this, create a vacuum for approximately 10 minutes using a water-jet pump or a membrane pump. Alternatively, use an ultrasonic bath or work with an eluent degasser.

- switching the valves. Using the pulsation absorber (6.2620.150) already built into the Metrohm ion chromatographs provides this protection.
- Regenerating separation columns* If separation columns are operated with clean eluents and filled with samples free of particles, you can expect the column to have a long service life. This means regenerating the column is not required and, after a multitude of injections, no longer possible.
- If the pressure in the column increases despite this or if the separating efficiency decreases, carry out the specified regeneration steps. Perform the regeneration outside the analytical line. For regeneration, connect the separation column to the pump directly. Route the regeneration solution through the column directly into a waste container. Rinse the separation column properly with fresh eluent. Then, reinstall the separation column.
- Shutting down the ion chromatograph* If the ion chromatograph is not used for a prolonged period (> 1 week), remove the separation column and seal it with the stoppers provided. Rinse the ion chromatograph, including all 3 suppressor chambers, with a methanol-water mixture (1:4). Store the separation column in the medium indicated on the column leaflet. Unless otherwise stated on the column leaflet, store the column at 4 to 30 °C.
- Prior to start-up, rinse the ion chromatograph with ultrapure water and then with fresh eluent. Bring the separation column back to ambient temperature before you install it. Then increase the temperature if necessary.
- Fun* Ion chromatography should be fun and should not stress you out. Metrohm puts all its work into ensuring you can work reliably with your IC systems with minimal maintenance, servicing and costs. Metrosep separation columns embody the attributes of quality, long service life and excellent results.
- Environmental protection* A significant advantage of ion chromatography is that most of the work involves aqueous media. As a result, the chemicals used in ion chromatography are largely non-toxic and do not impact the environment. When working with acids, bases, organic solvents or heavy metal standards, dispose of them properly after use.
- Guard columns* Guard columns are used to protect separation columns. Metrohm strongly recommends using guard columns. Guard columns normally contain the same stationary phase as separation columns. However, the quantity is significantly reduced to avoid impacting the chromatography. Guard columns remove critical contaminants that can react with column material. Guard columns also remove particles and bacterial contaminants. Replace the guard column in the following cases:
- If the backpressure in the system increases



- If the chromatography results deteriorate

Replace the guard column 3 to 4 times during the service life of the analytical column. Guard columns are available for all Metrosep separation columns.

Water quality

Aqueous media are mostly used in work involving ion chromatography. This means that water quality is a critical factor for good chromatography. If the water quality is inadequate, the results will be inadequate as well. Water with inadequate quality can damage instruments and separation columns. The ultrapure water being used must have a specific resistance greater than 18.2 M Ω -cm and should be free of particles. Therefore, filter the water using a 0.45- μ m filter and treat it with UV light. Modern ultrapure water systems for laboratory use ensure this level of water quality (Type I).

3 Eluent production

Metrohm recommends selecting a high degree of purity for chemicals for standard production and eluent production.

3.1 Production of standard eluent

Proceed as follows to produce 2 L of standard eluent with 3.1 mmol/L of sodium hydrogen carbonate and 5.6 mmol/L of sodium carbonate:

Producing 2 L of standard eluent

Accessories

- Eluent bottle (6.1608.120)
- Lid (6.1602.200) equipped with CO₂ adsorber
- Ultrapure water
- Sodium carbonate
- Sodium hydrogen carbonate

1 Pre-rinse the eluent bottle with ultrapure water several times.

2 Fill 2 L of ultrapure water into the eluent bottle.

3 Degas the ultrapure water.

Use the eluent degasser.

If no eluent degasser is available, degas the ultrapure water for 5–10 minutes using a vacuum pump. Degassing prevents problems with air bubbles in the high-pressure pump.

4

- Weigh 520.9 mg of sodium hydrogen carbonate.
- Weigh 1,187.1 mg of sodium carbonate.
- Add the weighed amounts of sodium hydrogen carbonate and sodium carbonate to the ultrapure water.

5 Rinse the column with eluent for 1 h.

This eluent (3.1 mmol/L of sodium hydrogen carbonate and 5.6 mmol/L of sodium carbonate) and chemical suppression can be used to achieve background conductivity of less than 23 $\mu\text{S}/\text{cm}$. The noise is typically less than 0.2 nS/cm.

4 Start-up

4.1 Connecting and rinsing the guard column

Guard columns protect separation columns and significantly increase their service life. The guard columns available from Metrohm are either actual guard columns or guard column cartridges used together with a cartridge holder. The process of installing a guard column cartridge into the corresponding holder is described in the cartridge leaflet.



NOTE

Metrohm recommends always working with guard columns. Guard columns protect the separation columns and can be replaced regularly as needed.



NOTE

Information regarding which guard column is suitable for your separation column can be found in the **Metrohm Column Program** (which is available from your regional Metrohm representative), the column leaflet and the product information or in consultation with your regional Metrohm representative.

You can find product information for your separation column at <http://www.metrohm.com> in the Ion Chromatography product area.



CAUTION

New guard columns are filled with solution and sealed with stoppers or caps on both sides.

Before inserting the guard column, ensure that this solution can be mixed with the eluent being used (follow the manufacturer specification).



NOTE


Do not connect the guard column until after the initial start-up of the instrument. Until then, replace the guard column and the separation column with couplings (6.2744.040).

Accessories

For this step, you need the following accessories:

- Guard column (suitable for separation column)

Connecting the guard column



1 Removing the coupling

Remove the coupling (6.2744.040) installed between the column inlet capillary and the column outlet capillary for the initial start-up.

2 Preparing the guard column

- Remove the stoppers or the stopper and the sealing cap from the guard column.

3 Connecting the guard column



CAUTION

When inserting the guard column, ensure that it is inserted correctly based on the marked flow direction (if specified).

- Fasten the inlet of the guard column to the column inlet capillary using a short pressure screw (6.2744.070).
- If the guard column is connected to the separation column using a connection capillary, fasten this connection capillary to the guard column outlet with a pressure screw.


Rinsing the guard column

1 Rinsing the guard column

- Place a beaker under the guard column's outlet.

**NOTE**

Connect the separation column only after the initial start-up of the instrument. Until that point, insert a coupling (6.2744.040) instead of the guard column and separation column.

**Connecting the separation column****1 Removing the stoppers**

- Remove the stoppers from the separation column.

4.3 Conditioning

In the following cases, the system must be conditioned with eluent until a stable baseline has been reached:

- After installation
- After each time the instrument is switched on
- After each eluent change



NOTE

The conditioning time can lengthen considerably if the composition of the eluent is modified.

Conditioning the system

1 Preparing the software



CAUTION

Ensure that the configured flow rate is not higher than the flow rate permitted for the corresponding column (refer to the column leaflet and chip data record).

- Start the **MagIC Net** computer program.
- Open the **Equilibration** tab in MagIC Net: **Workplace ► Run ► Equilibration**.
- Select (or create) a suitable method.
See also: *MagIC Net tutorial* and online help.

2 Preparing the instrument

- Check whether the column is inserted correctly in accordance with the flow direction marked on the sticker (arrow has to point in the flow direction).
- Check whether the eluent aspiration tubing is immersed in the eluent and that there is enough eluent in the eluent bottle.

3 Starting the equilibration

- Start the equilibration in MagIC Net: **Workplace ► Run ► Equilibration ► Start HW**.



- Visually inspect whether all capillaries and their connections from the high-pressure pump to the detector are leak-tight. If eluent is leaking out anywhere, tighten the corresponding pressure screw further, or loosen the pressure screw, check the end of the capillary and shorten it using the capillary cutter if necessary and retighten the pressure screw.

4 Conditioning the system


Continue rinsing the system with eluent until the desired stability level for the baseline has been attained .

The instrument is now ready for measuring samples.

5 Applications

5.1 Standard chromatogram

- Columns:*
- Metrosep A Supp 20 - 150/4.0
 - Metrosep A Supp 20 - 250/4.0
- Sample preparation:* –
- Detection:* Conductivity
- Suppression:*
- 150/4.0: Sequential suppression with MSM A and MCS
 - 250/4.0: Sequential suppression with MSM-HC A and MCS
- Temperature:* 30 °C
- Loop:* 50 µL
- Flow rate:*
- 150/4.0: 0.85 mL/min
 - 250/4.0: 0.75 mL/min
- Eluent:* 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃



Metrosep A Supp 20 - xxx/4.0		mg/L
1	Fluoride	0.1
2	Chlorite	0.1




	Metrosep A Supp 20 - xxx/4.0	mg/L
3	Bromate	0.1
4	Chloride	1
5	Nitrite	0.1
6	Bromide	0.1
7	Chlorate	0.1
8	Nitrate	1
9	Phosphate	1
10	Dichloroacetate	1
11	Sulfate	1

5.2 Temperature influence

US EPA 300.1, Part A and Part B

<i>Column:</i>	Metrosep A Supp 20 - 150/4.0
<i>Sample preparation:</i>	–
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM A and MCS
<i>Temperature:</i>	10–70 °C
<i>Loop:</i>	50 µL
<i>Flow rate:</i>	0.85 mL/min
<i>Eluent:</i>	3.1 mmol/L NaHCO ₃ , 5.6 mmol/L Na ₂ CO ₃



Metrosep A Supp 20 - 150/4.0		mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5


The Metrosep A Supp 20 can be used at temperatures of 10–70 °C. When the temperature increases, the retention times of the monovalent anions decrease. At higher temperatures, the separation between chlorite and bromate deteriorates while the separation between bromide and chlorate improves. Therefore, 30 °C is the optimum temperature for this application. The retention times of phosphate and sulfate increase significantly as the temperature increases. At 10 °C, phosphate and chlorate are not completely separated. Dichloroacetate moves to the beginning of the chromatogram with increasing temperature: At 10 °C it elutes just before sulfate, while at 70 °C it elutes between nitrate and phosphate and the peak is significantly higher.



Increasing the temperature also causes the column backpressure to decrease considerably. At 70 °C, the backpressure is only approx. 7.3 MPa, whereas the column backpressure at 10 °C is more than twice as high (17.8 MPa).

DIN ISO 10304, Part 1 and Part 4.

- Column:* Metrosep A Supp 20 - 150/4.0
- Sample preparation:* –
- Detection:* Conductivity
- Suppression:* Sequential suppression with MSM A and MCS
- Temperature:* 10–70 °C
- Loop:* 50 µL
- Flow rate:* 0.85 mL/min
- Eluent:* 1.6 mmol/L NaHCO₃, 6.5 mmol/L Na₂CO₃



Metrosep A Supp 20 - 150/4.0		mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1


	Metrosep A Supp 20 - 150/4.0	mg/L
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Sulfate	5

The behavior of the temperature change under ISO 10304, Part 1 and Part 4 is very similar to that under EPA 300.1, Part A and Part B. At low temperatures, nitrate and phosphate coelute, while at high temperatures chlorate is very close to nitrate. Therefore, 25–30 °C is the optimum temperature for this application. At higher temperatures, the carbonate peak moves towards bromide and chlorate, which can interfere with the quantification of chlorate.

5.3 Variation of the eluent flow rate

Metrosep A Supp 20 - 150/4.0

<i>Column:</i>	Metrosep A Supp 20 - 150/4.0
<i>Sample preparation:</i>	–
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM A and MCS
<i>Temperature:</i>	30 °C
<i>Loop:</i>	50 µL
<i>Flow rate:</i>	0.3–1.5 mL/min
<i>Eluent:</i>	1.6 mmol/L NaHCO ₃ , 6.5 mmol/L Na ₂ CO ₃



Metrosep A Supp 20 - 150/4.0		mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Sulfate	5

The Metrosep A Supp 20 - 150/4.0 can be operated with a flow rate of up to 1.5 mL/min. As the flow rate increases, all ions are accelerated evenly, with sulfate eluting at 1.5 mL/min in less than 15 minutes. The pressure increases almost proportionally to the flow rate. Due to the higher flow rate, the dwell time of the analytes in the detector is shortened, resulting in smaller peak areas. The Metrosep A Supp 20 - 250/4.0 can be operated with a maximum flow of 1.2 mL/min.

5.4 Variation of the eluent

5.4.1 Constant $\text{Na}_2\text{CO}_3/\text{NaHCO}_3$ ratio under EPA 300.1, Part A and Part B

Column: Metrosep A Supp 20 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM A and MCS


Temperature: 30 °C

Loop: 50 μL

Flow rate: 0.85 mL/min

Eluent:

- 2.214 mmol/L NaHCO_3 , 4.0 mmol/L Na_2CO_3
- 2.657 mmol/L NaHCO_3 , 4.8 mmol/L Na_2CO_3
- 3.100 mmol/L NaHCO_3 , 5.6 mmol/L Na_2CO_3
- 3.543 mmol/L NaHCO_3 , 6.4 mmol/L Na_2CO_3
- 3.986 mmol/L NaHCO_3 , 7.2 mmol/L Na_2CO_3






	Metrosep A Supp 20 - 150/4.0	mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5

With increasing eluent concentration, all anions are significantly accelerated. The monovalent anions react less strongly than the polyvalent anions such as phosphate and sulfate. Consequently, phosphate elutes at high eluent concentrations (3.986 mmol/L NaHCO₃, 7.2 mmol/L Na₂CO₃) shortly before nitrate, while sulfate elutes together with dichloroacetate. For the monovalent anions, the selectivity does not change significantly with the change in the eluent strength. An eluent composition of 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃ is considered to be the ideal condition for the EPA 300.1, Part A and Part B application.

5.4.2 Constant Na₂CO₃/NaHCO₃ ratio under ISO 10304, Part 1 and Part 4 – Conditions

<i>Column:</i>	Metrosep A Supp 20 - 150/4.0
<i>Sample preparation:</i>	–
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM A and MCS
<i>Temperature:</i>	30 °C
<i>Loop:</i>	50 µL
<i>Flow rate:</i>	0.85 mL/min
<i>Eluent:</i>	A) 0.985 mmol/L NaHCO ₃ , 4.00 mmol/L Na ₂ CO ₃ B) 1.292 mmol/L NaHCO ₃ , 5.25 mmol/L Na ₂ CO ₃

- C) 1.600 mmol/L NaHCO₃, 6.50 mmol/L Na₂CO₃
 D) 1.908 mmol/L NaHCO₃, 7.75 mmol/L Na₂CO₃
 E) 2.215 mmol/L NaHCO₃, 9.00 mmol/L Na₂CO₃



Metrosep A Supp 20 - 150/4.0		mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Sulfate	5


An increasing eluent concentration leads to shorter chromatograms, as all anions are eluted faster. The polyvalent anions phosphate and sulfate react more strongly than the monovalent anions. This can lead to co-elutions: In cases of strong eluent, phosphate elutes between chlorate and nitrate; sulfate elutes immediately after nitrate. At low eluent concentrations, however, the peaks are well resolved and the chromatogram duration increases to almost 40 minutes. An eluent composition of 1.6 mmol/L



NaHCO₃, 6.5 mmol/L Na₂CO₃ is considered to be the optimum condition for the application according to ISO 10304, Part 1 and Part 4.

5.4.3 NaHCO₃ variation at constant Na₂CO₃

- Column:* Metrosep A Supp 20 - 150/4.0
- Sample preparation:* –
- Detection:* Conductivity
- Suppression:* Sequential suppression with MSM A and MCS
- Temperature:* 30 °C
- Loop:* 50 µL
- Flow rate:* 0.85 mL/min
- Eluent:*
 - A) 0.775 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃
 - B) 1.550 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃
 - C) 3.100 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃
 - D) 6.200 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃



Metrosep A Supp 20 - 150/4.0		mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2

	Metrosep A Supp 20 - 150/4.0	mg/L
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5


The effect of sodium hydrogen carbonate on the chromatogram depends mainly on the pH value of the eluent. The elution strength of sodium hydrogen carbonate is much lower than the elution strength of sodium carbonate. With increasing sodium hydrogen carbonate content, the pH value of the eluent shifts into the acidic range. This affects the polyvalent anions phosphate and sulfate. Phosphate in particular is shifted to the beginning of the chromatogram, as its effective charge is reduced at a lower pH value of the eluent. At 6.2 mmol/L NaHCO₃, phosphate elutes under nitrate, while at 0.775 mmol/L Na₂CO₃ it elutes between dichloroacetate and sulfate. Dichloroacetate elutes with a low pH value somewhat later.

5.4.4 Na₂CO₃ variation at constant NaHCO₃

<i>Column:</i>	Metrosep A Supp 20 - 150/4.0
<i>Sample preparation:</i>	–
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM A and MCS
<i>Temperature:</i>	30 °C
<i>Loop:</i>	50 µL
<i>Flow rate:</i>	0.85 mL/min
<i>Eluent:</i>	A) 3.1 mmol/L NaHCO ₃ , 4.0 mmol/L Na ₂ CO ₃ B) 3.1 mmol/L NaHCO ₃ , 4.8 mmol/L Na ₂ CO ₃ C) 3.1 mmol/L NaHCO ₃ , 5.6 mmol/L Na ₂ CO ₃ D) 3.1 mmol/L NaHCO ₃ , 6.4 mmol/L Na ₂ CO ₃



E) 3.1 mmol/L NaHCO₃, 7.2 mmol/L Na₂CO₃



	Metrosep A Supp 20 - 150/4.0	mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5

Sodium carbonate has a strong elution strength and has a considerable effect on the retention times of all anions. Phosphate and sulfate react most strongly. The retention times of all anions decrease with increasing sodium carbonate concentration. Due to the strong elution strength, the behavior in this setting is very similar to that in chapter 5.4.1 (see chapter 5.4.1, page 21): In high sodium carbonate concentrations, phosphate is under nitrate and dichloroacetate coelutes with sulfate, while phosphate coelutes with dichloroacetate in low sodium concentrations.

5.4.5 NaOH variation at constant Na₂CO₃

Column: Metrosep A Supp 20 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM A and MCS


Temperature: 30 °C

Loop: 50 µL

Flow rate: 0.85 mL/min

Eluent:

- A) 0.00 mmol/L NaOH, 5.6 mmol/L Na₂CO₃
- B) 0.25 mmol/L NaOH, 5.6 mmol/L Na₂CO₃
- C) 0.50 mmol/L NaOH, 5.6 mmol/L Na₂CO₃
- D) 1.00 mmol/L NaOH, 5.6 mmol/L Na₂CO₃
- E) 2.00 mmol/L NaOH, 5.6 mmol/L Na₂CO₃



Metrosep A Supp 20 - 150/4.0		mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2



	Metrosep A Supp 20 - 150/4.0	mg/L
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5

Sodium hydroxide can also be used in the eluent as an elution component instead of sodium hydrogen carbonate. This means that the pH value of the eluent shifts to a more alkaline range. With increasing sodium hydroxide concentration, the retention times of all anions are shortened due to the stronger eluent. The only exception is phosphate, the effective charge of which increases at a higher pH value of the eluent. The resolution between bromide and chlorate is reduced at a higher pH value of the eluent.

5.5 Variation with organic modifier

5.5.1 Variation of the acetone concentration

Column: Metrosep A Supp 20 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM A and MCS

Temperature: 30 °C


Loop: 50 µL

Flow rate: 0.85 mL/min

Eluent:

- A) 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃, 0% acetone
- B) 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃, 10% acetone
- C) 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃, 20% acetone





	Metrosep A Supp 20 - 150/4.0	mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5

In some cases, the use of an organic modifier is useful or even necessary. Use cases for an organic modifier can be:

- To increase the stability of the eluent against bacterial contamination.
- To rinse organic parts of a sample out of the separation column and to avoid contamination of the column.
- To change in the selectivity of the column to optimize the separation.
- To increase the ionizability of the ions during IC-MS couplings.


The separation selectivity is influenced by the addition of an organic modifier. In the case of acetone, all monovalent anions with shorter retention



times will elute when acetone is present. With an acetone content of 10% in the eluent, the dichloroacetate is greatly accelerated and coelutes with nitrate. Furthermore, the separation between chlorate and nitrate is not ideal. With an acetone content of 20%, all three anions (chlorate, nitrate, and dichloroacetate) coelute, and the separation between bromate and chloride is not ideal. In contrast to the monovalent anions, the polyvalent anions phosphate and sulfate react hardly at all to the addition of acetone to the eluent. The addition of an organic modifier generally increases the eluent viscosity, which leads to a higher backpressure in the column. When 20% acetone is added to the eluent, the pressure increases by approx. 50%.

5.5.2 Variation of the methanol concentration

<i>Column:</i>	Metrosep A Supp 20 - 150/4.0
<i>Sample preparation:</i>	–
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM A and MCS
<i>Temperature:</i>	30 °C
<i>Loop:</i>	50 µL
<i>Flow rate:</i>	0.85 mL/min
<i>Eluent:</i>	A) 3.1 mmol/L NaHCO ₃ , 5.6 mmol/L Na ₂ CO ₃ , 0% methanol B) 3.1 mmol/L NaHCO ₃ , 5.6 mmol/L Na ₂ CO ₃ , 10% methanol C) 3.1 mmol/L NaHCO ₃ , 5.6 mmol/L Na ₂ CO ₃ , 20% methanol



Metrosep A Supp 20 - 150/4.0		mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5

The addition of methanol to the eluent has an effect very similar to the addition of acetone. All monovalent anions elute earlier, while the polyvalent anions are hardly affected. At 10% methanol, dichloroacetate is in the tail end of the nitrate peak. If methanol is increased to 20%, then dichloroacetate will coelute with chlorate and bromate will shift into the chloride peak. In terms of pressure, adding 20% methanol to the eluent is expected to result in an increase of approx. 50%.



5.5.3 Variation of the acetonitrile concentration

Column: Metrosep A Supp 20 - 150/4.0

Sample preparation: –

Detection: Conductivity


Suppression: Sequential suppression with MSM A and MCS

Temperature: 30 °C

Loop: 50 µL

Flow rate: 0.85 mL/min

Eluent:
 A) 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃, 0% acetonitrile
 B) 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃, 10% acetonitrile
 C) 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃, 20% acetonitrile



	Metrosep A Supp 20 - 150/4.0	mg/L
1	Fluoride	1
2	Chlorite	1
3	Bromate	2
4	Chloride	1
5	Nitrite	2

	Metrosep A Supp 20 - 150/4.0	mg/L
6	Bromide	5
7	Chlorate	4
8	Nitrate	5
9	Phosphate	5
10	Dichloroacetate	10
11	Sulfate	5


Acetonitrile is also often used as an organic modifier in the eluent. Here, too, the effect is very similar: The monovalent anions are accelerated with increasing acetonitrile concentration. At 10% acetonitrile in the eluent, dichloroacetate shifts to the end of the nitrate peak, while at 20% acetonitrile it completely coelutes with nitrate. At this concentration, chlorate also coelutes with nitrate. With acetonitrile, the pressure increase is less pronounced than with acetone and methanol: For eluents with 20% acetonitrile, the increase is only approx. 25%.

5.6 Determination of standard anions and oxyhalides (chlorite, bromate, chlorate) in mineral water samples according to ISO 10304, Part 1 and Part 4.

5.6.1 Fast version

<i>Column:</i>	Metrosep A Supp 20 - 150/4.0
<i>Sample preparation:</i>	–
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM A and MCS
<i>Temperature:</i>	25 °C
<i>Loop:</i>	50 µL
<i>Flow rate:</i>	0.85 mL/min
<i>Eluent:</i>	3.1 mmol/L NaHCO ₃ , 5.6 mmol/L Na ₂ CO ₃


5.6 Determination of standard anions and oxyhalides (chlorite, bromate, chlorate) in mineral water samples according to ISO 10304, Part 1 and Part 4.



The gray chromatogram corresponds to a directly injected mineral water sample. The black chromatogram shows the same sample, enriched with 10 mg/L chloride, nitrate, and sulfate, 1 mg/L phosphate, 100 $\mu\text{g}/\text{L}$ fluoride, nitrite, and bromide, and 5 $\mu\text{g}/\text{L}$ chlorite, bromate, and chlorate. Under the given conditions, sample matrices with up to 250 mg/L chloride, 50 mg/L nitrate, 250 mg/L sulfate, and 300 mg/L of hydrogen carbonate can be analyzed.

5.6.2 Improved matrix compatibility

Column:	Metrosep A Supp 20 - 250/4.0
Sample preparation:	–
Detection:	Conductivity
Suppression:	Sequential suppression with MSM A and MCS
Temperature:	30 °C
Loop:	50 μL
Flow rate:	0.75 mL/min
Eluent:	1.6 mmol/L NaHCO_3 , 6.5 mmol/L Na_2CO_3



The gray chromatogram corresponds to a directly injected mineral water sample. The black chromatogram shows the same sample with 10 mg/L chloride, nitrate, and sulfate, 1 mg/L phosphate, 100 µg/L fluoride, nitrite, and bromide, and 5 µg/L chlorite, bromate, and chlorate. Under the given conditions, sample matrices can be analyzed with up to 500 mg/L chloride, 100 mg/L nitrate, 500 mg/L sulfate, and 300 mg/L hydrogen carbonate.

5.7 Determination of standard anions, oxyhalides (chlorite, bromate, chlorate), and dichloroacetates in mineral water samples according to US EPA 300.1, part A and part B

5.7.1 Fast version

Column: Metrosep A Supp 20 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM A and MCS


Temperature: 30 °C

Loop: 50 µL

Flow rate: 0.85 mL/min

5.7 Determination of standard anions, oxyhalides (chlorite, bromate, chlorate), and dichloroacetates in mineral water samples according to US EPA 300.1, part A and part B

Eluent: 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃



The gray chromatogram corresponds to a directly injected mineral water sample. The black chromatogram shows the same sample with 10 mg/L chloride, nitrate, and sulfate, 1 mg/L phosphate and dichloroacetate, 100 $\mu\text{g}/\text{L}$ fluoride, nitrite, and bromide, and 5 $\mu\text{g}/\text{L}$ chlorite, bromate, and chlorate. Under the given conditions, sample matrices can be analyzed with up to 250 mg/L chloride, 50 mg/L nitrate, 250 mg/L sulfate, and 300 mg/L hydrogen carbonate.

5.7.2 Improved matrix compatibility

Column: Metrosep A Supp 20 - 250/4.0

Sample preparation: –

Detection: Conductivity


Suppression: Sequential suppression with MSM A and MCS

Temperature: 30 °C

Loop: 50 μL

Flow rate: 0.75 mL/min

Eluent: 3.1 mmol/L NaHCO₃, 5.6 mmol/L Na₂CO₃



The gray chromatogram corresponds to a directly injected mineral water sample. The black chromatogram shows the same sample with 10 mg/L chloride, nitrate, and sulfate, 1 mg/L phosphate and dichloroacetate, 100 $\mu\text{g/L}$ fluoride, nitrite, and bromide, and 5 $\mu\text{g/L}$ chlorite, bromate, and chlorate. Under the given conditions, sample matrices can be analyzed with up to 500 mg/L chloride, 100 mg/L nitrate, 500 mg/L sulfate, and 300 mg/L hydrogen carbonate.

6 Troubleshooting

6.1 Regeneration



CAUTION

Do not regenerate the column as a preventive measure!

Each regeneration process causes stress on the separation column and reduces its service life *see "Regenerating separation columns", page 5.*

Problem

- The backpressure increases.
- Double peaks occur.
- Tailing effects occur.
- The retention times become shorter.
- The resolution deteriorates.

Correction

Regenerating the separation column

Start by replacing the guard column if the above problems occur. Regenerate the separation column as described below if this measure does not help.

1 Disconnecting the separation column from the IC system

Disconnect the separation column outlet from downstream functional units such as the suppressor or the detector.

Collect the flow of liquid in a beaker.

2 Regenerating the separation column



NOTE

Check whether the maximum pressure is not exceeded during regeneration. If the pressure becomes too high, then reduce the flow rate.

Depending on the type of contamination, regenerate the separation column as follows:

- Contamination with organic components (*see table 3, page 39*)

- Contamination with inorganic components (see table 4, page 39)

When using organic modifiers for the regeneration, pay attention to the maximum backpressure.

Table 3 Contamination with organic components

	Rinse with	Duration [h]	Flow rate [mL/min]	Flow direction
1	Ultrapure water	1	0.4	Direction against the flow
2	Acetonitrile-water mixture (50:50)	2	0.4	Direction against the flow
3	Ultrapure water	1	0.4	Direction against the flow
4	Eluent	0.5	0.6	Regular

Table 4 Contamination with inorganic components

	Rinse with	Duration [h]	Flow rate [mL/min]	Flow direction
1	56 mmol/L Na ₂ CO ₃ , 31 mmol/L NaHCO ₃	2	0.4	Regular
2	Eluent	0.5	0.6	Regular

6.2 Decreasing resolution and asymmetrical peaks

Problem

The resolution of peaks deteriorates or peak shapes are asymmetrical.

Causes and prevention

Causes

The separation column has been overloaded.

Prevention or correction

The separation column can be overloaded by factors such as a high salt content in the sample matrix.

- Dilute the sample.
- Inject less sample.



Causes	Prevention or correction
There are dead volumes in the IC system.	<ul style="list-style-type: none"> Check that all of the capillaries have an inner diameter of ≤ 0.25 mm (6.1831.010). If not, use capillaries with a smaller inner diameter. Check that all of the capillaries are correctly installed. The IC Maintenance multimedia guide shows the installation process step-by-step.

6.3 Unstable retention times

Problem

The retention times are unstable.

Causes and prevention

Causes	Prevention or correction
Carbonate in the eluent	<p>Carbon dioxide from the air affects the carbonate / hydrogen carbonate balance in the eluent. The eluent becomes weaker over time.</p> <ul style="list-style-type: none"> Always keep the eluent bottle and bottles with eluent concentrates well sealed. Always use a CO₂ adsorber.
Air bubbles in the eluent	<p>Air bubbles make the eluent flow rate unstable. Backpressure is one indicator of an unstable flow rate. The backpressure must remain stable within ± 0.1 MPa.</p> <ul style="list-style-type: none"> Purge the high-pressure pump. Use the eluent degasser.

6.4 Unknown peaks

Problem

The chromatogram contains wide, unknown peaks.

Causes and prevention

Causes	Prevention or correction
Analytes eluting late	<p>Somewhat wider, unknown peaks can be the result of sample components eluting late. They are the result of the previous injection.</p> <ul style="list-style-type: none"> Extend the chromatogram duration.



6.5 Increasing backpressure

Problem

The backpressure increases.

Causes and prevention

Causes	Prevention or correction
Particles on the guard column	<ul style="list-style-type: none"> ▪ Replace the guard column.
Particles on the separation column	<p>Rinse the separation column at a reduced flow rate in the direction opposite to the flow direction.</p> <ul style="list-style-type: none"> ▪ Hold the column outlet in a beaker. ▪ Rinse the separation column for approximately 1 h. ▪ Install the separation column back in the flow direction.
Particles in the sample	<ul style="list-style-type: none"> ▪ Sample preparation, e.g. removing particles through Inline Ultrafiltration.

7 Literature

Metrohm recommends the following literature for more detailed information:

- Application Note S-236: Drinking water quality by EPA 300.1
- Column catalog, 8.000.5347

Index

B

Baseline
 Condition 14

C

Capacity 2
 Column
 see "Separation column" 10
 Conditioning 14

E

Eluent 7
 Equilibration 13

F

Flow rate 2

G

Guard column
 Installation 8
 Rinse 8, 9

I

IC column
 see "Separation column" 10
 Installation
 Guard column 8
 Separation column 10

O

Order number 1

R

Rinse
 Guard column 8, 9
 Separation column 12
 Rinsing
 Separation column 10

S

Separation column
 Installation 10
 Rinse 12
 Rinsing 10
 Specification 1