
6.2804.060 Deuteriumlampe 6.2804.060 Deuterium lamp

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6.2804.060 Deuterium lamp

WARNING

This product is classified as a special lamp in accordance with the regulations EC 244/2009, EU 1194/2012 and EU 2015/1428.

The lamp may only be used with the intended instruments. The lamp is not intended to be used to light rooms at home.

Intended use

The deuterium lamp 6.2804.060 is used in the following instruments:

- 2.844.0xxx UV/VIS Compact IC
- 2.887.0010 Professional UV/VIS Detector
- 2.944.0010 Professional UV/VIS Detector Vario

EU conformity

As this product is an accessory part, the EU declarations of conformity of the corresponding instruments are applicable. These instruments are listed under "Intended use".

WARNING

Harmful UV radiation

Damage to the eyes and skin by UV radiation if the lamp is removed while it is switched on.

- Switch off the instrument before removing the lamp.

CAUTION

Damage to the lamp

Reduced light transmission due to residues on the lamp surface. Damage to the lamp due to burnt-on residues on the lamp surface.

- Do not touch the lamp surface with your fingers!
- If the lamp surface is contaminated, clean the lamp with alcohol and a lint-free cloth.
- Make sure the lamp surface is dry before using the instrument.

Removing the old deuterium lamp

After a prolonged burning time, the radiation of the deuterium lamp starts to decrease. This can be noticed by an increased noise on the baseline. The service life of the lamp is approximately 1,000 operating hours.

Replace the deuterium lamp if the noise on the baseline becomes too strong and the burning time of the lamp is approximately 1,000 hours. Replace the deuterium lamp if the lamp does not light up anymore.

Prerequisites

- The instrument is switched off.
- The lamp has cooled down.

1. Disconnecting the deuterium lamp

Loosen the retaining ring (4) on the plug of the deuterium lamp. Remove the plug from the connection socket *UV*.

2. Loosening the fastening ring

Loosen and remove the fastening ring (6) of the deuterium lamp.

3. Removing the deuterium lamp

Carefully remove the deuterium lamp (5) from the housing.

Inserting the new deuterium lamp

Prerequisites

- The instrument is switched off.

1. Inserting the deuterium lamp

Insert the new deuterium lamp (6.2804.060) into the opening for the deuterium lamp on the lamp cooling element (1).

Align the groove on the lamp socket with the positioning pin in the lamp cooling element.

2. Fastening the deuterium lamp

Slide the fastening ring (6) over the cable of the deuterium lamp and screw it to the lamp cooling element.

3. Connecting the deuterium lamp

Plug the plug of the deuterium lamp into the connection socket *UV* of the instrument. Tighten the retaining ring (4).

Disposal

WEEE: This product is covered by European Directive 2012/19/EU, WEEE – Waste Electrical and Electronic Equipment.

The correct disposal of your old instrument will help to prevent negative effects on the environment and public health. More details about the disposal of your old instrument can be obtained from your local authorities, from waste disposal companies or from your local dealer.

RoHS: This product is covered by the directive 2011/65/EU, RoHS – Restriction of the use of certain hazardous substances in electrical and electronic equipment as well as its amendments 2015/863/EU.

Technical specifications

- Power consumption: 30 W
- Wavelength: 190 - 370 nm
- Service life: approx. 1,000 h
- Decrease in intensity after 1,000 h: ≤ 50%