## MSM-HC Rotor A, Hydroxid (6.2842.100)

<b>IC-Anionensuppressor MSM Hydroxid</b> .....	<b>1</b>
<b>IC Anion Suppressor MSM Hydroxide</b> .....	<b>5</b>
<b>Supresor de aniones Cl, MSM Hydroxid</b> .....	<b>9</b>
<b>Suppresseur d'anions Cl MSM Hydroxyde</b> .....	<b>12</b>

### DE

#### IC-Anionensuppressor MSM Hydroxid

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**EN**

## IC Anion Suppressor MSM Hydroxide

This leaflet contains information on cleaning and regenerating the Metrohm Suppressor Module "MSM" as well as its technical specifications and information on replacing the following spare parts:

- MSM-HC Rotor A, Hydroxide (6.2842.100)
- Connecting piece (6.2835.010)

The MSM-HC Rotor A, Hydroxide is suitable for applications with hydroxide eluents. The MSM-HC Rotor A, Hydroxide is used in systems with continuous eluent production.

### Standard conditions for MSM-HC Rotor A, Hydroxide

- With Metrosep A Supp 21: Gradients with 18–80 mmol/L KOH
- With Metrosep A Supp 19: Isocratic operation with 38 mmol/L KOH
- With Metrosep A Supp 18: Isocratic operation with 23 mmol/L KOH

### Technical specifications

Rotor	MSM-HC Rotor A, Hydroxide
Design	Micro packed-bed suppressor, robust, chemically resistant
Capacity	2.3 meq (2,300 µmol)*
Dead volume	< 120 µL
Regeneration	Chemical regeneration with STREAM or external rinsing with ultrapure water
Maximum flow rate	No irreversible damage possible
Backpressure limitations	No irreversible damage possible, no specific backpressure required
Solvent stability	100% solvent stability For applications with continuous eluent production, max. 20% methanol
Temperature range	No limit
Application	High-capacity MSM rotor for hydroxide gradient applications or isocratic hydroxide applications with very strong eluents
Storage	15 mmol/L KOH

\*Standard conditions: Metrosep A Supp 19 – 150/4.0, eluent: 38 mmol/L KOH, flow rate: 0.7 mL/min

## Regenerating

### Standard regeneration solution

- 0.2 mol/L H<sub>2</sub>SO<sub>4</sub>

If the suppressor units are loaded with certain heavy metals (such as iron) or organic contamination for long periods, then the standard regeneration solution may no longer be able to completely remove them. This constantly reduces the capacity of the suppressor units, which results in reduced phosphate sensitivity in mild cases and a large increase in the baseline in severe cases.

If these kinds of capacity problems occur at one or more positions, then all suppressor units must be regenerated at short notice with one of the following solutions:

### Regeneration solutions for restoration of suppressor capacity

- Contamination with heavy metals or elevated backpressure:**  
1 mol/L H<sub>2</sub>SO<sub>4</sub> + 0.1 mol/L oxalic acid
- Contamination with organic cationic complexing agents:**  
0.2 mol/L H<sub>2</sub>SO<sub>4</sub> / 0.1 mol/L oxalic acid / acetone 5%
- Heavy contamination with organic substances:**  
0.2 mol/L H<sub>2</sub>SO<sub>4</sub> / acetone ≥ 20%
- Contamination by certain environmental samples**  
1 mol/L H<sub>3</sub>PO<sub>4</sub>

#### NOTE

If phosphoric acid has been used once as regeneration solution in an IC system, then you will need to continue using phosphoric acid for regeneration. Regenerating it with sulfuric acid again can lead to problems in the baseline.

The MSM can be regenerated either with the peristaltic pump or with a Dosino. To restore the suppressor capacity, use the same technique as in continuous operation.

#### CAUTION

The pump tubing has only limited resistance to organic solvents. Insert a new pump tubing after regeneration.

### Procedure: Regeneration to restore suppressor capacity with peristaltic pump


- Replace the regeneration solution for continuous operation with the regeneration solution for restoring the suppressor capacity.
- Regenerate all 3 suppressor units, one after the other, for approx. 30 minutes using 1 of the solutions mentioned above.
  - In the software, set the peristaltic pump to Level 5.
  - In the software, use the **Step** command to switch to the next suppressor unit.
  - As soon as all three suppressor units have been regenerated, switch off the peristaltic pump.
- After regeneration, each of the three suppressor units must be rinsed with degassed ultrapure water for about 30 minutes.
  - Replace the regeneration solution with ultrapure water.
  - In the software, set the peristaltic pump to Level 5.
  - In the software, use the **Step** command to switch to the next suppressor unit.
  - As soon as all 3 suppressor units have been rinsed, switch off the peristaltic pump.
- Insert a new pump tubing.
- Replace the ultrapure water with regeneration solution for continuous operation.

### Procedure: Dosino regeneration to restore suppressor capacity

- Replace the regeneration solution for continuous operation with the regeneration solution for restoring the suppressor capacity.

2. Regenerate all 3 suppressor units, one after the other, for approx. 30 minutes using 1 of the solutions mentioned above.
  - a. Dose with a flow rate of 0.5 mL/min in the software.
  - b. In the software, use the **Step** command to switch to the next suppressor unit.
3. After regeneration, each of the 3 suppressor units must be rinsed with degassed ultrapure water for about 30 minutes.
  - a. Replace the regeneration solution with ultrapure water.
  - b. Dose with a flow rate of 0.5 mL/min in the software.
  - c. In the software, use the **Step** command to switch to the next suppressor unit.
4. Replace the ultrapure water with regeneration solution for continuous operation.

### Cleaning

*Suppressor – Parts*

- 1 Coupling nut**
- 2 Connecting piece (6.2835.010)**
- 3 MSM-HC Rotor A, Hydroxide**
- 4 Transparent suppressor housing**
- 5 Inspection window**

In the following cases, it may be necessary to clean the MSM:

- Increased backpressure at the MSM's connection tubing.
- Irremediable blockage of the MSM (solutions can no longer be pumped through the MSM).
- Irremediable jamming of the MSM (MSM can no longer be switched over).

### Procedure

1. Disconnect the MSM from the IC system.
  - a. In the software, stop the hardware and wait until the pressure in the high-pressure pump has been released.
  - b. Disconnect all capillaries of the MSM from the IC system.
2. Remove the MSM.
  - a. Unscrew the coupling nut **(1)** from the housing **(4)**.
  - b. Pull the connecting piece **(2)** out of the housing together with the rotor **(3)**.  
**NOTICE:** If the rotor gets stuck in the housing, you can push it out as follows:  
 Insert a pointed object into the inspection window **(5)** in the housing **(4)** and use it to push out the rotor.
  - c. Loosen the connecting piece from the rotor.
3. Rinse the capillaries.
  - a. Connect each of the 6 PTFE capillaries connected to the connecting piece **(2)** in sequential order to the outlet of the high-pressure pump using a coupling (6.2744.040) and pump ultrapure water through them.
  - b. Check whether water comes out at the connecting piece.

If one of the capillaries remains blocked, then the connecting piece **(2)** needs to be replaced (order number 6.2835.010 or 6.2832.010).
4. Clean the rotor.
  - a. Clean the sealing surface of the rotor **(3)** with ethanol using a lint-free cloth.

5. Reinsert the rotor into the housing.  
**CAUTION:** The rotor may be destroyed during start-up if it is not inserted correctly.
  - a. Insert the rotor **(3)** into the housing **(4)** in such a way that the tubing connections on the rear of the rotor fit into the corresponding recesses inside the housing and 1 of the 3 holes of the rotor is visible from below in the inspection window **(5)** of the housing **(4)**.  
**NOTICE:** The rotor's sealing surface is located approx. 4 mm deep inside the housing when the rotor is inserted correctly. If this is not the case, then the rotor needs to be moved into the correct position from below by means of a pointed object.
6. Clean the connecting piece.
  - a. Clean the sealing surface of the connecting piece **(2)** with ethanol using a lint-free cloth.
7. Insert the connecting piece.
  - a. Insert the connecting piece **(2)** into the housing **(4)** in such a way that connector 1 is on top and the 3 pins of the connecting piece fit into the corresponding recesses on the housing.
  - b. Reattach the coupling nut **(1)** and tighten it by hand (do not use a tool).
8. Connect and condition the MSM.
  - a. Reconnect all capillaries of the MSM to the IC system.
  - b. Before stepping the MSM for the first time, rinse the 3 suppressor units with solution for 5 minutes.

### Replacing parts

Parts of the MSM may need to be replaced in the following cases:

- Irremediable loss of suppressor capacity (reduced phosphate sensitivity and/or significant rise in the baseline).
- Irremediable blockage of the MSM (solutions can no longer be pumped through the MSM).

Both the rotor and the connecting piece can be replaced.

### Procedure

1. Disconnect the MSM from the IC system.
  - a. In the software, stop the hardware and wait until the pressure in the high-pressure pump has been released.
  - b. Disconnect all capillaries of the MSM from the IC system.
2. Remove the MSM.
  - a. Unscrew the coupling nut **(1)** from the housing **(4)**.
  - b. Pull the connecting piece **(2)** out of the housing together with the rotor **(3)**.  
**NOTICE:** If the rotor gets stuck in the housing, you can push it out as follows:  
Insert a pointed object into the inspection window **(5)** in the housing **(4)** and use it to push out the rotor.
  - c. Loosen the connecting piece from the rotor.
3. Clean the new rotor.
  - a. Clean the sealing surface of the rotor **(3)** with ethanol using a lint-free cloth.
4. Insert the rotor into the housing.  
**CAUTION:** The rotor may be destroyed during start-up if it is not inserted correctly.
  - a. Insert the rotor **(3)** into the housing **(4)** in such a way that the tubing connections on the rear of the rotor fit into the corresponding recesses inside the housing and 1 of the 3 holes of the rotor is visible from below in the inspection window **(5)** of the housing **(4)**.  
**NOTICE:** The rotor's sealing surface is located approx. 4 mm deep inside the housing when the rotor is inserted correctly. If this is not the case, then the rotor needs to be moved into the correct position from below by means of a pointed object.
5. Clean the new connecting piece.
  - a. Clean the sealing surface of the connecting piece **(2)** with ethanol using a lint-free cloth.
6. Insert the new connecting piece.
  - a. Insert the connecting piece **(2)** into the housing **(4)** in such a way that connector 1 is on top and the 3 pins of the connecting piece fit into the corresponding recesses on the housing.
  - b. Reattach the coupling nut **(1)** and tighten it by hand (do not use a tool).
7. Connect and condition the MSM.
  - a. Reconnect all capillaries of the MSM to the IC system.
  - b. Before stepping the MSM for the first time, rinse the 3 suppressor units with solution for 5 minutes.

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