## Thermal paper/Thermopapier 60 mm / 25 m

**EN**

### Properties

Table 1 Paper properties

Stability properties	Unit	Value	Test method
<b>Heat resistance:</b>			
Image loss (optical density)	OD	Max. 10%	24 h / 60 °C
<b>Humidity resistance:</b>			
Image loss (optical density)	OD	Max. 10%	24 h / 40 °C / 90% RH
<b>Light resistance:</b>			
Image loss (optical density, Xenotest)	OD	Max. 10%	5,000 lux / 100 h
<b>Stability of image &amp; printability</b>			
Correct storage for thermal paper		25 years	

**EN**

### General conditions

Thermoreactive coated paper can be used for many years. The printing also remains legible for several years if the paper is stored properly and some particularities in the handling of this paper are considered. These particularities are briefly described below.

#### Contrast

The contrast of the printout can be decreased if the paper has become yellow. All paper yellows more or less strongly if it is exposed to ultraviolet radiation. Thermal paper must be stored in a dark place protected from UV radiation, just like other archivable paper.

#### Storage temperature

The thermoreactive coating must be protected from strong heat. In case of storage in extremely warm conditions, reactions may occur, which make the coating turn gray and finally black. This blackening can be prevented if the storage temperature is not above normal ambient temperature.

Ideal: 24 °C (±5 °C).

#### Relative humidity

As for other papers too, the humidity in the warehouse should not be too high. Ideally, the relative humidity should be 50%. In any case, a relative humidity over 75% must be avoided.

Ideal RH: 50% (±10%).

#### Solvents and softeners

Thermoreactive coated paper is very sensitive to the effects of organic solvents and softeners, compared to uncoated paper of the same type. These substances can attack the coloring agents and lead to fading of the printout as well as to color reactions at room temperatures. Thermal paper must be protected against any contact with plastics containing softeners, as for example plasticized PVC. Carbonless copy paper reacts similarly. Be careful when using ink containing solvents to print on thermal paper.

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