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## 6.1204.XXX Electrode tips

- The electrode has to be protected against mechanical load!
- The tips of the electrodes may not be exposed to concentrated mineral acids as they will otherwise be irreversibly damaged!
- Please avoid extensive use of the electrode tips in organic solvents.
- It is recommended to use an individual electrode kit, consisting of working electrode (WE), reference electrode (RE), and auxiliary electrode (AE) for each application.
- Do not use the electrode outside the recommended temperature range.

### Electrode preparation / regeneration

The electrode should be rinsed with ethanol and dist. H<sub>2</sub>O before use.

The RDE is a solid electrode with a stationary surface. With increasing use, this surface becomes contaminated with the products of the electrode redox processes. The surface of the 6.1204.XXX electrode tips must therefore be regenerated from time to time by mechanical cleaning with extremely fine aluminum oxide powder. Proceed as follows:

1. Affix polishing cloth (part of 6.2802.000 Polishing kit) to a flat surface (e.g. to the bottom of a Petri dish) then add a little Al<sub>2</sub>O<sub>3</sub> powder (part of 6.2802.000 Polishing kit) and dist. H<sub>2</sub>O to form a thick paste.
2. Polish electrode tip by hand in the paste for ca. 10 s with small, ∞(eight)-shaped movements.
3. Rinse electrode tip with dist. H<sub>2</sub>O.

With unsatisfactory baselines and bad reproducibility of the measurement results, the electrode must be repolished and rinsed again.