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6.1204.1X0 Electrode tips

- *Please remove the protection cap before using the electrodes.*
- *The electrode has to be protected against mechanical load!*
- *PEEK electrode tips may not be exposed to concentrated mineral acids as they will otherwise be irreversibly damaged!*
- *Please avoid extensive use of the electrode tips in organic solvents.*
- *It is recommended to use an individual electrode kit, consisting of working electrode (WE), reference electrode (RE), and auxiliary electrode (AE) for each application.*
- *The electrodes are polished extensively ex works. Do **not** repolish the electrodes mechanically. Repolishing the electrode declines the signal-to-noise ratio and damages the electrode irreversibly!*
- *Do not use the electrode outside the recommended temperature range.*

Electrode preparation / regeneration

The electrode should be rinsed with ethanol and dist. H₂O before use.

The RDE is a solid electrode with a stationary surface. With increasing use this surface becomes contaminated with the products of the electrode redox processes. To clean the electrode, proceed as follows:

1. Put the electrode tip with the Pt surface into a 0.5 molar NaOH solution (at least p.A.-quality) for ca. 10 minutes.
2. Rinse electrode tip with dist. H₂O.

For electroplating bath applications (CVS), the electrode can (alternatively or additionally to cleaning with NaOH) be cleaned in pure Virgin Make-up Solution (VMS) by applying cyclic voltage sweeps.