6.1204.180 Ultra-Trace-Graphit-Electrode

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6.1204.180: Together with 6.1204.2x0 drive shaft makes up a rotating disk electrode (RDE).

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Attention

Do only polish 6.1204.180 Ultra Trace Graphite Electrode tips with the Polishing kit for Ultra Trace Electrodes 6.2802.020! All other electrode tips 6.1204.1XX get destroyed mechanically by this treatment!

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6.1204.180 Electrode tips

CAUTION

- The electrode tip must not be exposed to concentrated mineral acid, because it will be de-destroyed irreversibly otherwise!
- Please avoid extensive use of the electrode tip in organic solvents.
- It is recommended to use an individual electrode kit consisting of working electrode (WE), reference electrode (RE) and auxiliary electrode (AE), as well as and an individual measuring vessel for each application.
- Always clean and dry the electrode tip before long time storage.
- Do not use the electrode outside the recommended temperature range.
- Do not treat electrode tips in ultrasonic baths!

Electrode preparation / regeneration

The electrode should be rinsed with ethanol and ultrapure H₂O before use.

The RDE is a solid electrode with a stationary surface. With increasing use, this surface becomes contaminated with the products of the electrode redox processes, which can affect reproducibility and baseline stability. The surface of the 6.1204.180 electrode tips must therefore be regenerated from time to time by mechanical cleaning with the aluminum oxide Polishing strap 6.2802.030. Proceed as follows:

1. Insert a Polishing strap 6.2802.030 (part of the Polishing kit for Ultra Trace Electrodes 6.2802.020) with the yellow side up into the lower slot of the Polishing cylinder (part of the Polishing kit for Ultra Trace Electrodes 6.2802.020).
2. Remove the reference electrode (RE) and the auxiliary electrode (AE) from your VA stand.
3. Put the working electrode (WE) into the measuring head of you VA stand and let the electrode rotate at 2000 rpm.
4. Put the Polishing cylinder with the Polishing strap bottom-up over the electrode. The working surface of the electrode should touch the Polishing strap.
5. Press on the Polishing cylinder gently for 5 s, without stopping the driving motor (s. Fig.2). Then, move back and forth the Polishing strap through the slot for 5 s (simultaneously, press the Polishing strap gently on the electrode). It is important, that not only material from the Electrode shaft is ablated, but also from the graphite zone. The rubbed-off parts leave distinct marks on the Polishing strap.
6. Remove the Polishing cylinder from the electrode tip and turn off the driving motor of the VA stand.
7. To remove dirt and coal particles, rinse the electrode tip with ultra pure water thoroughly and rub it with a lint-free paper towel or filter paper.
8. Reinstall the reference electrode (RE) and the auxiliary electrode (AE) into the VA stand. Make sure, that all electrodes are connected correctly to the corresponding electrode cable.

If baselines or reproducibility are still not satisfying, repeat this procedure.

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