6.0728.xx0 Ag/AgCl Referenzelektrode 6.0728.xx0 Ag/AgCl Reference Electrode

DE

6.0728.xx0 Ag/AgCl Referenzelektrode

Allgemeines

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6.0728.xx0 Ag/AgCl Reference Electrode

General

- Check electrodes for transport damage when unpacking.
- Rinse off crystallized electrolyte from diaphragm.
- For special applications, it is recommended to use an own set of electrodes, consisting of a Working Electrode (WE), a Reference Electrode (WE) and an Auxiliary Electrode (AE).
- During the electrolyte exchanging process, crystallization problems can occur. A solution that produces a stable potential should be used as internal electrolyte for all reference electrodes. For example:
 - Aqueous systems: $c(\text{KCl}) = 3 \text{ mol/L}$
 - Nonaqueous systems: $c(\text{LiCl}) = 2 \text{ mol/L}$ in EtOH $c(\text{TEABr}^*) = 0.4 \text{ mol/L}$ in ethylene glycol
* tetraethylammonium bromide

Voltages of the Ag/AgCl reference electrode vs. a standard hydrogen electrode (concentration data at 25 °C):

Electrolyte	20 °C U [mV]	25 °C U [mV]
c(KCl) = sat.	201.9	197.0
c(KCl) = 3 mol/L	210.0	207.0
c(LiCl) = sat. in EtOH	148	143
c(LiCl) = 2 mol/L in EtOH	164	157
c(TEABr) = 0.4 mol/L in ethylene glycol	217	213

When using double junction electrodes (in combination with the electrolyte vessel 6.1245.010) the bridge electrolyte must be compatible with the internal electrolyte.

Setting up for operation

Fill in internal electrolyte

The filling of the internal system is only necessary for the (dryly delivered) electrodes 6.0728.000 and 6.0728.010. Refill all other electrode as required up to the deaerating holes:

1. Hold the internal system horizontally.
2. Screw off the diaphragm mount.
3. Fill the electrolyte area completely with the desired internal electrolyte.
4. Screw on the diaphragm mount.
5. Now, hold the internal system in a way that the diaphragm points downwards, and check whether the internal system is filled completely. Remove air bubbles by tapping slightly on the electrode shaft.

Do not use the electrode until 24 h after the filling process!

Fill in bridge electrolyte

1. Fill the electrolyte vessel 6.1245.010 with a suitable bridge electrolyte, normally $c(\text{KCl}) = 3 \text{ mol/L}$. The composition of the bridge electrolyte should match the analytical requirements (aqueous / nonaqueous solutions).

Metrohm recommends the following bridge electrolytes:

- For solutions containing perchlorate: $c(\text{NaCl}) = 3 \text{ mol/L}$
- For applications in the electroplating bath analysis: $c(\text{KNO}_3) = 1 \text{ mol/L}$

Screw together reference electrode

1. Put the internal system into the filled bridge electrolyte vessel (6.1245.010) and screw it tight.

A part of the bridge electrolyte is pressed out of the deaerating hole. Metrohm recommends to wait 20 minutes until the diaphragm is soaked with bridge electrolyte.

Measurement and storage

- Before each measurement, make sure that the internal electrolyte and bridge electrolyte are refilled up to their deaerating holes.
- If the reference electrode is not used for a longer period of time, store the internal system in the corresponding internal electrolyte solution.
If internal and bridge electrolyte are identical, the complete reference electrode can be stored in the corresponding electrolyte solution.

Maintenance and care

6.0728.000, 6.0728.010, 6.0728.x20

For trace analysis, periodically (weekly) exchanging the bridge electrolyte should be sufficient (use the same electrolyte type!).

In case your internal system dried out, rinse it thoroughly with internal electrolyte (Procedure *see "Fill in internal electrolyte"*).

6.0728.x30

The «Long Life» reference system must be filled necessarily with $c(\text{KCl}) = 3 \text{ mol/L}$.

For electroplating bath analysis, Metrohm recommends to renew the internal electrolyte weekly and the bridge electrolyte daily (Procedure *see "Setting up for operation"*; use the same electrolyte type!).

In case your Long Life Electrode dried out, fill the internal system as follows:

1. Screw off the diaphragm mount.
2. Put the internal system into a beaker filled with $c(\text{KCl}) = 3 \text{ mol/L}$ (deaerating holes must be covered with KCl) and evacuate it in an exsiccator ($< 50 \text{ mbar}$).
3. Afterwards, let air in slowly, so that the internal system is filled completely with KCl.
4. Screw on the diaphragm mount.
5. Now, hold the internal system in a way that the diaphragm points downwards, and check whether the internal system is filled completely. Air bubbles should be removed by tapping slightly on the electrode shaft

Do not use the electrode until 24 h after the filling process!