## Electrode tips for rotating disk electrode (RDE) Elektrodentips zu rotierender Scheibenelektrode (RDE)

6.1204.120 / 6.1204.130 / 6.1204.600

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### EN

## Electrode tips for rotating disk electrode (RDE)

An electrode tip (6.1204.1x0 or 6.1204.600) together with a driving axle (6.1204.210 or 6.1204.220) constitute the rotating disk electrode for VA Stands.

### CAUTION

- Remove the protective cap before using the electrode.
- Protect the electrode against mechanical load!
- PEEK electrode tips may not be exposed to concentrated mineral acids as they will otherwise be irreversibly damaged!
- Avoid extensive use of the PEEK electrode tips in organic solvents.
- PEEK electrode tips may be used within a temperature range of 0 °C to 40 °C.
- Glass electrode tips may be used within a temperature range of 0 °C to 50 °C.
- It is recommended to use an individual electrode kit, consisting of working electrode (WE), reference electrode (RE), auxiliary electrode (AE), and measuring vessel for each application.
- Do not treat electrode tips with ultrasound.

### Preparation

Rinse the electrode with ethanol and distilled water before use.

### Regeneration

The RDE is a solid-state electrode with a stationary surface. During use, this surface becomes increasingly contaminated with the products of the electrode redox processes.

Therefore, the surface of the electrode tip must be regenerated every now and then by mechani-

cally cleaning it with finest aluminum oxide powder.

For regeneration, the polishing set (6.2802.000) is required.

Proceed as follows:

1. Stick the polishing cloth onto a flat surface (e.g. the bottom of a petri dish).
2. Mix a small amount of aluminum oxide powder with distilled water to form a thick pulp.
3. Polish the electrode tip by hand forming little eights with the electrode in the pulp for approx. 10 s.
4. Thoroughly rinse the electrode tip with distilled water.

If the baseline is not satisfactory or if the reproducibility of the measuring results is poor, then the electrode must be polished and rinsed once again.

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