Pd-Referenzelektrode Pd reference electrode

6.1257.740

DE

Pd-Referenzelektrode

Allgemein

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Pd reference electrode

General


Upon receipt, please check the electrodes for damages during shipment.

Electrodes are consumables!

Usage

The Pd reference electrode (6.1257.740) has been conceived as a long-lasting and low-maintenance reference electrode for the use in the measuring cell (6.1257.xxx). Its potential is pH-dependent. It can be applied especially well with alkaline eluents (pH 8 to 14).

Potential-determining half-reaction

$E \approx +0.14 \text{ V}$ in 1 M NaOH at 25 °C, 0.1 MPa, vs. SHE

At pH 12.8 and 0.1 mol/L sodium hydroxide, the potential corresponds approximately the potential of the Ag/AgCl/KCl(sat.) reference electrode. If the pH value of the eluate is known, a potential given for a Ag/AgCl/KCl(sat.) reference electrode E_{vsAgCl} can be calculated using the following formula:

$$E_{\text{vsPd}} \approx E_{\text{vsAgCl}} + 0.05 \text{ V} \times (\text{pH} - 12.8)$$

Example 1:

$E_{\text{vsAgCl}} = 0.1 \text{ V}$; pH = 10.3 (carbonate/bicarbonate-eluent):

$$E_{\text{vsPd}} \approx 0.1 \text{ V} + 0.05 \text{ V} \times (10.3 - 12.8) = -0.025 \text{ V}$$

These values provide a basis for optimizations.

Care

Cleaning: rinse the Pd reference electrode or softly wipe it clean with a moist paper towel.

CAUTION

- Do not polish or rub the reference electrode.
- Avoid contact of the Pd reference electrode with concentrated acids.

- Avoid contact of the Pd reference electrode with strong acidic reducing agents (e. g. citric acid), which can reduce the PdO coating.
- Do not expose the Pd reference electrode to any ultrasonic fields.
- Connect the Pd reference electrode as a reference electrode only, never as a cathode!

Performance

After an eluate exchange the potential takes up to 24 h to stabilize (deviation up to 0.07 V).

If you return to the original eluate after an eluate exchange, the original potential stability of 0.03 V is reached only after approx. 48 h, at constant pH. Check standards must be performed at least once a day.

Disposal

The Pd reference electrode can either be recycled with noble metals or disposed of as metal containing laboratory waste.