



## 1 Temperature sensor

### 1.1 General

Immediately after receiving the temperature sensor, check to make sure that it works properly. Temperature sensors that do not work properly must be sent back for warranty processing within two months (starting from the day of delivery). If the defect is proven to be due to a material or manufacturing defect, the temperature sensor will be replaced at no charge. The transport costs are to the customer's account.

### 1.2 Selecting the temperature sensor in accordance with the application

It is imperative that you select the appropriate temperature sensor for the application:

Temperature sensors made of glass are suitable for all types of liquid samples. For semisolid materials, such as cheese, fruits, etc., the temperature sensor made of stainless steel can be used. This sensor, however, is not resistant to corrosive and strongly acidic media.

### 1.3 Cleaning



#### CAUTION

Do not use the ultrasonic bath for temperature sensors, as they may be damaged by such a treatment.

The temperature sensor can be cleaned by briefly immersing it in a suitable solvent or wiping it with a dampened cloth.

## **1.4 Storage**

The temperature sensor must be stored dry.