6.5405.000

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1. Preparation of the cell

Dry all parts well!

1. Mount titration vessel (1) in the titration vessel holder and fix it onto rod by means of clamping ring.
2. Put stirring bar (10) into the titration vessel.
3. Cut ground joint sleeves into right length and provide all ground joints with them.
4. Put indicator electrode (5) into titration vessel.
5. Take the electrode cable with the smaller plug (plug F), screw it onto indicator electrode and connect it to the KF Coulometer on the socket for the indicator electrode.
6. Put generator electrode (2) into titration vessel.
7. Fill drying tube (9) with molecular sieve and put it into the generator electrode.
8. Screw the electrode cable with the bigger plug (plug H) onto generator electrode and connect it to the KF Coulometer on the socket for the generator electrode.


Note: Avoid mixing up the two electrode cables! In order not to confound them, mark one of the cables and the corresponding electrode socket accordingly.

9. Introduce a septum into the screw cap (8) and screw it onto the glass connection piece.
10. Fill the cell with reagent up to the lower margin of the article number (app. 80–90 ml). Use only reagents which are especially designated for applications with cells without diaphragm.
11. Close the last ground joint aperture with stopper (7) .

Alternatively, you may use the inlet/outlet tube (6) instead of stopper (7) for the connection of a **pump**, to add fresh solvent and aspirate used solution:

12. Put stopper (7) onto inlet/outlet tube (6) and place it in titration vessel.
13. Screw on outlet and inlet tubing.
14. If the pump is not connected during a certain period, the stoppers with thread M8 (6.1446.080) are used to close the tubing connection of the inlet/outlet tube.

For work with the **KF Oven** the inlet/outlet tube (6) is transformed as follows:

1. Remove the screw connector for tubing with thread M8 on the inlet tube (longer tube) and replace it by screw connector for tubing with thread M6 (6.1811.000).
2. Screw tubing from KF Oven onto inlet tube.
3. If the KF Oven is not connected for a certain period, the tube is closed with a stopper with thread M6 (6.1446.040) and M8 (6.1446.080) each.

2. Scope of delivery and ordering designations

Cell without diaphragm including the following accessories:	6.5405.000
1 Double platinum electrode	6.0341.100
1 Generator electrode	6.0342.100
1 Drying tube	6.1403.030
1 Glass stopper SGJ 14	6.1437.000
1 Inlet/outlet tube	6.1439.010
1 Stopper with thread M6 for work on the KF Oven	6.1446.040
2 Stoppers with thread M8 for work with a Pump	6.1446.080
2 Sets of septa (5 items each)	6.1448.020
1 Titration vessel, 200 ml	6.1455.313
1 Tubing connection, 50 cm for work with a Pump	6.1805.200
1 Screw connector for work with a KF Oven	6.1811.000
1 Stirring bar, 25 mm	6.1903.030
1 Titration vessel holder	6.2047.000
1 Electrode cable, plug F, for indicator electrode	6.2104.020
1 Electrode cable, plug H, for generator electrode	6.2104.120
1 Screw cap	6.2701.040
3 PTFE ground sleeves SGJ 14	6.2713.000
1 PTFE ground sleeve SGJ 29	6.2713.010
1 PTFE ground sleeve SGJ 19	6.2713.020