Ion exchanger cartridges

IC-RP 6.1012.000 (50) / 6.1012.100 (10)

Sample preparation cartridge with C18 (octyl) for non-polar solid phase extraction.

Application

- When organic substances which could interfere with the chromatogram or harm the column have to be removed.
- For the enrichment of heavy metals as a complex and subsequent determination (also extraction from an interfering matrix; e.g. heavy metals in brine).

Activation

Before IC-RP cartridges are used, they must be activated with 10 mL methanol p.a. and flushed with 10 mL deionised water.

Regeneration

Regeneration is performed in the same manner, except that the amount of methanol is increased to ca. 20...50 mL.

IC-H 6.1012.010 (50) / 6.1012.110 (10)

Cation exchanger cartridge in the H⁺ form for anion analysis.

Application

- When the front peak is too broad owing to cations (e.g. Ca²⁺, Mg²⁺) and hence anions which are eluted early are masked.
- For the removal of carbonate/bicarbonate (possibly also treated briefly in an ultrasonic bath or with nitrogen).
- With alkaline samples (NaOH + R-SO₃H⁺ → H₂O + R-SO₃Na⁺).
- For the removal of cations from samples for ion exclusion chromatography.

Preparation

Before usage rinse with 10 mL deionised water. If possible rinse afterwards with sample to avoid dilution of the sample caused by dead volume of the cartridge.

Regeneration

- HCl p.a. 10% (10mL) is forced slowly through the cartridge using a syringe or sucked through by vacuum. The filter is then flushed free from chloride with ultra pure water (control injection).
- If the regenerated exchanger cartridge is not used for quite some time (2 days), it should be flushed again with 5 mL ultra pure water before use.

Another acid (e.g. HNO₃ p.a. 15%) can be used instead of HCl p.a. 10%, depending on the ions which have to be determined.

IC-Ag 6.1012.020 (50) / 6.1012.120 (10)

Cation exchanger cartridge in the Ag⁺ form for anion analysis.

Application

- Removal of the halides chloride, bromide, iodide (the content of fluoride, nitrite, nitrate and sulphate is influenced only to a minor extent or not at all).

Preparation

Before usage rinse with 10 mL deionised water. If possible rinse afterwards with sample to avoid dilution of the sample caused by dead volume of the cartridge.

- To retain excess silver, the sample must be injected via the cation exchanger cartridge IC-H 6.1012.010.

Regeneration

Cannot be regenerated.

IC-OH 6.1012.030 (50) / 6.1012.130 (10)

Anion exchanger cartridge in the OH⁻ form for cation analysis.

Application

- With highly acidic samples (pH < 2), the pH value can be raised and then set to pH=3 with HNO₃ (HCl + R-NH₃⁺ + OH⁻ → H₂O + R-NH₂ + Cl⁻).

Preparation

Before usage rinse with 10 mL deionised water. If possible rinse afterwards with sample to avoid dilution of the sample caused by dead volume of the cartridge.

Regeneration

Anion exchangers are regenerated in the same way as cation exchangers. A base, chiefly NaOH p.a. or KOH p.a., is used in place of acid.

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