6.2245.010 Datenkarte (SRAM)

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
6.2245.010 Data Card (SRAM)

Technical Specifications

Type	SRAM IC Memory Card according to JEIDA ICMC 4.0 / PCMCIA (68 pin)
Storage capacity	128 kB
Battery type	CR2025 (3 V)
Battery life time	at least 5 years (with use at $\leq 25^{\circ}\text{C}$)

Precautions

- Do not touch the terminals of the data card.
- Do not bend, drop, or expose to impact the data card.
- Avoid heat, moisture and direct sunlight.

Caution


- Danger of explosion if battery is incorrectly replaced!
- Use only CR2025 batteries for replacement.
- Discard used batteries according to manufacturer's instructions.

Battery replacement

In order to retain stored data and avoid data loss with a used data card, the battery should be replaced while the card is inserted in the instrument, supplying voltage to the card. Proceed as follows:

1. Place the battery holder lock switch in "OFF" position.
2. Pull the battery holder out of the card.
3. Remove the old battery and put the new battery (type CR2025) on the battery holder with "+" face up.
4. Insert the battery holder firmly into the card.
5. Place the battery holder lock switch in "ON" position.
6. Enter new date for battery change in the instrument.

Write protect

The data card has a write protect (WP) switch to protect data in the card from accidental overwriting.

Position "ON":

Prohibit data-writing.

Position "OFF":

Enable data-writing.