
6.2832.000 / 6.2842.000 / 6.2844.000 / 6.2835.000

IC-Anionensuppressoren .....	1
IC anion suppressors .....	4
Supresores de aniones Cl .....	8
Suppresseurs d'anions Cl .....	11

**DE**

**IC-Anionensuppressoren**

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## EN

### IC anion suppressors

This leaflet contains information on cleaning and regenerating the Metrohm Suppressor Module "MSM" and the Sample Preparation Module "SPM", as well as their technical specifications and information on replacing the following spare parts:

- MSM Rotor A (6.2832.000) with adapter (6.2842.020)
- MSM-LC Rotor A (6.2844.000) with adapter (6.2842.020)
- MSM-HC Rotor A (6.2842.000)

- SPM Rotor A (6.2835.000)
- Connecting piece (6.2835.010 or 6.2832.010)

#### NOTICE

The term "rotor" is used in this leaflet as a generic term for all three MSM rotors and the SPM rotor.

#### Regenerating the suppressor/SPM

If the suppressor units are loaded with certain heavy metals (such as iron) or organic contamination for long periods, then the standard regeneration solution may no longer be able to completely remove them. This constantly reduces the capacity of the suppressor units, which results in reduced phosphate sensitivity in mild cases and a large increase in the baseline in severe cases.

If such capacity problems occur at one or more positions, all anion suppressor units must be regenerated with one of the following solutions:

#### Regeneration solutions

- **Contamination with heavy metals or elevated backpressure:**  
1 mol/L H<sub>2</sub>SO<sub>4</sub> + 0.1 mol/L oxalic acid
- **Contamination with organic cationic complexing agents:**  
0.1 mol/L H<sub>2</sub>SO<sub>4</sub> / 0.1 mol/L oxalic acid / acetone 5%
- **Heavy contamination with organic substances:**  
0.2 mol/L H<sub>2</sub>SO<sub>4</sub> / acetone ≥ 20%
- **Contamination by certain environmental samples**  
1 mol/L H<sub>3</sub>PO<sub>4</sub>

#### NOTICE

If phosphoric acid has been used once as regeneration solution in an IC system, you will have to continue using phosphoric acid for regeneration. Regenerating it with sulfuric acid again can lead to problems in the baseline.

#### CAUTION

The pump tubing made of PVC has only limited resistance to organic solvents.

If a regeneration solution **with organic solvents** is used, then a separate pump tubing is used exclusively for regeneration. As an alternative, use the high-pressure pump for regeneration.

If a regeneration solution **without organic solvents** is used, then use the pump tubing for standard operation for regeneration as well.

#### Procedure

1. If a regeneration solution with organic solvents is used, then insert the dedicated pump tubing for regeneration into the peristaltic pump.
2. Replace the regeneration solution for continuous operation with the regeneration solution for suppressor regeneration.
3. Regenerate all three suppressor units, one after the other, for approx. 30 minutes using one of the solutions mentioned above.
  - a. In the software, set the peristaltic pump to Level 5.
  - b. In the software, use the **Step** command to switch to the next suppressor unit.
  - c. As soon as all three suppressor units have been regenerated, switch off the peristaltic pump.

4. After regeneration, each of the three suppressor units must be rinsed with degassed ultrapure water for about 30 minutes.
  - a. Replace the regeneration solution with ultrapure water.
  - b. In the software, set the peristaltic pump to Level 5.
  - c. In the software, use the **Step** command to switch to the next suppressor unit.
  - d. As soon as all three suppressor units have been rinsed, switch off the peristaltic pump.
5. If a dedicated pump tubing was used for the regeneration, then use the pump tubing again for standard operation.
6. Replace the ultrapure water with regeneration solution for continuous operation.

### Cleaning the suppressor/SPM

In the following cases, it may be necessary to clean the suppressor:

- Increased backpressure at the suppressor's connection tubing.
- Irremediable blockage of the suppressor (solutions can no longer be pumped through the suppressor).
- Irremediable jamming of the suppressor (suppressor can no longer be switched over).

### Procedure

1. Disconnect the suppressor from the IC system.
  - a. In the software, stop the hardware and wait until the pressure in the high-pressure pump has been released.
  - b. Disconnect all capillaries of the suppressor from the IC system.
2. Dismantle the suppressor.
  - a. Unscrew the union nut **(1)** from the housing **(4)**.
  - b. Pull the connecting piece **(2)** out of the housing together with the rotor **(3)**.  
**NOTICE:** If the rotor gets stuck in the housing, you can push it out as follows:  
 Put a pointed object into the slot **(5)** in the housing **(4)** and use it to push out the rotor.
  - c. Loosen the connecting piece from the rotor.
3. Rinse the capillaries.
  - a. Connect each of the 6 PTFE capillaries connected to the connecting piece **(2)** in sequential order to the outlet of the high-pressure pump using a coupling (6.2744.040) and pump ultrapure water through them.
  - b. Check whether water comes out at the connecting piece.

If one of the capillaries remains blocked, the connecting piece **(2)** has to be replaced (order number 6.2835.010 or 6.2832.010).
4. Clean the rotor.
  - a. Clean the sealing surface of the rotor **(3)** with ethanol using a lint-free cloth.
5. Reinsert the rotor into the housing.
 

**CAUTION:** The rotor may be destroyed during start-up if it is not inserted correctly.

  - a. Insert the rotor **(3)** into the housing **(4)** in such a way that the tubing connections on the rear of the rotor fit into the corresponding recesses inside the housing and one of the three holes of the rotor is visible from below in the slot of the housing **(5)**.  
**NOTICE:** The rotor's sealing surface is located approx. 4 mm deep inside the housing if the rotor is inserted correctly. If this is not the case, then the rotor has to be moved into the correct position from below by means of a pointed object.
6. Clean the connecting piece.
  - a. Clean the sealing surface of the connecting piece **(2)** with ethanol using a lint-free cloth.
7. Insert the connecting piece.
  - a. Insert the connecting piece **(2)** into the housing **(4)** in such a way that connector 1 is on top and the 3 pins of the connecting piece fit into the corresponding recesses on the housing.
  - b. Reattach the union nut **(1)** and tighten by hand (do not use a tool).
8. Connect and condition the suppressor.
  - a. Reconnect all capillaries of the suppressor to the IC system.
  - b. Before switching the suppressor over for the first time, rinse the 3 suppressor units with solution for 5 minutes.

## Replacing parts of the suppressor/SPM

Parts of the suppressor may have to be replaced in the following cases:

- Irremediable loss of suppressor capacity (reduced phosphate sensitivity and/or significant rise in the baseline).
- Irremediable blockage of the suppressor (solutions can no longer be pumped through the suppressor).

Both the rotor and the connecting piece can be replaced.

### Procedure

1. Disconnect the suppressor from the IC system.
  - a. In the software, stop the hardware and wait until the pressure in the high-pressure pump has been released.
  - b. Disconnect all capillaries of the suppressor from the IC system.
2. Dismantle the suppressor.
  - a. Unscrew the union nut **(1)** from the housing **(4)**.
  - b. Pull the connecting piece **(2)** out of the housing together with the rotor **(3)**.  
**NOTICE:** If the rotor gets stuck in the housing, you can push it out as follows:  
Put a pointed object into the slot **(5)** in the housing **(4)** and use it to push out the rotor.
  - c. Loosen the connecting piece from the rotor.
3. Clean the new rotor.
  - a. Clean the sealing surface of the rotor **(3)** with ethanol using a lint-free cloth.
4. Insert the rotor into the housing.  
**CAUTION:** The rotor may be destroyed during start-up if it is not inserted correctly.
  - a. Insert the rotor **(3)** into the housing **(4)** in such a way that the tubing connections on the rear of the rotor fit into the corresponding recesses inside the housing and one of the three holes of the rotor is visible from below in the slot of the housing **(5)**.  
**NOTICE:** The rotor's sealing surface is located approx. 4 mm deep inside the housing if the rotor is inserted correctly. If this is not the case, then the rotor has to be moved into the correct position from below by means of a pointed object.
5. Clean the new connecting piece.
  - a. Clean the sealing surface of the connecting piece **(2)** with ethanol using a lint-free cloth.
6. Insert the new connecting piece.
  - a. Insert the connecting piece **(2)** into the housing **(4)** in such a way that connector 1 is on top and the 3 pins of the connecting piece fit into the corresponding recesses on the housing.
  - b. Reattach the union nut **(1)** and tighten by hand (do not use a tool).
7. Connect and condition the suppressor.
  - a. Reconnect all capillaries of the suppressor to the IC system.
  - b. Before switching the suppressor over for the first time, rinse the 3 suppressor units with solution for 5 minutes.

### Technical specifications

Rotor	MSM	MSM-HC	MSM-LC
Structure	Micro packed-bed suppressor, robust, chemically resistant		
Capacity	approx. 90 min* (0.45 meq)	approx. 350 min* (1.8 meq)	approx. 45 min** (0.13 meq)
Dead volume	< 50 µL	< 250 µL	< 10 µL
Regeneration	Chemical regeneration with STREAM		
Maximum flow rate	No irreversible damage possible		
Backpressure limitations	No irreversible damage possible, no specific backpressure required		
Solvent stability	100% solvent stability		
Temperature range	No limit		

Rotor	MSM	MSM-HC	MSM-LC
Application	Universal use	High-capacity suppressor for gradient applications or isocratic applications with very strong eluents	Low-capacity suppressor with minimum dead volume for applications with 2 mm columns

\* Standard conditions: Metrosep A Supp 5 - 100/4.0, eluent: 3.2 mmol/L Na<sub>2</sub>CO<sub>3</sub> / 1.0 mmol/L NaHCO<sub>3</sub>, flow rate: 0.7 mL/min

\*\* Standard conditions: Metrosep A Supp 10 - 100/2.0, eluent: 5.0 mmol/L Na<sub>2</sub>CO<sub>3</sub> / 5.0 mmol/L NaHCO<sub>3</sub>, flow rate: 0.2 mL/min

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