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
| | |
|--|-----------|
| IC-Kationensuppressor | 2 |
| Aufbau | 2 |
| Funktionsweise und Anwendung | 2 |
| Regenerieren | 3 |
| Reinigen | 3 |
| Teile ersetzen | 4 |
| Technische Daten | 5 |
| | |
| IC cation suppressor | 5 |
| Structure | 5 |
| Mode of operation and application | 6 |
| Regenerating | 6 |
| Cleaning | 7 |
| Replacing parts | 8 |
| Technical specifications | 8 |
| | |
| Suppresseur de cations pour CI | 9 |
| Structure | 9 |
| Fonctionnement et application | 9 |
| Régénération | 10 |
| Nettoyage | 11 |
| Remplacement des pièces | 11 |
| Caractéristiques techniques | 12 |
| | |
| Supresor de cationes CI | 13 |
| Estructura | 13 |
| Funcionamiento y aplicación | 13 |
| Regeneración | 14 |
| Limpieza | 14 |
| Sustitución de piezas | 15 |

DE

IC-Kationensuppressor

Aufbau

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
4 Transparent suppressor housing

5 Inspection window

This leaflet contains information on the cleaning and regeneration of the cation suppressor and on the replacement of the following spare parts:

- Rotor MSM-HC C (6.2842.200)
- Connecting piece (6.2835.010)

Mode of operation and application

The rotor (MSM-HC C) is comprised of three suppressor units filled with anion exchanger material.

The following three processes can run simultaneously as a result:

- Suppression takes place in the first suppressor unit.
- The second suppressor unit is regenerated with an aqueous solution of 70 mmol/L sodium carbonate (Na_2CO_3) and 70 mmol/L sodium bicarbonate (NaHCO_3).
- The third suppressor unit is rinsed with ultrapure water or suppressed eluate.

Before each analysis, the suppressor is rotated 120° so that a freshly regenerated and rinsed suppressor unit is always available for suppression.

The regeneration solution can be conveyed two different ways:

- by peristaltic pump
Speed rate 1, 6.1826.420 pump tubing (PharMed® pump tubing (orange/yellow), 3-stopper)
- by Dosino Regeneration
2 mL dosing unit, dosing rate 0.5 mL/min, volume 5 mL

The conductivity value in a suppressed cation system is typically less than 0.3 $\mu\text{S}/\text{cm}$.

CAUTION

Use of unsuitable separation columns

If an unsuitable separation column is used, irreversible damage could occur to the ion chromatograph and the detector.

Use **polymer-based cation columns exclusively** in combination with the IC cation suppressor.

CAUTION

Conversion of an IC system to a suppressed cation system

Precipitations could occur after the conversion if the IC system was previously used

- with silica gel-based separation columns.
- for measuring anions.

Therefore use a **conductivity detector** for the suppressed cation system that is used **exclusively for this application**.

Regenerating

If the suppressor units are exposed to certain contaminations for prolonged periods, then these contaminations can no longer be removed completely with the standard regeneration solution. The performance of the suppressor units declines continuously, as indicated by a rise of the baseline or asymmetrical peaks.

If such problems occur at one or more positions, then treat all of the suppressor units as follows:

Option A

1. Stop the regeneration solution feed.
2. Rinse the system with eluent until the suppressor unit is completely exhausted (indicated by a significant increase in conductivity). This can take up to 240 minutes under standard conditions.
3. In the **Manual control** program part of the software, use the **Step** command to switch to the next suppressor unit. Repeat step 2.

4. Repeat step 3 for the third suppressor unit.
5. As soon as all three suppressor units are completely exhausted, restore the regeneration solution feed.
6. Equilibrate the system as usual (*see chapter "Conditioning" in the manual for the ion chromatograph*).

Option B

1. Rinse the three suppressor units with approx. 150 mL 2 mol/L nitric acid at a maximum flow rate of 5 mL/min (2 mL dosing unit or peristaltic pump).
At the same time, use the **Step** command in the **Manual control** program part in the software to switch to the next suppressor unit every five minutes.
2. Rinse the three suppressor units with approx. 150 mL standard regeneration solution.
At the same time, use the **Step** command in the **Manual control** program part of the software to switch to the next suppressor unit every five minutes.
3. Equilibrate the system as usual (*see chapter "Conditioning" in the manual for the ion chromatograph*).

Cleaning

The IC cation suppressor must be cleaned under the following circumstances:

- The backpressure on the connection capillaries rises.
- The suppressor is clogged (solutions can no longer be conveyed through the suppressor).
- The suppressor is jammed (suppressor can no longer be switched over).

Procedure

1. Disconnect the suppressor from the IC system.
 - a. In the software, stop the hardware and wait until the pressure in the high-pressure pump has been released.
 - b. Disconnect all suppressor capillaries from the IC system.
2. Dismantle the suppressor.
 - a. Unscrew the union nut (2-**1**) from the housing (2-**4**).
 - b. Pull the connecting piece (2-**2**) out of the housing together with the rotor (2-**3**).
NOTE: If the rotor gets stuck in the housing, you can push it out as follows:
Insert a pointed object into the inspection window (2-**5**) at the bottom of the housing. Use it to push out the rotor.
 - c. Loosen the connecting piece from the rotor.
3. Rinse the capillaries.
 - a. Connect each of the six PTFE capillaries connected to the connecting piece in sequential order to the outlet of the high-pressure pump using a coupling (6.2744.040). Pump ultrapure water through them.
 - b. Check whether water comes out at the connecting piece.
If one of the capillaries remains blocked, the connecting piece must be replaced (order number 6.2835.010).
4. Clean the rotor.
 - a. Clean the sealing surface of the rotor with ethanol using a lint-free cloth.
5. Reinsert the rotor into the housing.
CAUTION: The rotor may be destroyed during start-up if it is not inserted correctly.
 - a. Insert the rotor into the housing in such a way that the tubing connections on the rear of the rotor fit into the corresponding recesses inside the housing and one of the three holes of the rotor is visible in the inspection window at the bottom of the housing.
NOTE: The rotor's sealing surface is located approx. 4 mm deep inside the housing if the rotor is inserted correctly. If this is not the case, then the rotor has to be moved into the correct position from below by means of a pointed object.
6. Clean the connecting piece.
 - a. Clean the sealing surface of the connecting piece (2-**2**) with ethanol using a lint-free cloth.
7. Insert the connecting piece.
 - a. Insert the connecting piece into the housing in such a way that connector 1 is on top and the three knobs of the connecting piece fit into the corresponding recesses on the housing.
 - b. Reattach the union nut (2-**1**). At first, use only the finger to rotate it shut. Afterwards, tighten by applying more force until the full stop is reached (approx. 3/4 rotation) (do not use any tools).

8. Connect and condition the suppressor.
 - a. Reconnect all suppressor capillaries to the IC system (*see manual for the ion chromatograph*).
 - b. Before switching the suppressor over for the first time, rinse the three suppressor units with regeneration solution for five minutes.

Replacing parts

Parts of the suppressor must be replaced under the following circumstances:

- The suppressor capacity diminishes (indicated by an increased baseline $> 0.3 \mu\text{S}/\text{cm}$).
- The suppressor remains clogged (solutions can no longer be conveyed through the suppressor).

Procedure

1. Disconnect the suppressor from the IC system.
 - a. In the software, stop the hardware and wait until the pressure in the high-pressure pump has been released.
 - b. Disconnect all suppressor capillaries from the IC system.
2. Dismantle the suppressor.
 - a. Unscrew the union nut (2-**1**) from the housing (2-**4**).
 - b. Pull the connecting piece (2-**2**) out of the housing together with the rotor (2-**3**).
NOTE: If the rotor gets stuck in the housing, you can push it out as follows:
 Insert a pointed object into the inspection window (2-**5**) at the bottom of the housing. Use it to push out the rotor.
 - c. Loosen the connecting piece from the rotor.
3. Clean the new rotor.
 - a. Clean the sealing surface of the rotor with ethanol using a lint-free cloth.
4. Reinsert the rotor into the housing.

CAUTION: The rotor may be destroyed during start-up if it is not inserted correctly.

 - a. Insert the rotor into the housing in such a way that the tubing connections on the rear of the rotor fit into the corresponding recesses inside the housing and one of the three holes of the rotor is visible in the inspection window at the bottom of the housing.
NOTE: The rotor's sealing surface is located approx. 4 mm deep inside the housing if the rotor is inserted correctly. If this is not the case, then the rotor has to be moved into the correct position from below by means of a pointed object.
5. Clean the new connecting piece.
 - a. Clean the sealing surface of the connecting piece (2-**2**) with ethanol using a lint-free cloth.
6. Insert the connecting piece.
 - a. Insert the connecting piece into the housing in such a way that connector 1 is on top and the three knobs of the connecting piece fit into the corresponding recesses on the housing.
 - b. Reattach the union nut (2-**1**). At first, use only the finger to rotate it shut. Afterwards, tighten by applying more force until the full stop is reached (approx. 3/4 rotation) (do not use any tools).
7. Connect and condition the suppressor.
 - a. Reconnect all suppressor capillaries to the IC system (*see manual for the ion chromatograph*).
 - b. Before switching the suppressor over for the first time, rinse the three suppressor units with solution for five minutes.

Technical specifications


| Rotor | MSM-HC C |
|--------------|---|
| Structure | Micro packed-bed suppressor, robust, chemically resistant |
| Capacity | approx. 50 min* (0.25 meq) *with standard conditions: Metrosep C Supp 1 - 250/4.0; eluent 5 mmol/L HNO ₃ / 50 µg/L Rb ⁺ , flow rate 1.0 mL/min |
| Dead volume | < 250 µL |
| Regeneration | Chemical regeneration with STREAM |

| Rotor | MSM-HC C |
|--------------------------|---|
| Maximum flow rate | No irreversible damage possible |
| Backpressure limitations | No irreversible damage possible, no specific backpressure required |
| Solvent stability | 100% solvent stability |
| Temperature range | No limit |

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