IC-Probenvorbereitungskartuschen IC sample preparation cartridges Cartouches de préparation d'échantillons CI Cartuchos de preparación de muestras CI

Typ / Type / Type / Tipo	Anzahl / Quantity / Quantité / Cantidad	Nummer / Number / Numéro / Número
IC-RP	50 10	6.1012.000 6.1012.100
IC-H	50 10 25	6.1012.010 6.1012.110 6.1012.210
IC-Ag	50 10 25	6.1012.020 6.1012.120 6.1012.220
IC-OH	50 10	6.1012.030 6.1012.130
IC-Na	50	6.1012.040
IC-C18	50	6.1012.050

DE

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IC sample preparation cartridges

IC-RP

Sample preparation cartridge with C18 (octadecyl) for nonpolar solid-phase extraction.

Use

- If organic substances, which could interfere with the chromatogram or damage the separation column, have to be removed.
- For enriching heavy metals as a complex and subsequent determination (also for extraction

from an interfering matrix, e.g. heavy metal in soil).

Activation

- Before using the IC-RP cartridges, they have to be activated with 10 mL of methanol p.a. and rinsed with 10 mL of methanol.

NOTE

The sample preparation cartridges are supplied **WITHOUT** methanol.

Regeneration

- The regeneration is carried out analogous to the activation; however, the methanol quantity is increased to approx. 20 to 50 mL.

IC-H

Cation exchange cartridge in H⁺ form for anion analysis.

Use

- If the injection peak becomes too wide due to cations (e.g. Ca²⁺, Mg²⁺), thereby overlapping early eluting anions.
- For removing carbonate/hydrogen carbonate (in addition, poss. briefly treat in an ultrasonic bath or with nitrogen).
- For alkaline samples (NaOH + R-SO₃-H⁺ → H₂O + R-SO₃-Na⁺).
- For removing cations from samples for ion chromatography.

Preparation

- Prior to use, rinse with 10 mL deionized water. If possible, rinse with the sample afterwards to avoid diluting the sample by the dead volume of the cartridge.

Regeneration

- Slowly push 10 mL HCl p.a. 10% through the cartridge using a syringe or suction it through via vacuum. Afterwards, rinse the filter chloride-free using ultrapure water (control injection).

- If the regenerated exchange cartridge will not be used over an extended time period (2 days), it should be rinsed again with 5 mL ultrapure water prior to use.
- Depending on the ions to be determined, a different acid than HCl p.a. 10% can be used (e.g. HNO₃ p.a. 15%).

IC-Ag

Cation exchange cartridge in Ag⁺ form for anion analysis.

Use

- Removing the halogenides chloride, bromide, iodide (the content of fluoride, nitrite, nitrate and sulfate is hardly affected or not at all).

Preparation

- Prior to use, rinse with 10 mL deionized water. If possible, rinse with sample afterwards to avoid diluting the sample by the dead volume of the cartridge.
- To retain excess silver, the sample must be injected via the cation exchange cartridge IC-H (6.1012.010).

Regeneration

- This sample preparation cartridge cannot be regenerated.

IC-OH

Anion exchange cartridge in OH⁻ form for the cation analysis.

Use

- For heavily acidic samples (pH < 2), the pH value can be raised and subsequently be set to pH = 3 with HNO₃ (HCl + R-NH₃⁺ + OH⁻ → H₂O + R-NH₃ + Cl⁻).

Preparation

- Prior to use, rinse with 10 mL deionized water. If possible, rinse with sample afterwards to avoid diluting the sample by the dead volume of the cartridge.

Regeneration

- The regeneration of the anion exchangers is done analogous to that of the cation exchangers (see above). Instead of acid, leach (primarily NaOH p.a. or KOH p.a.) is used.

IC-Na

Cation exchange cartridge in Na⁺ form for anion analysis.

Use

- To exchange cations with Na⁺. It can also be used to remove or retain cations from the sample without changing the pH value of the sample.
- Useful for enriching cations or when removing cations from an interfering matrix.
- To remove cations from the sample in the ion-exclusion chromatography.

Preparation

- Prior to use, rinse with 10 mL deionized water. If possible, rinse with sample afterwards to avoid diluting the sample by the dead volume of the cartridge.

Regeneration

- Regeneration is not recommended.

IC-C18

Reversed phase cartridge C18 for manual sample preparation to remove lipophilic components.

Use

For the nonpolar solid-phase extraction. The cartridge removes organic substances.

Conditioning

- Rinse with 3 to 5 mL methanol.
- Rinse with 3 to 5 mL water or buffer solution. Do not allow the cartridge to dry out prior to adding the sample.

Application

- Add the sample at the top and push or pull through the cartridge at a rate of 1 to 5 mL/min.

If the desired component was not retained, provide the sample for the analysis.

If the desired component was retained:

- Rinse out retained interfering components with a polar solution.
- Elute the desired component with 1 to 2 mL nonpolar solution and provide it for the analysis.

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