DE

Dis-Cover-Schwenkarm (6.1462.080) für 786 Swing Head

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EN

Dis-Cover robotic arm (6.1462.080) for 786 Swing Head

This leaflet describes the installation and the mode of operation of the Dis-Cover robotic arm (6.1462.080) for the 786 Swing Head. The Dis-Cover robotic arm is used in the Sample Processor 814 / 815 (2.XXX.0110 / 2.XXX.0130).

Prerequisite:

- The 786 Swing Head is mounted.

Configuring the Sample Processor


Enter the following configuration data in the Sample Processor or in the control software prior to the assembly of the Dis-Cover robotic arm:

- Swing direction: +
- Robotic arm offset: 8°
- Maximum swing range: 117°
- Swing radius: 112 mm

Mounting the Dis-Cover robotic arm

The Dis-Cover robotic arm (6.1462.080) is designed for left-swinging use. Therefore, it is to be mounted on tower 2 of the Sample Processor.

- Switch off the instrument.
- Switch on the instrument.
- Mount the Dis-Cover robotic arm with the 3 supplied screws. Hold the robotic arm in such a way that it reaches the maximum swing range.
- Guide the connection cable of the Dis-Cover robotic arm through the guide chain of the Sample Processor.
- Connect the plug to the *ext. pump* socket on the rear of the tower.

EN

Mode of operation

The head of the Dis-Cover robotic arm consists of a permanent magnet and an electromagnet. When switched off, vessel covers with metal inserts are attracted by the permanent magnet. In order to remove a vessel cover (6.2037.050 or 6.2037.060) from the Dis-Cover robotic arm, the electromagnet is switched on with the PUMP command of the Sample Processor. The electrical field of the electromagnet neutralizes the magnetic field of the permanent magnet. The vessel cover is dropped.

CAUTION

Heating up of permanent magnet

High temperatures weaken the permanent magnet over time.

- To avoid the permanent magnet from heating up too much:
 - Do not keep the Dis-Cover robotic arm switched on for longer than 3 seconds.
 - Switch off the Dis-Cover robotic arm at the end of the sequence.

Switch on the Dis-Cover robotic arm with a PUMP command immediately before moving to the sample vessel to ensure that vessel covers are picked up in a controlled manner. The magnetic effect is thus suspended.

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