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New thermal switch

Thermal switch

In addition to the temperature control system of the heater, the instruments 873 Biodiesel Rancimat and 743 Rancimat are equipped with a safety switch that protects the instrument against overheating. The temperature control system serves to achieve and control the temperature defined in the Rancimat software. The safety switch is only triggered in case of a failure of the temperature control system. In such a case, a thermal switch disrupts the power supply to the heater at a certain temperature and thus prevents the instrument from heating up in an uncontrolled manner.

Biodiesel and biodiesel blends

The ignition temperature of biodiesel (FAME = fatty acid methyl ester) is approximately 250 °C, whereas the ignition temperature of diesel is approximately 220 °C. The ignition temperature (or autoignition temperature) describes the temperature at which a flammable substance can ignite in the presence of oxygen without a source of ignition and only due to the temperature. Autoignition can only occur if, in addition to a sufficiently high temperature, a specific mixture of flammable substance and oxygen is present in the gaseous phase. The value for diesel is between 0.6% and 7.5% at a natural ambient oxygen concentration of 21%.

873 Biodiesel Rancimat up to series 44

873 Biodiesel Rancimats up to serial number 1873001444xxx or 1873001544xxx were built and delivered ex works with a thermal switch which is activated at 260 °C (± 11 °C). As the ignition temperatures of biodiesel and diesel are below this temperature, spontaneous autoignition of the sample cannot be completely ruled out in case of a failure of the temperature control system. For safety reasons, we recommend that you replace the original thermal switch with the new thermal switch, which is activated already at a temperature of 180 °C (± 8 °C).

873 Biodiesel Rancimat series 45 and higher

873 Biodiesel Rancimats with the serial number 1873001445xxx or 1873001545xxx or higher were built and delivered ex works with the new thermal switch, which is activated at 180 °C (± 8 °C). (These instruments can also be identified by the label "Maximum operation temperature: 170 °C" attached to the back of the instrument.) No further measures are necessary for these instruments.

743 Rancimat used for biodiesel applications

743 Rancimats used for biodiesel applications are delivered ex works with a thermal switch which is activated at 260 °C (± 11 °C). As the ignition temperatures of biodiesel and diesel are below this temperature, spontaneous autoignition of the sample cannot be completely ruled out in case of a failure of the temperature control system. For safety reasons, we recommend that you replace the original thermal switch with the new thermal switch, which is activated already at a temperature of 180 °C (± 8 °C).

Note

If a Rancimat is equipped with a new thermal switch with a shutdown temperature of 180 °C, no methods involving temperatures higher than 170 °C must be used. In the Rancimat software, the setting of higher temperatures is not disabled; however, such a setting would prompt the thermal switch to shut down the heater.