

Column manual



Metrosep A Supp 19 (6.01034.4x0)

Manual

8.0107.8013EN / v1 / 2023-03-08



Metrohm AG
CH-9100 Herisau
Switzerland
+41 71 353 85 85
info@metrohm.com
www.metrohm.com

Column manual

Metrosep A Supp 19 (6.01034.4x0)

Manual

This documentation is protected by copyright. All rights reserved.

This documentation has been prepared with great care. However, errors can never be entirely ruled out. Please send comments regarding possible errors to the address above.

Disclaimer

Deficiencies arising from circumstances that are not the responsibility of Metrohm, such as improper storage or improper use, etc., are expressly excluded from the warranty. Unauthorized modifications to the product (e.g. conversions or attachments) exclude any liability on the part of the manufacturer for resulting damage and its consequences. Instructions and notes in the Metrohm product documentation must be strictly followed. Otherwise, Metrohm's liability is excluded.

Table of contents

1	General information	1
1.1	Ordering information	1
1.2	Technical specifications	1
2	Key aspects of working with separation columns	3
3	Eluent production	7
3.1	Production of standard eluent	7
4	Start-up	8
4.1	Connecting and rinsing the guard column	8
4.2	Connecting and rinsing the separation column	10
4.3	Conditioning	14
5	Applications	16
5.1	Standard chromatogram	16
5.2	Effects of temperature	17
5.3	Variation of the eluent flow rate	21
5.4	Variation of the eluent	22
5.4.1	Constant Na ₂ CO ₃ /NaHCO ₃ ratio	22
5.4.2	NaHCO ₃ variation at constant Na ₂ CO ₃	25
5.4.3	Na ₂ CO ₃ variation at constant NaHCO ₃	29
5.5	Variation with organic modifier	31
5.5.1	Variation of the acetone concentration	31
5.5.2	Variation of the methanol concentration	34
5.5.3	Variation of the acetonitrile concentration	38
5.5.4	Variation of the ethanol concentration	41
5.6	Determination of standard anions in mineral water samples	45
5.7	Determination of standard anions and organic acids in boiler feed water of power plants	46
5.8	Direct determination of standard anions in bioethanol ..	47
5.9	Determination of fluoride in dental gel according to USP	48
5.10	Direct determination of standard anions in lactose-free milk	49
5.11	Determination of standard anions and 13 organic acids in food samples with conductivity	50

1 General information

The Metrosep A Supp 19 is a high-performance separation column and is especially suitable for the determination of inorganic anions and low-molecular organic anions with chemical and sequential suppression. Due to its high capacity, the Metrosep A Supp 19 can readily handle samples with high ionic strength and with large fluctuations in concentration. The outstanding peak symmetries and high number of theoretical plates allow universal use in ion chromatography.

1.1 Ordering information

Table 1 Separation columns

Order number	Designation
6.01034.410	Metrosep A Supp 19 - 100/4.0
6.01034.420	Metrosep A Supp 19 - 150/4.0
6.01034.430	Metrosep A Supp 19 - 250/4.0

Table 2 Guard column

Order number	Designation
6.01034.500	Metrosep A Supp 19 Guard/4.0

1.2 Technical specifications

Column material Hydrophilic polystyrene/divinylbenzene copolymer with quaternary ammonium groups

Particle size 4.6 µm

Measurements

Order number	Measurements
6.01034.410	100 x 4.0 mm
6.01034.420	150 x 4.0 mm
6.01034.430	250 x 4.0 mm

pH range eluent 0–14

ph range sample 0–14

Temperature range 10–70 °C

25 °C

Maximum pres-
sure

Order number	Maximum pressure
6.01034.410	20 MPa (200 bar)
6.01034.420	25 MPa (250 bar)
6.01034.430	25 MPa (250 bar)

Flow rate

Order number	Recommended flow rate	Maximum flow rate
6.01034.410	0.7 mL/min	1.3 mL/min
6.01034.420	0.7 mL/min	1.2 mL/min
6.01034.430	0.7 mL/min	1.0 mL/min

Standard eluent

8.0 mmol/L sodium carbonate (Na_2CO_3) and 0.25 mmol/L sodiumhydrogencarbonate (NaHCO_3)

Permitted organic additives

0–100% acetonitrile, acetone and methanol

Capacity

Order number	Capacity
6.01034.410	94 μmol (Cl ⁻)
6.01034.420	140 μmol (Cl ⁻)
6.01034.430	234 μmol (Cl ⁻)

Preparation

Use a flow gradient to set the column to the standard flow within 5 minutes. Then wait until the baseline is given.

Storage

Store the column in the standard eluent and at 4 to 30 °C.

Typical pressure

For columns with a guard column under standard conditions with chemical suppression:

Order number	Typical pressure
6.01034.410	11 MPa
6.01034.420	14 MPa
6.01034.430	18 MPa

Column housing

Smart column with a chip, called an iColumn, made of PEEK

Application

Determination of inorganic anions and low-molecular organic anions with chemical and sequential suppression

2 Key aspects of working with separation columns

<i>Storage</i>	Once the backpressure in the ion chromatograph has dissipated, remove the column at ambient temperature. Seal the column at both ends using the original stoppers (6.2744.060). Store the column in the standard eluent and at 4 to 30 °C.
<i>Bacterial growth</i>	<p>Bacterial growth significantly affects chromatography and ruins separation columns. A vast array of problems in chromatography are caused by the growth of algae, bacteria and fungi.</p> <p>In order to prevent bacterial growth, always use fresh eluents, rinsing solutions and regeneration solutions. Do not use any solutions that have not been used for a prolonged period. Metrohm recommends cleaning all vessels as follows before filling them:</p> <ol style="list-style-type: none"> 1. Thoroughly rinse with ultrapure, UV-treated water (> 18.2 MΩ). 2. Swirl an acetonitrile-water mixture around in the vessel. 3. Rinse again with ultrapure water. <p>If you notice the growth of bacteria or algae despite these precautionary measures, add 5% methanol, acetonitrile or acetone to the eluent. This is only possible if you are <i>not using membrane suppressors</i>. Organic solvents can destroy membrane suppressors. The Metrohm Suppressor Modules (MSM, MSM-HC and MSM-LC) are 100% resistant to solvents.</p>
<i>Chemical quality</i>	All chemicals must have at least a quality of p.a. or puriss. Standard solutions must be intended specifically for ion chromatography.
<i>Chemical stress</i>	Even though specifications may indicate that separation phases do cover a large pH range, this does not mean they are chemically inert. Separation columns last longest when subjected to constant chemical conditions. Never allow a column to dry out and ensure the column is sealed well at all times.
<i>CO₂</i>	Carbon dioxide from the air affects the carbonate / hydrogen carbonate balance in the eluent. The eluent becomes weaker over time. In order to prevent this, always outfit the eluent bottle with a CO ₂ adsorber with the adsorber material soda lime.
<i>Degassing the eluent</i>	In order to prevent bubbles from forming, degas the produced eluent before using it in the IC system. To do this, create a vacuum for approximately 10 minutes using a water-jet pump or a membrane pump. Alternatively, use an ultrasonic bath or work with an eluent degasser.

Regenerating separation columns

switching the valves. Using the pulsation absorber (6.2620.150) already built into the Metrohm ion chromatographs provides this protection.

If separation columns are operated with clean eluents and filled with samples free of particles, you can expect the column to have a long service life. This means regenerating the column is not required and, after a multitude of injections, no longer possible.

If the pressure in the column increases despite this or if the separating efficiency decreases, carry out the specified regeneration steps. Perform the regeneration outside the analytical line. For regeneration, connect the separation column to the pump directly. Route the regeneration solution through the column directly into a waste container. Rinse the separation column properly with fresh eluent. Then, reinstall the separation column.

Shutting down the ion chromatograph

If the ion chromatograph is not used for a prolonged period (> 1 week), remove the separation column and seal it with the stoppers provided. Rinse the ion chromatograph, including all 3 suppressor chambers, with a methanol-water mixture (1:4). Store the separation column in the medium indicated on the column leaflet. Unless otherwise stated on the column leaflet, store the column at 4 to 30 °C.

Prior to start-up, rinse the ion chromatograph with ultrapure water and then with fresh eluent. Bring the separation column back to ambient temperature before you install it. Then increase the temperature if necessary.

Fun

Ion chromatography should be fun and should not stress you out. Metrohm puts all its work into ensuring you can work reliably with your IC systems with minimal maintenance, servicing and costs. Metrosep separation columns embody the attributes of quality, long service life and excellent results.

Environmental protection

A significant advantage of ion chromatography is that most of the work involves aqueous media. As a result, the chemicals used in ion chromatography are largely non-toxic and do not impact the environment. When working with acids, bases, organic solvents or heavy metal standards, dispose of them properly after use.

Guard columns

Guard columns are used to protect separation columns. Metrohm strongly recommends using guard columns. Guard columns normally contain the same stationary phase as separation columns. However, the quantity is significantly reduced to avoid impacting the chromatography. Guard columns remove critical contaminants that can react with column material. Guard columns also remove particles and bacterial contaminants. Replace the guard column in the following cases:

- If the backpressure in the system increases

3 Eluent production

Metrohm recommends selecting a high degree of purity for chemicals for standard production and eluent production.

3.1 Production of standard eluent

Proceed as follows to produce 2 L of standard eluent with 0.25 mmol/L of sodium hydrogen carbonate and 8.0 mmol/L of sodium carbonate:

Producing 2 L of standard eluent

Accessories

- Eluent bottle (6.1608.120)
- Cover (6.1602.200) equipped with CO₂ adsorber
- Ultrapure water
- Sodium carbonate
- Sodium hydrogen carbonate

1 Pre-rinse the eluent bottle with ultrapure water several times.

2 Fill 2 L of ultrapure water into the eluent bottle.

3 Degas the ultrapure water.

Use the eluent degasser.

If no eluent degasser is available, degas the ultrapure water for 5 to 10 minutes using a vacuum pump. Degassing avoids problems with air bubbles in the high-pressure pump.

- 4**
- Weigh 42.0 mg of sodium hydrogen carbonate.
 - Weigh 1,695.8 mg of sodium carbonate.
 - Add the weighed amounts of sodium hydrogen carbonate and sodium carbonate to the ultrapure water.

5 Rinse the column with eluent for 1 h.

This eluent (0.25 mmol/L of sodium hydrogen carbonate and 8.0 mmol/L of sodium carbonate) and chemical suppression can be used to achieve background conductivity of less than 23 µS/cm. The noise is typically less than 0.2 nS/cm.

4 Start-up

4.1 Connecting and rinsing the guard column

Guard columns protect separation columns and significantly increase their service life. The guard columns available from Metrohm are either actual guard columns or guard column cartridges used together with a cartridge holder. The process of installing a guard column cartridge into the corresponding holder is described in the cartridge leaflet.



NOTICE

Metrohm recommends always working with guard columns. Guard columns protect the separation columns and can be replaced regularly as needed.



NOTICE

Information regarding which guard column is suitable for your separation column can be found in the **Metrohm Column Program** (which is available from your regional Metrohm representative), the column leaflet and the product information or in consultation with your regional Metrohm representative.

You can find product information for your separation column at <http://www.metrohm.com> in the Ion Chromatography product area.



CAUTION

New guard columns are filled with solution and sealed with stoppers or caps on both sides.

Before inserting the guard column, ensure that this solution can be mixed with the eluent being used (follow the manufacturer specification).


**NOTICE**

Do not connect the guard column until after the initial start-up of the instrument. Until then, replace the guard column and the separation column with couplings (6.2744.040).

Accessories

For this step, you need the following accessories:

- Guard column (suitable for separation column)

Connecting the guard column**1 Removing the coupling**

Remove the coupling (6.2744.040) installed between the column inlet capillary and the column outlet capillary for the initial start-up.

2 Preparing the guard column

- Remove the stoppers or the stopper and the sealing cap from the guard column.

3 Connecting the guard column**CAUTION**

When inserting the guard column, ensure that it is inserted correctly based on the marked flow direction (if specified).

- Fasten the inlet of the guard column to the column inlet capillary using a short pressure screw (6.2744.070).
- If the guard column is connected to the separation column using a connection capillary, fasten this connection capillary to the guard column outlet with a pressure screw.

Rinsing the guard column

1 Rinsing the guard column

- Place a beaker under the guard column's outlet.
- Start manual control in MagIC Net and select the high-pressure pump: **Manual ► Manual control ► Pump**
 - **Flow:** in accordance with column leaflet
 - **On**
- Rinse the guard column with eluent for approx. 5 minutes.
- Stop the high-pressure pump in the manual control in MagIC Net again: **Off**.

4.2 Connecting and rinsing the separation column

The smart separation column (iColumn) is the heart of ion chromatographic analysis. It separates the different components according to their interactions with the column. Metrohm separation columns are equipped with a chip where their technical specifications and history (start-up, operating hours, injections etc) are stored.



NOTICE

Information regarding which separation column is suitable for your application can be found in the **Metrohm Column Program**, the product information for the separation column or it can be obtained from your regional Metrohm representative.

You can find product information for your separation column at <http://www.metrohm.com> in the Ion Chromatography product area.




A test chromatogram accompanies every column. The column leaflet can be found online at <http://www.metrohm.com> with the corresponding article. Detailed information on special IC applications can be found in the corresponding **Application Bulletins** or **Application Notes**. You can find these online at <http://www.metrohm.com> in the Applications area or request them free of charge from your responsible regional Metrohm representative.

**CAUTION**

New IC Columns are filled with solution and sealed with stoppers on both sides. Before inserting the column, ensure that this solution can be mixed with the eluent being used (follow the information provided by the manufacturer).

**NOTICE**

Connect the separation column only after the initial start-up of the instrument. Until that point, insert a coupling (6.2744.040) instead of the guard column and separation column.



1 Removing the stoppers

- 12 ■■■■■■

2 Installing the inlet of the separation column



CAUTION

When inserting the column, ensure that it is inserted correctly based on the marked flow direction.

There are 3 possibilities:

- Attach the column inlet directly onto the guard column or,
- if the guard column is connected to the separation column using a connection capillary: Connect the column inlet to the guard column outlet capillary using a PEEK pressure screw (6.2744.070) or,
- if no guard column is used (not recommended): Connect the column inlet capillary to the inlet of the separation column using a short pressure screw (6.2744.070).

3 Rinsing the separation column

- Place a beaker under the outlet of the separation column.
- Start manual control in MagIC Net and select the high-pressure pump: **Manual ► Manual control ► Pump**
 - **Flow:** Increase gradually up to the flow rate recommended in the column leaflet.
 - **On**
- Rinse the separation column with eluent for approx. 10 minutes.
- Stop the high-pressure pump in the manual control in MagIC Net again: **Off**.

4 Removing the coupling

- Remove the coupling (6.2744.040) from the column outlet capillary.

5 Installing the outlet of the separation column

- Fasten the column outlet capillary to the column outlet using a short PEEK pressure screw (6.2744.070).

6 Inserting the separation column

- Insert the separation column with the chip into the column holder until you hear it snap in place.

The separation column is now detected by MagIC Net.

4.3 Conditioning

In the following cases, the system must be conditioned with eluent until a stable baseline has been reached:

- After installation
- After each time the instrument is switched on
- After each eluent change



NOTICE

The conditioning time can lengthen considerably if the composition of the eluent is modified.

Conditioning the system

1 Preparing the software



CAUTION

Ensure that the configured flow rate is not higher than the flow rate permitted for the corresponding column (refer to the column leaflet and chip data record).

- Start the **MagIC Net** computer program.
 - Open the **Equilibration** tab in MagIC Net: **Workplace ► Run ► Equilibration**.
 - Select (or create) a suitable method.
- See also: *MagIC Net tutorial* and online help.

2 Preparing the instrument

- Check whether the column is inserted correctly in accordance with the flow direction marked on the sticker (arrow has to point in the flow direction).
- Check whether the eluent aspiration tubing is immersed in the eluent and that there is enough eluent in the eluent bottle.

3 Starting the equilibration

- Start the equilibration in MagIC Net: **Workplace ► Run ► Equilibration ► Start HW.**

- Visually inspect whether all capillaries and their connections from the high-pressure pump to the detector are leak-tight. If eluent is leaking out anywhere, tighten the corresponding pressure screw further, or loosen the pressure screw, check the end of the capillary and shorten it using the capillary cutter if necessary and retighten the pressure screw.

4 Conditioning the system

Continue rinsing the system with eluent until the desired stability level for the baseline has been attained .

The instrument is now ready for measuring samples.

	Metrosep A Supp 19 - XXX/4.0	mg/L
1	Fluoride	2
2	Chloride	2
3	Nitrite	5
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

5.2 Effects of temperature

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity


Suppression: Sequential suppression with MSM and MCS

Temperature: 10–70 °C

Loop: 20 µL

Flow rate: 0.7 mL/min

Eluent: 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃




Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

The Metrosep A Supp 19 can be used at temperatures from 10 to 70 °C. The retention times of the monovalent anions decrease with increases in temperature. Especially the polarizable ions nitrite, bromide and nitrate are strongly accelerated at higher temperature. At the same time, the peaks become somewhat sharper at higher temperatures. The retention time of phosphate increases slightly when the temperature increases. The retention time of sulfate increases significantly with increasing temperatures. At 10 °C, phosphate and sulfate are not completely separated.

Increasing the temperature also causes the column backpressure to decrease considerably. At 70 °C, the backpressure is only approx. 7.2 MPa, whereas the column backpressure at 10 °C is more than twice as high (about 17.8 MPa). Due to the high backpressure at low temperatures, standard flow measurements at 10 °C are not possible on the long Metrosep A Supp 19 - 250/4.0 column.

Organic acids




	Metrosep A Supp 19 - 150/4.0	mg/L
1	Glycolate	10
2	Acetate	10
3	Formate	10
4	Pyruvate	10
5	Tartrate	10
6	Succinate	10
7	Oxalate	10

The monovalent organic acids such as glycolate, acetate, formate and pyruvate elute between fluoride and chloride. Glycolate, acetate and formate are largely separated from one another at 30 °C. Their retention times are affected only slightly by temperature. The retention time of pyruvate, on

the other hand, increases with increasing temperature. At 10 °C, formate and pyruvate coelute. At 70 °C, pyruvate is separated from formate and elutes after formate. The divalent organic acids, such as tartrate, succinate and oxalate, behave similarly to the divalent inorganic anions, such as sulfate: With increasing temperature, the retention times also increase.

Haloacetic acids



Metrosep A Supp 19 - 150/4.0		mg/L
1	Monochloroacetate	10
2	Monobromoacetate	10
3	Dichloroacetate	10
4	Dibromoacetate	10
5	Trichloroacetate	10

The haloacetic acids are another class of components that are often analyzed in ion chromatography. Although haloacetic acids are monovalent ions, their reaction to temperature influences differs from that of inorganic anions. When the temperature increases, the retention times of all ions increase. It should be noted that haloacetic acids degrade at higher temperatures, making them difficult to analyze.

5.3 Variation of the eluent flow rate

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity


Suppression: Sequential suppression with MSM and MCS

Temperature: 30 °C

Loop: 20 µL

Flow rate: 0.3–0.9 mL/min

Eluent: 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃



Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

The Metrosep A Supp 19 - 150/4.0 can be operated with a flow of up to 1.2 mL/min. As the flow rate increases, all ions are accelerated uniformly, which means that sulfate elutes in less than 14 minutes at 0.9 mL/min. The pressure increases almost proportionally to the flow. Due to the higher flow rate, the dwell time of the analytes in the detector is reduced, resulting in smaller peak areas. The long Metrosep A Supp 19 - 250/4.0 can be operated with a maximum flow of 1.0 mL/min. The temperature must be increased in the event of high flow rates on the long column. If not, there will be an overpressure on the column.

5.4 Variation of the eluent

5.4.1 Constant $\text{Na}_2\text{CO}_3/\text{NaHCO}_3$ ratio

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: —

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS


Temperature: 30 °C

Loop: 20 μ L

Flow rate: 0.7 mL/min


Eluent:

- A) 0.125 mmol/L NaHCO₃, 4.0 mmol/L Na₂CO₃
- B) 0.188 mmol/L NaHCO₃, 6.0 mmol/L Na₂CO₃
- C) 0.250 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃
- D) 0.313 mmol/L NaHCO₃, 10.0 mmol/L Na₂CO₃
- E) 0.375 mmol/L NaHCO₃, 12.0 mmol/L Na₂CO₃




Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

With increasing eluent concentration, all anions are significantly accelerated. In this process, the multivalent anions phosphate and sulfate are accelerated faster than the monovalent anions. With a stronger eluent, the peaks are sharper and correspondingly higher. For a strong eluent (0.313 mmol/L NaHCO_3 , 10 mmol/L Na_2CO_3), nitrate coelutes with phosphate. For an even stronger eluent (0.375 mmol/L NaHCO_3 , 12 mmol/L Na_2CO_3), the elution order is phosphate, nitrate, sulfate. The separation between nitrate and sulfate is not optimal.



	Metrosep A Supp 19 - 150/4.0	mg/L
1	Glycolate	10
2	Acetate	10
3	Formate	10
4	Tartrate	10
5	Succinate	10
6	Oxalate	10

The organic acids are greatly accelerated with increasing eluent concentration. The monovalent organic acids such as glycolate, acetate and formate are accelerated to an extent similar to that of fluoride or chloride. The divalent organic acids such as tartrate, succinate and oxalate, on the other hand, behave like sulfate and are accelerated to a much greater extent. As the eluent concentration increases, the peaks become sharper and higher. This is particularly visible with divalent organic acids.



Metrosep A Supp 19 - 150/4.0		mg/L
1	Monochloroacetate	10
2	Monobromoacetate	10
3	Dichloroacetate	10
4	Dibromacetate	10
5	Trichloroacetate	10

The haloacetic acids react strongly on the Metrosep A Supp 19 to the eluent strength. With increasing eluent strength, all haloacetic acids are greatly accelerated.

5.4.2 NaHCO_3 variation at constant Na_2CO_3

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS

[illegible]

Eluent:

- A) 0.250 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3
- B) 1.0 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3
- C) 4.125 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3
- D) 8.0 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3

Figure 1 displays four stacked conductivity chromatograms showing the separation of Na_2CO_3 and NaHCO_3 mixtures. The y-axis represents Conductivity [$\mu\text{S}/\text{cm}$] from 0 to 100, and the x-axis represents Time [min] from 0 to 25. The traces are labeled with their respective concentrations:

- Top trace (red): 8.00 mmol/L Na_2CO_3 and 8.00 mmol/L NaHCO_3 . Peaks are labeled 1, 2, 3, 4, 5, and 7.
- Second trace (black): 8.00 mmol/L Na_2CO_3 and 4.13 mmol/L NaHCO_3 . Peak 6 is labeled.
- Third trace (blue): 8.00 mmol/L Na_2CO_3 and 1.00 mmol/L NaHCO_3 . Peak 7 is labeled.
- Bottom trace (orange): 8.00 mmol/L Na_2CO_3 and 0.25 mmol/L NaHCO_3 . No peaks are labeled.


Dashed lines connect the peaks across the traces, indicating the elution time of each component. The baseline conductivity decreases from top to bottom trace due to the decreasing concentration of NaHCO_3 .

Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10

Metrosep A Supp 19 - 150/4.0		mg/L
7	Sulfate	10

Sodium hydrogen carbonate influences the eluent strength much less than sodium carbonate. The retention times of the anions are therefore only slightly shortened by increasing the content of sodium hydrogen carbonate in the eluent. Only the retention time of phosphate is shortened to a significantly greater extent, because the pH value of the eluent changes significantly, which means that the effective charge of the phosphate ion is also reduced. The retention time of sulfate also shortens as the hydrogen carbonate content increases. No significant change in peak heights was observed in the tested range of sodium hydrogen carbonate. In an eluent composition of 4.125–8.000 mmol/L NaHCO₃ and 8.0 mmol/L Na₂CO₃, bromide and phosphate are not optimally separated.

Organic acids



Metrosep A Supp 19 - 150/4.0		mg/L
1	Glycolate	10
2	Acetate	10
3	Formate	10

The hydrogen carbonate content in the eluent does not affect the haloacetic acids. The retention times and shapes of the peaks of these analytes remain constant and are not affected by hydrogen carbonate.

5.4.3 **Na₂CO₃ variation at constant NaHCO₃**

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS


Temperature: 30 °C

Loop: 20 µL

Flow rate: 0.7 mL/min

Eluent:

- A) 0.25 mmol/L NaHCO₃, 4.0 mmol/L Na₂CO₃
- B) 0.25 mmol/L NaHCO₃, 6.0 mmol/L Na₂CO₃
- C) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃
- D) 0.25 mmol/L NaHCO₃, 10.0 mmol/L Na₂CO₃
- E) 0.25 mmol/L NaHCO₃, 12.0 mmol/L Na₂CO₃



Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

Due to the strong elution strength of sodium carbonate, the influence of the sodium carbonate content in the eluent is much stronger than the influence of the sodium hydrogen carbonate content. Increasing the sodium carbonate content in the eluent significantly reduces the retention times of all anions. The multivalent anions phosphate and sulfate are accelerated the most here. In a weak eluent (0.25 mmol/L NaHCO_3 and 4.0 mmol/L Na_2CO_3), phosphate and sulfate elute late after nitrate. In a concentrated eluent (0.25 mmol/L NaHCO_3 and 10.0 mmol/L Na_2CO_3), the baseline separation between nitrate and sulfate is no longer present. In a very strong eluent (0.25 mmol/L NaHCO_3 and 12.0 mmol/L Na_2CO_3), phos-

phate elutes before nitrate. The acceleration of the anions due to the use of a stronger eluent also causes higher peaks.

5.5 Variation with organic modifier

5.5.1 Variation of the acetone concentration

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS


Temperature: 30 °C

Loop: 20 µL

Flow rate: 0.7 mL/min

Eluent:

- A) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 0% acetone
- B) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 10% acetone
- C) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 20% acetone




Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

In some cases, the use of an organic modifier is useful or even necessary. The eluent can be made more stable against bacterial contamination by adding a modifier, or the modifier can help improve the rinsing-out of the organic parts of a sample from the separation column. With the addition of an organic modifier, the backpressure of the column and the selectivity of all anions change. An increase of the acetone content in the eluent accelerates the elution of chloride, nitrite, bromide and nitrate. Phosphate and sulfate, on the other hand, are retarded. Adding acetone to the eluent tends to cause deterioration of the shapes of the peaks and thus also the peak heights.


The column backpressure increases with increasing acetone content in the eluent due to the increasing viscosity of the eluent mixture. Without acetone in the eluent, the pressure is approx. 12 MPa. With 20% acetone in the eluent, it is already approx. 20 MPa.

Organic acids



Metrosep A Supp 19 - 150/4.0		mg/L
1	Glycolate	10
2	Acetate	10
3	Formate	10
4	Tartrate	10
5	Succinate	10
6	Oxalate	10

The reactions of the different organic acids vary greatly when acetone is added to the eluent. Glycolate, formate, tartrate and oxalate remain almost unchanged with 0 to 20% acetone in the eluent. Acetate and succinate, on the other hand, are accelerated by the addition of acetone. Glycolate and acetate as well as tartrate and succinate coelute with 10% acetone. With 20% acetone, succinate overtakes tartrate.



Metrosep A Supp 19 - 150/4.0		mg/L
1	Monochloroacetate	10
2	Monobromoacetate	10
3	Dichloroacetate	10
4	Dibromoacetate	10
5	Trichloroacetate	10

Acetone in the eluent accelerates the haloacetic acids. The retention times of all haloacetic acids shorten uniformly and linearly with increasing acetone content.

5.5.2 Variation of the methanol concentration

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: —

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS


Temperature: 30 °C

Loop: 20 µL

Flow rate: 0.7 mL/min

Eluent:
 A) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 0% methanol
 B) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 10% methanol
 C) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 20% methanol

Standard anions




Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

Methanol is also commonly used as an organic modifier. The retention times of the anions fluoride, chloride, nitrite, bromide and nitrate do not change with the addition of methanol. In contrast, the retention times of phosphate and sulfate increase with increasing methanol content. The elution order of the standard anions does not change due to methanol in the eluent. In contrast, the peak areas and peak heights of nitrite and phosphate in particular become significantly smaller with increasing methanol content.

Similar to acetone, methanol contributes to an increase in eluent viscosity. As a result, the column backpressure increases significantly, from approx. 12 MPa with 0% methanol in the eluent to approx. 20 MPa with 20% methanol.

Organic acids




Metrosep A Supp 19 - 150/4.0		mg/L
1	Glycolate	10
2	Acetate	10
3	Formate	10
4	Tartrate	10
5	Succinate	10

Metrosep A Supp 19 - 150/4.0		mg/L
6	Oxalate	10

The monovalent organic acids, similar to the monovalent inorganic anions, react hardly at all to the addition of methanol to the eluent. Only the shapes of the peak are affected by the addition of methanol to the eluent. Glycolate and acetate coelute starting with 10% methanol in the eluent. In the case of multivalent organic acids, succinate behaves differently than tartrate and oxalate. Succinate undergoes minimal additional retardation with increasing methanol content. Tartrate and oxalate, by contrast, exhibit a significant shift in retention time. This reverses the elution order between succinate and tartrate.

Haloacetic acids



Metrosep A Supp 19 - 150/4.0		mg/L
1	Monochloroacetate	10
2	Monobromoacetate	10
3	Dichloroacetate	10
4	Dibromoacetate	10
5	Trichloroacetate	10

The smaller haloacetic acids react hardly at all to the addition of methanol to the eluent. Only the retention time of trichloroacetate decreases with increasing methanol content.

5.5.3 Variation of the acetonitrile concentration

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: —

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS

Temperature: 30 °C


Loop: 20 μ L

Flow rate: 0.7 mL/min

Eluent:

- A) 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3 , 0% acetonitrile
- B) 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3 , 10% acetonitrile
- C) 0.25 mmol/L NaHCO_3 , 8.0 mmol/L Na_2CO_3 , 20% acetonitrile

Standard anions




Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

The addition of acetonitrile as an organic modifier contributes significantly less to the pressure increase than does the addition of acetone or methanol. Even with 20% acetonitrile in the eluent, the column backpressure increases by only approx. 33%.

The effect of acetonitrile on retention times is very similar to the effect of acetone: Fluoride, chloride, nitrite, bromide and nitrate elute earlier as acetonitrile content increases. The retention times of phosphate and sulfate are prolonged. In the process, the peak heights and peak areas become only slightly smaller. With 20% acetonitrile, phosphate and sulfate coelute.

Organic acids



	Metrosep A Supp 19 - 150/4.0	mg/L
1	Monochloroacetate	10
2	Monobromoacetate	10
3	Dichloroacetate	10
4	Dibromoacetate	10
5	Trichloroacetate	10

Acetonitrile has only a minimal effect on the retention times of monochloroacetate and monobromoacetate. In contrast to this, the haloacetic acids dichloroacetate, dibromoacetate, and trichloroacetate are significantly accelerated with the addition of acetonitrile.

5.5.4 Variation of the ethanol concentration

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity

Suppression: Sequential suppression with MSM and MCS


Temperature: 30 °C

Loop: 20 µL

Flow rate: 0.7 mL/min

Eluent:

- A) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 0 % ethanol
- B) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 10% ethanol
- C) 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃, 20% ethanol




Metrosep A Supp 19 - 150/4.0		mg/L
1	Fluoride	10
2	Chloride	10
3	Nitrite	10
4	Bromide	10
5	Nitrate	10
6	phosphate	10
7	Sulfate	10

Ethanol is also used as an organic modifier in isolated cases. Ethanol increases the eluent viscosity by far the most, which is then reflected in the column backpressure. Adding 20% ethanol almost doubles the pressure.

The retention times of the monovalent anions are only slightly affected by the addition of ethanol. The peak height and peak area of nitrite and phosphate decrease with increasing ethanol content. The multivalent anions are retarded with the addition of ethanol. This improves the separation.


ration between phosphate and sulfate. The effect of ethanol on the separation behavior is comparable to that of methanol.

Organic acids



Metrosep A Supp 19 - 150/4.0		mg/L
1	Glycolate	10
2	Acetate	10
3	Formate	10
4	Tartrate	10
5	Succinate	10
6	Oxalate	10

The retention times of the monovalent organic acids are not affected by the addition of ethanol to the eluent. Due to peak widening, glycolate coelutes with acetate starting with 10% ethanol. Succinate is slightly accelerated with the addition of ethanol, whereas tartrate and oxalate become retarded with increasing ethanol content.



	Metrosep A Supp 19 - 150/4.0	mg/L
1	Monochloroacetate	10
2	Monobromoacetate	10
3	Dichloroacetate	10
4	Dibromoacetate	10
5	Trichloroacetate	10

Except for trichloroacetate, the haloacetic acids show no change following addition of ethanol to the eluent. Their retention time does not change, regardless of ethanol content. Trichloroacetate is slightly accelerated with increasing ethanol content.

5.6 Determination of standard anions in mineral water samples

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: –

Detection: Conductivity


Suppression: Sequential suppression with MSM and MCS

Temperature: 25 °C

Loop: 20 µL


Flow rate: 0.7 mL/min

Eluent: 0.25 mmol/L NaHCO₃, 8.0 mmol/L Na₂CO₃



5.8 Direct determination of standard anions in bioethanol

Column:	Metrosep A Supp 19 - 150/4.0
Sample preparation:	–
Detection:	Conductivity
Suppression:	Sequential suppression with MSM and MCS
Temperature:	30 °C
Loop:	20 µL
Flow rate:	0.7 mL/min
Eluent:	0.25 mmol/L NaHCO ₃ , 8.0 mmol/L Na ₂ CO ₃



The red chromatogram corresponds to the direct injection of the ethanol sample onto the column. The Metrosep A Supp 19 - 150/4.0 is inert when injecting organic solvents, such as ethanol. The black chromatogram

shows the same sample to which 2 µg/L of the standard anions were added.

5.9 Determination of fluoride in dental gel according to USP

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: —

Detection: Conductivity


Suppression: Sequential suppression with MSM and MCS

Temperature: 40 °C

Loop: 20 μ L

Flow rate: 1.0 mL/min

Eluent: 15.0 mmol/L KOH




In this method, the current monograph was repeated with the Metrosep A Supp 19. The red chromatogram corresponds to the dental

gel that contains a nominal 1 µg/mL NaF. The black chromatogram shows the same sample to which 1 µg/mL NaF was added. Thus, it shows that the Metrosep A Supp 19 is also suitable for USP methods.

5.10 Direct determination of standard anions in lactose-free milk


<i>Column:</i>	Metrosep A Supp 19 - 150/4.0
<i>Sample preparation:</i>	Dilution 1:50
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM and MCS
<i>Temperature:</i>	30 °C
<i>Loop:</i>	20 µL
<i>Flow rate:</i>	0.75 mL/min
<i>Eluent:</i>	0.25 mmol/L NaHCO ₃ , 8.0 mmol/L Na ₂ CO ₃




Lactose-free milk was analyzed with Metrosep A Supp 19 for standard anions (fluoride, chloride, nitrite, bromide, nitrate, phosphate and sulfate). For this purpose, the sample was diluted 1:50 in ultrapure water. The red chromatogram shows the original sample. The black chromatogram shows the same sample to which 1 µg/L of the standard anions was added. The standard anions can be readily determined under these conditions. An organic acid elutes between fluoride and chloride.

5.11 Determination of standard anions and 13 organic acids in food samples with conductivity

<i>Column:</i>	Metrosep A Supp 19 - 250/4.0
<i>Sample preparation:</i>	-
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM and MCS
<i>Temperature:</i>	25 °C
<i>Loop:</i>	20 µL
<i>Flow rate:</i>	0.7 mL/min
<i>Eluent:</i>	A) 1.0 mmol/L Na ₂ CO ₃ B) 50.0 mmol/L Na ₂ CO ₃ 0–18 min: 100% A 18–20 min: 100–90% A 20–70 min: 90% A 70–75 min: 90–40% A 75–107 min: 40% A 107–120 min: 100% A



In addition to the determination of standard anions, the Metrosep A Supp 19 - 250/4.0 is particularly suitable for the determination of low-molecular-weight organic acids. Quinate, galacturonate, fluoride, lactate, shikimate, acetate and formate can be reliably resolved in the front part of the chromatogram with gradient elution. The column can thus be of excellent use for determinations in food samples.



In the back part of the chromatogram, higher eluent concentrations are used to elute the multivalent analytes. Organic acids are also dissolved after sulfate.

Metrosep A Supp 19 - 250/4.0		mg/L
1	Quinate	2.5
2	Galacturonate	2.5
3	Fluoride	2.5
4	Lactate	2.5
5	Shikimate	2.5
6	Acetate	2.5
7	Formate	2.5
8	Chloride	1
9	Nitrite	1
10	Bromide	1
11	Nitrate	1
12	phosphate	1
13	Sulfate	1

	Metrosep A Supp 19 - 250/4.0	mg/L
14	Malate	5
15	Tartrate	5
16	Succinate	5
17	Oxalate	5
18	Fumarate	5
19	Citrate	5
20	Isocitrate	10

5.12 Determination of standard anions and 16 organic acids with IC-MS

Column: Metrosep A Supp 19 - 150/4.0

Sample preparation: -

Detection: Mass spectrometry

Suppression: Sequential suppression with MSM and MCS

Temperature: 60 °C

Loop: 10 µL

Flow rate: 0.75 mL/min

Eluent: High-pressure gradient

A) 8.0 mmol/L Na₂CO₃, 0.25 mmol/L NaHCO₃, 10% MeOH (v/v)

B) 80.0 mmol/L Na₂CO₃, 2.5 mmol/L NaHCO₃, 10% MeOH (v/v)

0–8 min: 99% A



8–15 min: 99–80% A

15–20.9 min: 80–75% A

20.9–21 min: 75–20% A

21–28 min: 20% A



28–29.9 min: 99% A



Mass spectrometry can be used to determine analytes that are not completely separated on the analytical column. The analysis time of 30 minutes can thus be significantly reduced. In the application shown here, 200 µg/L of a wide variety of organic acids were added to a tea. In the front part of the chromatogram, 9 analytes can be determined: Quinate, fluoride, shikimate, lactate, glycolate, propionate, acetate, formate and pyruvate. For multivalent organic acids, malate, malonate, tartrate, succinate, glutarate, oxalate, fumarate and citrate are determined.

5.13 27 anions with one column

<i>Column:</i>	Metrosep A Supp 19 - 150/4.0
<i>Sample preparation:</i>	-
<i>Detection:</i>	Conductivity
<i>Suppression:</i>	Sequential suppression with MSM and MCS
<i>Temperature:</i>	30 °C
<i>Loop:</i>	20 µL
<i>Flow rate:</i>	0.7 mL/min
<i>Eluent:</i>	0.25 mmol/L NaHCO ₃ , 8.0 mmol/L Na ₂ CO ₃



A large number of different anions can be separated isocratically with the Metrosep A Supp 19 - 150/4.0.

6 Troubleshooting

6.1 Regeneration



CAUTION

Do not regenerate the column as a preventive measure!

Each regeneration process causes stress on the separation column and reduces its service life see "*Regenerating separation columns*", page 5.

Problem

- The backpressure increases.
- Double peaks occur.
- Tailing effects occur.
- The retention times become shorter.
- The resolution deteriorates.

Correction

Regenerating the separation column

Start by replacing the guard column if the above problems occur. Regenerate the separation column as described below if this measure does not help.

1 Disconnecting the separation column from the IC system

Disconnect the column outlet from the downstream functional units such as suppressor or detector.

Collect the flow of liquid in a beaker.

2 Regenerating the separation column



NOTICE

Ensure that the maximum pressure is never exceeded during regeneration. If the pressure becomes too high, reduce the flow rate.

Depending on the type of contamination, regenerate the separation column as follows:

- Contamination with organic components (*see table 3, page 59*)

- Contamination with inorganic components (see table 4, page 59).

When using organic modifiers for the regeneration, pay attention to the maximum backpressure.

Table 3 Contamination with organic components

	Rinse with	Duration [h]	Flow rate [mL/min]	Flow direction
1	Ultrapure water	1	0.4	Direction against the flow
2	Acetonitrile-water mixture (50:50)	2	0.4	Direction against the flow
3	Ultrapure water	1	0.4	Direction against the flow
4	Eluent	2	0.7	Regular

Table 4 Contamination with inorganic components

	Rinse with	Duration [h]	Flow rate [mL/min]	Flow direction
1	80 mmol/L Na ₂ CO ₃ , 2.5 mmol/L NaHCO ₃	2	0.4	Regular
2	Eluent	2	0.7	Regular

6.2 Decreasing resolution and asymmetrical peaks

Problem

The resolution of peaks deteriorates or peak shapes are asymmetrical.

Causes and prevention

Causes	Prevention or correction
The separation column has been overloaded.	<p>The separation column can be overloaded by factors such as a high salt content in the sample matrix.</p> <ul style="list-style-type: none"> ▪ Dilute the sample. ▪ Inject less sample.

6.5 Increasing backpressure

Problem

The backpressure increases.

Causes and prevention

Causes	Prevention or correction
Particles on the guard column	<ul style="list-style-type: none"> ▪ Replace the guard column.
Particles on the separation column	<p>Rinse the separation column at a reduced flow rate in the direction opposite to the flow direction.</p> <ul style="list-style-type: none"> ▪ Hold the column outlet in a beaker. ▪ Rinse the separation column for approximately 1 h. ▪ Install the separation column back in the flow direction.
Particles in the sample	<ul style="list-style-type: none"> ▪ Sample preparation, e.g. removing particles through Inline Ultrafiltration.

Index

B

Baseline	
Condition	15

C

Capacity	2
Column	
see "Separation column"	10
Conditioning	15

E

Eluent	7
Equilibration	14

F

Flow rate	2
-----------------	---

G

Guard column	
Installation	8
Rinse	8, 10

I

IC column	
see "Separation column"	10
Installation	
Guard column	8
Separation column	10

O

Order number	1
--------------------	---

R

Rinse	
Guard column	8, 10
Separation column	13
Rinsing	
Separation column	10

S

Separation column	
Installation	10
Rinse	13
Rinsing	10
Specification	1